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
INDICATIONS FOR USE

- Indication requiring immobilization and/or range of motion restriction. This may include:
 - Stable fractures of the elbow, distal humerus, proximal radius or ulna.
 - Tendon and ligament injuries or repairs.
 - Collateral ligament reconstructions.
 - Elbow hyperextension.
 - Triceps tendon repair.
 - Chronic elbow injuries.
 - Tennis elbow release.

KONTRAINDIKATIONEN

None known.

INTENDED USE

Indicated range-of-motion for post-traumatic, post-surgical repair of the elbow. This device is for single patient use. This device should be fitted by a medical professional.

GENERAL SAFETY INSTRUCTIONS

Inform your patient about everything in this document that is required for safe use of this product.

WARNING: If the device shows signs of damage or wear hindering normal function, the patient should stop using the device and contact a healthcare professional.

WARNING: In case of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

CAUTION: Before adjusting the device, the patient should contact a healthcare professional.

PRODUCT DESCRIPTION (FIGURE 2)

- Length and Cuff Position Adjustment Buttons (4)
- Adjustable Straps & Buckles (4)
- Drop-Lock Stop
- Range-of-Motion (ROM) Hinge with Flexion and Extension Stops
- Malleable cuffs (4)

FITTING

- Locate the hinge pivot center over the lateral humeral epicondyle (outer elbow) (Figure 3).
- To adjust the length and cuff position:
 - Push down on length adjustment button to unlock paddle (Figure 4).
 - Slide paddle to desired length or location.
 - Repeat adjustment for all paddles. Ensure all length adjustment buttons are in locked position.
- Important: Brace length should be adjusted to provide optimal elbow stabilization.
 - Lower extension (ulnar side of arm) should extend up to the styloid process (bump on the medial side of wrist).
 - Upper extension (humeral side of arm) should extend up to the most comfortable position at the axilla or arm pit.
- Beginning at cuffs proximal to elbow, bend cuffs to securely wrap around circumference of arm (Figure 5).
- Adjusting strap length:
 - Beginning with buckles closest to elbow, secure buckles and fasten and tighten straps (Figure 6a).
 - Straps may be trimmed if necessary (Figure 6b).
- To adjust range-of-motion, pull the flexion or extension button and rotate until center of button is pointing at desired position and release (Figure 7).
- (Optional) Place range-of-motion (ROM) stop-locks over flexion and extension buttons to limit ROM adjustment (Figure 8).

Neck Strap Application

- Attach neck strap clasp with d-ring to the loop in forearm padding closest to wrist (Figure 9).
- Place neck strap over shoulder opposite injured arm, across back, and under injured arm.
- Use second clasp to attach free end of neck strap back to d-ring on strap (Figure 10).
- Adjust strap length to support brace in a comfortable position (Figure 11).
- Slide neck strap pad over shoulder so neck strap feels comfortable.

Drop-Lock (Figure 12)

- Determine desired lock angle.
- Align blue setting arrow on proximal strut to appropriate angle as marked by blue lines on hinge.
- Slide drop-lock button toward center of hinge to lock hinge angle.

NOTE: The side bars may be contoured to fit the anatomy of the individual. Side bars must be contoured with final length.

Care and cleaning

- Do not wash in washing machine or place in dryer.
- Use mild detergent and cold water (30°C/86°F).
- Remove pads and straps and wash separately.
- Do not immerse hinge & uprights in water.
- Do not use bleach or strong detergents.

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe: The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations.

LIABILITY

- Ossur does not assume liability for the following:
 - Device not maintained as instructed by the instructions for use.
 - Device assembled with components from other manufacturers.
 - Device used outside of recommended use condition, application, or environment.

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