
3R15, 3R49

DE Gebrauchsanweisung	5
EN Instructions for use	11
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	22
ES Instrucciones de uso	28
PT Manual de utilização	34
NL Gebruiksaanwijzing	40
SV Bruksanvisning	46
DA Brugsanvisning	51
PL Instrukcja użytkowania	56
HU Használati utasítás	62
CS Návod k použití	67
RO Instrucțiuni de utilizare	72
HR Upute za uporabu	78
TR Kullanma talimatı	83
EL Οδηγίες χρήσης	88
JA 取扱説明書	94
ZH 使用说明书	100


● 4D1



Einzelteile-Pack
Single-Component Pack
Kit de pièces de rechange
Imballo dei singoli componenti
Kit componentes
Conjunto de peças de reposto
Onderdelenpakket
Service-Set
Komponentpakke
Zestaw naprawczy
Alkatrész-csomag 4D1-hez
Sada jednotlivých dílů
Set componente individuale
Pakiranje pojedinačnih dijelova
Münferit parça ambalajı
Σετ μεμονωμένων εξαρτημάτων
セット販売
零件包

▲
Mindestmenge
Minimum order quantity
Quantité minimum
Quantità minima
Cantidad mínima
Quantidade minima
Minimikvantitet
Minimum aantal
Mindste mængde
Ilość minimalna
Minimum mennyiség
Minimální množství pro objednání
Cantitate minimă
Minimalna količina
Minimum miktar
Ελάχιστη ποσότητα
最低発注数
最小起订量

■
Einzelteile
Single components
Pièces à l'unité
Singoli componenti
Componentes
Onderdelen
Separata delar
Komponenter
Podzespoły
Alkatrészek
Jednotlivé díly
Componente individuale
Pojedinačni dijelovi
Münferit parçalar
Μεμονωμένα εξαρτήματα
構成部品
零件

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Last update: 2017-02-13

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1 Components (fig. 1)

(Components of Titanium in brackets)

- | | |
|---|--------------------------------------|
| (1) 4G30=N (4G143) Upper Joint Section with Pyramid | (15) 4B74=20.4x24.5-K Brake Bushing |
| (2) 4G97 (4G144) Lower Joint Section | (16) 4G31 Bushing Safety Device |
| (3) 4G20=N (4G145) Swing Block | (17) 4Z47 Axis Screw |
| (4) 4G28=SET Knee Axis | (18) 4Z60=25 Brass Washer |
| (5) 4A62=N Swing Axis | (19) 507S16=10x1 External Lock Ring |
| (6) 4G26 Pressure Block | (20) 21Y113=2 Extension Assist Cable |
| (7) 506G3=M6x12 Set Screw | (21) 21Y66 Distal Anchor Shell |
| (8) 4G38 Plastic Block | (22) 21Y67 Adjustment Ring |
| (9) 4L7=2.0 Compression Spring | (23) 506G1=M8x25 Adjustment Screw |
| (10) 4R36 Plastic Guide Pin | (24) 513D28 Compression Spring |
| (11) 4Z32 Stance Control Adjustment Screw | (25) 21Y70=N Flexible Plastic Cover |
| (12) 4G37 Extension Stop | (26) 4Z60=17 Brass Washer |
| (13) 4Z46 Bumper | (27) 4Z59=4x6 Extension Stop |
| (14) 501T8=M4x6 Cap Screw | |

1.1 Single-Component Pack


4D1 Single-Component Pack for 3R49 and 3R15 (fig. 1, ●)

Consisting of: 1 brake bushing, 1 extension stop, 1 bumper, 2 large and 2 small washers, 2 lock rings, 2 extension stops, 1 bushing safety device, 1 axis screw.

2 Description

2.1 Application and recommendations

The monocentric 3R15 (stainless steel) and 3R49 (titanium) Modular Knee Joints with a stance phase friction brake activated by weight bearing were designed to be used **exclusively** for the prosthetic fitting of lower extremities.

Field of application according to the **Ottobock Mobility System MOBIS**: Recommendation for **mobility degrees 1 and 2** (indoor walkers and walkers confined to a limited radius outdoors). Approved for a patient with a weight of up to **100 kg (220 lbs)**.

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

The joint is especially suited for above-knee amputees who require a high degree of stability. The brake increases the stability of the knee and, for example, compensates for weak musculature.

**Attention!**

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.
- Inform the patient.

**Attention!**

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

2.2 Construction and function

The Swing Block (3) is connected to the Upper Joint Section (1) through the Swing Axis (5) and with the Lower Joint Section (2) through the Knee Axis (4). The slotted Swing Block (3) accepts the slotted Brake Bushing (15) which surrounds the knee axis. During weight bearing, resistance to movement between the knee axis and brake bushing increases the stability of the knee.

In the unloaded condition the Compression Spring (9) and the elasticity of the Swing Block (3) release the brake effect.

Extension of the knee joint is provided by the Extension Assist (20–25). The Extension Assist Cable (20) runs round the Distal Anchor Shell (21) and is hooked to the Swing Axis (5) (fig. 2) with the end of the Extension Assist Cable pointing to the rear (fig 3, arrow).

2.3 Service life

INFORMATION

As a basic principle, all Ottobock knee joints are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

3 Technical data

Article number	3R15	3R49
Proximal connection	Adjustment pyramid	Adjustment pyramid
Distal connection	Adjustment pyramid	Adjustment pyramid
Knee flexion angle	150°	150°
Material	Stainless Steel	Titanium
Weight	490 g	360 g
System height	9 mm	9 mm
Proximal system height up to alignment reference point	8 mm	8 mm
Distal system height up to alignment reference point	1 mm	1 mm
Max. user weight	100 kg	100 kg
Mobility Grade	1 + 2	1 + 2

4 Handling

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis.

The advantages of the 3R15/3R49 can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimisation.

4.1.1 Bench alignment with alignment tool

(the following steps refer to Fig. 11)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run 10 mm in front of the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter.
- ⑤ Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

Attention!

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. If required, use the 4R118 Adapter Plate for optimal positioning of the prosthetic knee joint.

- ⑦ Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static alignment optimization using L.A.S.A.R. Posture 743L100

(the following steps refer to Fig. 126)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).

- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **40 mm anterior to the upper front knee axis** (see Fig. 12).
- ③ After step 2, dynamic optimisation can take place between parallel bars.

4.2 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

4.3 Adjustment and final assembly

4.3.1 Adjusting swing and stance phase

This Modular Knee must be adjusted to the individual needs of each patient. Swing and stance phase are individually adjusted through the degree of extension assist force, friction and the braking action which is dependent upon weight bearing.

Before adjusting the knee joint, read the following instructions:

4.3.2 Adjusting swing phase with extension assist

Adjust the tension of the extension spring by turning the adjustment screw. Use a 4 mm hexhead wrench (fig. 3).

Clockwise rotation:	increases tension
Counter-clockwise rotation:	decreases tension

4.3.3 Setting swing phase speed with the friction adjustment

Resistance to movement between the knee axis and the brake bushing may be adjusted with the Set Screw (7). Use a 3 mm hex head wrench (fig. 4).

Clockwise rotation:	movement slower
Counter-clockwise rotation:	movement faster

Adjustments of as little as 15 to 20° provide a perceptible change. Settings which are too loose may cause noise and lead to premature wear. Those which are too tight will cause the brake to be activated even when not under load.

4.3.4 Adjusting the stance phase control

The responsiveness of the brake effect may be regulated by the stance control adjustment screw. Use a 6 mm hex head wrench (fig. 5).

Clockwise rotation: **braking is less responsive**
Counterclockwise rotation: **braking is more responsive**

The knee joint functions well only when the extension assist, friction and brake effect are harmoniously adjusted. Take care in making individual adjustments.

4.4 Foam cosmetic cover

Use foam cover 3R6 or 3S106 for joints 3R15/3R49 .

Attention!

To eliminate noise in the cosmetic foam cover, use 519L5 Silicone Spray. Do not use talcum powder! Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.

Useful information

To optimize sliding and for eliminating noises, please apply 519L5 Silicone Spray directly on the friction surfaces of the foam cosmetic covers.

5 Service directions

5.1 Lubrication

Apply 633F16=0.5 Lubricant to all moving metal surfaces (except the 4B74=20.4x24.5-K Brake Bushing) in the modular knee joint.

5.2 Replacing the brake bushing

If the brake bushing must be replaced, proceed as follows:

Unhook the Extension Assist Cable (20) (fig. 2). Remove the Axis Screw (17) from the joint section using a 6 mm hex head wrench (fig. 6). Push the guide surface of the knee axis out of the lower joint section to the rear and up (fig. 7). Remove the Brass Washers (18). Turn the stance control adjustment screw clockwise to open the axis hole (fig. 8). Push out the knee axis and remove the brake bushing.

Before installing the new brake bushing, be sure the Bushing Safety Device (16) is seated properly (fig. 9). Install the new brake bushing (fig. 10). Clean the knee axis with **634A58 (isopropyl alcohol)** and insert it into the bushing.

Reassemble in the reverse order.

Tighten the Axis Screw (17) using a torque wrench 710D4.

Torque: 25 Nm/18.4 ft.lbf./221 lbf.in.

Adjust joint functions as described under section 4.

6 Maintenance information

Ottobock recommends the readjustment of the knee joint's settings after the patient has spent a period of time getting used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the wear and functionality of the knee joint at least once a year and make follow-up adjustments if need be. Pay special attention to resistance to movement and the development of unusual sounds. Complete flexion and extension must be guaranteed.

Attention! Please inform your patients.

Due to different environmental conditions and conditions of application, the function of the knee joint can be impaired. To avoid the risk of accident and possible injury to the patient, the knee joint must no longer be used after a noticeable loss of function has occurred. This noticeable loss of function can be for example stiffness, non-attainment of the extension stop, decreasing swing phase or stance phase stability, abnormal noises, etc.

Measure

Have the knee joint examined by a specialist workshop where service parts or, if necessary, the entire knee joint will be exchanged.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH