3R36


3R20


3R20, 3R36

DE Gebrauchsanweisung	5
EN Instructions for use	11
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	23
ES Instrucciones de uso	29
PT Manual de utilização	35
NL Gebruiksaanwijzing	41
SV Bruksanvisning	47
DA Brugsanvisning	52
PL Instrukcja użytkowania	57
HU Használati utasítás	63
CS Návod k použití	68
RO Instrucțiuni de utilizare	74
EL Οδηγίες χρήσης	80
TR Kullanma talimatı	86
RU Руководство по применению	91
ZH 使用说明书	97

■	● 4D13	▲
Einzelteile	Einzelteile-Pack	Mindestmenge
Single components	Single-Component Pack	Minimum order quantity
Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
Singoli componenti	Imballo dei singoli componenti	Quantità minima
Componentes	Kit componentes	Cantidad mínima
Componentes	Conjunto de peças de reposo	Quantidade mínima
Onderdelen	Onderdelenpakket	Minimikvantitet
Separata delar	Service-Set	Minimum aantal
Komponenter	Komponentpakke	Mindeste mængde
Podzespoły	Zestaw naprawczy	Ilość minimalna
Alkatrészek	Alkatrész-csomag	Minimum mennyiség
Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
Componente individuale	Set componente individuale	Canitate minimă
Münferit parçalar	Münferit parça ambalajı	Minimum miktar
Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
Отдельная деталь	Комплект деталей	Минимальное количество
单个部件	维修组件	最低起订量

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Last update: 2017-02-10

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1 Components (Fig. 1)

1.1 Single components

(for 3R36 in brackets)

- | | |
|--|---|
| (1) 4G60 (4G111) Upper Joint Section | (12a) 513D34=L Spring for Extension Assist |
| (2) 4G61 (4G112) Lower Joint Section | (13) 4Z41 Adjustment Screw |
| (3) 4G117 Plastic Knee Cap | (14) 4V100 Spring Housing |
| (4) 501Z2=M4x12 Cap Screw with Hexagon Socket | (15) 509Y1=8.0 Bearing Ball, stainless steel |
| (5) 4H58 Extension Stop, complete | (16) 506G3=M6x6 Set Screw with Hexagon Socket |
| (6) 4G62=N (4G114=N) Posterior Linkage Bar | (17) 4B90=12-K Bushing |
| (7) 4Z40 Stop Bumper | (18) 502S19=M6 Lock Nut |
| (8) 4G53 (4G113=N) Anterior Linkage Bar | (19) 4A50 Posterior Axis Pin |
| (9) 501T15=M6x40 Cap Screw with Hexagon Socket | (20) 4Z49=20 Brass Washer |
| (10) 4B86 Knee Axis Bushing | (21) 4B91 Friction Stop Bumper |
| (11) 4V99 Pin for Extension Assist | (22) 507U16=6.4 Washer |
| (12) 513D34 Spring for Extension Assist | (23) 709S18 Wrench |

1.2 Single-Component Packs (Fig. 1, ●)

4D13 Single-Component Pack for 3R20 and 3R36

Consisting of: 1 Plastic Knee Cap, 1 Cap Screw, 1 Extension Stop, complete, 1 Stop Bumper, 1 Pin for Extension Assist, 1 Spring for Extension Assist, 1 Adjustment Screw, 1 Bearing Ball, 2 Lock Nuts, 2 Posterior Axis Pins, 4 Brass Washers, 4 Washers.


2 Description

2.1 Indications for use

The Ottobock Habermann Modular Knee Joints 3R20 (stainless steel) and 3R36 (titanium) are to be used **exclusively** for the prosthetic fitting of amputations of the lower limb.

2.2 Field of application

Field of application according to the **Ottobock Mobility System MOBIS**:

Recommendation for **mobility grades 1 and 2** (indoor walkers, restricted outdoor walkers). Approved for **patients with a weight of up to 100 kg/220 lbs.**

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

2.3 Service life

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.
- Inform the patient.

Attention!

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

INFORMATION

As a basic principle, all Ottobock knee joints are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

2.4 Construction and function

Upper Joint Section (1) and Lower Joint Section (2) are bound together through four rotation axes according to the engineering principle of a closed kinematic chain. The two Anterior Linkage Bars (8) with cone shaped axis pins and the Posterior Linkage bar (6) provide the moveable connection. The Extension Assist (11–15) accommodated in the Lower Joint Section extends the joint until this is limited by the adjustable, elastic Extension Stop (5).

The joint provides a rotational and translational movement due to its polycentric action. During flexion the centre of rotation changes its position, dependent on the degree of flexion (instantaneous centre of movement/rotation). In the extended position the instantaneous centre of rotation is located above the joint and posterior to the weight bearing line, thus providing knee stability in the stance phase (see section 4).

3 Technical data

Article number	3R20	3R20
Proximal connection	Pyramid adapter	
Distal connection	Pyramid adapter	
Knee flexion angle	110°	
Weight	690 g	445 g
System height	41 mm	
Proximal system height up to alignment reference point	–3 mm	
Distal system height up to alignment reference point	44 mm	
Max. user weight	100 kg / 220 lbs	
Mobility Grade	1 + 2	

4 Handling

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position of the axes. In extension, the instantaneous centre of rotation is located above the pyramid adapter and behind the alignment reference line, thus providing knee stability in the stance phase (Fig. 11).

The advantages of the 3R20 and 3R60 Knee Joints can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector.

Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimisation.

4.1.1 Bench Alignment with Alignment Tool (the following steps refer to Fig. 11)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run through the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee-ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter.
- ⑤ Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.
- ⑦ Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static alignment optimisation using L.A.S.A.R. Posture 743L100

(the following steps refer to Fig. 12)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **35 mm anterior to the upper knee axis** (see Fig. 12).
- ③ After step 2, dynamic optimisation can take place between parallel bars.

4.2 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

4.3 Adjustments and final assembly

Contrary to single-axis knee joints the 3R20/3R36 Modular Knee Joints are stable from heel strike to mid-stance phase. During roll-over of the foot knee flexion is initiated. The position of the anterior and posterior linkage bars to each other and thus the location of the instantaneous centre of rotation determine knee stability and initiation of flexion.

By varying the point of contact of the Extension Stop (5) stability of the knee joint can be adjusted according to the individual requirements of the patient. To provide swing phase adjustment the axis friction and extension assist tension may be regulated separately.

Before adjusting the knee joint, read the following instructions:

4.3.1 Adjusting stance phase stability with extension stop setting (Fig. 2)

To vary the position of the knee extension stop, turn the Adjustment Screw (5) using a 4 mm hex head wrench.

Clockwise rotation = decreases stability of the joint
(initiation of flexion becomes easier).

Counter-clockwise = increases stability.
rotation

4.3.2 Setting swing phase speed with Axis Friction adjustment (Fig. 3)

The axis pins are located in maintenance-free Knee Axis Bushings (10) and are attached by a Cap Screw (9) each to Washer (22) and Lock Nut (18). The bearing allows adjustment of the axis friction (resistance to motion) and readjustment of tolerance variation.

Fix Lock Nut (18) using enclosed Wrench (23), and turn Cap Screw (9) using a 5 mm hex head wrench.

Clockwise rotation = increases friction

Counter-clockwise = decreases friction
rotation

Adjustments of as little as 5° to 10° provide a perceptible change.

4.3.3 Regulating the swing phase with Extension Assist adjustment (Fig. 4)

Adjust the tension of the Spring for Extension Assist (12) by turning the Adjustment Screw (13). Use a 4 mm hex head wrench. The maximum depth of the adjustment screw in the Spring Housing (14) should not exceed 8 mm.

To have a lower effect from the extension assist, the 513D34=L (12a) Extension Assist Spring can be ordered.

5 Service directions

5.1 Lubrication

Apply 633F16=0.250 Lubricant to all moving metal surfaces in the modular knee joint.

Sparingly apply 633G6 Lubricant to the axis pins in the cone shaped bushings (10).

5.2 Replacing the Extension Stop (5) (Figs. 5–7)

Flex the knee joint. Loosen the Cap Screw (4) using a 3 mm hex head wrench. Remove the Plastic Knee Cap (3). Turn the Adjustment Screw (5) until the Extension Stop drops out. Use a 4 mm hex head wrench.

When replacing the Extension Stop first turn the Adjustment Screw (5) into the threaded hole by hand, two to three threads. At the same time, place the guide groove in the Extension Stop opposite the guide piece on the Upper Joint Section. Hold the Extension Stop with the left hand while turning the Adjustment Screw using a hexagonal wrench until the guide groove engages the guide piece.

5.3 Replacing the Stop Bumper (7) (Figs. 8–10)

Remove the Spring for Extension Assist (12) according to the instructions presented under section 4.3. Loosen the Set Screw (16) which holds the upper axis pin of the Posterior Linkage Bar (6) using a 3 mm hex head wrench.

Push out the Posterior Axis Pin (19) using a blunt object. Rotate the Posterior Linkage Bar distally and replace the Stop Bumper (7).

Assemble in the reverse order and lubricate the metal surfaces with 633F16 Lubricant. When inserting the Posterior Axis Pin (19), check to make sure the flat area in the middle of the pin is opposite the Set Screw (16).

Tighten Set Screw (16) using a torque wrench. Torque: 7 Nm/5.2 ft.lbf./62 lbf.in.

To secure the screws use 636K13 Loctite.

6 Cosmetic foam cover

Use foam cover 3R24, 3S124, 3R57 for joints 3R20/3R36.

Attention!

To eliminate noise in the cosmetic foam cover, use 519L5 Silicone Spray. **Do not use talcum powder!** Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. **Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.**

Notice:

To optimise sliding and for eliminating noises, please apply 519L5 Silicone Spray directly on the friction surfaces of the cosmetic foam covers.

7 Maintenance instructions

Attention!

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the joint to Ottobock.

**Attention! Please inform your patients.**

Due to different environmental conditions and conditions of application, the function of the knee joint can be impaired. To avoid the risk of accident and possible injury to the patient, the knee joint must no longer be used after a noticeable change of function has occurred. This noticeable change of function can be for example stiffness, non-attainment of the extension stop, decreasing swing phase control or stance phase stability, abnormal noises, etc.

Measure

Have the knee joint examined by a specialist workshop.

Ottobock recommends the readjustment of the knee joint's settings after the patient has spent a period of time getting used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the wear and functionality of the knee joint at least once a year and make readjustments if need be. Pay special attention to resistance to movement and the development of unusual noises. Complete flexion and extension must be guaranteed.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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A series of horizontal lines spanning the width of the page, providing a template for writing. The lines are evenly spaced and extend from the left margin to the right margin.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt/Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com