



3R80, 3R80=ST

DE	Gebrauchsanweisung (Fachpersonal)	18
EN	Instructions for use (qualified personnel)	28
FR	Instructions d'utilisation (personnel spécialisé)	37
IT	Istruzioni per l'uso (personale tecnico specializzato)	47
ES	Instrucciones de uso (personal técnico especializado)	58
PT	Instruções de utilização (pessoal técnico)	68
NL	Gebruiksaanwijzing (deskundig personeel)	78
SV	Bruksanvisning (fackpersonal)	88
PL	Instrukcja używania (personel fachowy)	97
HU	Használati útmutató (szakszemélyzet)	107
CS	Návod k použití (odborný personál)	117
RO	Instrucțiuni de utilizare (personal de specialitate)	127
HR	Upute za uporabu (stručno osoblje)	137
TR	Kullanım kılavuzu (uzman personel)	146
RU	Руководство по применению (для специалистов)	156
UK	Інструкція щодо використання (для кваліфікованого персоналу)	167
JA	取扱説明書(有資格担当者)	177
ZH	使用说明书(专业人员)	186


1 Recommended combinations | Unallowable combinations


2 Scope of delivery


3 Technical data

1 max. 150 kg
max. 330 lbs






3 3R80
1196 g


4 3R80=ST
1191 g


4 Technical data



5 Bench alignment of the prosthesis

6 Bench alignment of the prosthesis


7 Adapter inserts


8 Installing the tube adapter


1



2


3


4


5




9 Static alignment of the prosthesis | 3D Mode

10 Static alignment of the prosthesis | 2D Mode




11 Setting options | Swing phase flexion resistance



- 1 Setting range: 1,5 revolutions
- 2 Decreasing swing phase flexion resistance
- 3 Increasing swing phase flexion resistance

12 Setting options | Swing phase extension resistance




- 1 Setting range: 1,5 revolutions
- 2 Decreasing swing phase extension resistance
- 3 Increasing swing phase extension resistance

13 Setting options | Switching threshold

- 1 Setting range: 5x A – F
- 2 Decreasing the switching threshold
- 3 Increasing the switching threshold


14 Setting options | Stance phase flexion resistance


- 1 Setting range: 1 – 6
- 2 Decreasing the stance phase flexion resistance
- 3 Increasing the stance phase flexion resistance

15 Using the lock

1 Activating the lock


2 Deactivating the lock


16 Checking the settings and resetting to factory settings

1 Swing phase flexion resistance


2 Swing phase extension resistance




17 Checking the settings and resetting to factory settings

1 Stance phase flexion resistance


2 Threshold

3 Bumper without play, pre-stressed


18 Checking the settings and resetting to factory settings | Lock

19 Change between stance phase and swing phase

1 Stance phase flexion resistance at heel strike activated

2 Stance phase flexion resistance deactivated at heel relief

20 Exercises and adjustments

1

2


3


4

5

6



7


21 Exercises and adjustments | Sitting down



22 Exercises and adjustments | Walking

23 Exercises and adjustments | Walking


24 Exercises and adjustments | Walking down ramps

25 Exercises and adjustments | Cycling


26 Finishing the prosthesis


4 mm
10 Nm


4 mm
15 Nm
636K13- LOCTITE


- 506G3=M8X10
- 506G3=M8X12-V
- 506G3=M8X14
- 506G3=M8X16




27 Finishing the prosthesis




28 Checking the sacrificial anode

29 Maintenance | 4G1226=* Disassembling the stop bumper

30 Maintenance | 4G1226=* Assembly of stop bumper

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INFORMATION

Date of last update: 2025-02-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The "3R80/3R80=ST" prosthetic knee joint is referred to as the product in the following.

This document is intended for qualified personnel with manufacturing knowledge of lower limb prostheses.

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 30).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Observe the operating conditions of the product (see page 29) and do not expose the product to prohibited environmental conditions (see page 29).
- ▶ To avoid the risk of pinching, do not reach into the joint mechanism.
- ▶ In the event of leaks, stop using the product and have it repaired.
- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Only use the approved cleaning agents and disinfectants.
- ▶ Observe the instructions for cleaning and care (see page 35).
- ▶ Follow the maintenance instructions (see page 35).

Signs of changes in or loss of functionality during use

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

3 Product description

The product has the following key features:

- Monocentric prosthetic knee joint with rotation hydraulics
- Product components for stance phase stability:
 - Adjustable stance phase flexion resistance (hydraulic damping)
 - Adjustable trigger threshold for switching between swing phase resistance to stance phase flexion resistance
 - Lock (activated and deactivated by the user)
- Product components for swing phase control:
 - Adjustable swing phase flexion resistance (hydraulic damping)
 - Adjustable swing phase extension resistance (hydraulic damping)
 - Hydraulic unit with extension assist function (the extension assist function can be enhanced by installing an additional spring)

4 Intended use

4.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

4.2 Operating conditions

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

Allowable field of application	
Recommended mobility grade: 3 + 4 (see fig. 3 item ①)	Everyday prosthesis
Allowable body weight: ≤ 150 kg (see fig. 3 item ②)	

4.3 Environmental conditions

Storage and transport
- Storage temperature: -20 °C to +60 °C - Relative humidity: 20 % to 90 % - No mechanical vibrations or impacts

INFORMATION

Clean the product after contact with humidity/chemicals/solids in order to avoid increased wear and tear and damage (see page 35).

Allowable environmental conditions
- Temperature range: -10 °C to +45 °C - Relative humidity: 20 % to 90 %
- Contact with soap solution (e.g. when taking a shower) - Contact with perspiration - Submersion in fresh water/chlorinated water (e.g. swimming pool): maximum immersion depth 2 m
- Submersion in salt water with a low salt content (14 days/year, 1 hour/day) - Permissible salt content: max. 3.5 % (e.g. ocean) - Maximum immersion depth swimming pool: 1 m - Maximum immersion depth sea: 0.5 m
Contact with dust, foam cover particles, drifting sand (e.g. walking on the beach)
Contact with salty, condensing air (e.g. in a 40 °C steam sauna)
UV resistant

Unallowable environmental conditions
- Contact with liquid-binding particles (e.g. talcum) - Contact with large amounts of sand and dust (e.g. digging, kneeling in sand, construction site work) - Contact with acids - Contact with urine
- Submersion in salt water with excessive salt content - Salt content: > 3.5 % (e.g. brine bath, Dead Sea)
Cleaning agents and disinfectants with solvents, chlorine and phosphate
High water pressure (e.g. diving, jumping into water)

4.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

The combination possibilities with prosthetic hip joints are shown in the following illustration: see fig. 1

	7E7	7E9	7E10
3R80			
3R80=ST			

4.5 Lifetime

This prosthetic component has been load tested by the manufacturer according to ISO 10328. The maximum lifetime is 5 years.

5 Scope of delivery

The scope of delivery is shown in the following illustration: see fig. 2

- | | |
|------------------------------|--|
| ① Prosthetic knee joint | ⑥ Tube adapter |
| ② Adjustment wrench | ⑦ Quick start guide |
| ③ Compression spring | ⑧ Instructions for use qualified personnel |
| ④ Adapter protector 2Z11=KIT | ⑨ Instructions for use – users |
| ⑤ Stop bumper set | |

6 Preparing the product for use

6.1 Bench alignment of the prosthesis

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

CAUTION

Initial use of the prosthesis by the user

Falling due to lack of experience of the user, incorrect alignment or adjustment of the prosthesis

- For the safety of the user when standing and walking for the first time, use a suitable aid (e.g. parallel bars, handrail and anterior walker).

The prosthesis must be aligned using the alignment images as well as the instructions for use of all prosthetic components used (see fig. 5; see fig. 6). The prosthesis should be aligned so as to enable the user to stand safely.

Optional: Using the adapter protector 2Z11=KIT

The connecting section of the prosthetic joint can be protected against scratches when the adapter protector is used during alignment and fitting.

- 1) Follow the instructions in the accompanying document for the adapter protector.
- 2) Remove the adapter protector before leaving the fitting area.

- 3) Readjust the prosthetic alignment.

Alignment device 743A220 PROS.A.* Assembly: Installing adapter inserts

- Mount the adapter inserts 743Y580=* (see fig. 7) when using an alignment device.

6.2 Shortening the tube adapter

⚠ CAUTION

Incorrect processing of tube

Fall due to damage to the tube

- Do not clamp the tube into a vice.
- To shorten the tube, use only a tube cutter or a cutting device.

Shorten the tube adapter to the desired length (see fig. 8).

> **Recommended tools:** cutting device 704Y14=* or tube cutter 719R4, tube deburrer 718R1

- 1) Shorten the tube. Use a cutting device or tube cutter (item ❶) for this purpose.
- 2) Deburr the cut edge with the tube deburrer (item ❷).
- 3) Make sure that the tube is evenly deburred on the inside and outside (item ❸).

6.3 Installing the tube adapter

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- Clean the threads before every installation.
- Apply the specified torque values.
- Follow the instructions for thread lock.

⚠ CAUTION

Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- Clean the contact surfaces of the tube and the tube adapter using a degreasing cleaning agent.
- Position the tube in the tube adapter so the minimum insertion depth is met.

Securely install the tube adapter in the joint (see fig. 8).

> **Required tools:** Torque wrench 710D4, 4 mm hexagon socket adapter

- 1) Slide the tube into the product up to the stop (**0 mm**) (item ❹).
- 2) Tighten the cap screw **without** additional Loctite to **10 Nm** (item ❺).

6.4 Static Alignment

Ottobock recommends optimising the alignment with the 743L500 3D L.A.S.A.R. Posture measuring device. The 743L100 L.A.S.A.R. Posture can be used as an option.

The following optimisation methods are available:

- 3D mode
- 2D mode

Optimise the alignment of the prosthesis so that the user can stand comfortably and the reference points in the following illustrations are observed.

3D mode (see fig. 9)

- ❶ **PREREQUISITE:** The prosthetic knee joint is extended. The lock is activated.
- ❷ The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- ❸ The extension line of the auxiliary line points to the reference point (anterior superior iliac spine).

- 4 The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).

2D mode (see fig. 10)

- 1 **PREREQUISITE:** The prosthetic knee joint is extended. The lock is activated.
- 2 The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- 3 The extension line of the auxiliary line points to the reference point (anterior superior iliac spine).
- 4 The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).
- 5 The extension line of the auxiliary line points to the reference point (pyramid).

6.5 Dynamic Trial Fitting

CAUTION

Adjusting the settings

Falling due to improper or unfamiliar settings

- Only adapt the settings to the user gradually.
- Explain to the user the effects of the adjustments on the use of the prosthesis.

6.5.1 Overview of adjustment possibilities

Deactivate the lock for dynamic fitting.

Swing phase flexion resistance (see fig. 11)

- 1 Setting range
- 2 Decrease swing phase flexion resistance
- 3 Increase swing phase flexion resistance

Swing phase extension resistance (see fig. 12)

- 1 Setting range
- 2 Decrease swing phase extension resistance
- 3 Increase swing phase extension resistance

Switching threshold (see fig. 13)

- 1 Setting range
- 2 Decrease the switching threshold (lower weight required to initiate the stance phase)
- 3 Increase the switching threshold (higher weight required to initiate the stance phase)

Stance phase flexion resistance (see fig. 14)

- 1 Setting range
- 2 Decrease stance phase flexion resistance
- 3 Increase stance phase flexion resistance

Use lock (see fig. 15)

CAUTION

Locking function

Risk of injury due to unintentional activation or deactivation of the locking function

- Check that the joint is functioning correctly each time you activate or deactivate the lock.
- Make sure that cosmetic covers and clothing cannot activate the lock.

To activate or deactivate the lock, put weight on the prosthetic knee joint, push against the extension stop and activate the lock lever.

- ❶ Activate: Push the lock lever down (e.g. for relaxed standing)
- ❷ Deactivate: Pull the lock lever up (e.g. for walking)

Use in wet areas

CAUTION

Use in wet areas

Risk of injury due to unlocked joint or walking too fast on wet surfaces

- Lock the joint before using it in wet areas.
- Walk carefully on wet or slippery surfaces.

- 1) Observe the environmental conditions prior to use in wet areas (see page 29).
- 2) Activate the lock.
- 3) Walk carefully on wet and slippery surfaces.
- 4) Deactivate the lock after leaving wet areas.
- 5) Dry the joint.

6.5.2 Checking settings and restoring factory settings

The product is set to factory settings on delivery.

- 1) Check the settings of the product prior to the dynamic fitting.
- 2) In case of deviations, reset the product to the factory settings.
→ See factory settings table (see page)

Function	Adjustment position	Factory settings	Meaning
Swing phase flexion resistance	Adjustment valve "F"	Stop "–"	Minimum resistance
Swing phase extension resistance	Adjustment valve "E"	Stop "–"	Minimum resistance
Switching threshold to initiate the stance phase	Adjustment ring "I"	Bumper without play, pre-stressed	Low trigger threshold (Stance phase flexion resistance is activated with low load)
Stance phase flexion resistance	Adjustment ring "S"	Stop "6"	Maximum resistance

Swing phase flexion resistance (see fig. 16)

- ❶ Swing phase flexion resistance
- ❷ Swing phase extension resistance

Stance phase flexion resistance (see fig. 17)

- ❶ Stance phase flexion resistance
- ❷ Switching threshold to initiate the stance phase
- ❸ Bumper without play, pre-stressed

Lock (see fig. 18)

- ❶ Lock (deactivated)

6.5.3 Switching between the swing phase and stance phase

- ❶ Stance phase flexion resistance activated at heel strike (see fig. 19)
- ❷ Stance phase flexion resistance deactivated when heel is relieved

6.5.4 Exercises and adjustments

- 1) Take safety into account first if the settings cannot be fully adjusted to the needs of the user.
- 2) The full extension position has to be reached with every step and at any walking speed.
- 3) Adapt the prosthetic knee joint through fine adjustments and exercises.
- 4) Check the prosthesis settings during annual inspections and consultations.
- 5) Instruct the user to have the prosthesis inspected in case of changes in functionality.
- 6) Observe the recommended sequence of the exercises and settings (see page 11).

Sitting down (see fig. 21)

- Practise sitting down with the user.

Walking (see fig. 22)

- 1) **CAUTION!** Allow the user to walk only with safeguards in place, as the product may switch to swing phase flexion resistance too late with the factory settings.
 - 2) Reduce the switching threshold if stance phase resistance cannot be activated at the beginning of the stance phase, until it works for all walking speeds and stride lengths.
 - 3) Check the effect on the gait pattern after each change.
 - 4) Adjust the settings first at normal walking speed, then with short, quick steps and finally with long, fast steps.
 - 5) Adjust the swing phase flexion resistance so the lower leg of the prosthesis does not swing too far towards the back (dorsal position $\leq 70^\circ$) and reaches full extension in time for the next heel strike.
 - 6) Adjust the swing phase extension resistance so the product does not swing too hard against the extension stop and reaches full extension in time for the next heel strike.
- Install the additional spring to enhance the extension assist function and reach full extension at the next heel strike (see fig. 23).

Walking down ramps and stairs (see fig. 24)

- Adjust the stance phase flexion resistance so it provides the user with adequate safety yet does not generate excessive resistance.

Walking (repeat for readjustment) (see page 12)

- Check the swing phase flexion resistance again. Adjust it so the lower leg of the prosthesis does not swing too far towards the back (dorsal position $\leq 70^\circ$) and reaches full extension in a timely manner for the next heel strike.

Use lock (see fig. 15)

- Practise using the lock.

Cycling (see fig. 25)

- 1) Perform physical exercises such as cycling (e.g. on an exercise bike).
- 2) Deactivate the lock when doing so.

6.6 Attaching the cosmetic cover

CAUTION

Use of liquid-binding substances (e.g. talcum)

Risk of injury, damage to the product due to extraction of fluid from lubricant

- Prevent contact of the product with liquid-binding substances.

The recommended foam cover for the 3R80* is **3S107**.

- When attaching the cosmetic covers, note the settings from the dynamic fitting (see page 32).
- Apply silicone spray 519L5 directly onto the contact surfaces of the cosmetic foam cover to reduce friction and to eliminate noise.

6.7 Finishing the prosthesis

- 1) Replace set screws that are too short and could lead to a risk of falling (e.g. because a set screw slides in and does not secure the components) with longer ones.

- 2) Replace set screws that are too long and can damage the product (e.g. due to collision or friction) with shorter ones.
- 3) **NOTICE! Align the prosthesis so the prosthetic socket and joint do not come into contact with each other. Otherwise, the protector may partially open or become damaged when kneeling (see fig. 27). If contact is unavoidable when kneeling, show the user how to close the open protector.**
- 4) Clean the threads of the screw connections on the prosthetic socket or foot using a degreasing cleaner.
- 5) Only use Loctite 636K13 on the screw connections of the prosthetic socket or foot(see fig. 26).
- 6) Finish work on the prosthesis by tightening all screw connections. Observe the assembly tightening torques and thread locks for all prosthetic components (see fig. 26).

7 Cleaning

Permissible cleaning agents and disinfectants

Requirements:

Free of solvents, chlorine and phosphates

- > **Required materials:** 2 microfibre cloths, isopropyl alcohol 634A58, tap water
- 1) **In case of light soiling (e.g. small deposits):** Moisten the microfibre cloth with water.
- 2) **In case of heavy soiling (e.g. tough deposits):** Moisten the microfibre cloth with isopropyl alcohol.
- 3) **NOTICE! Verify the material compatibility when using isopropyl alcohol. Do not remove lubricant from bearing points.**
- 4) Clean the product with the moistened cloth.
- 5) Dry the product with a dry microfibre cloth.
- 6) Allow to air dry in order to remove residual moisture.

Cleaning after contact with salt water

- > **Required materials:** Microfibre cloth, tap water
- 1) **NOTICE! Remove salt crystals from the prosthetic knee joint to prevent corrosion.**
- 2) **Cleaning the prosthetic knee joint:** Immerse in tap water.
- 3) **Rinse out dirt and salt crystals:** Move the joint in the water (flex and extend several times).
- 4) **Remove residual water:** Remove the joint from the water and allow it to drain.
- 5) **Drying:** Dry with a microfibre cloth and allow to air dry.
- 6) **Perform a functional test:** Test the joint for proper function (flexing, extending, locking).
- 7) **Checking by qualified personnel:** After maximum use in salt water (14 days/year, 1 hour/day), have it checked by qualified personnel.

Cleaning with Disinfectants

- > **Required materials:** 2 microfibre cloths, colourless disinfectant
- 1) **NOTICE! Verify material compatibility when using disinfectants. Do not remove lubricant from bearing points.**
- 2) Apply disinfectant to the microfibre cloth.
- 3) Disinfect the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow to air dry in order to remove residual moisture.

8 Maintenance

- **NOTICE! Do not lubricate and grease the prosthetic joint.**
- **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- Schedule regular maintenance appointments with the user based on usage.

- ▶ After the user has familiarised themselves with the prosthetic joint, check the settings and adjust them if necessary.
- ▶ Check (inspect) the prosthetic components after 30 days of initial use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ Check the prosthetic joint for wear and functionality as part of the safety checks.
- ▶ Pay particular attention to the effortless ground compliance of the moveable joint components and unusual noises.
- ▶ Check the function of the set damping resistance.
- ▶ Check that full flexion and extension are always ensured.
- ▶ When used in salt water, check the product thoroughly for salt residues and traces of corrosion.
- ▶ Check all joint characteristics and adjustment options of the product for functionality.
- ▶ In case of functional limitations and damage, send the product to the manufacturer's service department.

8.1 Checking the sacrificial anode

The sacrificial anode ❶ protects the product against electrochemical corrosion, as it is destroyed rather than the product (see fig. 28).

- ▶ Instruct the user to have the prosthesis checked as soon as possible if corrosion is visible.
Example: In case of rust on the set screws that attach the proximal adapter to the prosthetic knee joint.
- ▶ Check the sacrificial anode during the annual safety check.
- ▶ Remove any dirt and crystals from the area around the sacrificial anode prior to checking.
 - Continue using the product if the sacrificial anode is present. (Entire sacrificial anode present; part of sacrificial anode present – dimension: **< 5 mm**)
 - **NOTICE!** If the sacrificial anode has been fully or partially destroyed down to the specified dimension, have a new sacrificial anode installed by the manufacturer's service department. (Sacrificial anode no longer present; part of sacrificial anode present – dimension: **≥ 5 mm**)

8.2 Replacing the stop bumper

Replace the stop bumper with every prosthesis service (see fig. 29; see fig. 30).

- > **Required tools and materials:** Flat screwdriver; flat nose pliers; compressed air; stop bumper set 4G1226=*
- 1) Deactivate the lock.
- 2) Move the product to maximum flexion.
- 3) Carefully push out the stop bumpers using the flat screwdriver.
- 4) Remove dirt and residues with compressed air.
- 5) Insert the new stop bumpers with the flat nose pliers.
- 6) Bring the product to maximum extension.
- 7) Activate the lock.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

11 Technical data

The total weight is shown in the figure below: see fig. 3

③ Total weight 3R80: **1196 g**

④ Total weight 3R80=ST: **1191 g**

The technical data are shown in the following figure: see fig. 4

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