






ottobock.

13E202=50 (50 Hz) / 13E202=60 (60 Hz)

DE	Gebruiksaanweisung (Fachpersonal)	3
EN	Instructions for use (qualified personnel)	12
FR	Instructions d'utilisation (Personnel spécialisé)	20
IT	Istruzioni per l'uso (Personale tecnico specializzato)	30
ES	Instrucciones de uso (Personal técnico especializado)	39
PT	Manual de utilização (Pessoal técnico)	48
NL	Gebruiksaanwijzing (Vakmensen)	57
SV	Bruksanvisning (Fackpersonal)	66
DA	Brugsanvisning (Faguddannet personale)	75
PL	Instrukcja użytkowania (Personel fachowy)	83
HU	Használati utasítás (szakszemélyzet)	93
CS	Návod k použití (Odborný personál)	101
RO	Instrucțiuni de utilizare (Personal de specialitate)	110
HR	Upute za uporabu (Stručno osoblje)	119
TR	Kullanma talimatı (Uzman personel)	127
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)	136
RU	Руководство по применению (Квалифицированный персонал) ...	146
JA	取扱説明書 (有資格担当者)	156
ZH	使用说明书 (专业人员)	164

INFORMATION

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INFORMATION

Last update: 2016-05-23

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The product "13E202 electrode" is referred to simply as the product/electrode below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

The product samples muscle action potential produced by muscle activity on the surface of the skin. Unwanted interference signals are filtered out by the product.

3 Application

3.1 Indications for use

The product is to be used **exclusively** for the control of myoelectric upper limb prosthetic components.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

3.3 Qualification

The fitting of a patient with the product may only be carried out by an orthopaedic technician.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 CAUTION
The heading describes the source and/or the type of hazard The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible: > E.g.: Consequence 1 in case of failure to observe the hazard > E.g.: Consequence 2 in case of failure to observe the hazard ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Fitting of children Injury due to swallowing small parts. ▶ Never leave children unattended during the fitting.
 CAUTION
Insufficient skin contact of the electrodes Injury due to faulty control or malfunction of the product. ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.

- ▶ In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary.
- ▶ If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

⚠ CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

⚠ CAUTION

Damage to the electrode cable due to kinking or abrasion

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Replace the damaged electrode cable immediately.

⚠ CAUTION

Turning the adjustment regulator too far

Injury due to faulty control or malfunction of the product.

- ▶ Do not turn the adjustment regulator further than the noticeable stop.

⚠ CAUTION

Slipping electrode

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ Inform the patient that uncontrolled movement of the prosthesis component may occur if the electrode slips.

⚠ CAUTION

Penetration of liquids

Injury due to faulty control or malfunction of the product.

- ▶ Make sure that no liquids can penetrate into the product.
- ▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Establishing a connection between the skin and metal components of the prosthesis with the use of carbon fibre

Injury due to unexpected product behaviour as a result of interference conducted by carbon fibres.

- ▶ During fabrication, make sure that no connection is formed between the skin and metal prosthesis components by carbon fibres.

NOTICE

Improper care of the housing

Damage to the casing through the use of acetone, white spirit or similar solvents.

- ▶ Only clean the housing with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

NOTICE

Thermal overloading due to external heat supply

Destruction of the product.

- ▶ Do not expose the electrode to hot air or other sources of heat.

4.4 Patient information

⚠ CAUTION

General use of electrodes

Skin irritation may occur under the electrodes.

- ▶ After removing the prosthetic components, check the skin beneath them for pressure sores and irritation.
- ▶ Observe safety information in the instructions for use.

⚠ CAUTION

Damaged electrode cable

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ Power down the system and contact your orthopaedic technician.

⚠ CAUTION

Use of a damaged product

Injury due to loss of product functionality.

- ▶ Prior to use, conduct an external visual inspection to verify that all parts of the product are undamaged.
- ▶ In case of damage, have the product repaired promptly.

⚠ CAUTION

Penetration of product with dirt and humidity

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see Page 20).

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points,...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Do not bring the surface of the electrodes into contact with a special cleaner for Ottobock cosmetic gloves, e.g 640F12.

5 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. 13E202=50/13E202=60 suction socket electrode
- 1 pc. 13E80 adjustment tool
- 1 pc. 647G334 instructions for use (qualified personnel)

Additional scope of delivery for socket fabrication

- 1 pc. 13E203 casting template for inner socket
- 1 pc. 13E204 lamination template for outer socket
- 2 pc. 503F3 socket screw with Allen head
- 2 pc. 507S15 lamination plate with bore

The following products are also required:

- 757M11=X-Change MyoBoy
- 13E129=* electrode cable
- 1 pc. 617H46 bonding agent

6 Preparation for use

Requirements for using the product

The following points must be observed for using the product:

The ideal electrode position was determined with the 757M11=X-Change MyoBoy.

Please note the 647G265 instructions for use.

The socket was fabricated according to the specifications in the 646T3=3.2 technical information.

6.1 Connecting the product

- 1) Shorten the electrode cable as needed.
- 2) Push the electrode cable with the grey side up into the plug connection all the way to the stop (see fig. 2).
- 3) Grease the connector receptacle of the product and the pin contacts with 633F11 silicone grease.
- 4) Apply 633F11 silicone grease to the electrode cable receptacle.
- 5) Insert the plug connection with connected electrode cable as far as possible into the product (see fig. 3).
- 6) Remove any silicone grease that seeps out.

Electrode cable replacement

- 1) Pull off the electrode cable at right angles to the product (see fig. 4).
- 2) Cut off at least 5 mm of the electrode cable and push it back into the plug connection.

6.2 Inserting the product

- 1) Apply 617H46 bonding agent into the groove of the electrode pad (see fig. 5).
- 2) Insert the product into the socket from the inside.

6.3 Adjusting the product

- 1) Use the adjustment tool to set the product sensitivity (see fig. 6).
- 2) After a pause of 10 to 15 minutes, readjust the sensitivity.

INFORMATION

The patient must be able to hold the maximal value of the control unit for at least 2 seconds. If this is not observed, the product setting is overly sensitive which makes it more difficult for the patient to control the prosthesis.

7 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

8 Legal information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

8.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9 Appendices

9.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Type BF applied part

LOT/PPPP YYYYY WW

Lot number

9.2 Technical data

Ambient conditions	
Storage and transport in original packaging	-15 °C/5 °F to +60 °C/140 °F
Storage and transport without packaging	-15 °C/5 °F to +60 °C/140 °F max. 80 % relative humidity, non-condensing
Operation	-15 °C/5 °F to +60 °C/140 °F max. 80 % relative humidity, non-condensing

General information	13E202=50	13E202=60
Frequency	50 Hz	60 Hz
Product service life	5 years	
Power supply	All batteries compatible with the Ottobock Myosystem	
Frequency bandwidth	90 - 450 Hz	
Sensitivity range	2,000 - 100,000 x	

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Otto Bock Healthcare Products GmbH