1K10, 1K30, 1S30

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	24
PT	Manual de utilização	29
NL	Gebruiksaanwijzing	35
SV	Bruksanvisning	40
DA	Brugsanvisning	45
NO	Bruksanvisning	50
FI	Käyttöohje	55
PL	Instrukcja użytkowania	60
HU	Használati utasítás	65
CS	Návod k použití	71
RO	Instrucțiuni de utilizare	76
HR	Upute za uporabu	81
SL	Navodila za uporabo	86
SK	Návod na používanie	91
BG	Инструкция за употреба	96
TR	Kullanma talimatı	102
EL	Οδηγίες χρήσης	107
RU	Руководство по применению	113
JA	取扱説明書	119
ZH	使用说明书	124
KO	사용 설명서	128

1

1K10

1K30

1S30

2

1 ■

2 ■

INFORMATION

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INFORMATION

Last update: 2016-02-08

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1K10, 1K30 and 1S30 prosthetic feet are suitable for use in both modular and exoskeletal prostheses. The functional properties are achieved through the combination of a contoured core and functional foam.

1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

2.2 Area of Application

1K10 Size 14 – 17, 1K30 Size 14 – 17, 1S30 Size 12 – 13

Approved for a body weight of up to **35 kg**.

1K10 Size 18 – 21, 1K30 Size 18 – 21

Approved for a body weight of up to **45 kg**.

2.3 Environmental Conditions

Allowable Environmental Conditions	
Temperature range for use: -10 °C to +60 °C (14 °F to 140 °F)	
Allowable relative humidity 0 % to 90 %, non-condensing	
Unallowable environmental conditions	
Mechanical vibrations or impacts	
Perspiration, urine, fresh water, salt water, acids	
Dust, sand, highly hygroscopic particles (e. g. talcum)	
Other Environmental Conditions	
UV light	Not resistant
Hydrolysis	Not resistant

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Excessive strain on the product Risk of injury due to breakage of load-bearing components ► Use the prosthetic components according to the classification (see section "Area of application").
 CAUTION
Reuse on another patient Risk of injury due to loss of functionality as well as damage to the product ► Only use the product for a single patient.
 CAUTION
Unallowable combination of prosthetic components Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover characteristics) is an indication of loss of functionality.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

1K10, 1K30, 1S30					
Fig.	Item no.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647G390
For prosthetic feet 1K10 Size 14 – 17, 1K30 Size 14 – 17, 1S30 Size 12 – 13					
2	1	■	0	Foot Adapter with screw connection	2R40=2
For prosthetic feet 1K10 Size 18 – 21, 1K30 Size 18 – 21					
2	2	■	0	Foot Adapter with screw connection	2R40=1

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- Observe the alignment and assembly instructions.

5.1 Installing the Foot Adapter (Modular Design)

- Assemble the foot adapter on the prosthetic foot according to its instructions for use.

5.2 Alignment

NOTICE

Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- Do not grind the prosthetic foot.

5.2.1 Bench Alignment

Foot size	Mid-foot anterior to the alignment reference line	Heel height
12 – 16 cm	15 mm	see Technical Data
17 – 21 cm	20 mm	

5.2.2 Static Alignment

- Ottobock recommends fitting the prosthesis using the L.A.S.A.R. Posture.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.2.3 Dynamic Trial Fitting

- Adapt the setup of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) so that correct heel contact and optimal roll-over are ensured.
- **TT fittings:** Make sure that physiological knee movement is achieved after the heel makes impact and the leg begins to bear weight.

6 Cleaning and Care

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct quarterly safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive.

The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH