50S20-1 Malleo Immobil Night Splint

DE	Gebruiksaanweisung	4
EN	Instructions for Use	7
FR	Instructions d'utilisation	10
IT	Istruzioni per l'uso	14
ES	Instrucciones de uso	17
PT	Manual de utilização.....	21
NL	Gebruiksaanwijzing	24
SV	Bruksanvisning	28
DA	Brugsanvisning.....	31
NO	Bruksanvisning	34
FI	Käyttöohje	37
PL	Instrukcja użytkowania	41
HU	Használati utasítás	44
CS	Návod k použití	47
RO	Instrucțiuni de utilizare.....	51
HR	Upute za uporabu	55
BS	Uputstvo za upotrebu	58
TR	Kullanma talimatı.....	61
EL	Οδηγίες χρήσης	64
RU	Руководство по применению	68
JA	取扱説明書.....	72

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INFORMATION


Date of the last update: 2013-10-21

- Please read this document carefully.
- Follow the safety instructions.

Explanation of Symbols

 Warnings regarding possible risks of accident or injury.

 Warnings regarding possible technical damage.

 Additional information on the fitting / use.

1 Description

1.1 Application

The 50S20-1 Malleo Immobil Night Splint is to be used **solely** for the orthotic treatment of the lower limbs and is **only** suitable for contact with healthy skin.

1.2 Indications and Effects

The 50S20-1 Malleo Immobil Night Splint is used for treatment following operative tendon lengthening in the calf muscle region and problems

resulting from excessive stress to the lower leg-foot region. Indication must be determined by the physician.

1.3 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. Parts to be fitted and those parts that come directly into contact with the skin can cause functional and hygienic risks if the orthosis is used by another person.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

2 Handling

2.1 Size Selection

This orthosis is universally applicable. Its size is determined based on shoe size (see sizing chart).

2.2 Adapting and Applying the Product

The initial fitting and application of the orthosis must be carried out by qualified personnel.

The orthosis is delivered ready for trial fitting.

1. To apply the orthosis, completely open it and position the foot in the middle of the orthosis with the knee flexed (Fig. 1).
2. Close the ankle strap and set the hook and loop closure to a comfortable length (Fig. 2). Close the foot and calf strap similarly. (Fig. 3).
3. Extend the knee to ensure that the ankle and foot straps are firmly in place. Readjust the hook and loop closure, if necessary (Fig. 4).
4. The hook and loop closures that support dorsal flexion can be incrementally adjusted to increase or reduce the angle of the foot (Fig. 5).

2.3 Instructions for Use and Care

NOTICE

Damage caused by wear and tear. Before each use, the orthosis must be checked for functional reliability and for possible wear or damage. If a closure or any other part of the orthosis shows any signs of wear, the orthosis should no longer be used.

Improper changes to the product are not permitted.

The orthosis must be used in strict conformity with the indication. The physician's instructions must be strictly observed.

Material: Foot bucket: Polypropylene; liner: polyester, foam, velour; hook and loop straps: nylon; closure pad: polyester, foam, Velour; closures: polypropylene; rivets: nylon; sole: Evazote.

Cleaning: Use a slightly moistened cloth to clean the outside of the orthosis. Hand-wash the textile elements in 30 °C warm water with a standard mild detergent. Soap residues can cause skin irritations; therefore, rinse the orthosis thoroughly so that all soap residue is removed. Let air dry. Do not expose to the sun and do not use oven or radiator heat for drying.

Disposal after usage: The material used in the product is recyclable and should be disposed of accordingly.

3 Further Usage Restrictions

The product is designed for use on one patient only.

The daily duration of use and period of application are dependent on the medical indication.

4 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH