6A20=10, 6A20=20

DE	Gebrauchsanweisung	3
EN	Instructions for use	9
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	21
ES	Instrucciones de uso	26
PT	Manual de utilização	32
NL	Gebruiksaanwijzing	38
SV	Bruksanvisning	44
NO	Bruksanvisning	50
PL	Instrukcja użytkowania	56
HU	Használati utasítás	62
CS	Návod k použití	68
RO	Instrucțiuni de utilizare	74
TR	Kullanma talimatı	80
EL	Οδηγίες χρήσης	85
RU	Руководство по применению	92
JA	取扱説明書	98
ZH	使用说明书	104

6A20=10

6A20=20

INFORMATION

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INFORMATION

Last update: 2015-12-03

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place

Overview of components

6A20=10 (Fig. 1)

Com- ponent Group	Description	Article no.	Torque	Description see section
1	Lamination Anchor <i>Clamping Screw</i>	4R111=N	10 Nm	4
2	Pin	6Y13=1		
3	Lamination Protection for Pin	5X440		
4	Ratchet Unit	6A52	20 Nm	4
5	Shuttle Lock Housing with Ad- justment Pyramid	6A51=10		4
6	Dummy Set with Screw	5X55		

6A20=20 (Fig. 2)

Com- ponent group	Description	Article no.	Torque	Description see section
1	Lamination Anchor <i>Clamping Screw</i>	4R111=N	10 Nm	4
2	Pin	6Y13=2		
3	Lamination Protection for Pin	5X440		
4	Ratchet Unit	6A52	20 Nm	4
5	Shuttle Lock Housing	6A51=20		4

Component group	Description	Article no.	Torque	Description see section
6	Dummy Set with Screw	5X55		
7	Set Screws (4 pieces)	506G3	15 Nm	4

1 Indications for use

The 6A20=10 and 6A20=20 Shuttle Lock Kits are to be **exclusively** used for the prosthetic fitting of the lower limb. They are not suitable for being used in bathing prostheses.

Caution!

If using unsuitable liners you can cause skin irritations. For further information, please see Liner Selection Tool 646D250=GB.

Field of application according to the **Ottobock Mobility System MOBIS**:

Approved for a patient with a **weight of up to 125 kg/275 lbs.**

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

2 Description

The shuttle lock is used to connect a liner into the socket of a modular prosthesis thus providing the connection to the distal end. A release button on the ratchet unit can be used to lock and unlock the 6Y13=1/6Y13=2 Pin to/from the shuttle housing.


Caution!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

3 Recommendation for socket reinforcement

For more details about the casting technique and incorporation of the 6A20=10/=20 Shuttle Lock please refer to the Technical Information 2.1.5 and 2.1.6.

	Caution! Sometimes the alignment line (solid line) differs from the long axis of the stump (broken line). The shuttle-lock must follow the broken line. In individual cases with the 6A20=10 this can require the use of a sliding adapter, e.g. 4R103 (observe weight limits!).
--	--

3.1 Preparation of socket adapter

Lubricate the threads of the adapter and clamping screw with silicone grease. Fill the openings near the clamping screw with Plasta Tape.

3.2 Reinforcement

- Affix the Domed Dummy with enclosed screw to the plaster positive extension.
- Pull on the soaked PVA bag
- Pull on a double layer of Nylglass Stockinette and tie it on in the groove of the Domed Dummy.
- Screw the enclosed Threaded Dummy onto the Socket Adapter and place it on the Domed Dummy.

 Attention! The anchor arms should be aligned in the AP and ML direction. The lamination anchor clamp must not be positioned so that it is aligned in the walking direction. If the clamp is in the wrong position, the lamination anchor cannot be supported. This will overload the shuttle lock and it will fail at the threaded connection.
--

**Attention!**

Prior to definitive use, clean the threads of the lamination anchor and of the screwed insert or connection adapter with acetone.

- Place the 616G4 Fiberglass Matting under the anchor arms
- Place 1 layer of 616G12 Carbon-Fiber Cloth under the anchor arms
- Thread the 699B1 Fiberglass Roving through the bore holes of the anchor arms and let it hang in loops
- Place the Carbon-Fiber Cloth over the anchor arms
- Place 1 layer of Carbon-Fiber Cloth medially and laterally up to the condylar area
- Place 1 layer of Carbon-Fiber Cloth (about 8 cm = 3.2 in. wide) circularly below the patella.
- Tie a double layer of the Nylglass Stockinette onto the adapter and reflect
- Tie 1½ layers of 616G13 Fiberglass Stockinette onto the Adapter and reflect.
- Tie a double layer of the Nylglass Stockinette onto the adapter and reflect.
- Pull over the soaked PVA bag and tie on proximally.

4 Assembly

**Caution!**

The shuttle lock housing may fracture if screwed out too much from the lamination anchor! This may cause the user of the prosthesis to fall and be injured.

To position the ratchet unit, never screw out the housing by more than 270°.

Fully screw the shuttle lock housing into the lamination anchor.

Then tighten the 501Z2=M5×22 Screw using the 710D4 Torque Wrench and secure it with 636K13 Loctite®.

Torque: 10 Nm.

**Caution!**

Improper tools may damage the ratchet unit or impair its function!

Never use an open-end wrench for unscrewing and retightening.
Use only the tool indicated in the Instructions for Use.

Secure the ratchet unit (4) with 636K13 Loctite® and screw it into the shuttle lock housing.

Torque: 20 Nm

Use the following tools:

- 710D4 Torque Wrench
- the suitable 709S26=1/4"×1/2" adapter piece and
- the S6A20-2 Socket Wrench Kit.

Caution!

Improper tools may damage the push button!

To remove the push button from the ratchet unit, please use a screw driver with 2 mm tip or a coin of comparable thickness.

Whenever grinding inside the socket, protect the shuttle lock from dirt by sealing the lock opening with tape or a piece of clean cloth. Please find further processing instructions in the Technical Information 2.1.5.

Caution!

Installation of the pin in a liner that is not approved

Falling due to loosening of the screw connection

1. Only install the pin in liners with a metal thread and observe the specified tightening torque.
2. **Inform the patient.**

Secure the Pin with Loctite® 241 (636K13), screw it into the liner's thread and tighten to 3 Nm (2.2 lbf.ft., 27 lbf.in.) using the 710D1 Torque Wrench or equivalent.

4.1 Shuttle Lock 6A20=20: Adjustment screws

Static corrections are easily made during alignment, trial fitting, or after finishing the prosthesis, using the four 506G3 Adjustment Screws (7). To exchange a module or during disassembly the previous position may be retained by loosening the two most deeply set, adjacent adjustment screws. Adjustment screws that appear to be too long at the desired adjustment may be replaced with shorter screws to avoid damage to the foam cover.

Use the 710D4 Torque Wrench to tighten the adjustment screws.

Torque moment: 15 Nm (11 lbf. ft., 133 lbf.in.).

Secure the adjustment screws in their final adjustment position of the prosthesis using 636K13 Loctite®.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH