4R23, 4R54, 4R73=A, 4R73=D, 4R74, 4R77

DE Gebrauchsanweisung	3
EN Instructions for Use	9
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	22
ES Instrucciones de uso	28
PT Manual de utilização	34
NL Gebruiksaanwijzing	41
SV Bruksanvisning	47
DA Brugsanvisning	53
NO Bruksanvisning	59
FI Käyttöohje	65
PL Instrukcja użytkowania	71
HU Használati utasítás	77
CS Návod k použití	83
RO Instrucțiuni de utilizare	90
HR Upute za uporabu	96
SL Navodila za uporabo	102
SK Návod na používanie	108
BG Инструкция за употреба	114
TR Kullanma talimatı	121
EL Οδηγίες χρήσης	127
RU Руководство по применению	134
JA 取扱説明書	141
ZH 使用说明书	146
KO 사용 설명서	152

4R54, 4R74

4R23

4R73=A

4R77


4R73=D

1

4R73=A, 4R73=D

2

4R77

3

1 Produktbeschreibung

INFORMATION

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INFORMATION

Date of the last update: 2012-03-19

- Please read this document carefully.
- Follow the safety instructions.

1.1 Construction and Function

The 4R23, 4R54, 4R73=A, 4R73=D, 4R74 and 4R77 Socket Adapters are designed for mounting on the socket of transfemoral and transtibial prostheses. They serve to connect with the distal prosthetic components. The 4R73=A and 4R73=D Socket Adapters each feature an eccentrically placed pyramid adapter. The adapters are equipped with different distal connections:

Reference number	Distal connection
4R23	Pyramid Adapter
4R54	Pyramid adapter with bore for Shuttle Lock
4R73=A	Pyramid adapter, 7 mm axial offset
4R73=D	Pyramid adapter, 5 mm diagonal offset

Reference number	Distal connection
4R74	Pyramid adapter with bore for Shuttle Lock
4R77	Pyramid adapter, rotating, with bore for Shuttle Lock

1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

2.2 Area of Application

Area of application according to the MOBIS® mobility system:

4R74

Approved for a body weight of up to **100 kg**.

4R23

Approved for a body weight of up to **125 kg**.

4R54, 4R73=A, 4R73=D and 4R77

Approved for a body weight of up to **150 kg**.

2.3 Environmental Conditions

Allowable Environmental Conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the user's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Non-observance of safety instructions

Injury and damage to the product

- Observe safety information in these instructions for use.
- Make the user aware of all safety instructions marked with "**Inform the user**".

CAUTION

Breakage of load-bearing components due to overloading

Injuries due to falling and damage to the prosthesis components

- Use the prosthesis components according to the MOBIS classification (see "Area of Application").
- **Inform the patient.**

CAUTION

Unallowable combination of prosthesis components

Injuries due to falling and damage to the prosthesis components

- Only combine the product with prosthesis components that are approved according to the chapter "Combination Possibilities".
- Based on the instructions for use, verify whether the prosthesis components may be combined with each other.

CAUTION

Use under environmental conditions that are not allowable

Injuries caused by falling and damage to the product

- ▶ Do not use the product under prohibited environmental conditions (see the chapter "Environmental Conditions").
- ▶ If the product has been exposed to prohibited environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the user.**

CAUTION

Exceeding the service life and reuse on other users

Injuries caused by falling, loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see the chapter "Service Life").
- ▶ Only use the product for a single user.
- ▶ **Inform the patient.**

CAUTION

Mechanical damage

Injuries and change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality" in this chapter).
- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the user.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

4R23, 4R54, 4R73=A, 4R73=D, 4R74, 4R77 Socket Adapter					
Fig.	Item No.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647H31
–	–	▲	–	Countersunk Head Screw	501S41=M6x12
–	–	▲	–	Countersunk Head Screw	501S41=M6x25
–	–	▲	–	Countersunk Head Screw	501S41=M6x30
The following parts are only for:				Socket Adapter	4R73=A, 4R73=D
2	1	▲	4	Countersunk Head Screw	501S41=M6x12x11.4
2	1	▲	4	Countersunk Head Screw	501S41=M6x25x11.4
2	1	▲	4	Countersunk Head Screw	501S41=M6x30x11.4
The following parts are only for:				Socket Adapter	4R77
3	1	▲	1	Pressure plate for 4R77	4Y19

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to falling and damage to the prosthesis components

- Observe the alignment and assembly instructions.

5.1 Installing the Socket Adapter

CAUTION

Risk of falling if care is not taken with the screw connections

Injuries due to falling and damage to the thread

- Clean the thread before every installation.
- Apply the installation torques provided.
- Pay attention to the instructions for securing the screw connections and the use of the correct length.

To install the 4R73=A* and 4R73=D* Socket Adapters, use the special screws according to the information provided in the following table.

Prosthetic component	Reference number (screw)
5R2 Lamination Disk	501S41=M6x12x11.4
4R101 Sliding Adapter	
5R1* Socket Attachment Block	501S41=M6x25x11.4
5R6=* Modular Adapter	501S41=M6x30x11.4

> **Recommended tools:**

710D4 Torque Wrench

- 1) Position the socket adapter and then use 2 countersunk head screws to fasten it on the posterior side (installation torque: **12 Nm**).
- 2) Tighten the 2 countersunk head screws on the anterior side (installation torque: **12 Nm**).

6 Maintenance

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

8 Legal Information

8.1 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Otto Bock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

8.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The

declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

8.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local manufacturer distribution company for more information.

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Otto Bock HealthCare GmbH