

4R136

DE	Gebrauchsanweisung	3
EN	Instructions for Use	6
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	11
ES	Instrucciones de uso	14
PT	Manual de utilização	17
NL	Gebruiksaanwijzing	20
SV	Bruksanvisning	23
DA	Brugsanvisning	29
NO	Bruksanvisning	29
FI	Käyttöohje	31
PL	Instrukcja użytkowania	34
HU	Használati utasítás	37
CS	Návod k použití	40
RO	Instrucțiuni de utilizare	43
HR	Upute za uporabu	46
BG	Инструкция за употреба	48
TR	Kullanma talimatı	51
EL	Οδηγίες χρήσης	54
RU	Руководство по применению	57
JA	取扱説明書	60
ZH	使用说明书	63

INFORMATION

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INFORMATION

Date of the last update: 2014-08-19

- Please read this document carefully.
- Follow the safety instructions.

Explanation of Symbols

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

The scope of delivery is shown on the front page.

1 Description**1.1 Intended use**

The Ottobock 4R136 exhaust valve is intended exclusively for use in prosthetic fittings of the lower limbs.

1.2 Field of application

In combination with a sealing sleeve, the 4R136 exhaust valve serves to create a vacuum in a prosthetic socket. Thanks to this technology, shearing and friction forces are effectively reduced.

1.3 Safety instructions

Please make your patients aware of the following safety instructions:

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to substances which could cause metal parts to corrode, such as fresh water, salt water and acids.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Failure to observe the installation instructions can result in complete loss of vacuum and therefore significantly compromise the adhesion.

1.4 Function

The 4R136 exhaust valve ejects excessive air out of the socket during the stance phase and closes during the swing phase to prevent further inflow of air. Thanks to this system, shearing and friction forces are effectively reduced.

2 Components

- V4 valve, 1 pc.
- Socket connector, angled, 1 pc.
- Socket connector, straight, 1 pc.
- Tube, 1 pc.
- Washer, 2 pcs.

3 Handling

3.1 Installation preparation

The following tools are needed in addition for installing the socket connector:

- Pliers
- 5, 6, 10 mm drill bit
- 2-component epoxy adhesive (5 min. curing time) or sealing resin-talcum mixture
- Insulating or polyethylene tape
- 1/4 UNF tap

3.2 Completion

Laminate the socket in the usual way. Drill a 10 mm bore hole into the posterior side of the socket, as close to the distal end as possible. Remove any burrs caused by drilling. Secure from inside with polyethylene adhesive tape. Fill the hole with lamination resin to achieve air tightness (Fig. 1). Drill a 5 mm hole in the centre of the sealing resin fill and cut the thread with a 1/4 UNF tap, or drill a 6 mm hole in the centre of the sealing resin fill and cut the thread in the hole using the connector. Verify that the thread of the

elbow connector is not longer than the laminate's wall thickness. If need be, use a washer (Fig. 2).

Mix the 2-component epoxy adhesive or the sealing resin-talcum mixture. Apply a thin layer of the mass to the thread of the connector (Fig. 3). When applying the adhesive, be careful not to let any adhesive get into the open end of the connector.

Screw the connector fully into the prepared hole using the pliers (Fig. 4). Slide the tube onto the connector.

Install the exhaust valve on the end of the hose. The arrow on the valve shows the direction of air flow (Fig. 5).

4 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Ottobock HealthCare GmbH