6A53, 6A54

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	24
PT	Manual de utilização	30
NL	Gebruiksaanwijzing	36
SV	Bruksanvisning	41
DA	Brugsanvisning	46
NO	Bruksanvisning	51
FI	Käyttöohje	56
PL	Instrukcja użytkowania	62
HU	Használati utasítás	67
CS	Návod k použití	72
RO	Instrucțiuni de utilizare	77
HR	Upute za uporabu	83
SL	Navodila za uporabo	88
SK	Návod na používanie	93
BG	Инструкция за употреба	98
TR	Kullanma talimatı	104
EL	Οδηγίες χρήσης	109
RU	Руководство по применению	115
JA	取扱説明書	121
ZH	使用说明书	126
KO	사용 설명서	130

1

2

3

4

INFORMATION

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INFORMATION

Last update: 2016-12-13

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Design and Function

The 6A53 and 6A54 Sliding Adapters permit the distal components of the modular prosthesis to be shifted in parallel, regardless of the angle adjustment by the pyramid adapters. The adjustments can be made either in the frontal plane (medial or lateral) or in the sagittal plane (anterior or posterior). The sliding adapters are equipped with different connections:

Product reference number	Connections
6A53	pyramid adapter/pyramid adapter
6A54	pyramid adapter/pyramid receiver

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

- The maximum approved body weight is specified in the technical data (see Page 13).

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Excessive strain on the product Risk of injury due to breakage of load-bearing components ► Use the product according to the specified area of application (see Page 9).
 CAUTION
Unallowable combination of prosthetic components Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Adapter	–
4	Set screw (only 6A54)	506G3=M8x12

5 Preparation for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Adjusting adapter

The adapter includes a pyramid adapter adjustable in 2 mm increments. The pyramid adapter can be aligned proximally or distally.

> **Required materials:** 710D4 torque wrench, 636K13 Loctite 241

- 1) Align the sliding adapter in the desired displacement direction (**a – p** or **m – l**).
- 2) Lift the ring on the adjustable pyramid adapter slightly in order to displace the pyramid adapter (see fig. 3).
- 3) Set the desired displacement and press the ring back down until it engages.

TIP: Mark the final position with a pen.

- 4) Connect the pyramid adapter with the applicable connecting component (see fig. 4). This fastens the pyramid adapter in its position.

5.2 Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> **Required materials:** 710D4 torque wrench, 636K13 Loctite® 241

1) Trial fitting:

Screw in the set screws.

Use the torque wrench to tighten the set screws (**10 Nm**).

2) Definitive mounting:

Use Loctite® to secure the set screws.

Screw in the set screws.

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).

- 3) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).

Selection table for set screws

Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X16	16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive.

The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH