


1M10 Adjust

DE	Gebrauchsanweisung	4
EN	Instructions for use	9
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	26
PT	Manual de utilização	32
NL	Gebruiksaanwijzing	37
SV	Bruksanvisning	42
DA	Brugsanvisning	48
NO	Bruksanvisning	53
FI	Käyttöohje	58
PL	Instrukcja użytkowania	63
HU	Használati utasítás	69
CS	Návod k použití	74
RO	Instrucțiuni de utilizare	80
HR	Upute za uporabu	85
SL	Navodila za uporabo	91
SK	Návod na používanie	96
BG	Инструкция за употреба	101
TR	Kullanma talimatı	107
EL	Οδηγίες χρήσης	112
RU	Руководство по применению	118
JA	取扱説明書	124
ZH	使用说明书	129
KO	사용 설명서	134

INFORMATION

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INFORMATION

Last update: 2018-03-26

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1.1 Construction and Function

The 1M10 Adjust prosthetic foot is suitable for use in modular prostheses. The prosthetic foot distinguishes itself with a high level of safety while standing and walking. This is ensured by a multi-axial joint in the load line. The heel characteristics can be individually adapted with a function module.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

Maximum body weight				
Size [cm]	22 to 23	24 to 25	26 to 27	28 to 30
Stiffness 1	52 kg	58 kg	72 kg	77 kg
Stiffness 2	68 kg	76 kg	95 kg	100 kg
Stiffness 3	80 kg	100 kg	125 kg	125 kg

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

Other environmental conditions	
UV light	Not resistant
Hydrolysis	Not resistant

2.4 Service life

This product has been tested by the manufacturer with 2 million load cycles according to ISO 22675. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION
Exceeding the service life and reuse on another patient
Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification.

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock	SL=Spectra-Sock-7

Spare parts/accessories (not included in the scope of delivery)

Fig.	Reference number	Reference number
1	Footshell	2C1
1	Insert for function module	2D11

5 Preparation for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

INFORMATION

- Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Alignment

5.1.1 Applying/removing the footshell

Inserting the foot into the footshell

- 1) Slide the prosthetic foot into the footshell (see fig. 3, see fig. 4).
- 2) Press the heel of the prosthetic foot into the footshell until it engages (see fig. 5).

Removing the footshell

- 1) Push the lock of the footshell to the rear and pull up the heel of the prosthetic foot.
- 2) Remove the prosthetic foot from the footshell.

5.1.2 Bench Alignment

TT bench alignment

Bench alignment process

Required materials: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:

Sagittal plane

①	Heel height: Effective heel height (shoe heel height - sole thickness of forefoot) + 5 mm
②	Exterior foot rotation: approx. 5°
③	Anterior placement, middle of the prosthetic foot to the alignment reference line: 30 mm
④	Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed.
⑤	Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: Individual residual limb flexion + 5°

Frontal plane

Bench alignment process	
6	Alignment reference line of prosthetic foot: between big toe and second toe Alignment reference line of prosthetic socket: along the lateral patella edge
7	Observe the abduction or adduction position.

TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

5.1.3 Static Alignment

- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.
- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

5.1.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, medialise the prosthetic foot. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation.

5.1.4.1 Optimising the heel characteristics

INFORMATION

Changing the configuration of the function module has no effect on the load line while standing.

The heel characteristics are adjusted with the function module. The following settings are possible:

Heel strike	Selected configuration (see fig. 2)	Position of insert
Soft	1	Insert removed, opening posterior
Moderate	2 (delivery condition)	Insert posterior
Firm	3	Insert anterior

- Turn the function module to the desired configuration (see fig. 6).

5.2 Optional: Installing the foam cover

A connecting element (such as a connection plate, connection cap or foam connection cap) provides a detachable connection between the foam cover and the prosthetic foot.

- > **Required materials:** Degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive
- 1) Measure the distance from the knee rotation point to the desired end of the foam cover and add the compression allowance: For PE foam **10 mm**, for PUR soft foam **30 mm**. For TF prostheses, add twice the compression allowance.
 - 2) Cut the pre-shaped foam block to length.
 - 3) Pull the foam cover over the prosthesis.
 - 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or sits on the foot adapter.
 - 5) Mount the prosthetic foot on the prosthesis.
 - 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
 - 7) Dismount the prosthetic foot and remove the connecting element.
 - 8) Clean the connecting element using a degreasing cleaner.

- 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
- 10) Let the glue dry (approx. **10 minutes**).
- 11) Mount the prosthetic foot and adapt the exterior cosmetic cover shape. Take into account compression due to cosmetic stockings or SuperSkin.

6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9 Technical data

Sizes [cm]	22	23	24	25	26	27	28	29	30
Heel height with 2C1=*N [mm]	10 ±5								
Heel height with 2C1=*S [mm]	20 ± 5					–			
System height with 2C1=*N [mm]	46	48	53		57	59	64		
System height with 2C1=*S [mm]	38		44		49	–			
Product weight without footshell [g]	255		320		385		465		
Max. body weight [kg]	80		100		125				
Mobility grade	1 + 2								
Colours	beige (4), light brown (15)								

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