

AxonArm Hybrid 12K500=*, AxonArm Ergo 12K501=*

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	21
FR Instructions d'utilisation (Personnel spécialisé)	39
IT Istruzioni per l'uso (Personale tecnico specializzato)	57
ES Instrucciones de uso (Personal técnico especializado)	75
PT Manual de utilização (Pessoal técnico)	93
NL Gebruiksaanwijzing (Vakmensen)	111
SV Bruksanvisning (Fackpersonal)	129
DA Brugsanvisning (Faguddannet personale)	147
NO Bruksanvisning (Fagpersonell)	165
FI Käyttöohje (Ammattihenkilöstö)	183
CS Návod k použití (Odborný personál)	201

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DE

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Table of contents

EN

1	Foreword	23
2	Product Description	23
2.1	Construction and Function	23
2.2	AFB (Automatic Forearm Balance) Flexion Aid	23
2.3	AxonMaster	24
2.4	EasyPlug	24
2.5	Forearm Connection Cables	24
2.6	Lock	24
2.6.1	Mechanical Release	24
2.6.2	Electronic Unlocking (AxonArm Ergo 12K501)	24
2.6.3	Controlling the Electronic Lock	24
2.7	Power Supply	24
2.8	Charging Receptacle	25
2.9	Rechargeable Battery	25
3	Application	25
3.1	Indications for Use	25
3.2	Orthotist / Prosthetist Qualifications	25
3.3	Conditions of Use	25
3.4	Application Area	25
4	Safety	25
4.1	Explanation of Warning Symbols	25
4.2	Structure of the safety instructions	26
4.3	General Safety Instructions	26
5	Scope of Delivery	27
6	Preparation for Use	27
6.1	Installing the AxonRotation Adapter	28
6.2	AxonMaster	28
6.2.1	Controlling the Electronic Lock	28
6.2.2	Connecting the Electrodes	29
6.2.3	Installing the AxonMaster in the Upper Arm Socket	29
6.2.4	Sealing the prosthesis	29
6.2.5	Connecting the AxonMaster to the EasyPlug	29
6.3	Installing the EasyPlug	29
6.4	Connecting the Elbow Joint to the Upper Arm	30
6.5	Setting the Preflexion	31
6.6	Adjusting the Humeral Rotation Feature	33
6.7	Installing the Pull Cable	33
6.8	Checking the Symmetry	33
7	Handling	33
7.1	Charging the Battery	33
7.2	Switching On the Axon-Bus Prosthetic System	33
7.3	Switching Off the Axon-Bus Prosthetic System	33
7.4	Determining the Battery Charge Level	33
7.5	Safety Shutoff	34
7.6	Open the Axon-Bus Gripping Component in an Emergency	34
7.7	Bluetooth Function	34
7.8	Beep Function	34
8	Settings	35
9	Maintenance	35
9.1	Service Information	35
10	Disposal	35

11	Legal Information	35
11.1	Liability	35
11.2	Trademarks	35
11.3	CE Conformity	35
12	Appendices	36
12.1	Fault Display and Fault Elimination	36
12.1.1	LED flashes red and beeper sounds for 3 seconds.....	36
12.1.2	The LED flashes three times after pressing the button in the colour of the current battery capacity	36
12.1.3	LED alternately flashes yellow/red and beeper sounds for 3 seconds	36
12.1.4	The movements of the Axon-Bus gripping component are becoming increasingly slower.....	36
12.2	Symbols Used	36
12.3	Technical Data.....	36
12.4	Accessories	37
12.5	Axon	37

1 Foreword

INFORMATION

Last update: 2014-11-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The 12K500 AxonArm Hybrid and the 12K501 AxonArm Ergo are hereinafter each referred to as the 'AxonArm' or the 'product'.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product Description

2.1 Construction and Function

12K500 AxonArm Hybrid

The 12K500 AxonArm Hybrid is a mechanically locked elbow joint.

Locking and unlocking is carried out mechanically, e.g. by means of body harnesses.

12K501 AxonArm Ergo

The 12K501 AxonArm Ergo is a myoelectrically locked elbow joint.

Locking and unlocking is controlled myoelectrically by means of electrode signals or with a switch.

12K500 AxonArm Hybrid/12K501 AxonArm Ergo

In combination with other Ottobock prosthetic components of the Axon-Bus prosthetic system such as the 8E500 Michelangelo Hand, it helps the wearer perform everyday tasks.


The key advantages of the AxonArm are as follows:

- High holding force
- Natural free swing behaviour
- Low noise level, no background noise during free swing
- Attractive appearance

The 757B501 AxonEnergy Integral supplies power to the Axon-Bus components.

If the Axon-Bus prosthetic system is off or the battery is empty, the elbow joint can still be locked and unlocked in any position, even under load, by means of the pull cable.

2.2 AFB (Automatic Forearm Balance) Flexion Aid

The flexion aid provides for harmonious movements and noiseless free swing.

When the prosthetic arm is extended, it stores the energy generated by gravity and uses it to support flexion. As a result, the patient only has to apply minimal force on the flexion cable or through body momentum. The level of flexion assistance can be adjusted to the individual weight of the prosthetic forearm with the Axon-Bus gripping component and different clothing by means of a dial.

Adjusting the flexion aid (see Page 31)

2.3 AxonMaster

Various control modes and switching methods as well as the option of implementing a fine-tuned patient-specific adjustment are available for optimum adaptation to the patient. For information please refer to the 647G590 Instructions for Use for the 13E500 AxonMaster.

Two MyoBock electrodes, e.g., 13E200/202=*, or if necessary a pull switch, e.g., 9X18, can be connected to the 13E500 AxonMaster.

2.4 EasyPlug

EasyPlug is a through-connection integrated into the elbow joint for connecting the Axon-Bus. Since there are no external cables, its overall durability is increased, reducing the risk of defects caused by broken cables and enhancing cosmetic appearance at the same time. With the AxonArm only the AxonMaster 13E500 is attached in the region of the elbow ball. All other slots have no function and are thus covered by caps.

2.5 Forearm Connection Cables

The forearm cables establish the electrical connection between the AxonArm and the Axon-Bus gripping component. They are already installed on the flexion aid and can be replaced if damaged.

2.6 Lock

The internal infinite-position lock on the AxonArm can be unlocked or locked under load in any position. In the locked state, a new AxonArm with a forearm length of 305 mm has a load capacity of 230 N (this value may decline slightly over the course of several years due to wear). The lock will slip under higher loads.

2.6.1 Mechanical Release

Mechanical unlocking and locking can be carried out by means of a pull cable when the AxonArm Ergo 12K501 is off or if the rechargeable battery is empty. The AxonArm Hybrid 12K500 can only be unlocked and locked using the pull cable because it has no electronic locking element.

The lock can also be released under load with a constant low applied force of approx. 10 N.

A strong pull on the cable (clicking sound) locks or unlocks the joint permanently.

2.6.2 Electronic Unlocking (AxonArm Ergo 12K501)

The locking and unlocking of the elbow is controlled myoelectrically by means of electrode signals or a switch. The control devices must be connected to the AxonMaster 13E500. For further information on control and the control programs please refer to the 647G590 Instructions for Use for the AxonMaster 13E500.

2.6.3 Controlling the Electronic Lock

The electronic lock can be controlled myoelectrically by the following means configured at the factory:

- A long co-contraction on channels 1 and 2 of the AxonMaster releases or locks the lock.
- A short co-contraction switches the gripping mode of the gripper.
- A pulse on channel 3 of the AxonMaster releases or locks the lock.

2.7 Power Supply

Power to the AxonArm and the Axon-Bus gripping component is supplied by the 757B501 AxonEnergy Integral positioned in the forearm. The capacity is sufficient for one day of normal everyday activities. Regular overnight recharging of the battery is recommended. LEDs on the battery charger provide information about the battery charge level.

INFORMATION

Switching off the Axon-Bus prosthetic system during extended passive pauses (e.g. air or rail travel, visit to a theatre or cinema, etc.) will make the battery last longer before requiring a recharge. Only the entire Axon-Bus prosthetic system with all connected components can be switched off. Individual Axon-Bus components cannot be switched off separately.

2.8 Charging Receptacle

The charging receptacle (Fig. 1, item 2) has the following functions when the AxonArm Ergo is connected with the AxonMaster:

- Contacts for charging the battery
- LED display to indicate the current battery charge level
- LED display for providing feedback on operating states
- Button to turn the Axon-Bus prosthetic system on and off, display the battery charge level and open the Axon-Bus gripping component in an emergency
- Beeper for providing feedback on operating states
- Button for activating the Bluetooth function

2.9 Rechargeable Battery

The rechargeable battery (Fig. 1, item 1) consists of Li-ion cells with integrated electronics. This protects the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range.

3 Application

3.1 Indications for Use

The AxonArm Hybrid 12K500 and the AxonArm Ergo 12K501 may only be used **exclusively** for exoprosthetic fittings of the upper limbs.

3.2 Orthotist / Prosthetist Qualifications

Fitting a patient with the AxonArm may only be carried out by prosthetists who have been authorised by Ottobock after completion of a corresponding training course.

3.3 Conditions of Use

The AxonArm was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist unit (pushups, downhill, mountain biking) or extreme sports (free climbing, paragliding, etc.).

The applicable environmental conditions are described in the technical data.

The AxonArm is intended for use by one user only. Use of the product by another person is not approved by the manufacturer.

3.4 Application Area

The AxonArm is intended solely for use with components of the Axon-Bus prosthetic system. It can be used for unilateral or bilateral amputees from a transhumeral amputation level to shoulder disarticulation.

4 Safety

4.1 Explanation of Warning Symbols

	Warning regarding possible serious risks of accident or injury.
	Warning regarding possible risks of accident or injury.
	Warnings regarding possible technical damage.

4.2 Structure of the safety instructions

CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General Safety Instructions

CAUTION

Signs of wear on system components

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ In the interest of user safety and in order to maintain the operational reliability and the validity of the warranty, the specified service intervals must be observed.

CAUTION

Operator errors during the adjustment procedure

Injury due to faulty control of the Axon-Bus prosthetic system.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During product training, you will receive a password which provides you with authorisation to access the 560X500=* AxonSoft software. To qualify for software updates, additional product training courses may be necessary.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

CAUTION

Unauthorised manipulation of Axon-Bus components

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Manipulations to the Axon-Bus prosthetic system other than the tasks described in these instructions for use are not permitted.
- ▶ Do not damage the battery or separate the connection cables between the battery packs.
- ▶ The Axon-Bus prosthetic system and any damaged Axon-Bus components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.

CAUTION

Incorrect electrode settings

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact Ottobock Myo-Service.
- ▶ Set the electrode gains as low as possible to reduce the risk of malfunctioning due to strong electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ Allow the patient to rest during the adjustment of the electrodes; otherwise, muscle fatigue will lead to inconsistent results. As a result, the therapist will tend to set the electrode gains too high.

CAUTION

Risk of pinching in elbow joint flexion region

- ▶ Injuries due to pinching of body parts.
- Ensure that fingers and other body parts are not in this area when bending the elbow joint.

⚠ CAUTION**Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

⚠ CAUTION**Penetration of dirt and moisture**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Do not let foreign particles or liquids get into the prosthetic arm.
- ▶ Do not expose the prosthetic arm, especially the elbow joint, to dripping or splashing water.
- ▶ Wear the prosthetic arm and especially the elbow joint under suitably resistant clothing in the rain.

NOTICE**Coating, gluing or painting the product**

Damage or fracture due to chemical processes.

- ▶ The product must not be coated, glued or painted.

5 Scope of Delivery

Piece (s)	Product name
1	12K500 AxonArm Hybrid or 12K501 AxonArm Ergo
1	Cable Lock
1	Lamination Ring (installed on AxonArm)
1	Lamination Protection Cover (hollow ball)
1	Lamination Protection Cover (disc with arrow)
1	13Z146 Lamination Ring for Axon-Bus gripping component
1	Spherical Cord Lock
2	Tongue Cap
4	Cable Ties 420mm x 4mm
1	Clamp Ring
4	Fixing Brackets
4	Self-Tapping Pan Head Screws
1	647G790 Instructions for Use
1	3.3.6 Technical Information (646T3=3.6-DE Casting and Lamination Instructions)
1	3.3.6 Technical Information (646T3=3.6-EN Casting and Lamination Instructions)

INFORMATION

The clamp ring, mounting brackets and pan head screws are also available as a set, "Alignment Aid for Elbow Components" under reference number 743A23.


INFORMATION

Only gripping components of the Axon-Bus prosthetic system, e.g. the 8E500 Michelangelo Hand, can be used with the 12K500 AxonArm Hybrid/12K501 AxonArm Ergo. Other gripping components may not be used.

6 Preparation for Use

For the fabrication of the patient's socket and cutting the forearm socket to length, consult the 646T3=3.6 Technical Information.

6.1 Installing the AxonRotation Adapter

1) Use silicone grease to seal the plug connection.

2) Connect the AxonRotation Adapter with the connecting cable.

3) Position the AxonRotation Adapter such that the **middle screw holes on the lamination ring** can be used.

INFORMATION: The other screw positions can be used for correcting the alignment of the AxonRotation Adapter after the lamination ring has been glued in.

4) Screw the AxonRotation Adapter to the lamination ring.

6.2 AxonMaster

Various control modes and switching methods as well as the option of implementing a fine-tuned patient-specific adjustment are available for optimum adaptation to the patient. For information please refer to the 647G590 Instructions for Use for the 13E500 AxonMaster.

Two MyoBock electrodes, e.g., 13E200/202=*, or if necessary a pull switch, e.g., 9X18, can be connected to the 13E500 AxonMaster.

Input on AxonMaster	Electrode
Input 1	Hand closing
Input 2	Hand opening
Input 3	Pull switch (optional)


► Remove the cover caps from the required inputs. The remaining inputs cannot be used at this time and therefore remain covered by the caps.

6.2.1 Controlling the Electronic Lock

The electronic lock can be controlled myoelectrically by the following means configured at the factory:

- A long co-contraction on channels 1 and 2 of the AxonMaster releases or locks the lock.
- A short co-contraction switches the gripping mode of the gripper.
- A pulse on channel 3 of the AxonMaster releases or locks the lock.

6.2.2 Connecting the Electrodes

- 1) Remove the cover caps from the required inputs.
- 2) Seal the slots on the AxonMaster with silicone grease.

- 3) Connect the electrode cables to the AxonMaster. Observe the slot assignments:
Slot 1: close
Slot 2: open
Slot 3: pull switch (optional)

6.2.3 Installing the AxonMaster in the Upper Arm Socket

The AxonMaster can be secured with adhesive tape or a hook-and-loop fastening strap. The AxonMaster can be covered with film to provide additional protection against moisture ingress.

Should the available space in the upper arm socket be too small an area for the AxonMaster (lamination dummy) must be reserved during lamination of the inner socket.

6.2.4 Sealing the prosthesis

- 1) After laminating the lamination ring, remove the grease and any dirt.
- 2) Lubricate the groove, inner edge and thread of the lamination ring with 633F30 Special Grease.

6.2.5 Connecting the AxonMaster to the EasyPlug

- 1) Grease the plug of the AxonMaster cable with 633F11 silicone grease.
- 2) Connect the cable.
- 3) Apply 633F11 silicone grease around the plug and blanking covers to seal.

6.3 Installing the EasyPlug

- 1) Lubricate the groove, inner edge and thread of the lamination ring with 633F30 special grease.
- 2) Push the EasyPlug with the attached cable into the lubricated lamination ring and press in firmly. The EasyPlug enclosure must engage completely.
→ **You must be able to hear a loud click!**

⚠ CAUTION

Misalignment of the EasyPlug during installation

- > Injury due to faulty control or malfunction of the product.
- > Reduction in mechanical strength due to thread misalignment.
- ▶ The circuit board and the edge of the lamination ring must be parallel to each other.


▶ **Left picture:**

The EasyPlug has been snapped in properly. The circuit board and the edge of the lamination ring are parallel to each other.

Right picture:

The EasyPlug has not been snapped in properly; it is crooked. If the EasyPlug is extremely crooked (as shown in the picture), you need to remove it and insert it again because otherwise it may become damaged.


6.4 Connecting the Elbow Joint to the Upper Arm


- 1) Unscrew the frictional clamp until the top thread is lined up with the fixed thread (right picture).


- 2) Position the lamination ring on the upper arm socket such that the notch in the lamination ring is across from the start of the thread.

- 3) Mark the base position for the flexion angle (preflexion) on the joint ball. Dimension $Y = 1.5$ cm.


- 4) Thread the belt under the axis.
- 5) Pull the belt up to the marking (dimension Y) with tweezers.

- 6) Twist the forearm to the upper arm socket until the notch on the lamination ring is visible above the guide of the strap clamp (left picture). Otherwise the strap clamp will not be able to be inserted (right picture).

INFORMATION: Ensure the strap is straight and not twisted.


- 7) Slide the strap clamp under the lamination ring.

- 8) Screw the strap clamp in.

6.5 Setting the Preflexion

- The level of AFB compensation (flexion assistance) can be adjusted and adapted to the weight of different clothing with the dial. Adjustments are easier to make when the forearm is flexed.
- The adjustment mechanism is protected by a slip coupling.
- The hydrodynamic damper of the flexion aid mechanism is designed to counterbalance the weight of the forearm plus Axon-Bus gripping component. Please be sure to stabilise the upper arm component of the AxonArm Ergo while testing the joint functions.

Region 1

Minimal compensation that allows natural free swing of the arm while walking.

Region 2

Compensation increases progressively with flexion of the arm and decreases during extension.

Region 3

Compensation remains constant. If set correctly, the weight of the forearm is balanced by the compensation and the elbow remains flexed.

Region 4

Reduced compensation until flexion stop.

⚠ CAUTION

Lock released during donning, doffing or storing of the extended arm


Injury due to sudden bending of the forearm or the upper arm socket due to released lock.

- The Axon-Bus prosthetic system must always be flexed for donning and doffing.
- The Axon-Bus prosthetic system must be stored in flexed condition.


⚠ CAUTION

Release of lock when arm horizontal or raised vertically

- Injury due to sudden bending of the arm due to released lock and lack of a load, especially with a high compensation force setting.
- Inform the patient of this danger.

- 1) Flex the AxonArm Ergo and minimise the compensation force using the dial.

- 2) Flex the AxonArm Ergo until a loop can be made in the strap using your finger (approx. 45°) and lock the joint using the pull cable.


- 3) Secure the strap (e.g. using an artery clamp).
INFORMATION: When using an artery clamp, verify that the jaws of the artery clamp have no sharp edges. Otherwise, tape over the jaws to protect the strap.
- 4) Loosen the screw on the strap clamp, but do not remove the strap clamp.
- 5) Shorten or lengthen the strap.
Lengthening the strap: preflexion is decreased.
Shortening the strap: preflexion is increased.
- 6) Tighten the screw on the strap clamp.
- 7) Remove the strap securing device (artery clamp).

A 3 mm change in the strap length results in an approx. 5° change in the flexion angle (preflexion).

The free swing range is changed by the same amount as the flexion angle. The start of the adjacent compensation region changes accordingly.

The factory-set preflexion is approx. 10°.

6.6 Adjusting the Humeral Rotation Feature

- The humeral rotation turntable has integral stops (+/- 80 degrees) to block extreme movements. Turntable friction is easily adjusted by turning of the external adjustment screw, using a coin or similar object.

6.7 Installing the Pull Cable

⚠ CAUTION

Removing the pull cable

Injury due to malfunction of the product.

- The pull cable on the 12K501 AxonArm Ergo is to be used for emergency operation only. For safety reasons, it must not be removed. When the Axon-Bus prosthetic system is switched off or defective and when the battery is empty, the AxonArm cannot be locked or unlocked without the pull cable.

- 1) Bring the forearm into neutral position (no outward or inward rotation).
- 2) Affix the cable guide to the socket with cap screws (centred in relation to the strap clamp and approx. 30 mm away from the lamination ring).
- 3) Test inward and outward rotation of the socket up to the stop. Ensure that the pull cable cannot be tensioned (loop).

6.8 Checking the Symmetry

- 1) Connect the Axon-Bus gripping component (e.g. Michelangelo Hand) to the AxonArm.
- 2) Check the symmetry with respect to the other hand (e.g. using the 743L20=230 Ottobock Laser-Line).

7 Handling

7.1 Charging the Battery

The 757B501 AxonEnergy Integral must be charged prior to initial use. The 757B501 AxonEnergy Integral may **only** be charged with the 757L500 AxonCharge Integral. The charging process can only be performed with the 13E500 AxonMaster connected.

- 1) Connect the charging plug of the 757L500 AxonCharge Integral to the charging receptacle.
 - The beeper sounds briefly 2x.
 - The Axon-Bus prosthetic system is deactivated and the charging process starts automatically.
- 2) To activate the Axon-Bus prosthetic system, disconnect the charging plug and switch on the Axon-Bus prosthetic system.

INFORMATION: The prosthesis cannot be used during the charging process.

7.2 Switching On the Axon-Bus Prosthetic System

- Press the charging receptacle button and hold for one second (Fig. 2, arrow).
- The beeper emits two short signals.

7.3 Switching Off the Axon-Bus Prosthetic System

- Press the charging receptacle button and hold for more than one second (Fig. 2, arrow).
- The beeper emits one long signal.

7.4 Determining the Battery Charge Level

The battery charge level can be queried at any time.

1. With the Axon-Bus prosthetic system switched on, press the charging receptacle button (Fig. 2, arrow) and hold for less than one second.
2. The LED is illuminated. The colour of the LED indicates the current battery charge level.

Battery empty

Orange LED is illuminated

Battery 50% charged

Yellow LED is illuminated

Battery fully charged

Green LED is illuminated


7.5 Safety Shutoff

The purpose of the safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit.

After a short circuit the battery charger must be connected and then disconnected in order for the electronics to be activated.

7.6 Open the Axon-Bus Gripping Component in an Emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open. A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

7.7 Bluetooth Function

Pressing the button in the charging receptacle and holding for more than four seconds with the Axon-Bus prosthetic system turned off activates the Bluetooth function of the connected AxonMaster 13E500. The blue LED lights up. Once a connection to the PC has been successfully established, the LED flashes blue.

7.8 Beep Function

Action	Beep 1x long	Beep 2x short	Beep 3x short	Beep 1x short, 1x long	Beep pulsating
Switch on the Axon-Bus prosthetic system Press and hold the button for one to four seconds.					
Switch off the Axon-Bus prosthetic system Press the button and hold for more than one second.					
Charge the Axon-Bus prosthetic system Connect charging plug, charging process starts.					
Stop charging Disconnect charging plug. Charging ends and the Axon-Bus prosthetic system is switched off.					
Empty battery indicator					
Switch on the Bluetooth function (if the Axon-Bus prosthetic system was previously switched					

Action	Beep 1x long	Beep 2x short	Beep 3x short	Beep 1x short, 1x long	Beep pulsating
off) Press the button and hold for more than four seconds.					
Open the Axon-Bus gripping component in an emergency Press the button and hold for approximately three seconds, until the Axon-Bus gripping component opens.					

8 Settings

Adjustments can be made to the AxonArm and the Axon-Bus gripping component via Bluetooth data transfer and the 560X500=* AxonSoft software. For this purpose, a wireless connection must be established between the 13E500 AxonMaster and the PC with the help of the 60X5 BionicLink PC.

For further functions, please see the instructions for use for the 13E500 AxonMaster and 60X5 BionicLink PC.

9 Maintenance

9.1 Service Information

Since all moving mechanical parts are subject to wear and tear, regular service is required within the guarantee period. In the course of service, Ottobock Service (Ottobock Myo-Service) checks the entire product. Wear parts are replaced as required. Omitting a scheduled service inspection voids the guarantee.

10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that does not comply with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

11 Legal Information

11.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

11.2 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The

declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

12 Appendices

12.1 Fault Display and Fault Elimination

12.1.1 LED flashes red and beeper sounds for 3 seconds

Probable cause: Error in the Axon-Bus prosthetic system.

Remedy:

1. Switch the Axon-Bus prosthetic system off and back on.
2. Check the plug connections in the socket.
3. Connect the charger. If the charger displays the service symbol (yellow spanner) or the system error symbol (red spanner) send the Axon-Bus prosthetic system to Ottobock Myo-Service.

12.1.2 The LED flashes three times after pressing the button in the colour of the current battery capacity

Probable cause: The connection from the AxonEnergy Integral to the Axon-Master has been interrupted.

Remedy:

1. Check the connection (e.g. plug connections, flex cable, EasyPlug etc.).
2. If no fault can be found, send the Axon-Bus prosthetic system to Ottobock Myo-Service.

12.1.3 LED alternately flashes yellow/red and beeper sounds for 3 seconds

Probable cause: The Axon-Bus prosthetic system is overheated.

Remedy:

1. Allow the Axon-Bus prosthetic system to cool for a few minutes.
2. Do not use the Axon-Bus prosthetic system at high ambient temperatures.

12.1.4 The movements of the Axon-Bus gripping component are becoming increasingly slower

The patient is alerted to the diminishing battery charge level as the Axon-Bus gripping component increasingly slows down and the gripping force decreases. Finally, when there is very little battery capacity remaining, the Axon-Bus prosthetic system switches off automatically to protect the battery against harmful deep discharge.

12.2 Symbols Used

Declaration of conformity according to the European Directive for Medical Devices 93/42/EEC.

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that does not comply with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions provided by the responsible authorities in your country regarding return and collection processes.

SN YYYY WW NNN

Serial number of device

Legal manufacturer

12.3 Technical Data

Description	12K500=* AxonArm Hybrid	12K501=* AxonArm Ergo
Operating temperature	0 °C to 60 °C/32 °F to 140 °F	0 °C to 60 °C/32 °F to 140 °F
Storage and transport in original packaging	-20 °C to 40 °C/-4 °F to 104 °F	-20 °C to 40 °C/-4 °F to 104 °F

Description	12K500=* AxonArm Hybrid	12K501=* AxonArm Ergo
Storage and transport without packaging	-20 °C to 40 °C/-4 °F to 104 °F	-20 °C to 40 °C/-4 °F to 104 °F
Relative humidity	max. 80 % relative humidity, non-condensing	max. 80 % relative humidity, non-condensing
Long-term storage temperature	-20 °C to 40 °C/-4 °F to 104 °F	-20 °C to 40 °C/-4 °F to 104 °F
Dimensions	approx. 345 mm x d55 mm to d80 mm / 13.6 inch x d2.16 inch to d3.15 inch (arm extended)	approx. 345 mm x d55 mm to d80 mm / 13.6 inch x d2.16 inch to d3.15 inch (arm extended)
Weight without Axon-Bus gripping component	approx. 840 g/29.6 oz	approx. 840 g/29.6 oz
Service life	5 years	5 years
Forearm length (as delivered)	305 mm / 12 inch	305 mm / 12 inch
Min. forearm length	187 mm / 7.36 inch	187 mm / 7.36 inch
Flexion angle	15° - 145°	15° - 145°
Maximum vertical load with locked elbow joint and a forearm length of 305 mm / 12 inch	23 kg/ 50.7 lb	23 kg/ 50.7 lb
Operating time when fully charged	approx. 2500 gripping and locking cycles or approx. 8 hours	approx. 2500 gripping and locking cycles or approx. 8 hours
Battery service life	2 years	2 years
Charging time	3.5 hours	3.5 hours
Charging temperature	0 °C to 40 °C/32 °F to 104 °F	0 °C to 40 °C/32 °F to 104 °F

12.4 Accessories

The following additional components are necessary for operating the AxonArm:

Piece (s)	Product name	Reference number
1	AxonCharge Integral	757L500
1	AxonMaster	13E500
1	AxonRotation Adapter	e.g. 9S501
1	AxonFlexion Adapter	e.g. 9S500=*
1	AxonSoft	560X500=*
1	Axon-Bus Gripping Component	e.g. 8E500 Michelangelo Hand
2	Electrodes	e.g. 13E200=*, 13E202=*
1	Pull Switch (optional)	e.g. 9X18
1	Prosthetic Glove	e.g. 8S500=R/L AxonSkin Natural when using the 8E500 Michelangelo Hand

Additional MyoBock components such as suction socket electrodes, electrode cables, pull-in tube and tube valve for suction socket if required.

12.5 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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Table of contents

EN

1	Foreword	16
2	Introduction	16
3	Purpose	16
3.1	Intended Use	16
3.2	Application	16
3.3	Conditions of Use	16
3.4	Qualification	16
4	Scope of Delivery, Tools and Instructions	17
4.1	Scope of Delivery	17
4.2	Accessories and Tools	17
4.3	Instructions	18
5	Preparation for Use	18
5.1	Adjustment, Fitting, Use	18
5.2	Disassembling the Patient's Prosthesis	18
5.3	Length Ranges	20
5.4	Adjusting the Length	21
5.5	Finishing	22
6	Service	24
7	Disposal	24
8	Legal Information	24
8.1	Liability	24
8.2	Trademarks	24
8.3	CE Conformity	24
9	Appendix	25
9.1	Technical data	25
9.2	Symbols Used	25

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The 12K501=M-S AxonArm Ergo (Service/Test) is hereinafter referred to as the AxonArm Ergo Service.

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2 Introduction

This document explains all relevant steps from the disassembly of the patient prosthesis to the assembly and delivery of the finished AxonArm Ergo Service. It is intended for qualified professionals and assumes that this qualified personnel is trained in the handling of the various materials, machines and tools.

- ▶ **INFORMATION: This document does not replace the 647G790 Instructions for use of the 12K500=* AxonArm Hybrid/12K501=* AxonArm Ergo containing detailed information on the handling and safety of your product. Please read the 647G790 Instructions for Use carefully!**

3 Purpose

3.1 Intended Use

The AxonArm Ergo Service is to be used solely for the duration of the service and for trial fittings for exoprosthetic fittings of the upper limbs.

The AxonArm Ergo Service is compatible with the upper arm sockets of the 12K500=* AxonArm Hybrid and 12K501=* AxonArm Ergo.

3.2 Application

The AxonArm Ergo Service is suitable for unilateral and bilateral above-elbow amputees.

The AxonArm Ergo Service is to be used **exclusively** for the fitting of adults.

3.3 Conditions of Use


The AxonArm Ergo Service was designed for everyday activities and must not be used for unusual activities such as extreme sports (i.e. free climbing, paragliding, etc.).

3.4 Qualification

The fitting of patients with the AxonArm Ergo Service may only be carried out by a prosthetist.

4 Scope of Delivery, Tools and Instructions

4.1 Scope of Delivery

Item	Piece(s)	Product name	Article number
1	1	AxonArm Ergo Service	12K501=M-S
2	1	Axon adapter, short	13Z152
2	1	Counter nut, large	13Z111
3	1	Axon adapter, long	13Z153
3	1	Counter nut, small	13Z115
4	1	Adapter, short	13Z112
5	1	Adapter, long	13Z113

4.2 Accessories and Tools

The following accessories and tools are supplied for the assembly of the AxonArm Ergo:

Item	Tool	Article number
1	Disposable foam cover	13Z82
2	Strap wrench	13Z116
3	Hook wrenches	13Z109
4	Case key	

4.3 Instructions

The following accompanying documents necessary for assembly can be found on the inner side of the case under the foam cover:

Item	Description	Product number
not illustrated	Instructions for Use – 12K500=* AxonArm Hybrid/12K501=* AxonArm Ergo	647G790
not illustrated	Alignment and Assembly Instructions for 12K501=M-S AxonArm Ergo (Service/Test)	647G976
not illustrated	Information slip for reusable package	647H324

5 Preparation for Use

5.1 Adjustment, Fitting, Use

Read the 647G790 Instructions for Use before fitting the patient with the AxonArm Ergo. This contains important additional information about the use of the AxonArm Ergo.

5.2 Disassembling the Patient's Prosthesis

Before adapting the product to the patient, the patient's elbow component must be removed from the socket.

Example: patient socket, AxonArm and Michelangelo Hand

Turn off the AxonArm.

- **INFORMATION:** The AxonArm can only be turned off when the 13E500 AxonMaster is connected.

Pull off the prosthetic component.

► **INFORMATION:** For pulling off the prosthetic component, consult the corresponding instructions for use for AxonBus gripping components, e.g. 447D115 Instructions for Use for the 8E500 Michelangelo Hand.

Bend the AxonArm.

Release the lock using the pull cable if required.


Unscrew the AxonRotation adapter.


Set the dial to minimum compensation.


Extend the AxonArm.


Lock the elbow joint in this position using the pull cable to fix the removal and installation position.

Remove the screw on the strap clamp and remove the clamp.


Turn the screw. The humeral rotation feature is loosened and the patient socket can be unscrewed more easily.

Determine the length. Measure the distance between the front edge of the socket and the rotation centre of the joint. Document value 'X'.

5.3 Length Ranges


Length ranges from 190 mm to 302 mm can be covered by combining long and short Axon adapters.

Length range 190 mm – 210 mm

Required accessories:


1 pc. Axon adapter, short with counter nut (item 2)

Length range 215 mm – 256 mm

Required accessories:

- 1 pc. Axon adapter, long with counter nut (item 3)
- 1 pc. adapter, short (item 4)

Length range 261 mm – 302 mm

Required accessories:

- 1 pc. Axon adapter, long with counter nut (item 3)
- 1 pc. adapter, long (item 5)


5.4 Adjusting the Length

Based on the established length, the corresponding version can be configured and assembled.

The length range of 215 mm – 256 mm is illustrated here as an example.


The adapter, short and the Axon adapter, short with counter nut were already threaded into the AxonArm Ergo Service.

Screw the Axon adapter, short with counter nut in or out as needed and set the documented measurement "X". See the section "Disassembling the Patient's Prosthesis".


Mark the radial position on the Axon adapter, short or Axon adapter, long. The position between the 2 mounting screw holes must lie in an imaginary line to the charging receptacle.

► **INFORMATION:** The threaded tube must be screwed in at least to the point where the marker bore holes (arrow) are covered by the counter nut. If this is not possible, the next longer version has to be chosen.

Hand-tighten the adapter and counter nut of the threaded tube clockwise.


Thread the AxonArm Ergo Service into the strap wrench and fix it. First apply the hook wrench to the adapter (optional) and finally to the counter nut of the threaded tube. Tighten the nut clockwise until the prescribed torque is reached. Ensure that the radial position of the Axon adapter is not twisted.

The tightening torque of 16 Nm is reached as soon as the insides of the handles touch each other at the hook wrench.


5.5 Finishing

The process of adapting the AxonArm Ergo Service to the patient is completed by pulling on the foam cover and installing the AxonRotation adapter.


Pull the disposable foam cover over the AxonArm Ergo Service.

Mark the corresponding length according to the measurement 'X' determined previously (see the section "Disassembling the Patient's Prosthesis").


Use scissors to cut the disposable foam cover to length.

Use silicone grease to seal the plug connection.


Connect the AxonRotation adapter with the connecting cable.

Mount the AxonRotation adapter to the Axon adapter, short or the Axon adapter, long. Observe the installation position when using the 9S503 AxonRotation adapter.

► **INFORMATION:** Do not twist the Axon adapter, short or Axon adapter, long after installing the AxonRotation adapter. The cables could be damaged.

Pull the disposable foam cover previously cut to length over the AxonArm Ergo Service.
Mount the prosthetic component.
Pull the cosmetic glove over the Michelangelo Hand if applicable.
In case of a fitting with the 9S503 AxonRotation Adapter, make sure that the rotational movement is working smoothly.

- 1) Remove the lamination ring of the AxonArm Ergo Service (see the section "Disassembling the Patient Prosthesis").
- 2) Connect the elbow joint to the upper arm and configure the settings. See the 647G790 Instructions for Use, sections "Connecting the Elbow Joint to the Upper Arm", "Setting the Preflexion", "Adjusting the Humeral Rotation Feature", "Installing the Pull Cable" and "Checking the Symmetry".

6 Service

After each fitting, the AxonArm Ergo Service has to be sent to the authorised Ottobock Service Centre for inspection. There, the complete product is examined and readjusted if required. Wear parts are also replaced if required.

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8.3 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

9 Appendix

9.1 Technical data

Description	Value
Reference number	12K501=M-S
Operating temperature	0 °C to 60 °C/32 °F to 140 °F
Storage and transport temperature	-20 °C to 40 °C/-4 °F to 104 °F
Relative humidity	max. 80 % relative humidity, non-condensing
Long-term storage temperature	-20 °C to 40 °C/-4 °F to 104 °F
Dimensions	approx. 345 mm x d55 mm to d80 mm (arm extended)
Weight without APS gripper	approx. 840 g/29.6 oz
Service life	5 years
Forearm length (as delivered)	185 mm
Maximum adjustable forearm length	302 mm
Minimum adjustable forearm length	190 mm
Flexion angle	15 ° - 145 °
Maximum holding force for forearm length of 305 mm	230 N
Protection class	IP20

9.2 Symbols Used

Symbol	Meaning
	This product must not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions provided by the responsible authorities in your country regarding return and collection processes.
	Declaration of conformity according to the applicable European directives
	Serial number of device
	Lot number
	Legal manufacturer

[illegible]

[illegible]

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.