
Instruction for Use

BALANCE KNEE OFM1/KFM1

4

EN		Instructions for Use	6
DE		Gebrauchsanweisung	14
FR		Notice d'utilisation	22
ES		Instrucciones para el uso	30
IT		Istruzioni per l'uso	38
DA		Brugsanvisning	47
NL		Gebruiksaanwijzing	55
PT		Instruções de Utilização	63
PL		Instrukcja obsługi	72
RU		Инструкция по использованию	80

ALIGNMENT

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


The Balance Knee OFM1 and Balance Knee KFM1 joints are polycentric knee joints with a stable aluminium frame construction. The joints are suitable for users in mobility classes 1-2, and are approved for a maximum body weight of 136 kg.

Due to its unique construction, Balance Knee OFM1 and Balance Knee KFM1 can either be used as a locking joint with manual release or as a mechanical, freely moving knee joint as the user becomes more mobile. The dual function lever enables the user to switch easily between these two functions.

INDICATIONS FOR USE

The Balance Knee OFM1/KFM1 is designed for patients with transfemoral and knee disarticulation amputation.

CONTRAINDICATIONS FOR USE

N/A.

Proximal adaptation can be achieved in a number of ways:

- Balance Knee OFM1 with a pyramid adapter
- Balance Knee OFM1 with an IKF-Adapter
- Balance Knee KFM1 with a loop adapter
- Balance Knee KFM1 with a 3-arm socket adapter (particularly suitable for interim sockets)

For all adaptation options, Balance Knee OFM1/KFM1 is infinitely adjustable by up to 5 mm in the AP/ML directions and can be rotated without limitation (exception: IKF-adapter).

In addition, the flexion/extension angle and/or the adduction/abduction angle of the Balance Knee OFM1 version with a loop adapter, can be corrected by up to max. of 10° by using the adjusting discs. All the alignment adjustments can be carried out during the fitting phase or at any other later point in time.

The Balance Knee OFM1 version with an IKF adapter emulates the anatomical angle of strike of about 5 - 11° without deactivating the safe knee joint geometry. Using the IKF adapter reduces stress on the joint chain and gives the user a secure and energy-saving gait. The adapter also makes walking on gradual slopes and uneven ground considerably easier and safer. A number of different progressive buffers are available that can be varied individually depending on the walking habits of the prosthesis user. The IKF adapter can be used for all mobility classes.

TECHNICAL DATA

Name	Balance Knee OFM1 with a pyramid adapter	Balance Knee OFM1 with an IKF adapter	Balance Knee KFM1 with a loop adapter	Balance Knee KFM1 KP5 with a 3-arm socket adapter
Art. No.	1.721.126	1.721.127	1.721.128	1.721.129
Material	Aluminium + Carbon	Aluminium + Carbon	Aluminium + Carbon	Aluminium + Carbon
Maximum body weight	136 kg	136 kg	136 kg	136 kg
Amputation level	Above-knee	Above-knee	Knee exarticulation, long thigh	Knee exarticulation, long thigh
Total fitted height	approx. 159 mm	approx. 174 mm	approx. 151 mm	approx. 151 mm
Effective fitted height	ca. 107 mm	approx. 122 mm	approx. 116 mm	approx. 116 mm
Fitted height	approx. 18 mm	approx. 33 mm	approx. 27 mm	approx. 27 mm
Gewicht	approx. 640 g	approx. 770 g	approx. 650 g	approx. 735 g
Ground clearance	approx. 13 mm	approx. 13 mm	approx. 13 mm	approx. 13 mm
Prox. connection	Pyramid	Pyramid	Loop adapter	3-arm socket adapter
Swing phase	Integrated extension spring	Integrated extension spring	Integrated extension spring	Integrated extension spring
Stance phase	Optional locking of the axle geometry	Optional locking of the axle geometry Stance phase flexion	Optional locking of the axle geometry	Optional locking of the axle geometry
Axes	Polycentric	Polycentric	Polycentric	Polycentric
Flexion angle without socket	150°	150°	150°	150°

ASSEMBLY OF THE PROXIMAL CONNECTION

Balance Knee OFM1 with a pyramid adapter

- Slide the washers (34 & 35) onto the 40 mm central screw (33), guide the screw through the central opening in the head of the joint (26), and tighten the pyramid adapter (36) to a torque of 35 Nm.
- Secure the 40 mm central screw (33) with Loctite and use a torque wrench (8 mm Allen).

Balance Knee OFM1 with IKF-Adapter

- Slide the washers (34 & 35) onto the 27 mm central screw (31), guide the screw through the central opening in the head of the joint (26) and tighten the IKF adapter (40) to a torque of 35 Nm. Make sure you assemble the IKF adapter (40) in the correct AP direction.
- Secure the 27 mm central screw (31) with Loctite and use a torque wrench (8mm Allen key).

Selecting and inserting the buffer elements

The set is supplied with four buffers

- 3 back (P) progressive buffers (with holes)
- 1 front (A) buffer (without a hole).

Back buffers (with holes)		
Navy	X-Soft	Available as an accessory
Yellow	Soft	Supplied with set
Red	medium	Premounted
Brown	Hard	Supplied with set
Blue	X-hard	Available as an accessory
Front buffers (without hole)		
Brown	Standard	Premounted
Blue	X-hard	Available as an accessory

The hardness of the buffers must be selected and inserted to suit the patient's requirements.

Note: if the user feels unsafe while walking, e.g. has a gentle wobble during the stance phase, use the next hardest buffer.

If the user feels that there is too little damping in the stance phase, use the next softest buffer.

Mounting the buffer elements

The buffer without a hole is the front stop buffer (A).

The back buffer (P) with holes for progressive damping, is mounted in the IKF adapter with the holes facing upwards (**Fig. 2**).

Balance Knee KFM1 with loop adapter

- Fill the central thread of the loop adapter (39) completely with Plastelin, so that no resin can run into the thread during lamination and the central screw can later be completely screwed in. Safety relevant!
- Make up a dummy disc and fix this to the adapter with double-sided adhesive tape.
- Laminate the loop adapter (39) with sufficient carbon fibre and fibre glass to suit the patient's needs.

Important:

Please make sure it binds well in the groove!

- After lamination, grind the distal end of the adapter free and remove the Plastelin.
- Slide the washers (34 & 35) over the 29 mm central screw (32) and guide the screw through the central opening in the head of the joint (26).
- Slide the adjustment discs (38) over the thread of the central screw (32) and tighten the loop adapter (39) to a torque of 35 Nm.
- When mounting the two adjustment discs, make sure that their fine teeth are facing each other on the inner surfaces.
- Secure the 29 mm central screw (32) with Loctite and use a torque wrench (8mm Allen key).

Balance Knee KFM1 with 3-arm socket adapter

- Before adjusting (moulding) the 3-arm socket adapter to the individual socket form, screw the adapter with an external thread (37) into the internal thread of the 3-arm socket adapter. Secure the M5 locking screw of the 3-arm socket adapter to a torque of 10 Nm.
- Mould the 3-arm socket adapter to the individual socket shape.
- Loosen the M5 locking screw on the 3-arm socket adapter and remove the attached adapter.
- Cover the internal thread of the 3-arm socket adapter with grease and screw in the dummy supplied. Secure the M5 locking screw of the 3-arm socket adapter again to a torque of 10 Nm. Fill the slit in the 3-arm socket adapter, the screw head, and any exposed thread with Plastelin.
- Laminate the 3-arm socket adapter with sufficient carbon fibre and fibre glass to suit the patient's needs.
- After lamination, loosen the M5 locking screw of the 3-arm socket adapter and remove the screwed-in dummy and any remaining resin or Plastelin.
- For the final assembly, now completely screw the adapter with an external thread into the internal thread of the 3-arm socket adapter.

Attention: for any later adjustments, do not unscrew the adapter with the external thread (37) more than one full turn.

Secure the M5 locking screw of the 3-arm socket adapter with Loctite and to a torque of 10 Nm. Use a torque wrench (5 mm Allen key).

Slide the washers (34 & 35) onto the 29 mm central screw (32), guide the screw through the central opening in the head of the joint (26), and tighten the adapter (37) to a torque of 35 Nm.

Secure the central screw (32) with Loctite and use a torque wrench (8 mm Allen key).

Attention: never mix up the central screws! The pyramid adapter (36) must be screwed tight with the 40 mm long central screw (33), the loop adapter (39) and the adapter with the external thread (37) must be screwed tight with the 29 mm long central screw (32), and the IKF-Adapter (40) must be screwed tight with the 27 mm long central screw (31). Safety relevant!

OPTIMISING THE PROSTHESIS ALIGNMENT

After the respective central screw has been loosened, unlimited rotation (exception: IKF-adapter) and AP/ML corrections of 5 mm are possible.

Furthermore, in the case of version Balance Knee KFM1 with a loop adapter, alignment corrections in flexion, extension, adduction and abduction of up to max. 10° can be made with the aid of the adjustment discs (38).

Attention: always make sure that the respective central screw (36, 37 or 38) is tightened to a torque of 35 Nm at the dynamic fitting session to prevent any unwanted rotation and thus risk to the user.

ADJUSTMENT

Stance phase stability

The joints are configured as a conventional locking knee at the factory and are delivered in this setting.

Swing phase control when using the joint as a freely mobile joint with integrated extension spring

- Adjusting the extension assist spring:
 - Use a screwdriver to turn the adjusting screw (1) to the right ►►► to increase the action of the extension spring.
 - Use a screwdriver to turn the adjusting screw (1) to the left ◄◄◄ to decrease the action of the extension spring.
- Adjusting axle friction with the help of the friction screw (13)
 - Use a 5 mm Allen key to turn the friction screw (13) to the right ►►► to increase friction at the axle.
 - Use a 5 mm Allen key to turn the friction screw (13) to the left ◄◄◄ to decrease friction at the axle.

CONVERTING THE JOINTS

Locking ►►► mechanical

- Push the release lever (27) upwards and hold it in this position.
- Use a 2.5 mm Allen key to turn the two release screws (24) adjacent to the release lever (27) finger-tight to the right ►►► to secure it in this position.
- Caution: Tighten the release screws (24) to 5 Nm and always secure them with Loctite when finishing the prosthesis. Use a torque wrench.

Attention: Failure to observe point 1 when tightening the 2.5 mm release screws (24) results in loss of function of the joint and/or makes it impossible to release the joint.

Mechanical ►►► Locking

- Conversion to a fixed knee is performed in the reverse order by turning the two 2.5 mm release screws anticlockwise ◄◄◄ until the release lever (27) moves freely.

- Secure the release screws (24) again with Loctite and check that the release lever is freely movable (27).

Attention: After every conversion, carry out a functional check before the dynamic fitting.

Important additional information:

Play may be found when using the optional lock on the Balance Knee OFM1/KFM1. This play can be eliminated with the help of the screw (15) which is turned right ►►► with a 2.5 mm Allen key. This moves the extension stop (20), which compensates for the play. If the release lever (27) does not click into place easily, correct the position of the extension stop with the help of the screw (15), which you turn to the left ◄◄◄.

MOUNTING THE DUAL FUNCTION LEVER

For laminated sockets

- Laminate the threaded couple in the desired position.
- Fix the lamination tube in the laminate.
- Mount the dual function lever on the threaded couple. Secure the screw with Loctite.
- Guide the nylon lanyard through the lamination channel.
- Set the dual function lever in the lower position.
- Snap the release lever (27) into the knee joint.
- Fix the nylon lanyard to the release lever (27) on the joint.
- Secure the retaining screw (28) with Loctite and tighten by hand. Use a 2 mm Allen key.

For plastic sockets

- Drill a hole in the socket in the desired position and fix the threaded couple from the inside of the socket.
- Mount the dual function lever. Secure the screw with Loctite.
- Use the guide clips to determine where the protective duct will run.
- Fix the protective duct in the guide clips.
- Pass the nylon lanyard through the protective duct.
- Set the dual function lever in the lower position.
- Snap the release lever (27) into the knee joint.
- Fix the nylon lanyard to the release lever (27) on the joint.
- Secure the retaining screw (28) with Loctite and tighten by hand. Use a 2 mm Allen key.

Note: If the user only wants to use the dual function lever for releasing the knee joint in certain situations (e.g. sitting down), the enclosed M4 Allen screws can be screwed into the slots in the dual function lever. This stops the lever from slipping up beyond the taper and thus permanently inactivating the release lever (Fig. 3).

SAFETY WARNINGS

- Adapters used at the upper connection must always be tightened to the prescribed torque for the respective modular component. Use a torque wrench.
- The adapters used for the lower connection (e.g. tube adapter) must be inserted down to the end stop in the tube receiver of the knee joint.
- The screw (6) on the tube receiver must be tightened to 16 Nm and secured with Loctite. Use a torque wrench (5 mm Allen key).
- The central screws (31, 32, 33) must be secured with Loctite and tightened to a torque of 35 Nm. Use a torque wrench (8 mm Allen key).
- If the joint is being used with the locking function, the release screws (24) must not impede the function of the release lever (27) and must be secured with Loctite to rule out any inadvertent loosening and incorrect function.
- If the joint is being used with the freely moving function, the release screws (24) must be tightened to 5 Nm and secured with Loctite. Use a torque wrench (2,5 mm Allen key).
- If the joint is being used with the locking function, neither the release lever (27) nor the nylon cord of the dual function lever (41) may be impeded by the cosmesis and they must be able to move freely.
- If readjustment work is required, secure the release lanyard of the dual function lever (41) so that it cannot catch in grinding tools. High potential hazard!
- When readjusting the extension stop (20) with the adjusting screw (15) on the back stub axle, make sure that it functions or locks properly and carry out a functional check before the dynamic fitting.

- Observe the alignment instructions and adjustment notes.
- Do not unscrew any screws other than those described here.
- The use of talc, particularly together with lubricants, can cause the joint to jam. Therefore, please do not use any talc!
- Do not use compressed air for cleaning. We recommend the use of a brush.

When using with the IKF adapter

- If there are any visible changes, wear-and-tear, or functional limitations, replace the IKF adapter or its buffer.
- Please make sure that the patient is competent in the handling of the IKF adapter and his prosthesis before leaving the premises. Failure to observe this warning may cause the patient to fall.
- Any manipulation whatsoever leads to premature wear-and-tear (fracture) of the adapter.
- Make sure you assemble the IKF adapter in the correct AP direction.
- Make sure you mount the buffer elements correctly.

When using with a loop adapter

- To guarantee the safety of the user, make sure you use sufficient carbon/glass lay-up material when laminating the loop adapter (39).
- After lamination, make sure that the thread of the loop adapter (39) is free of Plastiline and resin residues, otherwise the central screw (32) cannot be screwed in completely.
- The fine teeth on the two adjustment discs must face each other.

When using with a 3-arm socket adapter

- The M5 locking screw on the 3-arm socket adapter must be tightened to a torque of 10 Nm and secured with Loctite. Use a torque wrench (5 mm Allen key).
- For readjustments, do not unscrew the adapter with an external thread (37) by more than one turn.

Observe all the above points, otherwise the guarantee becomes null and void!

MAINTENANCE

The joint must be inspected, and repaired if necessary, at least every 6 months!

Inspect:

- the alignment
- the screw connections
- glued joint and torque of the central screw
- the suitability of the patient (e.g. max. body weight, degree of mobility)
- loss of lubricant
- for damage to the joint and the accessories (e.g. the dual function lever)
- for soiling of the bushes
- the lanyard function for locking joints
- AP + ML play
- the extension stop (20)
- visual inspection of the proximal adapter for damage
- for damage to the two adjustment discs (38) when using the loop adapter (39)
- the threaded connection when using the adapter with an external thread (37)
- when using the IKF adapter: check that it functions properly, especially the mounted buffers
- for soiling with talc

CARE

Clean the joint with a soft cloth moistened with a little mild benzene. Do not use any more aggressive cleaning agents because these can damage the seals and bushes.

Do not use compressed air for cleaning! Compressed air can blow existing dirt into seals and sleeves where it causes damage and premature wear-and-tear.

LIABILITY

In case of damage: Össur can only consider complaints accompanied by a copy of the delivery note or the Össur invoice together with a detailed description of the reasons for returning the product. A manufacturer can only be

held liable for the failure of his own modular components. The manufacturer can only be held liable beyond this, if it can be proved that his modular components were causally responsible for the damage to or loss of function of modular components from other manufacturers.

CE MARK, COMPLIANCE

All the modular components have been tested in compliance with EC Directive 93/42/EEC and bear the CE mark.

General conditions

This product is for use by one patient only. If it is used for treating more than one patient, product liability pursuant to the medical Devices Act becomes null and void.

Össur products are guaranteed to be compatible with modular components from others, if the following points are observed:

- The modular components are used only in compliance with their intended purpose (see relevant Instructions for Use).
- If modular components with different max. weight limits are used, the weight limit of the weakest component applies.
- The use of tested individual components with the CE mark does not release the prosthetist from his obligation to check, to the best of his capabilities, the combination of fittings for their suitability, correct assembly and safety.
- Should there be any indications that a given combination of fittings may not assure the required degree of safety, the fittings must not be combined.
- The prostheses must be manufactured in compliance with the generally acknowledged professional rules of the orthopaedic appliances trade.
- If the patient's suitability for wearing the prosthesis changes (e.g. body weight, mobility, etc.), the prosthesis must be checked without delay.
- If the prosthesis has been exposed to unusually high stress, e.g. a fall, it must be inspected in a prosthetic workshop immediately for possible damage.
- Safety-relevant regulations for individual fittings (see relevant Instructions for Use) must be observed

Combination of different modular components:

- The above general conditions apply and additionally:
- Only use modular components that comply with the requirements of ISO 10328
- Only use modular components that comply with the requirements of the medical Devices Act
- The modular components must be suitable for each other and fit the standard pyramid or 30 mm system.

USEFUL LIFE

The maximum useful life of the Össur knee joint Balance Knee OFM1/KFM1 is 36 months.

REPLACEMENT PARTS, ACCESSORIES

No.	Name	Article Number
1	Adjusting screw for extension spring	1.970.818
2	Extension assist spring	1.970.819
6	Screw for tube receiver	1.970.822
13	Friction screw	1.970.827
15	Adjusting screw for extension stop	1.970.830
19	Plate for extension stop	1.970.831
20	Extension stop	1.970.832
24	Release screw	1.970.842
28	Retaining screw	1.970.835
31	Central screw for IKF adapter, 27 mm	M.631.254
32	Central screw for loop adapter (39) and adapter with an external thread (37), 29 mm	1.970.844
33	Central screw for pyramid adapter, 40 mm	1.970.859
34	Washer	1.970.845
35	Washer, flat	1.970.846
36	Pyramid adapter incl. 33, 34, 35	1.970.863
37	Adapter with an external thread M 36 x 1.5 incl. 32, 34, 35	1.633.121
38	Adjustment discs for loop adapter (39)	1.970.847
39	Loop adapter incl. 32, 34, 35, 38	1.970.862
40	IKF adapter incl. 31, 34, 35, buffer set, standard	M.631.100
-	IKF Standard buffer replacement part set	M.631.250
-	IKF buffer X-soft	M.631.251
-	IKF buffer X-hard	M.631.252
-	IKF buffer front X-hard	M.631.253
41	Dual function lever	M.721.222

Only the listed replacement parts and accessories are available.

To have a joint inspected and repaired, please send the joint requiring repair to Össur.

We will provide you with a service joint for the duration of the repair.

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EN – Caution: Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Ossur products are used with other recommended Ossur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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