

MODEL 116

# Baldwin Belt

## LSO PULLEY BELT

The Baldwin Belt is a universal lumbar orthosis which provides sagittal control. The belt features a flexible pulley system for maximum comfort and fit. A double handle pulley system allows for the patient to receive equalized compression throughout the brace.

- One universal size: S-4XL (25-60" circumference)



*Contoured Padded Posterior Panel*



*The Baldwin Belt provides superior 5:1 compression and stabilization*

### FEATURES

- Lightweight low profile rigid anterior panel is built into the belt and provides additional abdominal comfort and support
- Flexible pulley system contours and conforms to the patient's body type which allows the patient to adjust easily for desired comfort when standing or sitting
- Easy to use dual handle pulley helps provide equal compression as well as aids patients with arthritic hand issues and muscle weakness

### INDICATIONS

- Low back pain, sprains and strains, and scoliosis
- Provides relief from mild lower back pain.
- Superior sagittal control and support from L-1 to L-5.
- Chronic & Mechanical low back pain



F 800-521-2192 | 248-588-7480  
P 800-923-2537 | 248-588-2960  
**BeckerOrthopedic.com** | f t in


©2020 Becker Orthopedic Appliance Co.  
All rights reserved.  
Revision 11/30/20

EC REP  
Acorn Regulatory Consultancy Services Limited  
Knockmorris Cahir Co. Tipperary Ireland, Postcode: E21 R766  
P 012 4626 8456  
F 012 4626 8648



# Fitting Instructions


## 1 | LOCATE SIZE



## 2 | CUT TO SIZE




Lift sleeve and remove extension panels to adjust to appropriate size.



Cut to desired length with scissors (prior to cutting, refer to sizing chart for proper length).




If size adjustment is minimal, simply fold the extension panel to desired length.



Once adjustment is made, (cut or folded) simply lift the sleeve and attach the extension to Velcro. Press down to secure.

## 3 | FIT TO SIZE



Initially align each pulley handle to the side of the brace prior to fitting. Adjust tension for sitting or standing.



Center back panel on the spine.



Wrap the left side to the front, positioning the front panel. Then wrap right side on top of the left panel to close and attach in front.



Pull each side pulley at the same time out and away to desired compression and support.



Place the handles on the fabric for a secure fit.

## PULLEY ADJUSTMENT

To Shorten the Length of the Nylon Cord on the Pulley Handle:



Detach the cord from the front clip



Wrap the cord around the 360° side groove clockwise



Reattach the cord to the front clip