Bio Leg

Instructions for Use (Qualified Personnel)

Ver. 1.10 last update 2024-8-13

Table of Contents

1	Introduction.....	4
2	Design.....	4
3	Function.....	5
4	Application.....	6
4.1	Indications for use	6
4.2	Conditions of use	6
4.3	Indications	6
5	Safety Precautions	7
5.1	Contraindications.....	8
5.2	Warning.....	8
5.3	Caution	10
5.4	Notice.....	13
6	Preparing the Product for Use.....	13
6.1	Scope of delivery.....	13
6.2	Alignment.....	14
6.2.1	Bench alignment (sagittal plane).....	14
6.2.2	Bench alignment (front face value).....	14
6.2.3	Confirmation at maximum flexion	15
6.2.4	Dynamic alignment.....	15
6.3	Charging Bio Leg	16
6.3.1	Charging the knee joint battery.....	16
6.3.2	Display of the current charge level.....	17
6.4	About outer covers	18
6.4.1	How to install the outer cover	18
7	Setting with the “MoPro” Adjustment Software.....	19
7.1	Operating environment.....	19
7.2	Preparing the product to connect to the MoPro	19

7.3	Initial launch of adjustment software (MoPro)	19
8	Use	20
8.1	Switching on the product	20
8.1.1	Night mode	21
8.1.2	Automatic power off function	21
8.1.2.1	Power off prevention function	21
8.1.2.2	Seated power saving function	21
8.2	Movement patterns	22
8.2.1	Standing	22
8.2.1.1	Standing flexion resistance (yielding)	22
8.2.1.2	Standing assistance (rest)	22
8.2.2	Walking	23
8.2.2.1	Stance phase	23
8.2.2.2	Swing phase	24
8.2.3	Sitting/standing up	24
8.2.3.1	Sitting	25
8.2.3.2	Standing up	26
8.2.4	Ascending and descending the stairs	27
8.2.4.1	Walking up stairs	27
8.2.4.2	Walking down stairs	29
8.2.5	Ramp	30
8.2.5.1	Walking up a ramp	30
8.2.5.2	Walking down a ramp	31
8.3	Adjustment of prosthetic knee parameters	32
9	Storage	33
10	Cleaning	33
11	Product Warranty and Maintenance	33
11.1	Additional services	33
11.2	Warranty details	34
11.2.1	Maintenance	34
11.2.2	Scheduled maintenance	34
11.3	Product structural strength	34

12	Disposal	35
13	Legal Information	35
13.1	Liability	35
13.2	Local legal information	35
14	Technical Data	36
14.1	Bio Leg.....	36
14.2	AC Adapter (Battery charger)	37
15	Environmental Conditions.....	37
16	Electromagnetic Compatibility.....	38
16.1	Warnings for Electromagnetic Compatibility	38
16.2	Electromagnetic Emissions Compliance Information	40
16.3	Electromagnetic Immunity Compliance Information.....	40
16.4	RF wireless communications Compliance Information.....	41
16.5	Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields.....	41
17	Frequently Asked Questions.....	42
18	Signals	43
18.1	Signal for operating states	43
18.2	Charging status signals	43
18.3	Warning/error signals.....	44
18.3.1	Warning signals in use	44
18.3.2	Error signals in use	45
18.3.3	Error signals during charging	46
19	Symbols	47
20	Trademark.....	48
21	Software License.....	49
22	Link	50

1 Introduction

- Please read this manual carefully before using this product.
- Instruct the user of the safe product usage.
- Please contact the manufacturer if you have any questions about the product or in case of problems.
- Report each serious product-related incident to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state of the user.
- Please keep this manual in a safe place for easy access.

2 Design

The product consists of the following components

- ① Power Button
- ② Top Outer Cove

- ③ Proximal Pyramid Adapte
- ④ Back Button
- ⑤ Charging Socket

3 Function

This product uses motors and sensors to control the movement of the knee joint. The build-in microprocessor analyzes the user's movements based on sensor measurements and assists the knee joint in flexion and extension at appropriate timings.

It provides safer and more comfortable walking in a variety of situations.

In addition, the Bio Leg exclusive adjustment software settings allow each user to walk in a way that is more natural and comfortable.

This product is intended to assist the ability of the prosthetic knee user and is not intended to enable the user to easily climb stairs or stand up from a chair without training.

4 Application

4.1 Indications for use

This product is to be used solely for lower limb exoprosthetic fittings.

4.2 Conditions of use

- This product was developed for everyday use and must not be used for any other purpose.

Do not use the product in strenuous exercise or other activities that place an excessive load on the product.

- The product is intended for the exclusive use by an individual user.

Use of the product by users other than the intended user is not approved by the manufacturer.

4.3 Indications

- For users with knee disarticulation, transfemoral amputation, or hip disarticulation
- For unilateral or bilateral amputation
- For users with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation

To use this product, the user must fulfill the physical and mental requirements for perceiving visual/acoustic signals and/or mechanical vibration signals.

5 Safety Precautions

Please read this "Safety Precautions" carefully before to use the product correctly.

The precautions shown here are intended to ensure the correct use of this product and to prevent hazards or damages to the patient and other people.

The degree of danger or damage that may occur if the product is used incorrectly and the contents of the display are ignored is classified into "Warning," "Caution," and "Notice" and explained.

All of this important information must be followed.

WARNING

Warning regarding possible serious risks of accident or injury.

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

Particularly important supplementary information and other information.

5.1 Contraindications.

- This product should not be used by individuals weighing more than 125 kg. (275 pounds)
- This product is intended only for patients who have undergone a thigh amputation or knee disarticulation.
- This product must not be used in conjunction with osseointegrated implant systems.

5.2 Warning.

WARNING

Adherence to Safety Instructions

Please adhere to the safety notes specified in this manual. Improper use of the product may result in injury to the user or damage to the product.

- > Follow the safety and handling instructions outlined in this manual.

WARNING

Cease Use if Damaged

If you notice any damage or defects in the product or accessories, such as unusual noises, rattling, or visible damage, stop using the product immediately. Continued use in these conditions may lead to further damage and possibly injury to the user.

- > Contact the manufacturer for assistance.

WARNING

Do Not Disassemble or Modify

Do not disassemble or modify the product or accessories, including the AC adapter.

Doing so may cause damage or malfunction, which could lead to injury.

**WARNING****Prevent Entanglement During Use**

During the movement of the device's knee joint, there is a risk of injury from getting fingers or hands caught in the mechanism.

➤ Please keep your fingers and hands away from the moving parts of the device.

**WARNING****Power Off During Adjustment**

Always turn off the power when making alignment adjustments to the product. If the product moves unexpectedly during adjustment, it could cause injury.

**WARNING**

Do not use in areas where flammable anesthetic gases or combustible gases are present, or in high-concentration oxygen environments.

This may cause explosions or fires.

**WARNING****Risks Associated with Ingress of Dust and Moisture into the Device**

The device may malfunction unexpectedly due to faults, potentially causing injury to the user through falls or other accidents. Ingress of particles and liquids can lead to electric leakage or fire hazards.

- Ensure that dust and liquids do not enter the device.
- Take care to prevent splashes of water from contacting the product.
- On rainy days, protect the device from direct exposure to moisture by covering it with clothing or a cover.
- If you suspect that dust or liquid has entered the device, stop using it immediately. Wipe the device with a dry cloth and allow it to air dry thoroughly. Regardless of whether there are any operational abnormalities, always contact the relevant prosthetic manufacturing facility and have the product inspected by our company.

5.3 Caution

CAUTION

Use the product only in the environment specified in this manual

If used outside the specified environment, the device may be damaged, and users may be at risk of injury.

> Please refer to “Technical Data” and “Electromagnetic Compatibility” in this document.

CAUTION

Do Not Use with Other Medical Devices

This product should not be used in conjunction with other medical devices.

Using it with other devices could damage both the product and the other devices, compromising safety.

CAUTION

Do not use in locations with other medical equipment such as electric surgical units (electrosurgical units), MRI scanners, or CT scanners.

This product should not be used in conjunction with other medical devices.

Using it with other devices could damage both the product and the other devices, compromising safety.

CAUTION

Attention to Warnings/Messages

Failure to notice and respond to warnings or messages from the product could lead to falls and injuries. The product might perform unexpectedly due to changes in assistance force.

> Follow the procedures specified in this manual whenever warnings or messages are issued.

CAUTION

Always remove the device from the body before charging.

Charging while attached to the body may result in electrocution.

CAUTION

Use Specified Charger Only

Do not charge the device with any AC adapter or power cable not specified by the manufacturer. Using an unsuitable charger could damage the device.

CAUTION

Software Operation in Specified Environment Only

Run the adjustment software only in the operating environments specified.

Installing the software in an unspecified environment may lead to operational issues.

CAUTION

Assembly by Certified Professionals Only

The assembly of the device and any components should only be performed by certified prosthetists as designated by our company. Improper assembly (such as loose screws or misalignment) could lead to malfunctions and injuries during use, as well as damage to the device.

CAUTION

Adjustments by Certified Professionals Only

Adjustments to the prosthetic (such as alignment and control parameter settings) should only be made by certified prosthetists authorized by our company. Incorrect settings could cause unexpected movements of the product, potentially resulting in falls and injuries.

➤ Ensure settings are adjusted according to the procedures outlined in this manual.

**CAUTION****Caution When Using on Stairs**

When using the device on stairs, ensure that it is in the appropriate mode to prevent falls. Always use handrails and other fall-prevention measures when using the device on stairs.

➤ Follow the instructions in this manual to switch to the appropriate mode before using the device on stairs. Use the device on stairs only after proper training.

**CAUTION****Avoid Continuous Use in Intense Activities**

Do not use the device continuously in intense activities as it may overheat and damage the device.

**CAUTION****Keep the Device Clean**

Keep the device clean. Dirt and debris can cause corrosion and damage to the device.

➤ Clean the device before storage as outlined in the 'Device Cleaning' section of this manual.

**CAUTION****Storage/Transport Conditions**

Store and transport the device only in the environments specified. Improper storage or transport conditions can damage the device, and in the case of Li-ion batteries, can lead to deterioration, fire, or leakage.

**CAUTION****Avoid Mechanical Stress During Storage/Transport**

Do not expose the device or accessories to unnecessary mechanical stress during storage or transport, as this could lead to damage.

5.4 Notice

NOTICE

Even when the device is not in use, regular charging is necessary

The device contains a Li-ion battery, which may deteriorate if not charged for an extended period.

6 Preparing the Product for Use


6.1 Scope of delivery

Name	pc
Bio Leg (SN007) Knee joint (with pyramid adapter)	1
Top Outer Cover (Carbon Black)	1
Lower Outer Cover (Carbon Black)	1
Instruction for Use (this manual: for prosthetists and prosthetic fabrication facilities)	1

6.2 Alignment

Only when biomechanically aligned correctly can the knee joint function safely.
Please note that bench alignment should be performed with the power turned off.

6.2.1 Bench alignment (sagittal plane)

Align the distance so that the alignment reference line passes 20 mm forward from the alignment reference point (knee axis).
Sockets should be aligned to the socket style and feet to the manufacturer's recommended alignment.

6.2.2 Bench alignment (front face value)

The distal pyramid should be assembled in the middle position.

6.2.3 Confirmation at maximum flexion

Make sure that the socket or part does not contact the knee joint when it is maximally flexed.

6.2.4 Dynamic alignment

After adjustment using the adjustment software, please optimize the dynamic alignment during trial walking.

If necessary, take the following into account when making adaptations:

- Socket flexion angle (sagittal plane) with confirmation of stride symmetry
- Socket internal rotation angle and M-L position setting of the socket adapter (forehead)
- Rotational position of the axis of the knee joint and external rotation of the foot (horizontal plane)

6.3 Charging Bio Leg

A fully charged battery can be used for approximately 22 hours of average activity. When charging observe the following

- For constant use of the prosthetic knee, it is recommended to charge the battery at least once a day.
- Do not disconnect the power cord from the product until just before using the product to allow it to operate for the maximum time on a single charge.
- Charge the battery for at least 4 hours before using it for the first time.

Charge the battery until the LED on the back button of the unit turns blue and the unit vibrates and beeps.

6.3.1 Charging the knee joint battery

- 1) Connect the 3-pin plug of the power cable to the AC adapter.
- 2) Plug the AC adapter into an outlet.
- 3) Connect the charging plug of the power cable to the charging socket of the knee joint.
- 4) Start charging.
- 5) When charging is complete, a vibration and beep are emitted and the LED on the back button illuminates.
- 6) When charging is complete, disconnect the power cable from the knee joint.

When charging is correctly initiated, a vibration and beep are emitted and the LED on the back button flashes.

6.3.2 Display of the current charge level

By pressing the rear button, the charge level can be checked by the color of the rear LED, vibration and beep sound.


Battery charge status	Back LED Color	Beep signal	Vibration Signal
61% or more	Blue	4 times shorter	4 times shorter
41% – 60	Green	3 times shorter	3 times shorter
21% – 40	Yellow	2 times shorter	2 times shorter
20% or less	Violet	1 time shorter	1 time shorter

The charge level can also be checked from the adjustment software.

6.4 About outer covers

The top outer cover should always be attached, even when test walking.

Select the size of the lower outer cover according to the length of the connecting pipe to the foot section and make sure it does not interfere with the foot section.

6.4.1 How to install the outer cover

- 1) Attach the top outer cover by inserting it from the bottom and spreading it to the left and right.
- 2) Tighten the screws on the left and right sides (red arrows in the figure).
- 3) Attach the lower outer cover by inserting it from the front and spreading it to the left and right.
- 4) Tighten the screws on the left and right to secure the top and lower outer covers together (green arrow sections in the figure).

7 Setting with the “MoPro” Adjustment Software

An adjustment software (“MoPro”) can be used to change various parameters, such as assist force, within certain limits (e.g., when becoming accustomed to using the knee joint).

In addition, battery charge level and other information can be checked.

The app can be downloaded from the following URL or QR code.

<https://apps.apple.com/jp/app/mopro-bionicm/id1602657293>

7.1 Operating environment

MoPro has been tested on the following operating systems:

- iOS13 or later
- iPad OS13 or later
- Bluetooth 4.0 or later

7.2 Preparing the product to connect to the MoPro

Before connecting the MoPro to the prosthesis, make sure that

- Bluetooth on the mobile device must be activated.
- The mobile device must not be in “Airplane Mode” (offline mode).
- The mobile device must be connected to the Internet.
- Prepare the serial number and pairing password of the prosthetic knee to be connected.

The serial number is affixed as a sticker on the bottom of the product.

It is the six-digit number indicated under the word “SN” Pairing passwords can be found on our website.

7.3 Initial launch of adjustment software (MoPro)

To use MoPro, you need to agree to the Terms of Use and have Bluetooth connection, etc.

Please refer to the MoPro instruction manual for details.

8.1 Switching on the product

Switch Off

Switch On

The power button is located on the top of the knee joint. Press the power button and the power LED glows green to turn the power on. Press the power button again, and when the power LED turns off, the power is off.

- When fitting a prosthetic knee, make sure that the power is turned off and the prosthesis is properly fitted before turning on the power.
- When removing the prosthetic knee, make sure the power is turned off.
- The power-off prevention function automatically turns off the power after a certain period without the prosthetic leg attached.

8.1.1 Night mode

When the night mode is set to "on" in MoPro, the power LED and rear LED turn on/off settings will be changed as follows during the period from midnight to 5:00 a.m. The default setting of the night mode is "off".

		Night mode: On	Night mode: Off
Normal operation	Power LED	Lights off 10 seconds after power on	Always on
Charging	Power LED	Lights off 10 seconds after start of charging	Always on
	Back LED	Lights off 10 seconds after start of charging	Always blinking

Even with the night mode on, the rear LED will flash while the rear button is pressed while charging, the remaining battery level is indicated by LED color, and the power LED also lights up for 10 seconds.

8.1.2 Automatic power off function

In the following two situations, the system is equipped with a function that reduces battery consumption by automatically turning off the power.

8.1.2.1 Power off prevention function

After 90 minutes with the prosthetic knee removed and the knee joint is fully extended, the power is automatically turned off. This function is designed to prevent battery consumption due to forgetting to turn off the power when the prosthetic knee is not worn, such as while sleeping.

This feature is turned on by default and cannot be turned off.


8.1.2.2 Seated power saving function

After a certain period in the sitting position, the power turns off automatically. This function is designed to prevent battery consumption due to standby mode when working while seated for long periods of time. This function is initially turned off and can be turned on with MoPro. The time until the power is automatically turned off can also be set in MoPro.

8.2 Movement patterns

8.2.1 Standing

8.2.1.1 Standing flexion resistance (yielding)

It prevents sudden knee breakage by creating resistance to knee flexion.

This flexion resistance also works during stair descent and sitting movements to support stable and smooth movement.

8.2.1.2 Standing assistance (rest)

After standing still with the knee joint in slight flexion for 1.5 seconds, the knee extension force begins to function, and the user enters the standing assistance mode. When the prosthetic foot is lifted off the floor, the standing assistance mode is deactivated.

With the prosthetic knee slightly bent, weight can be placed on the prosthetic side, allowing the user to take a resting posture or lean against a wall with ease.

The knee extension force, the number of seconds of rest before the rest function starts, and the range of the knee flexion angle when the function starts can be adjusted in the "Rest" setting in MoPro.

8.2.2 Walking

The motor assists in smoothly swinging the next step from the stable stance phase to the swing phase smoothly.

To switch to the swing phase, step forward from full foot contact.

When test walking for the first time after wearing a prosthetic knee, be sure to perform the prescribed training.

8.2.2.1 Stance phase

With a certain amount of load applied to the knee joint, standing flexion resistance occurs.

This reduces the risk of knee breakage during heel ground contact.

The strength of flexion resistance can be adjusted.

8.2.2.2 Swing phase

In the swing phase, the resistance and assist force for each of flexion and extension can be adjusted to allow smooth swinging of the prosthetic leg. It can prevent excessive knee bending and assist with prosthetic leg extension while following walking speed.

Assist force can be adjusted.
The stronger the assist force, the faster the prosthesis will extend. However, depending on the degree of adjustment, the impact at full extension may be stronger, so adjust with caution.

8.2.3 Sitting/standing up

It controls knee flexion resistance and assists in knee extension during sitting and standing up activities in a chair.

8.2.3.1 Sitting

During the sitting motion, the flexion resistance of the knee joint as the flexion resistance of the standing position allows the left and right knee to bend equally.

Sitting in a chair, with the knee joint bent about 90° and no load applied to the prosthetic foot for 5 seconds, the user enters the sitting mode.

Once in the sitting mode, the resistance of the knee joint to both flexion and extension is adjusted to facilitate movement of the prosthetic knee.

remarks

- In the sitting mode, resistance acts in the flexion direction. By moving in the extension direction, the resistance in the flexion direction is reduced. After that, if the knee joint angle remains unchanged for 2 seconds, flexion resistance is applied again.
- The sit-to-stand power saving function automatically turns off the power after a certain period in a sitting position (the sit-to-stand power saving function is turned off by default). When the power is turned off by the sit-down power saving function, the power switch must be pressed again to turn the power on when using the stand-up assist function or when starting to walk.
- Use of the product with the power turned off may cause the user to fall over.
- The on/off of the sit-down power saving function and the time until the sit-down power saving function is activated can be set in MoPro.

8.2.3.2 Standing up

During the standing up motion, an assist in the direction of extension is applied to the prosthetic knee, which enters a standing up mode to lift the body.

- 1)Open both feet to an equal width.
- 2)Tilt the upper body slightly forward.
- 3)When the weight is applied equally to both legs with the prosthetic knee weighted, the knee assists in the direction of extension.
- 4)Make sure that the knee is fully extended before moving on to the walking motion.

Support the elbow rest or other support with hands and try to stand up.


Practice enough to allow full extension of the knee.

8.2.4 Ascending and descending the stairs

Assists smooth ascent and descent of the stairs by controlling knee flexion resistance and assisting knee extension.

8.2.4.1 Walking up stairs

- 1) In standing or walking mode, place the non-disconnected side on the first step of the staircase.
- 2) With the knee joint extended, extend the knee on the non-cutting side and lift the knee joint.
- 3) When the thigh on the prosthetic side is raised to place the knee joint on the top step, the knee joint flexes to avoid hitting the stairs.
- 4) When the prosthetic knee is placed on the upper step and weight is placed on the prosthetic knee, the knee joint of is extended to assist in the stair climbing motion.
- 5) When the patient finishes climbing the stairs, the prosthetic knee is swung forward as a walking motion to release the stair ascending mode.

Remarks


Assist force and knee joint flexion angle can be adjusted in the MoPro. Instruct the user to deactivate the stair ascending assist function when finished climbing stairs.

The two ways to unlock the system

- When the prosthetic foot is lifted off the floor, the knee joint is intentionally bent to enter the swing phase of gait (when released, the knee joint vibrates to indicate to the user that it has been released).
- Assume a standing position with both legs extended and wait for 5 seconds.

Walking without deactivating the stair ascending assist function may cause the knee to bend significantly, delaying the lower leg from coming forward and causing the user to fall. Please have the user practice deactivating the stair ascending assist function.

8.2.4.2 Walking down stairs

- 1) In standing or walking mode, place the non-disconnected side on the first step of the staircase.
- 2) Using the flexion resistance of the prosthetic knee, the knee joint is bent with the weight on the prosthetic knee side.
- 3) Once the non-cutting side is grounded on the bottom step, lift the prosthetic side and place it on the next step.
- 4) Just before the prosthetic leg touches the ground, the knee joint stops in a slightly flexed position.
- 5) After completing the stair descent, the prosthetic knee is swung forward as a walking motion to release the stair descent mode.


This action should be practiced thoroughly and with care. Only when the sole is properly grounded does the system react correctly and the flexion movement is controlled.

<div data-bbox="137 387 238 488" data-label="Image"> </div>	<div data-bbox="297 108 395 135" data-label="Section-Header"> <p>Remarks</p> </div> <div data-bbox="297 150 1001 296" data-label="Text"> <p>The strength of the flexion resistance and the angle at which the knee joint stops just before ground contact can be adjusted in the MoPro. Instruct the user to deactivate the stair descent assist function when finished descending stairs.</p> </div> <div data-bbox="297 311 693 338" data-label="Text"> <p>The two ways to unlock the system</p> </div> <div data-bbox="297 352 981 579" data-label="List-Group"> <ul style="list-style-type: none"> • When the prosthetic foot is lifted off the floor, the knee joint is intentionally bent to enter the swing phase of gait (when released, the knee joint vibrates to indicate to the user that it has been released). • Assume a standing position with both legs extended and wait for 5 seconds. </div> <div data-bbox="297 593 986 699" data-label="Text"> <p>Walking without the stair descent assist function deactivated may cause a delay in the lower leg coming forward, which may cause the user to fall.</p> </div> <div data-bbox="297 713 925 778" data-label="Text"> <p>Have the user practice disabling the stair descent assist function.</p> </div>
---	--

8.2.5 Ramp

8.2.5.1 Walking up a ramp


When swinging the prosthetic knee out, the extension of the knee joint is controlled, making it easier to get the prosthetic leg up the hill.

8.2.5.2 Walking down a ramp

When the prosthetic knee side touches the ground and the knee joint breaks, flexion resistance is increased, flexion of the knee joint is controlled, and the body's center of gravity is low and stable.

The swing phase is not initiated even if the knee joint is flexed.

Remarks

If the angle of the slope is steep, it may switch to stair descent mode.

If the mode is switched to stair descent mode, please deactivate the stair assist function and walk.

8.2.4 Walking up the stairs.

8.3 Adjustment of prosthetic knee parameters

MoPro can be used to adjust parameters such as assist force for each operation/function.

For details, please refer to the help in MoPro and the MoPro instruction manual.

The app can be downloaded from the following URL or QR code.

<https://apps.apple.com/jp/app/mopro-bionicism/id1602657293>

9 Storage

- If not used for a long period of time, store away from high temperature and humidity.
- Products should be recharged periodically, even when not in use for long periods of time.

10 Cleaning

- Clean the product with a cloth dampened with fresh water when needed.
- Dry the product with a lint-free cloth such as a non-woven cloth and allow it to air-dry fully.

11 Product Warranty and Maintenance

Regular check at 12-month intervals and regular maintenance at 24-month intervals are mandatory in the interest of user safety, maintaining operating reliability, protecting the warranty, ensuring basic safety and the essential performance characteristics, as well as reinforcing EMC safety.

The manufacturer grants a grace period of no more than one month before, or two months after, the due date.

11.1 Additional services

During maintenance, additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

- Bio Leg main unit

11.2 Warranty details

- 3-year
- Performing repairs
- Periodic maintenance

The warranty may not cover the product if it is damaged or broken due to surface scratches, negligent or improper use, intentional or careless handling, or rough handling.

11.2.1 Maintenance

To maintain product performance and ensure safe use, please have the specified periodic maintenance performed.

Regular maintenance includes inspection of sensors and replacement of worn parts.

11.2.2 Scheduled maintenance

- The 3-year warranty includes regular free check at 12 months and regular free maintenance at 24 months.

Failure to perform scheduled maintenance will void all warranties.

11.3 Product structural strength

This product has a strength that can withstand 3 million walking cycles, equivalent to a 125kg wearer for three years. Deterioration or damage caused by long-term use is not covered by the warranty.

*Tested in accordance with ISO10328-P6-125kg

12 Disposal

Some jurisdictions do not permit the disposal of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your local government may have a detrimental impact on health and the environment. Please observe the instructions of your local authority pertaining to return and collection.

Note that a Li-ion battery is used in this product.

13 Legal Information

All legal terms are subject to and may vary according to the national laws of the country of use.

13.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

13.2 Local legal information

All legal requirements are subject to the national laws of your country and may vary.

14 Technical Data

14.1 Bio Leg

Product name	Robotic Prosthetic Knee
Reference number	Bio Leg SN007.
Connection (proximal)	pyramid
Dimensions	Overall length: 282.3mm (11.1 in) Overall length (effective dimension): 254mm (10.0 in) To outer cover (bottom) (M): 358.1mm (14.1 in) To outer cover (bottom) (L): 111.8mm (16.2 in) Width: 111.8mm (4.4 in)
Weight	3kg (6.6 lbs)
Max. body weight	Max. 125kg (275 lbs)
Maximum possible flexion angle	132°
Operating temperature	-10°C/+14°F to +40°C/+104°F If the device is stored below 0°C /+32°F for an extended time, device will not perform as expected.
Internal power supply	Li-ion battery
Product temperature for charging	0°C/+32°F to 40°C/+104°F *Charging at between 10°C /+50°Fto 40°C /+104°F is recommended.

Please note that the specifications are subject to change without notice in order to improve product performance.

14.2 AC Adapter (Battery charger)

Model	PW-M065A2-1Y190A
Input	AC100-240V, 2-1A, 50-60Hz
Output	DC19.0V 3.24V 65.0W Max
Mark	IEC 60601-1

15 Environmental Conditions

The device is not waterproof against rain or submersion.

If wet, wipe the knee joint and parts with a cloth and dry thoroughly and naturally.

No salt water or chlorinated water allowed.

If the device has been submitted to these conditions, it must be sent for service to BionicsM

	Use	Shipping/Storage
Temperature	-10°C ~ 40°C	-20°C ~ 60°C
Relative humidity	30% ~ 75% non-condensing	10% ~ 100% non-condensing
Atmospheric pressure	70kPa ~ 106kPa	50kPa ~ 106kPa

16 Electromagnetic Compatibility

16.1 Warnings for Electromagnetic Compatibility

This product is a device that complies with the EMC (Electromagnetic Compatibility) standard IEC 60601-1-2:2014+A1:2020, which is required for the safe use of medical electrical equipment. The EMC standard is a regulation that ensures the safe use of medical electrical equipment by keeping the noise generated by the equipment from affecting other devices and by minimizing the influence of electromagnetic waves emitted by other devices to a certain level. As detailed information regarding the EMC environment necessary for the safe functioning of the equipment is required to be provided to users according to IEC 60601-1-2:2014+A1:2020 (Section 5.2.1.1), we provide the technical information below.

◆ What is EMC (Electromagnetic Compatibility)?

EMC (Electromagnetic Compatibility) is the ability to meet the following two requirements:

- Do not emit noise that causes unacceptable interference to other electronic devices in the vicinity. (Emission)
- Withstand the electromagnetic environment, such as noise emitted by other electronic devices in the vicinity, and perform its functions normally. (Immunity)

◆ Technical explanation related to EMC (Electromagnetic Compatibility)

Medical electrical equipment requires special attention concerning EMC, and failure to follow the EMC information provided may result in serious accidents.

- This product requires special attention regarding electromagnetic compatibility (EMC) and must be used according to the EMC information provided in the EMC technical documentation.
- This product may be affected by portable and mobile RF communication equipment.
- Do not use this product in close proximity to or stacked with other equipment (except during communication).

**WARNING**

Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

**WARNING**

Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

**WARNING**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

16.2 Electromagnetic Emissions Compliance Information

Emissions test	Standard	Compliance level
Conducted and radiated RF emissions	CISPR11	Group 1 – Class B

16.3 Electromagnetic Immunity Compliance Information

Immunity Test	Standard	Compliance Level
Electrostatic discharge	IEC 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM fields	IEC 61000-4-3	10V/m 80MHz~2.7GHz 1kHz 80%AM
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 16.4
Electrical fast transient/burst	IEC 61000-4-4	±2kV for power supply lines 100kHz
Surge	IEC 61000-4-5	±0.5kV, ±1kV line(s) to line(s) ±0.5kV, ±1kV, ±2kV line(s) to earth
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V rms 0,15MHz – 80MHz 6 V in ISM and amateur bands between 0,15MHz and 80MHz 80 % AM at 1 kHz
Power frequency magnetic field	IEC 61000-4-8	30A/m 50Hz and 60Hz
Voltage dips	IEC 61000-4-11	0%U _T :0.5 cycle 0°, 40°, 90°, 135°, 180° ,255°, 270° and 315° 0%U _T :1 cycle and 70%U _T :25 cycle Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0%U _T :250 cycle
Proximity magnetic fields	IEC 61000-4-39	See table 16.5
Note: UT is the A.C. mains voltage prior to application of the test level.		

16.4 RF wireless communications Compliance Information

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380~390	TETRA400	Pulse modulation 18Hz	1.8	0.3	27
450	430~470	GMRS460, FRS460	FM ± 5 kHz deviation 1kHz sine	2	0.3	28
710	704~787	LTE Band 13,17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800~960	GSM800/900 TETRA800, iDEN820, CDMA850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700~1990	GSM1800; CDMA1900; GSM1900; DECT; LET Band 1,3, 4,25;UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						
2450	2400~2570	Bluetooth, WLAN, 802.11b/g/n, RFID2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100~5800	WLAN802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

16.5 Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields

Test frequency	Modulation	IMMUNITY TEST LEVEL(A/m)
30kHz	CW	8
134.2kHz	Pulse modulation 2.1kHz	65
13.56MHz	Pulse modulation 50kHz	7.5

17 Frequently Asked Questions

Q: Is the product waterproof?	A: The product is not waterproof against rain or submersion, so please handle with care. If the product gets wet, wipe the product with a cloth and dry it thoroughly and naturally. If there is any abnormality after drying, please contact us.
Q: How many hours does it take to fully charge the battery?	Depending on the remaining battery level, ambient temperature, etc., the battery will be fully charged in about 3 hours.
Q: How long will the battery last?	When fully charged, it can be used continuously for approximately 22 hours, or 1 day of average use. If the prosthetic knee is to be used constantly, it is recommended that it be recharged once a day.

Please also check the “Frequently Asked Questions” on our website.

18 Signals

18.1 Signal for operating states

Beeps and vibrations notify the patient when switching between walking and stair ascending/descending.

Signal notifications regarding operating status can be toggled on and off in MoPro.

Event	Beep Signal	Vibration Signal
Operation transitions other than the following	Once	Once
Standing → Ascending stairs	Twice	Twice
Walking (small steps) → Standing	Twice	Twice
Ascending stairs → Walking	None	Once
Descending stairs → Walking	None	Once

18.2 Charging status signals

Events	Back LED Color	Beep	Vibration Signal
Start charging	Flashes in colors corresponding to charge level	Twice	Twice
Charging complete	Lights up in different colors according to charge level	Twice	Twice

Battery charge status

How to check the charge level.

18.3 Warning/error signals

18.3.1 Warning signals in use

Event	Back LED Color	Beep Signal	Vibration Signal	Required Action
Low temperature conditions (Warning No. "0028 0040")	Flashing yellow 5 times	5 times	5 times	The temperature of the product is low. Move the product to a warmer environment such as indoors. Although the product can be used continuously, if the temperature is outside the permissible range, the joints will be put on brake and the motion assist function will stop.
High temperature conditions (Warning No. "0028 0080" or "0032 0001")	Flashing yellow 5 times	5 times	5 times	The product temperature is high. Please refrain from activity. You can continue to use the product, but if the temperature is outside the acceptable range, the joints will be put on brake and the motion assist function will stop.

18.3.2 Error signals in use

Event	Back LED Color	Beep Signal	Vibration Signal	Required Action
Hypercooled state (Error No. "0028 0400")	Red blinking 10 times	10 times	10 times	The temperature of the product is out of the acceptable range. Move the product to a safe place (where you can sit) and a warm environment such as a room, then turn off the power and discontinue use until the temperature is within the permissible range. *If the temperature is not back within the allowable range, an error will be notified after restarting.
Excessively high temperature conditions (Error No. "0028 0800" or "0032 0200")	Red blinking 10 times	10 times	10 times	The product temperature is out of the acceptable range. Relocate to a place safe for sitting down, then turn off the power and discontinue use until the temperature is within the permissible range. *If the temperature is not back within the allowable range, an error will be notified after restarting.
Possible malfunction (Error No.: other than above)	Red blinking 10 times	10 times	10 times	There is a possibility of malfunction. Relocate to a safe place for sitting down, then turn off the power and restart it. If the error message appears repeatedly after restarting, discontinue use and contact the manufacturer for support.

When an error occurs, the error code is displayed on the MoPro.

18.3.3 Error signals during charging

Event	Back LED Color	Beep Signal	Vibration Signal	Required Action
Charging error	Red blinking 10 times	10 times	10 times	Charging error. Turn off the power and restart it. Also, check the connections of the power cable and AC adapter, and that they are properly connected to the product and the power outlet. If the error repeatedly appears, contact the manufacturer for support.

Consult instructions for use

Manufacturer

Type B Applied Part

Non-ionizing electromagnetic radiation

Device contains electronic components and/or batteries that should not be disposed of in regular waste

Indicates the temperature limits to which the equipment can be safely exposed.

Indicates the range of humidity to which the equipment can be safely exposed.

Indicates the range of atmospheric pressures to which the equipment can be safely exposed.

Date of manufacture (YYMMDD)

Serial number

Model number

Positive center connection

Direct current

20 Trademark

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such this does not justify the conclusion that the denotation in questions is free of third-parties rights.

21 Software License

////////////////////////////////////
Open-source software

This product uses software under the following open-source software license.

Zephyr

Copyright: The Zephyr Contributors.

Licensed under the Apache License, Version 2.0 (the "License"); and

You may not use this file except in compliance with the License.

You may obtain a copy of the License at

<http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "as is" basis, without warranties or conditions of any kind, either expressed or implied.

See the license for the specific language governing permissions and limitations under the license.

////////////////////////////////////

22 Link

You can view the operation and training methods for Bio Leg™ and the dedicated app MoPro™ via the following URL and QR code.

MoPro™ Guidebook for Professionals – App designed for Bio Leg™

https://bionicm.com/file/mopro_beg_en

The Bio Leg™ & MoPro™ Practitioner's Guide

https://bionicm.com/file/training_en

Instruction Manual for Professionals https://bionicm.com/file/ifu_pro_en

MoPro™ Guidebook for Users App designed for Bio Leg™

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BionicM Inc.