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Primary and secondary varicose veins, varicose veins during pregnancy, to support sclerotherapy, following vein surgery procedures, thrombophlebitis (superficial) and status post healed phlebitis, deep vein thrombosis of the leg, status post thrombosis, post-thrombotic syndrome, prevention of thrombosis in mobile patients, CVI of stages I – III according to Widmer or C1 – C6 according to CEAP, ulcer prevention, deep venous insufficiency, angiodysplasia, lymphedema, edema during pregnancy, post-traumatic edema, postoperative edema, idiopathic cyclical edema, lipedema, lipo-lymphedema, phlebolymphe-dema, congestion conditions following immobility (arthrogenic congestion syndrome, paresis and partial paresis of the extremity), improvement of venous symptoms, occupational edema (sedentary and standing jobs).

Decompensated heart failure (NYHA III + IV), septic phlebitis, acute bacterial, viral or allergic inflammation.

Advanced peripheral arterial disease (if one of these parameters applies):
ABPI < 0.5 ankle artery pressure < 60 mmHg, Pronounced weeping
 dermatoses, intolerance to pressure or of any of the components of the
 product, severe sensory disorders of the extremities, primary chronic
 polyarthritis, peripheral arterial disease (PAD)stage I/II, malignant
 lymphedema.

If there are open wounds in the area to be treated, the affected area must be covered with an appropriate wound dressing or bandage before the compression garment can be applied.

Severe adverse events due to compression treatment such as skin necrosis, nerve damage or thromboembolic events are rarely encountered if compression is correctly used and contraindications are considered. Discomfort, dry skin and itching are the most frequently reported adverse events related to compression use. In the event any side effects occur from use of a garment, you should consult with a medical professional.

Juzo Silver Garments should be removed prior to MRI Scans & Cardiac Defibrillations.

The manufacturer does not accept any liability for damage / injuries resulting from improper use or misuse.

WISTOCKING-0622

Freedom in Motion

Freedom to Make the Most of Every Day

Riding a bike. Working in the garden. Taking a walk. Going to work. These are the keep moving benefits of Juzo compression garments, designed to improve your quality of life, not inhibit it. Family-owned since 1912, Juzo has led the compression therapy garment industry by creating high-quality products that offer you the freedom to maintain a more active lifestyle. Advancements like our FiberSoft technology, seasonal colors, attractive styles and fun print designs ensure you move in comfort without compromising fashion. Thank you for choosing Juzo.

Care and Maintenance

To maintain your Juzo stockings, we recommend the following:

- Hand or machine wash warm daily to maximize therapeutic effectiveness and prolong the life of your Juzo garment.
- When machine washing and drying, use the lowest permanent press settings. A garment bag is strongly recommended.
- Wash dark colors separately from other clothing.
- Use Juzo washing detergent or another mild soap.
- DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives. These may damage your garment.
- Please refer to the washing symbols on the label in your Juzo stocking.
- Store your stockings in a dry place at room temperature.
- Largely resistant to sun, ointments, lotions, and perspiration.

Return Policy

Should it become necessary to return your garment for evaluation, for hygienic reasons Juzo requires that the garment be laundered prior to returning to the original place of purchase.


Fabric

Details are provided on label sewn inside your compression stocking.

Instructions for Use

It is recommended that you consult with your medical professional and fitter for proper guidelines on the use and application of this product. As therapy progresses or new symptoms develop, it may be necessary to take new measurements for a new compression stocking.

It is also recommended that you put your compression stockings on first thing in the morning every day. Rough fingernails, jagged rings, etc. can damage fabric; use a Juzo latex free donning glove to protect garment.

Useful tips to help with your compression stocking: For knee-highs, thigh-highs and pantyhose:

1 Turn the upper half of the stocking inside out to the heel by gently sliding the stocking down over the foot portion.

2 Slip your foot into the foot portion.

3 Gently ease the remaining part of the stocking over the heel.

4 Use the palms of your hands to massage the stocking upwards towards the knee. NOTE: Never pull at the upper edge of the stocking.

5 Follow the same procedure for the opposite leg, taking care that the stitch rows always run vertically.

6 For thigh-highs and pantyhose, position both stockings, as described above, until just below the knee

7 Gently ease the stocking onto the thigh. Work the stocking upwards into place using the palms of your hands until the stocking reaches the groin area. NOTE: If putting on a thigh-high stocking with silicone border, never pull the stocking up from the silicone border.

8 For pantyhose, ease the body part of the panty hose to your waist. The seams should run vertically up the front of the garment.

Always use Juzo latex free donning gloves when applying compression garments.

Enclosed with your JUZO OPEN TOE compression stocking is a complimentary Slippie designed to help you put on your open toe compression garment.

Follow the below instructions.

1 Insert half the Juzo Slippie into the top of the compression stocking.

2 Place the compression stocking w/ the inserted Juzo Slippie lengthwise on the Easy Pad. (Easy Pad is not included)

3 Slip your foot into the Juzo Slippie while it is inside the compression stocking.

4 Put weight onto the entire foot. Slide it forward and lift your toes slightly during the movement. The compression stocking will slide over the heel.

5 Repeat the last step until your entire foot is in the stocking.

6 When the heel of the stocking is in place adjust the foot portion by sliding your foot sideways to the left or right on the Juzo Easy Pad.

7 Place the heel of your opposite foot on the tip of the Juzo Slippie. Pull our leg back removing Slippie from inside the stocking.

8 While wearing Juzo donning gloves, work the stocking further up the leg. Constantly correct the fit of the compression stocking and evenly distribute the material.

For more information about Juzo donning aids, go to: www.juzousa.com • www.juzo.ca

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