




Instructions for Use

CHEETAH KNEE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	21
IT	Istruzioni per l'uso	26
NO	Bruksanvisning	31
DA	Brugsanvisning	36
SV	Bruksanvisning	41
FI	Käyttöohjeet	46
NL	Gebruiksaanwijzing	51
PT	Instruções de Utilização	56
PL	Instrukcja użytkowania	61
CS	Návod k použití	67
TR	Kullanım Talimatları	72
RU	Инструкция по использованию	77
JA	取扱説明書	83
ZH	中文说明书	88
KO	사용 설명서	93

ENGLISH

The Cheetah Knee is referred to as the device in the following. This document provides you important information on the indication for use, alignment, setup and handling of the device. This manual is intended for use by a certified prosthetist and the user of the device. The device should only be fitted and aligned by a certified prosthetist. Read this document carefully and follow the safety instructions.

PRODUCT DESCRIPTION

The device is a polycentric (four bar) mechanical prosthetic knee. It is intended to extend user activity to sport specific application such as long distance running and sprinting. The device is suitable for multi-speed ambulation, high impact levels and heavy loading.

The sprinting version of the device (available as special order) offers faster flexion during swing phase, allowing the user to run faster.

Adjustment of valve H is not available for the sprinting version as it is preset for fast running or sprinting.

NUMBERING IN FIGURES

1. Alignment reference line
2. Socket
3. Pivot axis
4. Front link
5. Hydraulic housing
6. Back link
7. Chassis
8. Valve F (Flexion)
9. Valve H (Heel Rise / Initial Flexion)
10. Valve E (Extension)
11. Bearings (lubrication points)

SAFETY

Warning: Statement that alerts the user to the possibility of serious injury, or serious adverse reactions associated with the use or misuse of the device.

Caution: Statement that alerts the user to the possibility of a problem with the device associated with its use or misuse.

Precaution: Please follow the safety precautions as failure to do so can lead to malfunction of the device and risk injury to the user. Prosthetists should make sure their clients are aware of, and understand, the safety precautions.

General Safety Instructions

Make sure that the end user of the device knows how the device functions and how to operate with it. Please inform the end user about all safety instructions before he starts using the device.

Make sure that the user is instructed to contact his/her prosthetist if experiencing change in performance of the device

Inform the end user to never operate this product if it is not working properly, if it has been dropped or damaged. In such cases the product should be returned to the healthcare provider for examination.

Caution: Avoid placing hands or fingers near moving joints.

FUNCTIONS

The device provides stability in stance through alignment with four bar

geometry and provides adjustable extension and flexion swing rate through hydraulic damping system. The device allows for smooth transition from walking to running and will accommodate for rapid flexion and extension.

MEDICAL INDICATIONS

The product is a non-invasive, single patient, reusable system. The system is not used in direct contact with the body. The device is to be used exclusively for the exoprosthetic fitting of amputations of the lower limb.

- Unilateral transfemoral amputation
- Bilateral transfemoral amputation
- Unilateral knee disarticulation amputation
- Bilateral knee disarticulation amputation

INDICATIONS FOR USE

This product is recommended for use by above knee amputees with a high activity level that subject the knee to repeated high impacts or walking/jogging longer distances on a regular basis.

- Activity level: Moderate to High (K3-K4)
- Impact level: High
- Maximum weight limit: 100 kg (220 lbs)

CONTRAINDICATIONS FOR USE

None

DEVICE SETUP

Refer to the Össur catalog for recommended selection of adapters.

Caution: Always use original Össur adapters. Adapters from other manufacturer may not fully engage all threads. They may compromise strength and void the warranty. Never tighten adapters against the hydraulic head.

Recommended Alignment of the Device

Bench Alignment

For bench alignment refer to Figure 2.

- Assemble prosthetic foot, pylon and knee.
- Establish knee center height, measured from the pivot axis (3)
- Attach socket with selected adapter combination to knee. Tighten to torque as specified in adapter IFU. Ensure flexion/extension and abduction/adduction angles are as determined in patient assessment
- Establish vertical pylon for proper knee function. Alignment reference line (1) from bisection of socket (2) on ischial level should pass 0–10 mm in front of the pivot axis (3) of the knee and an alignment reference point of the foot (as specified in selected foot IFU)

If alignment reference line (1) is too far anterior to pivot axis (3), more effort may be required to initiate knee flexion.

Static Alignment

For safety, make initial adjustments with patient standing between parallel bars.

- Fit prosthesis and check for correct length.
- Make sure the knee is in neutral position and socket angles are correct.
- Evaluate sagittal and coronal alignment for correct position of the alignment reference line.

Dynamic Alignment

The initial dynamic assessment is important for getting the user familiar with relying on the prosthesis. The best way to acquire a dynamic alignment is by having a symmetric gait pattern. Any deviations will be more pronounced in running gait.

1. Evaluate the user's level of voluntary control, the user should be able to maintain knee stability.
2. Use hip extension to maintain stability into stance.
3. Train the user to maintain normal step length on the sound side.

Precaution: Össur recommends that the initial assessment is completed at the parallel bars.

Adjustments

Valve adjustment – F, H, E

For valve adjustment refer to Fig. 3. Use small flathead screwdriver for adjustment.

Valve	Factory setting	Function
F (8)	¾ open	Affects swing flexion resistance from 60°–110° (fast walking and running only)
H (9)	¾ open	Affects swing flexion resistance from 0°–60° (slower walking)
E (10)	1 ½ open	Affects whole swing extension resistance, 110°–0°

- Preset swing flexion resistance, Valve F
- Sit patient on chair. Extend prosthetic knee. Let it fall into flexion.
- Increase resistance of Valve F until a small bump observed at 60° flexion angle.
- If you observe excessive heel rise, increase resistance of Valve F until heel rise has normalized.
- Use Valve H to smooth flexion resistance changes between initial flexion and 60°.
- Increase swing extension resistance with Valve E until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension.
- Fine tune Valves F and E until a smooth and secure gait pattern is achieved for slow and fast walking/running speeds.

Caution: Excessive swing extension resistance will prevent knee from moving into full extension.

Caution: Over tightening valves or flexing knee with all valves completely closed may damage valves.

PRODUCT SPECIFIC TRAINING

Instructing new users is essential to achieve a successful rehabilitation. For safety, initial adjustments and gait training should be done with patient standing between parallel bars.

Level Ground Walking/Running

- Patient's center of mass must be close to pivot axis (3) of knee. In this position, very little hip flexor effort is required to initiate knee flexion. Knee flexion will be initiated through rotation of the hip as weight is transferred onto the ball of the foot during terminal stance.
- Allow patient to stand in the parallel bars with their weight over the ball of the foot. Initiate knee flexion with pelvic rotation and slight hip flexion. Repeat several times.
- Start walking within parallel bars. As confidence increases, start walking outside of the parallel bars.
- As the user has improved control of their prosthesis running can be initiated. The hydraulic mechanism adapts to changes in speed.

MAINTENANCE

The device is designed for low-maintenance. The device should be inspected every six months for signs of unusual wear, by a prosthetist familiar with the product.

If the knee is subjected to excessive moisture or corrosive environment, it's recommended to clean and lubricate the knee frequently.

Cleaning

Wipe knee with soft cloth moistened with small amount of kerosene. DO NOT use solvent stronger than kerosene. If kerosene is not available, use general purpose oil (3 in 1) or sewing machine oil.

Precaution: Do not dip knee or pour solvent over knee. Bearings and seals may be damaged.

Precaution: Do not use compressed air to clean knee. Air forces pollutants into bearings and may cause malfunctions and wear.

Lubrication

Two sets of roller bearings are visible, see Fig. 2 (11). Place a few drops of machine oil on the roller bearings (general purpose oil (3 in 1) or sewing machine oil can also be used). Move the knee several times and wipe off excessive oil with a soft cloth. No other parts need external lubrication.

WARRANTY

The device is warranted against defective materials and workmanship for 18 months.

SPECIFICATIONS

Structural Specifications

Product weight: 0.62 kg (1.37 lbs)

Dimension: **(Fig. 1)**

Build height: 108 mm (4.2")

Max Knee flexion available: 110 degrees

Usage Environment

- Environmental operating temperature: -10 °C (+15 °F) to 40 °C (104 °F)
- Operating humidity: 0%–100% Relative humidity.
- Shipping and storage temperature: -15 °C (+5 °F) to 50 °C (+122 °F).
- Storage humidity: 10%–95% Relative humidity (non-condensing)

Caution: The knee should not be used in a dusty environment. Exposure to sand, talcum, water or other substances that may interfere with knee functionality should be avoided.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This device has been tested according to EN ISO 10328:2006 to three million load cycles. Depending on the amputee's activity this corresponds to a duration of use of three to five years. We recommend carrying out regular bi-annual safety checks.

ISO 10328 – P5 – 100 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

Caution! Avoid placing hands or fingers near moving joints.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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