ENGLISH

DESCRIPTION OF SYMBOLS

MD Medical Device

GENERAL SAFETY INFORMATION

WARNING: In cases of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

WARNING: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

USES

Emergency and post operative knee support.

CAUTION

Patients should be instructed to inspect their knee and this product before and after each use. If compression is desired, cover the knee with an elastic wrap before applying the brace.

PHOTO CAPTIONS

- A. Femoral Cuff and Straps
- B. Locking Buttons and Sizing Indicators
- C. Adjustable Contact Closure
- D. Posterior Stay
- E. Popliteal Pad
- F. Tibial Cuff and Straps

FITTING INSTRUCTIONS

1. Undo all contact closure straps.
2. Determine length of posterior stay assembly required for the patient. Product is pre-set to a 20" length.
3. Adjust femoral cuff by pressing the blue locking button as shown. Repeat adjustment for the tibial cuff. For best results, place the tibial cuff at or above the calf muscle.
4. Apply posterior stay assembly to the leg. Center the posterior stay along the midline of the leg. This product may be applied over a bare leg (under clothing), over a compression dressing, or over clothing. If over clothing, ensure pant leg is pulled down taut.
5. Center the knee between the two cuffs. Ensure that the popliteal pad is centered behind the knee.
6. Comfortably secure all four straps in order from bottom to top. Ensure that middle adjustable straps are not directly over patella area.
7. If necessary, slide straps up or down for maximum patient comfort. Straps should be positioned at least 2 inches away from the patella.
8. The cuff straps are length adjustable. If one or both of the cuff straps are too long, unfasten the hook strip and move until desired length is achieved. Refasten the hook strip to the strap.

WASHING INSTRUCTIONS

1. This product may be cleaned using a mild detergent and cold water. Remove cuff pads from posterior stay splint before washing. Do not use bleach or strong detergents. Do not wash in washing machine.
2. Lay cuff pads out flat on a towel and air dry. Do not squeeze or wring out the cuff pads. Do not place in dryer or near heater.

Exoform® Cool Knee Immobilizer

Part #	Description	Length	U/M
222100	Cool	18" - 24"	1 ea

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY TEXT

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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