
4R180 Harmony P4, 4R181 Harmony P4 HD

DE Gebrauchsanweisung	5
EN Instructions for use	14
FR Instructions d'utilisation	23
IT Istruzioni per l'uso	32
ES Instrucciones de uso	42
PT Manual de utilização	51
NL Gebruiksaanwijzing	61
SV Bruksanvisning	70
DA Brugsanvisning	79
NO Bruksanvisning	87
PL Instrukcja użytkowania	96
CS Návod k použití	105
SK Návod na používanie	114
TR Kullanma talimatı	123
JA 取扱説明書	132

INFORMATION

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INFORMATION

Last update: 2017-08-23

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The Harmony P4 and Harmony P4 HD mechanical vacuum pumps serve to generate an elevated vacuum in the prosthetic socket. The vacuum pump is mounted on the prosthetic socket. While walking, the vacuum pump is compressed and released again, thus generating the vacuum. Thanks to its design which includes an integrated elastomer rod, the vacuum pump dampens the vertical impact load and permits slight torsion of the prosthetic socket.

The Harmony system consists of a vacuum pump, a total surface weight-bearing socket, a polyurethane liner and a sealing sleeve.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.


2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

- The maximum approved body weight is specified in the technical data (see Page 23). The product must be used in TT prostheses only.

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C (14 °F) to +60 °C (140 °F)
Storage and transport conditions: -10 °C (14 °F) to +60 °C (140 °F), relative humidity: no restrictions
Moisture: fresh water

Unallowable environmental conditions
Mechanical vibrations or impacts
Salt water, perspiration, urine, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

2.5 Qualification

The product may only be integrated into the prosthesis by O&P professionals who have completed training as certified Harmony experts.

- **Fabricating the total surface weight-bearing socket:** Fabricating a total surface weight-bearing socket is required for installation of the product.

3 Safety

3.1 Explanation of warning symbols

CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION
Excessive strain on the product Risk of injury due to breakage of load-bearing components ► Use the product according to the specified area of application (see Page 15).

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by reduced prosthetic suspension.

4 Scope of delivery

Fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Harmony P4* with threaded ring and lamination disc	–
2	1	Lamination dummy	–
3	1	Pin wrench for Harmony P4*	4X901
4	1	Vice jaws for Harmony P4*	–
5	1	Allen key 2 mm	–
6	1	Rubber connector cap	–
7	1	Vacuum connector	2R117-0

Fig. 1, item	Quantity	Designation	Reference number
8	1	Vacuum test connector	–
9	1	Connection set	755Z39
–	4	Set screw	506G3=M8x16-"NIRO"
–	1	For 4R181 only: elastomer rod (red)	4Y347

Spare parts	
Designation	Reference number
Elastomer rod (red)	4Y347
Elastomer rod (yellow)	4Y348
4 set screws M4x6	4X902
Lamination disc set (1 lamination disc, 1 lamination dummy)	4X903
Threaded-ring set for 4R180 (1 threaded ring for 4R180, 4 set screws M4x6)	4X904
Threaded-ring set for 4R181 (1 threaded ring for 4R181, 4 set screws M4x6)	4X446

5 Preparation for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- Clean the threads before every installation.
- Apply the specified torque values.
- Follow the instructions regarding the length of the screws and about how to secure the screws.

INFORMATION

The layup described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layup.

5.1 Fabricating the thermoplastic socket

A prosthetic socket made of thermoplastic synthetic material is fabricated either as a check socket or the inner layer of the definitive socket, sealing the definitive socket. The difference between the two prosthetic sockets is the choice of material.

- > **Required materials:** 616T283* ThermoLyn clear or 616T252* ThermoLyn rigid for a check socket, 616T183* ThermoLyn PETG clear for a definitive socket, wet sanding paper, degreasing cleaner (such as 634A58 isopropyl alcohol)
 - 1) Vacuum form the prosthetic socket using thermoplastic synthetic material.
 - 2) **Check socket only:** Roughen the distal end of the prosthetic socket from the outside with wet sanding paper.
 - 3) **Definitive socket only:** Roughen the entire prosthetic socket from the outside with wet sanding paper.
 - 4) Degrease the outside of the prosthetic socket.
 - 5) Glue the lamination disc to the prosthetic socket (see Page 18).

5.2 Gluing the lamination plate

- > **Required materials:** 627B4 polyethylene adhesive tape, suitable adhesive (such as 636K49 SuperGlue, sealing resin-talcum mixture, Fabtech +PLUSeries 60-second adhesive), 699G30 synthetic support bandage
- 1) Screw the lamination plate with the threaded ring onto the vacuum pump.
- 2) Cover the hole in the centre of the lamination plate with polyethylene adhesive tape (see fig. 2).
- 3) Install the distal prosthetic components as described in the respective instructions for use.
- 4) Perform the bench alignment and determine the required structural height and alignment of the prosthetic components relative to the prosthetic socket.
- 5) **INFORMATION: The glued joint between the lamination plate and prosthetic socket must be air tight.**
Glue the lamination plate to the prosthetic socket using suitable adhesive (see fig. 3). In doing so, ensure that air inclusions are avoided.
- 6) Allow the adhesive to cure fully.
- 7) Once the adhesive has hardened, unscrew the vacuum pump from the lamination plate.
- 8) Either complete the check socket (see Page 18) or fabricate the definitive socket (see Page 18).

5.3 Completing the check socket

- > **Required materials:** 627B4 polyethylene adhesive tape, 699G30 synthetic support bandage
- 1) Secure the connection between the prosthetic socket and lamination plate with synthetic support bandage (see fig. 4). In doing so, wrap the support bandage firmly around the prosthetic socket up to the groove in the lamination plate.
- 2) **INFORMATION: In order for the sealing sleeve to seal, it has to overlap an air-tight surface by about 5 cm.**
If required, wrap the support bandage which is permeable to air with polyethylene adhesive tape to create an air-tight surface.
- 3) Drill a hole ($\varnothing$ 4 mm) through the middle of the lamination plate into the plaster model (see fig. 5).
- 4) Remove the prosthetic socket from the plaster model.
- 5) Deburr the hole on the inside of the socket.
- 6) Clean the prosthetic socket, lamination plate and bore hole.
- 7) Install the vacuum pump (see Page 19).
INFORMATION: Seal the hole in the prosthetic socket with a piece of adhesive tape if the prosthetic socket requires modification after pump installation. This prevents debris from entering the vacuum pump.
- 8) **Optional:** Install the 2R117-0 socket connector as described in the applicable instructions for use. This allows the vacuum level to be determined.

5.4 Fabricating the definitive socket

5.4.1 Laminating the transtibial socket

- > **Required tools and materials:**
99B81* PVA bag, 623T3=8 or 623T3=10 Perlon stockinette, 81A1=8 or 81A1=10 nylon tube sock, 616G15 carbon fibre woven hose, Nylglass stockinette 623T9*, 617H119 Orthocryl lamination resin 80:20 PRO, string
- > **Requirement:** The thermoplastic socket has been fabricated and the lamination plate has been glued on.
- 1) Apply a thin layer of Vaseline to the inside of the lamination plate (see fig. 6).
- 2) Insert the lamination dummy into the lamination plate (see fig. 7).
- 3) Seal the distal surface of the lamination plate with polyethylene adhesive tape (see fig. 8).
- 4) Pull a nylon stockinette over the plaster model.

- 5) Pull a piece of carbon fibre woven hose (2 times the length of the plaster model) over the plaster model up to the socket brim. Tie off the excess carbon fibre woven hose and fold it over the plaster model.
- 6) Pull a piece of nylglass stockinette over the plaster model.
- 7) Pull a piece of carbon fibre woven hose (2 times the length of the plaster model) over the plaster model up to the socket brim. Tie off the excess carbon fibre woven hose and fold it over the plaster model.
- 8) Pull a piece of Perlon stockinette (2 times the length of the plaster model) over the plaster model up to the socket brim. Tie off the excess Perlon stockinette and fold it over the plaster model.
- 9) Soak a PVA bag and pull it over the plaster model.
- 10) Complete the lamination process with Orthocryl.
- 11) Finish the prosthetic socket (see Page 19).

5.4.2 Finishing the prosthetic socket

> **Required materials:** 636W46 Loctite 242

- 1) Remove the PVA bag.
- 2) Sand off the laminate distally on the lamination plate and remove the lamination dummy.
- 3) Mark the trim line and cut it using an oscillating saw.
- 4) Drill a hole (**Ø 4 mm**) through the middle of the lamination plate into the plaster model.
- 5) Remove the prosthetic socket from the plaster model.
- 6) Sand the trim line and round off the edges. Deburr the hole on the inside of the socket.
- 7) Use compressed air to clean the prosthetic socket, bore hole and lamination plate.
- 8) Clean the inside of the lamination plate.
- 9) Install the vacuum pump (see Page 19).

INFORMATION: Seal the hole in the prosthetic socket with a piece of adhesive tape if the prosthetic socket requires modification after pump installation. This prevents debris from entering the vacuum pump.

- 10) **Optional:** Install the 2R117-0 socket connector as described in the applicable instructions for use. This allows the vacuum level to be determined.

5.5 Installing the vacuum pump

> **Required materials:** 636W46 Loctite 242

- 1) Slide the threaded ring onto the vacuum pump (see fig. 9).
- 2) With the supplied vice jaws, secure the vacuum pump in a vice (see fig. 10).
- 3) Set the lamination plate onto the vacuum pump and align it.
- 4) **Fitting:**
Use the pin wrench to tighten the threaded ring (**torque: 50 Nm**) (see fig. 11).
Screw the two (2) set screws into the threaded ring (see fig. 12).
- 5) **Definitive mounting:**
Use the pin wrench to tighten the threaded ring (**torque: 50 Nm**) (see fig. 11).
- 6) Apply Loctite to the two (2) set screws and screw them into the threaded ring (see fig. 12).
Screw in the set screws until they touch the lamination plate; then tighten them by half a turn.
- 7) Install the distal prosthetic components as described in the respective instructions for use.

5.6 Adjusting the vacuum pump

The vacuum pump offers the following adjustment possibilities: vacuum generation and damping by elastomer rod, rotation between vacuum pump and lamination plate, plantar flexion and dorsiflexion, and pronation and supination via the pyramid receiver. The distal prosthetic components have to be removed to adjust the elastomer rod. If there is no other possibility to do so, the connection between the vacuum pump and lamination plate has to be loosened to adjust the rotation of the prosthetic foot.

Vacuum generation/damping adjustment

> **Required materials:** Allen wrench 4 mm

- 1) Disassemble the distal prosthetic components.
- 2) Screw in the adjustment screw in the distal end of the vacuum pump fully to the stop (see fig. 13).
- 3) Loosen the adjustment screw as described in the adjustment table below.
- 4) **Optional:** To replace the elastomer rod, fully unscrew the adjustment screw. Repeat the adjustment process after inserting a new elastomer rod.
- 5) Install the distal prosthetic components as described in the respective instructions for use.

Adjustment table (turning out the elastomer rod)			
Body weight [kg]	Product	Red rod	Yellow rod
50 (110 lbs)	4R180, 4R181	3	–
60 to 70 (132 lbs to 154 lbs)		2.5	–
80 to 90 (176 lbs to 198 lbs)		2	–
100 (220 lbs)		1.5	3
110 (242 lbs)	4R181	–	3
120 to 130 (264 lbs to 286 lbs)		–	2.5
140 to 150 (308 lbs to 330 lbs)		–	2
The optimal setting depends on the user's level of activity and requests. Loosen the adjustment screw of the elastomer rod by a maximum of four (4) turns.			

Rotation adjustment


> **Required materials:** 636W46 Loctite 242

- 1) Loosen the 2 set screws in the threaded ring.
- 2) Loosen the threaded ring using the pin wrench until the vacuum pump can be turned.
- 3) Adjust the rotation.
- 4) Use the pin wrench to tighten the threaded ring (torque: **50 Nm**).
- 5) Apply Loctite to the two (2) set screws and screw them into the threaded ring (see fig. 12). Screw in the set screws until they touch the lamination plate; then tighten them by half a turn.

5.7 Testing the vacuum system

First steps

- **If a socket connector is installed:** Connect a vacuum gauge to the socket connector.
- **Without socket connector and vacuum gauge:** The sealing sleeve is pulled to the socket brim by the vacuum. Then the socket brim is clearly visible under the sealing sleeve.
- Ask the user to put on the prosthesis. Check the process to exclude errors during donning.
- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Perform a systematic inspection of the prosthesis.

To perform the check for leaks, the prosthesis is divided into areas.

- **Area 1:** Sealing sleeve and optional socket connector
- **Area 2:** Connection to the vacuum pump
- **Area 3:** Vacuum pump

Check area 1 (with socket connector)

Check socket connector

- Seal the opening to the socket connector on the inside of the prosthetic socket with a small piece of paper and air-tight adhesive tape (e.g. PVC tape).
- Connect a manual pump to the socket connector and build up a vacuum (**approx. 500 hPa or 15 inHg**) (see fig. 14).
- Check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Glue the socket connector again and re-check it.

Check the sealing sleeve

- Remove the adhesive tape in the prosthetic socket.
- Seal the opening to the vacuum pump on the inside of the prosthetic socket with air-tight adhesive tape (e.g. PVC tape).
- Get the user to put on the prosthesis and roll the sealing sleeve over the thigh.
- With the manual pump, build up a vacuum on the socket connector (**approx. 500 hPa or 15 inHg**).
- Check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Replace the sealing sleeve and check again.

Check area 1 (without socket connector and vacuum gauge)

- Replace the sealing sleeve with a new one.
- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is not maintained:** The problem lies in area 2 or 3.

Check area 3

- Clean the vacuum pump to exclude the possibility of dirty valves (see Page 22).
- Remove the vacuum pump from the prosthetic socket and the distal prosthetic components.
- Connect the pump vac test fitting to the proximal suction connection of the vacuum pump (see fig. 15).
- Set the vacuum pump with the hose down onto a suitable surface. Make sure the hose is not kinked.
- Produce a vacuum by pumping and check whether the vacuum is maintained (see fig. 16).
- **If the vacuum is not maintained:** Send the vacuum pump to the manufacturer.
- **If the vacuum is maintained:** The problem lies in area 2.

Check area 2

- Seal the opening to the vacuum pump on the inside of the prosthetic socket with air-tight adhesive tape (e.g. PVC tape).
- Connect a manual pump to the lamination plate and build up a vacuum (**approx. 500 hPa or 15 inHg**) (see fig. 17).
- Check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Fill the bore hole between the lamination plate and prosthetic socket, starting at the lamination plate, with sealing resin or Fabtech +PLUSeries 60-second adhesive and allow the adhesive to set.
- Drill a new hole with a smaller diameter and perform the test again.

6 Use

CAUTION

Use of powder (e.g. baby powder or talcum)

Skin irritation on the residual limb and loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- Keep the product away from powders.

6.1 Donning the prosthesis

- 1) Don the liner on the residual limb. Avoid wrinkles, air pockets and soft tissue displacement in the process.
- 2) **For liners without textile coating:** Pull on a nylon sock (451F21).
- 3) Slide into the prosthetic socket and unroll the sealing sleeve over the prosthetic socket onto the thigh.

6.2 Doffing the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Doff the residual limb and liner out of the prosthetic socket.

7 Cleaning

7.1 Cleaning the prosthetic socket

- Clean the product with a damp, soft cloth and then let it dry.

7.2 Cleaning the liner

- Clean the liner as described in the liner instructions for use.

7.3 Cleaning the vacuum pump

For qualified personnel: cleaning with compressed air

- With compressed air, carefully blow through the hole in the prosthetic socket to remove dirt in the vacuum pump and one-way valve.

For qualified personnel: flushing with distilled water

- 1) Remove all distal prosthetic components.
- 2) Fill the prosthetic socket with approx. 60 ml distilled water.
- 3) Place a towel on the floor or a table and set the vacuum pump onto it.
- 4) Run the vacuum pump until the water has been pumped through.
- 5) Install the distal prosthetic components as described in their instructions for use.

For users: cleaning in the course of use

If the vacuum pump no longer generates a sufficient vacuum and no O&P professional can be reached, the vacuum pump can be rinsed with a small amount of water.

- 1) Moisten a nylon sheath with distilled water.
- 2) Position the nylon sheath in the prosthetic socket.
- 3) Don the prosthesis and use it as normal. The vacuum pump draws the water through the hole in the bottom of the socket.
- 4) Contact the O&P professional to set an inspection date.

8 Maintenance

- A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.
- Perform a visual inspection of the elastomer bar annually and replace it if needed.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria out-

lined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R180	4R181
Weight [g]	465	590
System height [mm]	132	
Body weight [kg]	50 to 100 (110 lbs to 220 lbs)	50 to 150 (110 lbs to 330 lbs)
Mobility grade	2 to 4	

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Otto Bock HealthCare GmbH