
## ENGLISH

### CONTENT

- L-292011 Lock body 4-hole with guide (w/o clutch mechanism).
- L-292118 Delrin tooling body.
- L-292125 Lam/therm ring.
- L-292127 Pivot, lanyard.
- L-292129 Delrin screw.
- L-292130 Bushing, lanyard.
- L-293041 Screw m6x20.
- L-293042 Screw m6x35.
- L-292010 Cleat with cover and rivets.
- L-293012 Cord cover 40 cm.
- L-292016 Pull handle, lanyard.
- L-292018 Lanyard cord with connecting screw.
- L-292117 Elbow, lanyard.

### LAMINATION

- Flatten (7). The distal end enough to centrally attach the delrin tooling body (2) with the four nails (3).
- Apply the inner pva bag. Pull through the central hole of the aluminum ring (4). Flare the bag away from the central hole and attach the ring using the delrin screw (5), sealing the pva bag in the process. Ensure that the arrow (6) on the aluminum ring is aligned in the ap or ml direction. Coat the four spacing screws (7) and the delrin screw with a release agent.
- Tie in your lay-up material. Tie all but the first layer beyond the aluminum ring. Ensure that the lamination material is at least 6 mm thick (8) at the distal end and the four spacing screws are embedded evenly in the material.
- Apply outer pva bag and finish lamination.
- Grind the distal end accurately to expose the spacer screw heads. Do not proceed any deeper!
- Remove the spacer screws and attach the lock body (9) with the four long countersunk screws (10). Use locktite 425 and torque to 12 nm (9 ft.Lbs).

### Thermoplastic

- Flatten (1) the distal end enough to centrally attach the delrin tooling body (2) with the four nails (3).
- Starting at the edge of the delrin tooling body, drill a small hole with a wire drill through the full length of your plaster model to its base. This assures maximum vacuum at the distal end.
- Attach the aluminum ring (4) to the delrin tooling body using the delrin screw (5). Ensure that the arrow (6) on the aluminum ring is aligned in the ap or ml direction.
- Vacuum form the thermoplastic. Ensure proper forming around the aluminum ring. It is important that the thermoplastic is at least 6 mm thick (8) at the distal end.
- Grind the distal end accurately to expose the spacer screw heads. Do not

proceed any deeper!

- Remove the spacer screws (7) and attach the lock body (9) with the four long countersunk screws (10). Use locktite 425 and torque to 12 nm (9ft.Lbs)

### Lanyard components

- Shorten the elbow (11) to the desired length and glue it into the bushing (12) with locitite 425.
- Find the correct positioning for the cleat (13), drill two 3 mm (1/86) holes, and rivet it to the socket.
- Trim the cord cover (14) to the desired length. Glue it into the elbow with locitite 425 and to the socket in the correct location.

### Cleat positioning (Figure C)


- (16) Released
- (17) Locked

**IMPORTANT:** tie a knot in (15) the lanyard cord to stop it from inadvertently being pulled into the channel.

**NOTE:** the pull handle (not shown in diagram) can be used for getting a more secure grip on the lanyard cord when the socket is donned.

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

- 1 Year limited warranty
- It is recommended that the lanyard cord be regularly checked for wear.
- Rated up to high activity level (k3)
- Weight limit: 120 kg (265 lbs)

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
**EN – Caution:** Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Ossur products are used with other recommended Ossur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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www.ossur.com ©Copyright Ossur 2018 IFU 0023 IFU 1234\_001 Rev. 5