

















Instructions for Use







ICELOCK[®] 850 HYBRID
FOR 4-HOLE SOCKET
















		3
EN	Instructions for Use	8
DE	Gebrauchsanweisung	16
FR	Notice d'utilisation	25
ES	Instrucciones para el uso	33
IT	Istruzioni per l'uso	41
NO	Bruksanvisning	49
DA	Brugsanvisning	56
SV	Bruksanvisning	64
EL	Οδηγίες Χρήσης	72
FI	Käyttöohjeet	81
NL	Gebruiksaanwijzing	89
PT	Instruções de Utilização	97
PL	Instrukcja użytkowania	105
TR	Kullanım Talimatları	113
RU	Инструкция по использованию	121
JA	取扱説明書	129
ZH	中文说明书	137
KO	사용 설명서	144

DESCRIPTION

The device offers a combination of mechanical and vacuum suspension. The mechanical suspension uses a lanyard, allowing the patient to pull themselves into the socket. The device can be used exclusively as mechanical suspension or in combination with an expulsion valve to provide either passive or active suspension.

Components

- A. Central Screw
- B. Fixation Dummy
- C. Distal Connector
- D. Nails
- E. Manufacturing Dummy 1
- F. Manufacturing Dummy 2
- G. Aluminum Dummy
- H. Spacers
- I. Flat-head Screws for 4-hole
- J. Attachment Pin for Icelock 850
- K. Sealing Insert
- L. Lock Housing for 4-hole
- M. Large Square Ring
- N. Large O-ring
- O. Small Square Ring
- P. Quick Connect
- Q. Strap
- R. Routing Insert
- S. Latch
- T. Slider Button with screw
- U. Latch Housing Frame
- V. Latch Housing
- W. Strap Lock with screws

All parts are reused between the check socket and definitive socket except for the Latch Housing

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

INDICATIONS FOR USE AND TARGET PATIENT POPULATION

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

The weight limit for the device is 166 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur replacement parts for the device for ensured durability.

The device is for single patient use.

REQUIRED COMPONENTS

The device should be used in combination with an Iceross Seal-In® X Locking liner, Iceross Seal-In® X Seal and Icelock® 552 Expulsion Valve or Unity® Valve.

ASSEMBLY INSTRUCTIONS

The following socket types can be fabricated with this lock version.

- Thermoplastic check socket
Follow sections Plaster Positive and Thermoplastic Sockets
- Definitive socket with a thermoplastic rigid inner material.
Follow sections Plaster Positive, Thermoplastic Sockets (Steps 1-8) and Definitive Socket.
- Definitive socket with a thermoplastic flexible inner material.
Follow sections Plaster Positive, Flexible Inner Socket and Definitive Socket.

Socket Construction

Prepare plaster negative as usual.

Plaster Positive

Note: If using a moist plaster, isolate properly.

1. Model the plaster positive to agree with the residual limb.
2. Flatten the distal end until the diameter is the same size as the Manufacturing Dummy 1+2 (**Fig. 1**). For liner size 22 and above, use the two Manufacturing Dummies. For liner sizes smaller than size 22, only use the Manufacturing Dummy 1.
3. Use one of the two following methods to ensure air flow to the Distal Connector at the distal end:
 - Drill a hole through the whole plaster from the center of the distal end to the proximal end before fixing the Manufacturing Dummy.

OR

- Apply a thin stocking over the plaster.
- 4. Use the Nails supplied to attach the correct Manufacturing Dummy/ Dummies to the plaster (**Fig. 2**).

If using the Icelock 552 Expulsion Valve, perform steps 5 to 7.

5. Determine the valve location on the plaster and flatten the area to ensure the valve Dummy sits flush. In the center hole of the Dummy place a mark on the plaster with a permanent marker. Use the red Dummy when using a rigid inner material and the green one when using a flexible inner material.
6. Drill a vacuum hole from the Dummy location to the proximal end of the plaster or the previously drilled hole.
7. Attach the Dummy to the plaster, using the Nails provided.

If using the Unity Valve, perform steps 8 to 10 to fabricate the Socket Insert into the socket.

8. Fix the Socket Insert to the plaster using the Set Screw. Ensure it sits flush to the plaster.
9. Plaster may have to be removed or added.
10. Apply the Foam Cap to the Set Screw.

Thermoplastic Sockets

1. Put the Small Square Ring onto the cut-out on the Distal Connector. Use the Central Screw and the Fixation Dummy to secure the Distal Connector on the Manufacturing Dummy/Dummies and the plaster (**Fig. 3**).
- Note:**
- Orientate the Distal Connector to the line of progression.
 - For liners below size 22, replace the Large Square Ring with the Large O-ring before draping the material (**Fig. 4**).
2. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 5**).
 3. Tightly wind a strap into the groove of the Distal Connector during draping. This will help the socket become airtight. Additionally, guide the plastic around the valve Dummy to ensure proper molding.
 4. After fully hardening, grind and remove the material at the distal end until completely exposing the head of the Central Screw (**Fig. 6**).
 5. Remove the Central Screw.
 6. Continue grinding, until fully exposing the Fixation Dummy.
 7. Remove the Fixation Dummy and the Small Square Ring.
 8. Grind carefully until the material is flush with the surface of the Distal Connector (**Fig. 7**).
 9. If applicable, expose the valve Dummy by grinding away material until fully flat area is exposed.
10. Put the Aluminum Dummy on the socket and align the four holes of the Aluminum Dummy with the four holes on the Distal Connector (**Fig. 8**).
 11. Secure the Aluminum Dummy with the Central Screw.
 12. Drill four $\varnothing 6,5$ mm holes through the Aluminum Dummy until the heads of the Set Screws can be reached (**Fig. 8**).
 13. Remove the Set Screws, Central Screw, and Aluminum Dummy (**Fig. 9**).
 14. Remove the socket from the plaster.
 15. Push the Distal Connector out of the socket (**Fig. 10**).
 16. Replace the Large Square Ring with the Large O-ring (**Fig. 4**).
 17. Use a $\varnothing 10$ mm drill to drill again into the $\varnothing 6,5$ mm holes in the socket (**Fig. 11**).
- Note:** Do not damage the cylindrical walls inside of the socket, otherwise the sealing won't work.
- The Latch lock can be placed either anteriorly or laterally on the socket, according to the patient's preference. Consider the patient's dominant hand and hand dexterity when choosing where to place the Latch lock. The selected side cannot be changed after fabrication.
18. Place the Latch Housing perpendicular to the 4-hole connection, so the Strap will be directed straight from the Lock Housing (**Fig. 12**).
 19. Place the Latch Housing so that it is at least 2 cm below the expected proximal trimline. Orient the Latch Housing so the arrow points proximal.
- Note:** It is recommended to place the Latch Housing as proximally as possible for easier donning.

20. Mark the spot on the socket where the Latch Housing will be placed.
21. Roughen the surface on the chosen spot and glue the Latch Housing using Composite 1 adhesive, follow instructions on the adhesive packaging (**Fig. 12**).

Note: Adhesive working time is 1 minute.

22. Grind to expose the valve Dummy.

Caution: Make sure not to manipulate the distal attachment area when adjusting the socket with heat during the fitting process. Fix this area with rigid tape before heating.

Flexible Inner Socket

1. Put the Small Square Ring onto the cut-out on the Distal Connector. Use the Central Screw and the Fixation Dummy to secure the Distal Connector on the Manufacturing Dummy/Dummies and the plaster (**Fig. 3**).

Note:

- Orientate the Distal Connector to the line of progression.
 - For liners below size 22, replace the Large Square Ring with the Large O-ring before draping the material (**Fig. 4**).
2. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 5**).
 3. Tightly wind a strap into the groove of the Distal Connector during draping. This will help the socket become airtight. Additionally, guide the plastic around the valve Dummy to ensure proper molding.
 4. After fully hardening, grind and remove the material at the distal end until completely exposing the head of the Central Screw (**Fig. 6**).
 5. Remove the Central Screw, Fixation Dummy, and Small Square Ring.
 6. Cut the flexible inner socket so it sits flush to the flat surface on the Distal Connector (**Fig. 13**).

Note: Do not cut the flexible inner socket lower. The chamfered edge ensures that the inner socket is securely fastened in the laminated socket and contact to the seal ring on the Distal Connector ensures vacuum in the socket

7. Tie a PVA bag on the Distal Connector and trim the excess PVA material (**Fig. 14**).

Definitive Socket

A two-step lamination is required for making the definitive socket.

1. If doing a rigid inner socket, perform steps 1 through 8 in the Thermoplastic Sockets section. If doing a flexible inner socket, follow the Flexible Inner Socket section.
2. Put the Small Square Ring back on again.
Note: It prevents resin leaks towards the Distal Connector.
3. Align the four holes of the Aluminum Dummy with the four holes on the Distal Connector.
4. Secure the Aluminum Dummy with the Central Screw (**Fig. 15**).
5. Seal the transition between the inner socket and plaster at the proximal end.
6. Roughen the outer surface if using a rigid inner material.
7. Fill the head of the Central Screw with plasticine or similar.
8. Attach reinforcement material between the Distal Connector and Aluminum Dummy (**Fig. 16**).
9. Drape unidirectional carbon tape between the Distal Connector and Aluminum Dummy.
10. Apply the reinforcement material layup to agree with the expected stress.

The Latch lock can be placed either anteriorly or laterally on the socket, according to the patient's preference. Consider the patient's dominant hand and hand dexterity when choosing where to place the Latch lock. The selected side cannot be changed after fabrication.

11. Fold the top 2-3 layers of fiber off the plaster (**Fig. 17**). Place the Latch Housing Frame so that it is at least 2 cm below the expected proximal trimline and in a direct line from the middle of the 4-hole connection, so the Strap will be routed straight. Orient the Frame so the round end points towards the proximal end. If it is placed incorrectly, the socket will be unusable.

Note: It is recommended to place the Latch Housing as proximally as possible for easier donning.

12. Fold the layers of fiber over the Latch Housing Frame (**Fig. 18**).

13. Push the Latch Housing into the Frame, with the fiber layers in between, until it snaps securely into place (**Fig. 18**)

Note: Make sure that the arrow on the Latch Housing points towards the proximal end.

14. Put the upper PVA bag over the plaster and Aluminum Dummy.

15. When the resin is distributed, maneuver resin by hand to make sure no air pockets form around the Latch Housing and resin fully saturates the fibers underneath the Latch Housing (**Fig. 19**).

16. Before the resin gels, fasten a cord around the PVA bag, distally to the Aluminum Dummy.

17. Once resin sets, remove PVA bag and excess resin on top of the Aluminum Dummy.

18. Drill four $\varnothing 6,5$ mm holes through the Aluminum Dummy until the heads of the Set Screws can be reached (**Fig. 20**).

19. Remove the Set Screws, Central Screw, and Aluminum Dummy (**Fig. 21**).

20. Grind to clean the edges if needed.

21. Remove the socket from the plaster.

22. Push the Distal Connector out of the socket (**Fig. 10**).

23. Replace the Large Square Ring with the Large O-ring (**Fig. 4**).

24. Grind to expose the valve Dummy.

If using the Icelock 552 Expulsion valve with a flexible inner socket, perform steps 25 to 27.

25. Remove the flexible inner socket from the laminated socket.

26. Drill the laminated outer socket with a 44 mm hole saw at the valve location and sand the edge.

27. Grind the flexible inner socket to expose the valve Dummy.

Lock Assembly

1. Remove all dust, dirt and other contamination from the socket and Lock Housing.

2. Remove plastic cap and apply pressurized air into the hole in the soft Dummy to loosen it from the socket and pull the soft Dummy out of the Latch Housing.

3. Apply medium strength threadlocker to the screw on the Latch lock and fasten the Latch lock by sliding it into the slot on the Latch Housing and hand-tightening the screw (**Fig. 22**).

4. Route the Quick Connect into the Lock Housing and the Sealing Insert (**Fig. 23**).

5. Slide the Routing Insert into the Lock Housing.

6. Assemble the Distal Connector and the lock to the socket, along with a 4-hole adapter of choice. Apply medium strength threadlocker to the Flat-Head Screws and torque to 10 Nm (**Fig. 24**).

Note: If doing a thermoplastic check socket only, put the four Spacers into the \varnothing 10 mm holes.

7. Refer to the instructions for use for the Icelock 552 Expulsion Valve or Unity Valve for the assembly process.

Note: If the Strap end gets pulled out of the Routing Insert during assembly, it can be reassembled and threaded (**Fig. 25**).

FITTING INSTRUCTIONS

Warning: Failure to adhere to the following instructions may cause the pin to detach from the liner resulting in loss of suspension, which may lead to injury.

Note: Setting the length of the Strap should be done before donning the system to ensure correct engagement.

1. Apply Loctite 410/411 or an equivalent compound onto the threads on the Attachment Pin and screw it into the liner. Torque the Attachment Pin to 4 Nm.
2. Pull the Quick Connect out of the socket and fasten to the Attachment Pin on the liner (**Fig. 26**).

Caution: Make sure to fully engage the pin into the Quick Connect.

3. Pull hard on the outermost Strap (a) so the liner is fully donned into the socket (**Fig. 27**).

Note: This ensures correct final position of the residual limb inside the socket once the lock has been fastened.

4. Adjust the length by using the inner two Straps. Pull the middle Strap (b) to elongate the system and the inner Strap (c) to shorten it (**Fig. 27**). Adjust the Straps so they become tight when the Latch is nearly closed (**Fig. 28**).

5. Fasten the **two** inner Straps (b) and (c) with the Strap Lock and hand-tighten the screws (**Fig. 29**).

Note: The screws have a pre-applied threadlocking patch.

6. Fit the device to the patient using the Donning and Doffing instructions described in the next section.

Note: Only minor length adjustments can be made on the Strap after fitting the device. The Strap should not be lengthened by more than a few millimeters to ensure the Quick Connect will still enter the Lock Housing and seal with vacuum.

7. When the fitting is finished, cut the excess Strap (c) 4-6 cm above the Routing Insert and heat seal the edge to prevent fraying (**Fig. 30**).

Note: If necessary, during the lifetime of the device, the strap should be retightened by a healthcare professional. Correct strap length ensures optimal function of device.

Caution: Make sure to fasten the two inner Straps with the Strap Lock, not only one, as it will not provide a secure enough fastening.

USAGE

Donning Instructions

1. Don the liner and the seal.
2. Secure the Quick Connect onto the liner pin (**Fig. 26**).
Caution: Make sure to fully engage the pin into the Quick Connect.
3. While sitting or standing, pull the latch while donning the socket to keep tension on the strap. This will help guide the Quick Connect into the distal end (**Fig. 31**).
4. Stand up to push remaining air out of the socket before securing the latch (**Fig. 32**).

Caution: The latch won't be secured if the socket fit is too tight. Do not put excessive force into pulling the residual limb into the socket as this may cause too much tension on the distal end.

Doffing Instructions

1. Slide the latch lock down to release the latch (**Fig. 33**). Remove the latch from the socket.
2. Push the valve to allow airflow inside the socket while pulling the residual limb out (**Fig. 34**).
3. Push the black button on the Quick Connect and remove from the liner pin (**Fig. 35**).

Store the device so the Quick Connect is pulled out of the socket and the latch is not locked (**Fig. 36**).

Note: If the latch is locked while no liner is attached, the Quick Connect will be unreachable by hand in the distal end. To retrieve the Quick Connect, take the liner and put it into the socket. Push it down so the liner pin engages the Quick Connect and drag the liner up (**Fig. 37**).

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P7 - 166 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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