











Instructions for Use











---











ICELOCK<sup>®</sup> 850 HYBRID FOR  
LAMINATED SOCKET

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	13
FR	Notice d'utilisation	21
ES	Instrucciones para el uso	28
IT	Istruzioni per l'uso	35
NO	Bruksanvisning	42
DA	Brugsanvisning	48
SV	Bruksanvisning	55
EL	Οδηγίες Χρήσης	62
FI	Käyttöohjeet	70
NL	Gebruiksaanwijzing	77
PT	Instruções de Utilização	84
PL	Instrukcja użytkowania	91
TR	Kullanım Talimatları	98
RU	Инструкция по использованию	105
JA	取扱説明書	112
ZH	中文说明书	119
KO	사용 설명서	125

### DESCRIPTION

The device offers a combination of mechanical and vacuum suspension. The mechanical suspension uses a lanyard, allowing the patient to pull themselves into the socket. The device can be used exclusively as mechanical suspension or in combination with an expulsion valve to provide either passive or active suspension.

### Components

- A. Funnel Dummy
- B. Funnel
- C. Sealing Insert
- D. Attachment Pin for Icelock 850
- E. Fabrication Screw
- F. Lock Housing for Lamination
- G. Routing Insert
- H. Quick Connect
- I. Strap
- J. Latch
- K. Latch Housing Frame
- L. Latch Housing
- M. Slider Button with screw
- N. Strap Lock with screws

All parts are reused between the check socket and definitive socket except for the Funnel Dummy, Lock Housing for Lamination, and the Latch Housing.

### INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### INDICATIONS FOR USE AND TARGET PATIENT POPULATION

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

The weight limit for the device is 166 kg.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Only use Össur replacement parts for the device for ensured durability.

The device is for single patient use.

## REQUIRED COMPONENTS

The device should be used in combination with an Iceross Seal-In® X Locking liner, Iceross Seal-In® X Seal and Icelock® 552 Expulsion Valve or Unity® Valve.

## ASSEMBLY INSTRUCTIONS

The following socket types can be fabricated with this lock version.

- Thermoplastic check socket  
Follow sections Plaster Positive and Thermoplastic Sockets
- Definitive socket with a thermoplastic rigid inner material.  
Follow sections Plaster Positive, Thermoplastic Sockets (steps 1-4) and Definitive Socket.
- Definitive socket with a thermoplastic flexible inner material.  
Follow sections Plaster Positive, Flexible Inner Socket and Definitive Socket.

### **Socket Construction**

Prepare plaster negative as usual.

### **Plaster Positive**

**Note:** If using a moist plaster, isolate properly.

1. Model the plaster positive to agree with the residual limb.
2. Flatten the distal end until the diameter is the same size as the Funnel Dummy (**Fig. 1**).
3. Use one of the two following methods to ensure air flow to the Distal Connector at the distal end:
  - Drill a hole through the whole plaster from the distal end to the proximal end before fixing the Funnel Dummy. Align the drill so the hole does not go through the center of the distal end, as that space will be used for the center screw.

OR

- Apply a thin stocking over the plaster.

If using the Icelock 552 Expulsion Valve, perform steps 4 to 6.

4. Determine the valve location on the plaster and flatten the area to ensure the valve Dummy sits flush. In the center hole of the Dummy place a mark on the plaster with a permanent marker. Use the red Dummy when using a rigid inner material and the green one when using a flexible inner material.
5. Drill a vacuum hole from the Dummy location to the proximal end of the plaster or the previously drilled hole.
6. Attach the Dummy to the plaster, using the Nails provided.

If using the Unity Valve, perform steps 7 to 9 to fabricate the Socket Insert into the socket.

7. Fix the Socket Insert to the plaster using the Set Screw. Ensure it sits flush to the plaster.
8. Plaster may have to be removed or added.
9. Apply the Foam Cap to the Set Screw.

### **Thermoplastic Sockets**

The Latch lock can be placed either anteriorly or laterally on the socket, according to the patient's preference. Consider the patient's dominant

and hand dexterity when choosing where to place the Latch lock. The selected side cannot be changed after fabrication.

1. Attach the Funnel Dummy on to the plaster with the Fabrication Screw and orient the arrow towards the selected side where the Latch lock will be placed (**Fig. 2**).
2. Screw the Lock Housing onto the Funnel Dummy and tighten by hand until reaching the end stop. Hold the Funnel Dummy in place when tightening. Make sure that the red surface faces the selected side after tightening (**Fig. 3**).
3. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 4**).
4. Tightly wind a strap into the groove of the Funnel Dummy and Lock Housing during draping. This will help the socket become airtight. Additionally, guide the plastic around the valve Dummy to ensure proper molding.
5. Apply buildup material for a socket adapter of choice. Make sure not to put it in front of the red surface on the Lock Housing (**Fig. 5**).
6. Remove the socket from the plaster and trim the proximal edge of the socket.
7. Grind the thermoplastic material on the side of the Lock Housing until the red surface has disappeared (**Fig. 6**).
8. Place the Latch Housing in a direct line from the opening on the Lock Housing, so the Strap will be directed in a straight line (**Fig. 7**).
9. Place the Latch Housing so that it is at least 2 cm below the expected proximal trimline. Orient the Latch Housing so the arrow points proximal.  
**Note:** It is recommended to place the Latch Housing as proximally as possible for easier donning.
10. Mark the spot on the socket where the Latch Housing will be placed.
11. Roughen the surface on the chosen spot and glue the Latch Housing using Composite 1 adhesive, follow instructions on the adhesive packaging (**Fig. 7**).  
**Note:** Adhesive working time is 1 minute.
12. Grind to expose the valve Dummy.

**Caution:** Make sure not to manipulate the Lock Housing area when adjusting the socket with heat during the fitting process. Fix this area with rigid tape before heating.

### Flexible Inner Socket

The Latch lock can be placed either anteriorly or laterally on the socket, according to the patient's preference. Consider the patient's dominant hand and hand dexterity when choosing where to place the Latch lock. The selected side cannot be changed after fabrication.

1. Attach the Funnel Dummy on to the plaster with the Fabrication Screw and orient the arrow towards the selected side where the Latch lock will be placed (**Fig. 2**).
2. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 8**).
3. Cut the flexible inner socket so it aligns with the top of the thin wall below the threads (**Fig. 9**).  
**Note:** Do not cut the flexible inner socket lower. The chamfered edge ensures that the inner socket is securely fastened in the laminated socket and contact to the seal ring on the Funnel ensures vacuum in the socket.
4. Screw the Lock Housing onto the Funnel Dummy and tighten by hand until reaching the end stop. Hold the Funnel Dummy in place when

tightening. Make sure that the red surface faces the selected side after tightening (**Fig. 3**).

5. Tie a PVA bag tightly onto the Lock Housing, just above the Funnel Dummy, and trim the excess PVA material (**Fig. 10**).

### **Definitive Socket**

A two-step lamination is required for making the definitive socket.

1. If doing a rigid inner socket, perform steps 1 through 4 in the Thermoplastic Sockets section. If doing a flexible inner socket, follow the Flexible Inner Socket section.
2. Roughen the outer surface if using a rigid inner material.
3. Laminate the first layer of structural material. When lamination is set, pull off outer PVA bag.
4. Apply buildup material for a socket adapter of choice. Make sure not to put it in front of the red surface on the Lock Housing (**Fig. 5**).
5. Attach reinforcement material to the socket adapter and drape over the plaster.
6. Fold the top 2-3 layers of fiber off the plaster (**Fig. 11**). Place the Latch Housing Frame so that it is at least 2 cm below the expected proximal trimline and in a direct line from the center of the red surface on the Lock Housing, so the Strap will be routed straight. Orient the Frame so the round end points towards the proximal end. If it is placed incorrectly, the socket will be unusable.

**Note:** It is recommended to place the Latch Housing as proximally as possible for easier donning.

7. Fold the layers of fiber over the Latch Housing Frame (**Fig. 12**).
  8. Push the Latch Housing into the Frame, with the fiber layers in between, until it snaps securely into place (**Fig. 12**)
- Note:** Make sure that the arrow on the Latch Housing points towards the proximal end.
9. Pull PVA bag over plaster and laminate.
  10. When the resin is distributed, maneuver resin by hand to make sure no air pockets form around the Latch Housing and resin fully saturates the fibers underneath the Latch Housing (**Fig. 13**).
  11. When lamination is set, pull off outer PVA bag and remove the socket from the plaster.

12. Grind the material on the side of Lock Housing until the red surface has disappeared (**Fig. 6**).

13. Trim socket proximally.

14. Grind to expose the valve Dummy.

If using the Icelock 552 Expulsion valve with a flexible inner socket, perform steps 15 to 17.

15. Remove the flexible inner socket from the laminated socket.
16. Drill the laminated outer socket with a 44 mm hole saw at the valve location and sand the edge.
17. Grind the flexible inner socket to expose the valve Dummy.

### **Lock Assembly**

1. Remove all dust, dirt and other contamination from the socket and Lock Housing.
2. Remove plastic cap and apply pressurized air into the hole in the soft Dummy to loosen it from the socket and pull the soft Dummy out of the Latch Housing.
3. Apply medium strength threadlocker to the screw on the Latch lock and fasten the Latch lock by sliding it into the slot on the Latch Housing and hand-tightening the screw (**Fig. 14**).

4. Unscrew Funnel Dummy from Lock Housing using a 19 mm socket wrench.
5. Route the Quick Connect into the Lock Housing and into the Sealing Insert and the Funnel (**Fig. 15**).
6. Slide the Routing Insert into the Lock Housing.
7. Slide the Sealing Insert to the bottom, into the Routing Insert and screw the Funnel on using a cross driver with a socket wrench extension (**Fig. 16**). The end position of the sticker on the Funnel should align with the opening on the Lock Housing (**Fig. 17**).
8. Remove the sticker on the Funnel.
9. Refer to the instructions for use for the Icelock 552 Expulsion Valve or Unity Valve for the assembly process.

**Note:** If the Strap end gets pulled out of the Routing Insert during assembly, it can be reassembled and threaded (**Fig. 18**).

## FITTING INSTRUCTIONS

**Warning:** Failure to adhere to the following instructions may cause the pin to detach from the liner resulting in loss of suspension, which may lead to injury.

**Note:** Setting the length of the Strap should be done before donning the system to ensure correct engagement.

1. Apply Loctite 410/411 or an equivalent compound onto the threads on the Attachment Pin and screw it into the liner. Torque the Attachment Pin to 4 Nm.
2. Pull the Quick Connect out of the socket and fasten to the Attachment Pin on the liner (**Fig. 19**).

**Caution:** Make sure to fully engage the pin into the Quick Connect.

3. Pull hard on the outermost Strap (a) so the liner is fully donned into the socket (**Fig. 20**).

**Note:** This ensures correct final position of the residual limb inside the socket once the lock has been fastened.

4. Adjust the length by using the inner two Straps. Pull the middle Strap (b) to elongate the system and the inner Strap (c) to shorten it (**Fig. 20**). Adjust the Straps so they become tight when the Latch is nearly closed (**Fig. 21**).
5. Fasten the **two** inner Straps (b) and (c) with the Strap Lock and hand-tighten the screws (**Fig. 22**).

**Note:** The screws have a pre-applied threadlocking patch.

6. Fit the device to the patient using the Donning and Doffing instructions described in the next section.

**Note:** Only minor length adjustments can be made on the Strap after fitting the device. The Strap should not be lengthened by more than a few millimeters to ensure the Quick Connect will still enter the Lock Housing and seal with vacuum.

7. When the fitting is finished, cut the excess Strap (c) 4-6 cm above the Routing Insert and heat seal the edge to prevent fraying (**Fig. 23**).

**Note:** If necessary, during the lifetime of the device, the strap should be retightened by a healthcare professional. Correct strap length ensures optimal function of device.

**Caution:** Make sure to fasten the two inner Straps with the Strap Lock, not only one, as it will not provide a secure enough fastening.

## USAGE

### **Donning Instructions**

1. Don the liner and the seal.
2. Secure the Quick Connect onto the liner pin (**Fig. 19**).  
**Caution:** Make sure to fully engage the pin into the Quick Connect.
3. While sitting or standing, pull the latch while donning the socket to keep tension on the strap. This will help guide the Quick Connect into the distal end (**Fig. 24**).
4. Stand up to push remaining air out of the socket before securing the latch (**Fig. 25**).  
**Caution:** The latch won't be secured if the socket fit is too tight. Do not put excessive force into pulling the residual limb into the socket as this may cause too much tension on the distal end.

### **Doffing Instructions**

1. Slide the latch lock down to release the latch (**Fig. 26**). Remove the latch from the socket.
2. Push the valve to allow airflow inside the socket while pulling the residual limb out (**Fig. 27**).
3. Push the black button on the Quick Connect and remove from the liner pin (**Fig. 28**).

Store the device so the Quick Connect is pulled out of the socket and the latch is not locked (**Fig. 29**).

**Note:** If the latch is locked while no liner is attached, the Quick Connect will be unreachable by hand in the distal end. To retrieve the Quick Connect, take the liner and put it into the socket. Push it down so the liner pin engages the Quick Connect and drag the liner up (**Fig. 30**).

### **Cleaning and care**

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

### **Environmental Conditions**

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

## MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

### **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

## Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P7 - 166 kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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