



















Instructions for Use

MIAMI TLSO 464


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	25
NO	Bruksanvisning	29
DA	Brugsanvisning	33
SV	Bruksanvisning	37
EL	Οδηγίες Χρήσης	41
FI	Käyttöohjeet	46
NL	Gebruiksaanwijzing	50
PT	Instruções de Utilização	55
PL	Instrukcja użytkowania	60
CS	Návod k použití	65
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	73
JA	取扱説明書	78
ZH	中文说明书	82
KO	사용 설명서	86

19

ENGLISH

PRODUCT DESCRIPTION

Miami TLSO 464 is a reusable rigid thoraco-lumbo-sacral spine orthosis implemented as a modular compression belt. The device is designed to support and provide motion restriction to a patient's spine during routine daily use.

The device is a non-invasive, adjustable brace for single patient use. Device length can be adjusted to patient morphology requirements. A pulley system is used to tighten the belt and belt panels to the patient.

List of parts (Fig. 1)

- a. Yoke
- b. Strut
- c. Overlap closure/ belt end
- d. Buckle
- e. PTE
- f. Belt
- g. Pectoral disc
- h. ATE

INDICATIONS FOR USE

Gross thoracic and lumbar spine immobilization in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative spinal pathologies
- Post-surgical stabilization
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

CONTRAINDICATIONS

Do not use with unstable, displaced fractures.

INTENDED USE

Patients requiring a gross thoracic and lumbar spine immobilization in all three planes. Sizing is appropriate for most adults; the product is non-gender specific.

GENERAL SAFETY INSTRUCTIONS

Warning:

The Miami TLSO with anterior and posterior thoracic extensions is a supportive brace only, and is not intended or guaranteed to prevent spinal injury. It should be applied by a licensed healthcare practitioner and as directed by a physician, or other qualified medical authority. Össur accepts no responsibility for any injury that might occur while using the orthosis.

Precaution:

Do not overtighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system. If brace causes discomfort, irritation, or is damaged, contact your

healthcare provider for further instructions and discontinue use.

Make sure the hinged adjustment mechanism is tightened sufficiently to lock it into position.

The brace is single patient use only.

The brace is not designed to be worn while bathing.

INSTRUCTIONS FOR FITTING

MIAMI TLSO BRACE (ATE+PTE) PREPARATION:

No preparation is needed for over-the-shoulder strapping configuration.

For axillary strapping configuration:

- Remove strap from the posterior panel (PTE)
- Reattach the free end of the strap to the strap anchor (**Fig. 2**)

For over-the-shoulder and axillary strapping configuration (requires purchase of additional kit):

- Remove pectoral disc assemblies from the ball joints on the yoke (**Fig. 3**)
- Attach secondary buckles from accessory kit over the ball joints into the groove (**Fig. 4**)
- Reattach the pectoral disc assemblies onto the ball joints
- Attach the free end of the strap from the accessory kit to the strap anchors (**Fig. 2**)

OVER-THE-SHOULDER MIAMI TLSO BRACE (ATE+PTE) APPLICATION

1. Log roll the patient onto their side. Omit this step if fitting in standing position. If standing, patient must stand with back straight and chest out.
2. Position the belt by centering the posterior panel (PTE) over the spine, with bottom of posterior panel approximately at the sacrococcygeal joint. (**Fig. 5**).
3. Log roll the patient back into the supine position. Omit this step if fitting in standing position.
4. Wrap the lumbar belt around the waist and through Össur Fit Tool (**Fig. 6**).
5. Adjust length of belt arms to position tool centrally. Secure the hook and loop tabs to desired sizing (**Fig. 7**).
6. Remove Össur Fit Tool from belt.
Note: If an Össur Fit Tool is not available, use sizing guide (**Fig. 8**) to estimate best patient fit based on waist circumference. Simply fold or trim the belt arms at the suggested location. If patient's waist exceeds 50" (127 cm), an extension belt can be added to each belt arm to reach a maximum belt circumference of 70" (178 cm).
7. Place belt ends into the overlap closures (**Fig. 9**).
8. To obtain the lowest profile fit, the belt arm ends can be trimmed if later adjustments are not required.

Note:

- Once belt ends are secured in overlap closures, the overlap closures will serve as belt ends.
 - Panels can be customized by heat-molding, grinding, or trimming, to achieve correct patient fit.
9. Place the anterior thoracic extension (ATE) onto the patient's front, centering the anterior panel on the abdomen with the bottom edge just above the symphysis pubis, while still allowing the patient to sit

comfortably (**Fig. 10**).

10. Adjust the height of the ATE by pushing the button on the anterior panel and sliding the strut to the appropriate length. The superior edge of the pectoral pads should be located approximately one inch below the sternal notch (**Fig. 11**).
11. Using the provided tool, loosen the screw by one and a half turn and then adjust the angle on the hinged adjustment mechanism to accommodate sternal contours (**Fig. 12**).
Note: Make sure the hinged adjustment mechanism is tightened enough to lock it into position.
12. Fasten the non-labeled belt end onto the ATE panel.
13. Fasten the retaining straps onto the ATE panel soft good, over the non-labeled belt end (**Fig. 13**). Close the belt.
14. If necessary, adjust the compression system cord length (**Fig. 14**).
15. Connect both quick-release buckles, tighten the shoulder straps, and position the shoulder pad for optimal comfort (**Fig. 15**).
16. If desired, remove alligator, trim the strap, and replace alligator.
17. Slide fingers through the holes in the two compression system handles and pull until the belt is at the appropriate tightness. (**Fig. 18**).
18. The brace may require fine tuning for comfort, fit and stabilization when sitting, standing, or ambulating. This may include adjusting or tightening straps, compression system, rigid panels, and/or ATE/PTE.

AXILLARY STRAP ONLY OPTION / ADJUSTMENTS

- Replace step 14 with the following: fasten the strap anchor onto a smooth rearward area of the belt, then tighten strap (**Fig. 16**)
- Resume steps 15 and 16.

ATE+PTE DUAL STRAP OPTION / ADJUSTMENTS

- Before step 15, attach axillary strap to secondary buckles (**Fig. 16**)
- Fasten the strap anchor onto the belt and then tighten strap (**Fig. 17**)
- Resume steps 15 and 16.

HINGE/LENGTH ADJUSTMENTS

If brace needs to be taller (**Fig. 19**):

1. Remove hinge screw carefully with provided tool; do not lose the pre-load spring.
2. Attach extension piece to lower stay.
3. Attach yoke to extension piece.
4. Using the provided tool, adjust the angle on the hinged adjustment mechanism to accommodate sternal contours.
Note: Make sure the hinged adjustment mechanism is tightened sufficiently to lock it into position.

BRACE REMOVAL

1. Detach strap buckle(s) on right shoulder.
2. Detach compression system handles and loosen. Fix them onto the belt after loosening.
3. Detach the overlap closure and remove brace.
4. Before reapplying the brace and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart. Reattach right strap buckle before reapplying the brace.

PANEL ADJUSTMENTS

The posterior and anterior rigid panels can be modified as necessary to optimize patient fit and comfort. Simply remove the panels from their fabric sleeve and adjust the rigid panel with a heat gun and/or trimming device.

MAINTENANCE & CLEANING

The Miami TLSO system is designed for minimum maintenance and care. The padded covers may be removed from the rigid panels for washing.

- Hand-wash using mild soap and rinse thoroughly
- Air dry flat
- Do not machine wash or place in tumble dryer

To remove pectoral pads from the ATE for cleaning, pull and twist the pectoral pads until the ball snaps out of the socket. To reassemble, press the ball into the socket until it snaps back in place.

PRODUCT MATERIAL LIST

Nylon, glass-filled nylon, polyamide, spandex, thermoplastic polyurethane, polyurethane, ester polyurethane, polyester, polyurethane foam, bamboo, acrylonitrile butadiene styrene, steel, stainless steel, ethylene vinyl acetate, glass, polyethylene, low density polyethylene, aluminum.

ENVIRONMENTAL CONDITIONS

The device can be used in normal conditions (home or healthcare), both indoors and outdoors.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

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EN – Caution: This product has been designed and tested as a single use product, and is not intended for multiple use. In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure. If any problems occur with the use of this product, immediately contact your medical professional.

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