










Instructions for Use



i - LIMB[®] HAND

3

EN		Instructions for Use	5
DE		Gebrauchsanweisung	23
FR		Notice d'utilisation	41
ES		Instrucciones para el uso	59
IT		Istruzioni per l'uso	77
NO		Bruksanvisning	95
DA		Bruksanvisning	113
SV		Bruksanvisning	129
EL		Οδηγίες Χρήσης	146
FI		Käyttöohjeet	164
NL		Gebruiksaanwijzing	182
PT		Instruções de Utilização	200
PL		Instrukcja użytkowania	218
CS		Návod k použití	236
TR		Kullanım Talimatları	254
RU		Инструкция по использованию	272
JA		取扱説明書	290
ZH		中文说明书	307
KO		사용 설명서	324
AR		تعليمات الاستخدام	341

REGULATORY MARKINGS

Consult instructions for use

Class II equipment - provides double isolation to protect against electric shock

IP22

Protected against solid foreign objects of 12.5 mm diameter and greater and protected against splashing water.

Serial Number

The unique serial number for i-limb™ quantum devices is a "M" with a 6 digit alpha/ numeric number.

The unique serial number for i-limb™ ultra devices is a "U" with a 6 digit alpha/ numeric number.

The unique serial number for i-limb™ access devices is a "A" with a 6 digit alpha/ numeric number.

WEEE Compliance: This marking on the product, packaging, accessories or literature indicates that the product contains electronic components and/or batteries that should not be disposed of in regular waste at the end of its usable life. To prevent possible harm to the environment or human health from uncontrolled waste disposal users are requested to separate these items from other types of waste and recycle them responsibly to support sustainable reuse of material resources. Users should contact their local government office for information on how these items can be recycled or disposed of in an environmentally sound manner.

To protect natural resources and to promote material reuse, please separate batteries and electrical components from other types of waste and recycle them through your local, free electronic parts return system.

Manufacturer - YYYY

Manufacturer and year of manufacturing (YYYY)

European Conformity

BF Applied Part

INSTRUCTIONS FOR USE

I-LIMB HAND

The i-Limb® hand is referred to as “the device” in the following document. This document provides information on the indications for use, and handling of the device. It is intended for the user of the device. The device may only be configured and fitted by a qualified practitioner authorized by Össur after completing the corresponding training. These “Instructions for Use” relate to: i-Limb Quantum, i-Limb Ultra, i-Limb Access, analog electrodes, magnetic charger port, mains and car charger.

PRODUCT DESCRIPTION

The i-Limb hands are a range of prosthetic hands consisting of individual motorized digits, stall detection and a proprietary control app for Apple iOS devices (Fig. 1a). The product label can be found on the side of the device. (Fig. 1b). The device is assembled as part of your prosthesis together with a socket custom made for you by your clinician.

Users can choose from a selection of automated grips and gestures to help complete their daily tasks. Depending on model, grips can be customized and automated for additional control. An overview of available control options can be found in the features comparison table below:

Features comparison			
	i-Limb Quantum	i-Limb Ultra	i-Limb Access
Control options			
Gesture control	Yes	-	-
App control	Yes	Yes	Yes
Muscle control	Yes	Yes	Yes
Proximity control	Yes	-	-
Grips available	24	18	12
My grips	12	-	-
Speed boost	Yes	Yes	Yes
Powered thumb rotation	Yes	Yes	-
Vari-grip	Yes	Yes	-
Natural hand mode	Yes	Yes	-

While digits are not controlled individually, automated grips allow i-Limb hands users to move certain digits for pinching an object or make a gesture. Training is required to fully utilize the i-Limb hand and fully understand all the functional benefits.

INDICATIONS FOR USE

Individuals with upper limb absence.

CONTRAINDICATIONS

None known.

INTENDED USE

Individuals with upper limb absence with the potential to use an upper limb prosthesis.

REQUIRED DEVICES

The My i-Limb and Biosim apps can be downloaded from the Apple Store. The two apps require an Apple iOS device supported by the manufacturer, e.g. an iPhone or iPad. See the Apple Store for device compatibility.

SAFETY INSTRUCTIONS

Warnings

i-Limb Hand:

- The end user is the intended operator of the device and is responsible for its use.
- The i-Limb hand does not provide sensation, heat and moisture cannot be felt. The i-Limb hand is for low to moderate impact activities only.
- Do not use without an approved cover.
- Do not use with a damaged cover.
- Do not disassemble componentry or modify in any way.
- Do not service or perform maintenance when in use.
- Do not carry objects using only the tips of the digits. Carry objects by evenly distributing weight across the digits, as close to the knuckles and palm of the hand.
- Do not use for heavy lifting.
- Do not use with machinery with moving parts that may cause personal injury or damage.
- Do not use for extreme activities that may cause injury to a natural hand.
- Do not expose to vibration.
- Do not expose to excessive or high forces, particularly on the fingertips and on the side of the digits.
- Do not expose to water.
- Do not expose to excessive moisture, liquids, dust, high temperatures, or shock.
- Do not use in hazardous environments.
- Do not expose to high temperatures.
- Do not expose to flames.
- Do not use in or expose to explosive atmospheres.
- The electrode, magnetic charger port, DC port and switch block are an APPLIED PART.
- The electrode may contain nickel

Batteries:

- Do not cut or modify battery wires.
- Do not bend or exert excessive pressure on the battery.
- Do not pierce the battery.
- Do not disassemble battery.
- Do not expose to high temperatures.
- Do not incinerate batteries.
- Do not alter battery terminal wires.
- Do not short circuit the battery.
- Do not store batteries inside a vehicle.
- Dispose of batteries in accordance with US, European or local regulations.

Precautions

i-Limb Hand:

- Users must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind and any other motorized vehicle or device. It is entirely the user's responsibility to seek confirmation that they are physically and legally able to drive using the i-Limb hand and to the fullest extent permitted by law.
- Only use with approved Össur accessories and tooling.
- Damaged covers must be replaced or repaired by a qualified Össur technician or technical partner.
- Maintenance, repairs and upgrades may only be performed by qualified Össur technicians and technical partners. Össur will provide upon request information to assist service personnel in repair of device.
- Do not use an i-Limb device to operate electronic devices connected to a mains outlet, as this may affect EMG signal.

- It is not recommended to use your device whilst undergoing specific investigations or treatments when in close proximity of other medical electrical equipment.
- Do not use while charging is in progress.
- Only use with gloves supplied by Össur.
- Always use with gloves to avoid risk of electrostatic build-up and discharge.
- Do not use oil-based lotions on the skin, e.g. Vaseline.
- Do not expose electrode to dirt or fluids.

Batteries:

- Only Össur batteries should be used with this device.
- Only use the Össur charger to charge Össur batteries.
- Internal batteries must not be replaced by the end user.
- Batteries are to be replaced annually by service personnel only.

Do not use your i-Limb hand if the battery has visibly ballooned or swelled. Contact your clinician and do as follows:

- discontinue the charging process immediately
- disconnect the battery
- remove to a safe area
- leave and observe for 15 minutes
- replace with new battery
- do not re-use
- dispose of any leaking batteries in an appropriate manner

If the device will not be used for a longer period of time, it is advisable to remove the battery from the prosthesis. Contact your clinician for this to be completed.

SAFE USE GUIDELINES

- Push the thumb to the side. Hold handles or similarly shaped objects securely in palm near the base of the fingers. Move the thumb to the side so that it doesn't block the index from closing (**Fig. 2**).
- Do not block the index finger from closing around thinner objects. The object will be less secure within the i-Limb hand if the digits are not able to conform around its shape (**Fig. 3**).
- Hold objects close to palm with all fingers fully closed. Ensure all digits are fully closed around objects (**Fig. 4**).
- Do not hold objects with fingertips or side of the fingers (**Fig. 5**).
- Position objects close to the knuckles when pulling/pushing (**Fig. 6**).
- Do not pull or push objects with fingertips (**Fig. 7**).
- Push up with fully closed fist with force at knuckles (**Fig. 8**).
- Do not push up on fingers (**Fig. 9**).

POWER

A fully charged device can be used for up to 16 hours, depending on use.

To charge the i-Limb hand, remove it first from the residual limb. Depending on model, turn either off or to standby by moving power switch to the left position (**Fig. 10, 11**).

Magnetic Charging Port

The magnetic charging port allows the battery to be charged, the power switched on/off and the remaining battery level to be monitored in one place (**Fig. 12**).

If your i-Limb has an internal battery and supports magnetic charger port, please follow these steps:

Turning the Device On and Off

- Press switch on charge port for 1 second and release to turn on/off.
- When turned on: the display will fully light up and then dim to a low level.
- When turned off: the display will fully light up and then switch off.

Checking the Battery Charge Level Status

- The battery charge level is shown by the bars when the power is on.
- Each bar of light shows a 20% charge level. The number of bars lit up shows how much power is left in the battery.
- When the battery level reaches 5%, the red low battery warning will light up (**Fig. 13**). The light will shine for 3 minutes, and then the device will turn itself off.

NOTE: an emergency reserve of power remains available, to turn on and release the device from an object, ensuring the safety of the user.

Charging the battery

WARNING: Do not wear an i-Limb hand when the battery is charging.

- Plug the mains charger into the power outlet. Attach the magnetic charger to the charger port on the prosthesis.
- When in standby, the charger unit will display a faint green light.
- When charging is in progress, a red light is displayed.
- When fully charged a green light is displayed.
- Charging time: 90 minutes to 3 hours.
- A car charger is also available.

External Battery

- If you have external batteries, remove the batteries from your prosthesis and insert into the charger base unit (**Fig. 14**). Connect the charger base unit with the power cable. Plug the power cable into the power outlet.
- The charging state is indicated on the back of the base unit:
 - Middle light on: Charger is plugged in
 - 2nd and 5th lights flashing green: batteries are charging
 - 2nd and 5th lights solid green: batteries are charged
 - 1st and 4th red lights on: battery fault, unplug and try again. If lights continue to illuminate, contact your clinician.

Power Supply

CAUTION: Only use the Össur power supply to charge Össur batteries.

Manufacturer:	PowerSolve
Model No.:	000311A
Input:	100-240Vac, 50-60 Hz, 0.3A Max
Output:	8.4VDC, 1A

CONTROLLING THE I-LIMB HAND

Identifying the Device Number:

The device number can be found proximal to the thumb base (**Fig. 15**). For extra small hands, or those with a flexion wrist, the number is located within the chassis.

When connecting your i-Limb hand with the My i-Limb app, the device number will be displayed on the connection screen. Selecting the number connects the app to your i-Limb hand. Alternatively, when connected to the app your device number is displayed in the 'about' section.

Gesture Control (available only on i-Limb Quantum hands)

Enables an automated grip to be accessed through a smooth motion of the prosthesis in one of four directions (forward, back, left or right). The grips programmed for each direction are customized to the user's requirements using the My i-Limb app.

To access gesture control:

1. Hold i-Limb hand parallel to the ground (elbow bent to 90°)
2. Maintain an open signal until index finger twitches
3. Move the i-Limb Quantum hand (within 1 second) in direction assigned to the desired grip
4. The i-Limb Quantum hand will adopt the grip
5. To exit the grip, hold an open signal until the hand exits the grip

Default setting for gesture control is hold open, however it can alternatively be accessed using co-contraction. This is changed using the My i-Limb app.

App Control (available on all i-Limb hands)

An automated grip can be accessed at the touch of an icon on the My i-Limb app. These are called quick grips. The i-Limb hand will exit the grip when the icon is tapped again or by selecting another grip icon to enter.

Muscle Control (available on all i-Limb hands)

Triggers are specific muscle signals that can be used to access an automated grip. There are 4 potential muscle triggers: hold open, double impulse, triple impulse and co-contraction.

You can use the app to activate and program muscle control.

Proximity Control (available only on i-Limb Quantum hands)

Grip Chips are small Bluetooth devices which change the function of the i-Limb Quantum hand when the hand comes close to them.

To use Grip Chips, make sure the i-Limb hand is not connected to the My i-Limb app.

Proximity is accessed by holding the i-Limb hand fully opened near (15 cm / 6") to a Grip Chip. Wait until the grip is accessed, which may take up to 3 seconds. To exit a grip accessed via proximity control, give a long open signal.

Double tapping a Grip Chip activates the grip. Quickly tap a Grip Chip twice, just like double-clicking a computer mouse. The LED on the Grip Chip flashes once when tapping is successful.

NOTE: a 3 second pause between each double tap is required. The pause prevents the grip chip from incorrectly detecting multiple taps in a very short period of time. This can cause the hand to enter a grip and then immediately exit it again.

The i-Limb must be fully open/digits stalled open for proximity or tapping to be successful.

The individual Grip Chips are programmed using the My i-Limb app and can be reprogrammed by the user at any time.

I-LIMB HAND COVERS

Donning a Cover

To don a cover, bring i-Limb hand into the position shown in figure 16 and turn the hand off. The position can also be achieved using the don/doff quick grip.

Donning the i-Limb Skin Active, i-Limb Skin Contour and i-Limb Skin Natural covers:

1. Line up cover with the fingers of the i-Limb hand, and slide cover downwards.
2. When the fingers are mostly donned, pull the thumb opening over the thumb (**Fig. 17**).
3. Slide the remainder of the cover over the i-Limb hand (**Fig. 18**).
4. Make sure each digit tip is fitted to the cover.
5. Do not pull the cover too tightly over the hand.

Donning the i-Limb Skin Match cover:

1. Spray the outer surface with isopropyl alcohol (**Fig. 19**).
2. Invert the cover to the level of the finger openings (**Fig. 20**).

3. Ensure the fingers are straight and not bunched.
4. Align the cover finger holes with the digits of the hand (**Fig. 21**).
5. Pull down onto digits of the i-Limb hand.
6. Place thumb opening over thumb digit.
7. Pull the cover carefully over the remainder of the i-Limb hand, without applying too much pressure on the thumb (**Fig. 22**).
8. Check the cover for bunching and make sure the cover tips are fully fitted to the digit tips.
9. Check the function of the i-Limb hand and ensure full opening and closing is possible and the digit tips align.

Doffing all Covers

1. Bring the i-Limb hand into the same position as for donning and switch off/ to standby.
2. Pull the cover upwards on each digit to release.
3. Carefully ease the full cover off, without applying too much pressure on the thumb.
4. Continue to pull cover upwards, until fully removed.

CLEANING

Cleaning instructions

The i-Limb hand covering, electrode surface and magnetic charger port can only be cleaned with a soft damp cloth and mild soap.

Clean the electrode surface regularly.

Clean the cover with isopropyl alcohol once a week to help with disinfection.

NOTE: do not submerge the i-Limb hand, cover, electrode surface or magnetic charger port in water. Do not use any strong chemicals.

TECHNICAL SPECIFICATIONS

Voltage	Range 6 – 8.4 V
Max. Current	6 A
Battery Capacity	Rechargeable lithium polymer 7.4. V (nominal); 2000 mAh capacity; 1,300 mAh capacity
Max hand load limit (static limit)	40 kg/ 88 lb. (extra small) 90 kg/ 198 lb. (small/ medium/ large)
Finger carry load (static limit)	20 kg/ 44 lb. (extra small) 32 kg/ 71 lb. (small/ medium/ large)
Time from full open to full close	0.8 seconds
Expected service life	5 years

Device weight (Quantum and Ultra) Note: Titanium digits add an additional 30g/ 0.07lb. per hand		Extra small	Small	Medium/ large
	QWD	472 g/ 1.04 lb.	512 g/ 1.13 lb.	528 g/ 1.16 lb.
	WD	432 g/ 0.95 lb.	472 g/ 1.04 lb.	488 g/ 1.08 lb.
	Flexion	572 g/ 1.26 lb.	612 g/ 1.35 lb.	628 g/ 1.38 lb.
	Friction	467 g/ 1.03 lb.	507 g/ 1.12 lb.	523 g/ 1.15 lb.

Device weight (Access) Note: Titanium digits add an additional 30 g/ 0.07 lb. per hand		Extra small	Small	Medium/ large
	QWD	432 g/ 0.95 lb.	468 g/ 1.03 lb.	478 g/ 1.05 lb.
	WD	392 g/ 0.86 lb.	428 g/ 0.94 lb.	438 g/ 0.97 lb.
	Flexion	532 g/ 1.17 lb.	568 g/ 1.25 lb.	578 g/ 1.27 lb.
	Friction	427 g/ 0.94 lb.	463 g/ 1.02 lb.	473 g/ 1.04 lb.

LIST OF PARTS

i-Limb Quantum	
Code	Description
PL700000A	i-Limb Quantum QWD RHXS
PL701000A	i-Limb Quantum QWD LHXS
PL702000A	i-Limb Quantum QWD RHS
PL703000A	i-Limb Quantum QWD LHS
PL704000A	i-Limb Quantum QWD RHM
PL705000A	i-Limb Quantum QWD LHM
PL706000A	i-Limb Quantum QWD RHL
PL707000A	i-Limb Quantum QWD LHL
PL708000A	i-Limb Quantum Flexion RHXS
PL709000A	i-Limb Quantum Flexion LHXS
PL710000A	i-Limb Quantum Flexion RHS
PL711000A	i-Limb Quantum Flexion LHS
PL712000A	i-Limb Quantum Flexion RHM
PL713000A	i-Limb Quantum Flexion LHM
PL714000A	i-Limb Quantum Flexion RHL
PL715000A	i-Limb Quantum Flexion LHL
PL716000A	i-Limb Quantum WD RHXS
PL717000A	i-Limb Quantum WD LHXS
PL718000A	i-Limb Quantum WD RHS
PL719000A	i-Limb Quantum WD LHS
PL720000A	i-Limb Quantum WD RHM
PL721000A	i-Limb Quantum WD LHM
PL722000A	i-Limb Quantum WD RHL
PL723000A	i-Limb Quantum WD LHL
PL724000A	i-Limb Quantum Friction RHXS
PL725000A	i-Limb Quantum Friction LHXS
PL726000A	i-Limb Quantum Friction RHS
PL727000A	i-Limb Quantum Friction LHS

PL728000A	i-Limb Quantum Friction RHM
PL729000A	i-Limb Quantum Friction LHM
PL730000A	i-Limb Quantum Friction RHL
PL731000A	i-Limb Quantum Friction LHL
PL792000A	i-Limb Quantum Ti QWD RHS
PL793000A	i-Limb Quantum Ti QWD LHS
PL794000A	i-Limb Quantum Ti QWD RHM
PL795000A	i-Limb Quantum Ti QWD LHM
PL796000A	i-Limb Quantum Ti QWD RHL
PL797000A	i-Limb Quantum Ti QWD LHL
PL798000A	i-Limb Quantum Ti Flexion RHS
PL799000A	i-Limb Quantum Ti Flexion LHS
PL800000A	i-Limb Quantum Ti Flexion RHM
PL801000A	i-Limb Quantum Ti Flexion LHM
PL802000A	i-Limb Quantum Ti Flexion RHL
PL803000A	i-Limb Quantum Ti Flexion LHL
PL804000A	i-Limb Quantum Ti WD RHS
PL805000A	i-Limb Quantum Ti WD LHS
PL806000A	i-Limb Quantum Ti WD RHM
PL807000A	i-Limb Quantum Ti WD LHM
PL808000A	i-Limb Quantum Ti WD RHL
PL809000A	i-Limb Quantum Ti WD LHL
PL810000A	i-Limb Quantum Ti Friction RHS
PL811000A	i-Limb Quantum Ti Friction LHS
PL812000A	i-Limb Quantum Ti Friction RHM
PL813000A	i-Limb Quantum Ti Friction LHM
PL814000A	i-Limb Quantum Ti Friction RHL
PL815000A	i-Limb Quantum Ti Friction LHL

i-Limb Ultra

Code	Description
PL732000A	i-Limb UltraPT QWD RHXS
PL733000A	i-Limb UltraPT QWD LHXS
PL734000A	i-Limb UltraPT QWD RHS
PL735000A	i-Limb UltraPT QWD LHS
PL736000A	i-Limb UltraPT QWD RHM
PL737000A	i-Limb UltraPT QWD LHM
PL738000A	i-Limb UltraPT QWD RHL
PL739000A	i-Limb UltraPT QWD LHL
PL740000A	i-Limb UltraPT Flexion RHXS

PL741000A	i-Limb UltraPT Flexion LHXS
PL742000A	i-Limb UltraPT Flexion RHS
PL743000A	i-Limb UltraPT Flexion LHS
PL744000A	i-Limb UltraPT Flexion RHM
PL745000A	i-Limb UltraPT Flexion LHM
PL746000A	i-Limb UltraPT Flexion RHL
PL747000A	i-Limb UltraPT Flexion LHL
PL748000A	i-Limb UltraPT WD RHXS
PL749000A	i-Limb UltraPT WD LHXS
PL750000A	i-Limb UltraPT WD RHS
PL751000A	i-Limb UltraPT WD LHS
PL752000A	i-Limb UltraPT WD RHM
PL753000A	i-Limb UltraPT WD LHM
PL754000A	i-Limb UltraPT WD RHL
PL755000A	i-Limb UltraPT WD LHL
PL756000A	i-Limb UltraPT Friction RHXS
PL757000A	i-Limb UltraPT Friction LHXS
PL758000A	i-Limb UltraPT Friction RHS
PL759000A	i-Limb UltraPT Friction LHS
PL760000A	i-Limb UltraPT Friction RHM
PL761000A	i-Limb UltraPT Friction LHM
PL762000A	i-Limb UltraPT Friction RHL
PL763000A	i-Limb UltraPT Friction LHL
PL816000A	i-Limb UltraPT Ti QWD RHS
PL817000A	i-Limb UltraPT Ti QWD LHS
PL818000A	i-Limb UltraPT Ti QWD RHM
PL819000A	i-Limb UltraPT Ti QWD LHM
PL820000A	i-Limb UltraPT Ti QWD RHL
PL821000A	i-Limb UltraPT Ti QWD LHL
PL822000A	i-Limb UltraPT Ti Flexion RHS
PL823000A	i-Limb UltraPT Ti Flexion LHS
PL824000A	i-Limb UltraPT Ti Flexion RHM
PL825000A	i-Limb UltraPT Ti Flexion LHM
PL826000A	i-Limb UltraPT Ti Flexion RHL
PL827000A	i-Limb UltraPT Ti Flexion LHL
PL828000A	i-Limb UltraPT Ti WD RHS
PL829000A	i-Limb UltraPT Ti WD LHS
PL830000A	i-Limb UltraPT Ti WD RHM
PL831000A	i-Limb UltraPT Ti WD LHM

PL832000A	i-Limb UltraPT Ti WD RHL
PL833000A	i-Limb UltraPT Ti WD LHL
PL834000A	i-Limb UltraPT Ti Friction RHS
PL835000A	i-Limb UltraPT Ti Friction LHS
PL836000A	i-Limb UltraPT Ti Friction RHM
PL837000A	i-Limb UltraPT Ti Friction LHM
PL838000A	i-Limb UltraPT Ti Friction RHL
PL839000A	i-Limb UltraPT Ti Friction LHL
i-Limb Access	
Code	Description
PL764000A	i-Limb Access QWD RHXS
PL765000A	i-Limb Access QWD LHXS
PL766000A	i-Limb Access QWD RHS
PL767000A	i-Limb Access QWD LHS
PL768000A	i-Limb Access QWD RHM
PL769000A	i-Limb Access QWD LHM
PL770000A	i-Limb Access QWD RHL
PL771000A	i-Limb Access QWD LHL
PL772000A	i-Limb Access Flexion RHXS
PL773000A	i-Limb Access Flexion LHXS
PL774000A	i-Limb Access Flexion RHM
PL775000A	i-Limb Access Flexion LHM
PL776000A	i-Limb Access Flexion RHL
PL777000A	i-Limb Access Flexion LHL
PL778000A	i-Limb Access WD RHXS
PL779000A	i-Limb Access WD LHXS
PL780000A	i-Limb Access WD RHS
PL781000A	i-Limb Access WD LHS
PL782000A	i-Limb Access WD RHM
PL783000A	i-Limb Access WD LHM
PL784000A	i-Limb Access WD RHL
PL785000A	i-Limb Access WD LHL
PL786000A	i-Limb Access Friction RHS
PL787000A	i-Limb Access Friction LHS
PL788000A	i-Limb Access Friction RHM
PL789000A	i-Limb Access Friction LHM
PL790000A	i-Limb Access Friction RHL
PL791000A	i-Limb Access Friction LHL
PL840000A	i-Limb Access Ti QWD RHS

PL841000A	i-Limb Access Ti QWD LHS
PL842000A	i-Limb Access Ti QWD RHM
PL843000A	i-Limb Access Ti QWD LHM
PL844000A	i-Limb Access Ti QWD RHL
PL845000A	i-Limb Access Ti QWD LHL
PL846000A	i-Limb Access Ti Flexion RHS
PL847000A	i-Limb Access Ti Flexion LHS
PL848000A	i-Limb Access Ti Flexion RHM
PL849000A	i-Limb Access Ti Flexion LHM
PL850000A	i-Limb Access Ti Flexion RHL
PL851000A	i-Limb Access Ti Flexion LHL
PL852000A	i-Limb Access Ti WD RHS
PL853000A	i-Limb Access Ti WD LHS
PL854000A	i-Limb Access Ti WD RHM
PL855000A	i-Limb Access Ti WD LHM
PL856000A	i-Limb Access Ti WD RHL
PL857000A	i-Limb Access Ti WD LHL
PL858000A	i-Limb Access Flexion RHS
PL859000A	i-Limb Access Flexion LHS
PL860000A	i-Limb Access Friction RHXS
PL861000A	i-Limb Access Friction LHXS

LIST OF MAIN COMPONENTS

Part Code	Description
PL000290A	Magnetic Charger Base
PL000291A	Magnetic Charger Plug
PL000292A	Magnetic Charger Car
SA000192A	Battery 1300 mAh
SA000235A	Battery 2000 mAh
PL091050A	Electrodes Compact 50 Hz 300 mm
PL091060A	Electrodes Compact 60 Hz 300 mm
PL091127A	Electrodes Compact 50 Hz 600 mm
PL091128A	Electrodes Compact 60 Hz 600 mm
PL069466B	Electrodes Remote 50 Hz 300 mm
PL069467B	Electrodes Remote 60 Hz 300 mm
PL069468B	Electrodes Remote 50 Hz 600 mm
PL069469B	Electrodes Remote 60 Hz 600 mm
PL544008A	Digital Electrode Remote

PL544009A	Digital Electrode Compact
092023A	i-limb skin active ultra Medium RH Clear
092024A	i-limb skin active ultra Medium LH Clear
092025A	i-limb skin active ultra Small RH Clear
092026A	i-limb skin active ultra Small LH Clear
092027A	i-limb skin active ultra Medium RH Black
092028A	i-limb skin active ultra Medium LH Black
092029A	i-limb skin active ultra Small RH Black
092030A	i-limb skin active ultra Small LH Black
PL353000A	i-limb skin active RH Extra Small Black
PL353001A	i-limb skin active LH Extra Small Black
PL353002A	i-limb skin active RH Extra Small Clear
PL353003A	i-limb skin active LH Extra Small Clear
PL353004A	i-limb skin contour RH Extra Small Clear
PL353005A	i-limb skin contour LH Extra Small Clear
PL353014A	i-limb skin contour RH Extra small Black
PL353015A	i-limb skin contour LH Extra small Black
PL354004A	i-limb skin contour RH Small Clear
PL354005A	i-limb skin contour LH Small Clear
PL354014A	i-limb skin contour RH small Black
PL354015A	i-limb skin contour LH small Black
PL355004A	i-limb skin contour RH Medium Clear
PL355005A	i-limb skin contour LH Medium Clear
PL355014A	i-limb skin contour RH medium Black
PL355015A	i-limb skin contour LH medium Black

ENVIRONMENTAL CONDITIONS

Do not use, transport or store i-Limb hand outside of the boundaries in the table below:

	Use	Shipping	Extended Storage
Temperature	-0 °C to +40 °C	-25 °C + 70 °C	-25 °C + 70 °C
Relative Humidity	10 % to 100 %	10 % to 100 %	10 % to 100 %
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Re	Type and Frequency Characteristics	Effective Radiated Power
Bluetooth Low Energy Dual Mode Module Model BR-LE4.0-D2A	FCC Contains FCC ID: XDULE40-D2 Canada Contains IC: 8456A-LE4D2 Japan Contains transmitter with certificate number  R 205-160268	(Dual Mode) Version V2.1 +ED (GFSK + $\pi/4$ DQPSK + 8DPSK) 2402 -2480 MHz Version V4.0 (GFSK) 2402 -2480 MHz	Adjustable Power (-23 dBm to 10.5 dBm) short to long range

ELECTROMAGNETIC COMPATIBILITY

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation

In order to regulate the requirements for Electromagnetic Compatibility (EMC) with the aim to prevent unsafe product situations, the BS EN 60601-1-2:2015/ IEC 60601-1-2: 2015 standard has been implemented. This standard defines levels of electromagnetic emissions for medical devices.

The Össur Myoelectric Prosthetic Devices are manufactured by Össur hf. conforms to this BS EN 60601-1- 2:2015/ IEC 60601-1-2: 2015 standard for both immunity and emissions.

The i-Limb hand is suitable for use in any environment except where immersion in water or any other fluid is possible, or where exposure to highly electrical and/or magnetic fields can occur (e. g. electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Refer to further guidance below regarding the EMC environment in which the device should be used:

Guidance and manufacturer's declaration – electromagnetic emissions		
Össur Myoelectric Prosthetic Devices is intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that it each are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Not applicable Battery Powered	Össur Myoelectric Prosthetic Devices use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Össur Myoelectric Prosthetic Devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low- voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable Battery Powered	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable Battery Powered	

Guidance and manufacturer's declaration – electromagnetic immunity			
Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >3 m
Surge IEC 61000-4-5	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered
Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >30 m
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m 50/ 60 Hz	Battery Powered

Guidance and manufacturer's declaration – electromagnetic immunity			
Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customers or the users of Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable	Not applicable Battery Powered No Cables >3 m	Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter. Recommended separation distance d = 1.2 √P d = 1.2 √P 80 MHz to 800 MHz d = 2.3 √P 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	12 V/m 26 MHz to 1000 MHz 10 V/m 1000 MHz to 2700 MHz	12 V/m 26 MHz to 1000 MHz 10 V/m 1000 MHz to 2700 MHz 1kHz 80% AM	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Össur Myoelectric Prosthetic Devices are used exceeds the applicable RF compliance level above, Össur Myoelectric Prosthetic Devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Össur Myoelectric Prosthetic Devices			

Recommended separation distance between portable and mobile RF communications equipment and the Össur Myoelectric Prosthetic Devices			
Össur Myoelectric Prosthetic Devices are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the user of Össur Myoelectric Prosthetic Devices can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Össur Myoelectric Prosthetic Devices as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter in Watt	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7GHz
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

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PL809000A	i-Limb Quantum Ti WD LHL
PL810000A	i-Limb Quantum Ti Friction RHS
PL811000A	i-Limb Quantum Ti Friction LHS
PL812000A	i-Limb Quantum Ti Friction RHM
PL813000A	i-Limb Quantum Ti Friction LHM
PL814000A	i-Limb Quantum Ti Friction RHL
PL815000A	i-Limb Quantum Ti Friction LHL
i-Limb Ultra	
Code	Description
PL732000A	i-Limb UltraPT QWD RHXS
PL733000A	i-Limb UltraPT QWD LHXS
PL734000A	i-Limb UltraPT QWD RHS
PL735000A	i-Limb UltraPT QWD LHS
PL736000A	i-Limb UltraPT QWD RHM
PL737000A	i-Limb UltraPT QWD LHM
PL738000A	i-Limb UltraPT QWD RHL
PL739000A	i-Limb UltraPT QWD LHL
PL740000A	i-Limb UltraPT Flexion RHXS
PL741000A	i-Limb UltraPT Flexion LHXS
PL742000A	i-Limb UltraPT Flexion RHS
PL743000A	i-Limb UltraPT Flexion LHS
PL744000A	i-Limb UltraPT Flexion RHM
PL745000A	i-Limb UltraPT Flexion LHM
PL746000A	i-Limb UltraPT Flexion RHL
PL747000A	i-Limb UltraPT Flexion LHL
PL748000A	i-Limb UltraPT WD RHXS
PL749000A	i-Limb UltraPT WD LHXS
PL750000A	i-Limb UltraPT WD RHS
PL751000A	i-Limb UltraPT WD LHS
PL752000A	i-Limb UltraPT WD RHM
PL753000A	i-Limb UltraPT WD LHM
PL754000A	i-Limb UltraPT WD RHL
PL755000A	i-Limb UltraPT WD LHL
PL756000A	i-Limb UltraPT Friction RHXS
PL757000A	i-Limb UltraPT Friction LHXS
PL758000A	i-Limb UltraPT Friction RHS
PL759000A	i-Limb UltraPT Friction LHS
PL760000A	i-Limb UltraPT Friction RHM
PL761000A	i-Limb UltraPT Friction LHM

PL762000A	i-Limb UltraPT Friction RHL
PL763000A	i-Limb UltraPT Friction LHL
PL816000A	i-Limb UltraPT Ti QWD RHS
PL817000A	i-Limb UltraPT Ti QWD LHS
PL818000A	i-Limb UltraPT Ti QWD RHM
PL819000A	i-Limb UltraPT Ti QWD LHM
PL820000A	i-Limb UltraPT Ti QWD RHL
PL821000A	i-Limb UltraPT Ti QWD LHL
PL822000A	i-Limb UltraPT Ti Flexion RHS
PL823000A	i-Limb UltraPT Ti Flexion LHS
PL824000A	i-Limb UltraPT Ti Flexion RHM
PL825000A	i-Limb UltraPT Ti Flexion LHM
PL826000A	i-Limb UltraPT Ti Flexion RHL
PL827000A	i-Limb UltraPT Ti Flexion LHL
PL828000A	i-Limb UltraPT Ti WD RHS
PL829000A	i-Limb UltraPT Ti WD LHS
PL830000A	i-Limb UltraPT Ti WD RHM
PL831000A	i-Limb UltraPT Ti WD LHM
PL832000A	i-Limb UltraPT Ti WD RHL
PL833000A	i-Limb UltraPT Ti WD LHL
PL834000A	i-Limb UltraPT Ti Friction RHS
PL835000A	i-Limb UltraPT Ti Friction LHS
PL836000A	i-Limb UltraPT Ti Friction RHM
PL837000A	i-Limb UltraPT Ti Friction LHM
PL838000A	i-Limb UltraPT Ti Friction RHL
PL839000A	i-Limb UltraPT Ti Friction LHL
i-Limb Access	
Code	Description
PL764000A	i-Limb Access QWD RHXS
PL765000A	i-Limb Access QWD LHXS
PL766000A	i-Limb Access QWD RHS
PL767000A	i-Limb Access QWD LHS
PL768000A	i-Limb Access QWD RHM
PL769000A	i-Limb Access QWD LHM
PL770000A	i-Limb Access QWD RHL
PL771000A	i-Limb Access QWD LHL
PL772000A	i-Limb Access Flexion RHXS
PL773000A	i-Limb Access Flexion LHXS
PL774000A	i-Limb Access Flexion RHM

PL775000A	i-Limb Access Flexion LHM
PL776000A	i-Limb Access Flexion RHL
PL777000A	i-Limb Access Flexion LHL
PL778000A	i-Limb Access WD RHXS
PL779000A	i-Limb Access WD LHXS
PL780000A	i-Limb Access WD RHS
PL781000A	i-Limb Access WD LHS
PL782000A	i-Limb Access WD RHM
PL783000A	i-Limb Access WD LHM
PL784000A	i-Limb Access WD RHL
PL785000A	i-Limb Access WD LHL
PL786000A	i-Limb Access Friction RHS
PL787000A	i-Limb Access Friction LHS
PL788000A	i-Limb Access Friction RHM
PL789000A	i-Limb Access Friction LHM
PL790000A	i-Limb Access Friction RHL
PL791000A	i-Limb Access Friction LHL
PL840000A	i-Limb Access Ti QWD RHS
PL841000A	i-Limb Access Ti QWD LHS
PL842000A	i-Limb Access Ti QWD RHM
PL843000A	i-Limb Access Ti QWD LHM
PL844000A	i-Limb Access Ti QWD RHL
PL845000A	i-Limb Access Ti QWD LHL
PL846000A	i-Limb Access Ti Flexion RHS
PL847000A	i-Limb Access Ti Flexion LHS
PL848000A	i-Limb Access Ti Flexion RHM
PL849000A	i-Limb Access Ti Flexion LHM
PL850000A	i-Limb Access Ti Flexion RHL
PL851000A	i-Limb Access Ti Flexion LHL
PL852000A	i-Limb Access Ti WD RHS
PL853000A	i-Limb Access Ti WD LHS
PL854000A	i-Limb Access Ti WD RHM
PL855000A	i-Limb Access Ti WD LHM
PL856000A	i-Limb Access Ti WD RHL
PL857000A	i-Limb Access Ti WD LHL
PL858000A	i-Limb Access Flexion RHS
PL859000A	i-Limb Access Flexion LHS
PL860000A	i-Limb Access Friction RHXS
PL861000A	i-Limb Access Friction LHXS

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PL809000A	i-Limb Quantum Ti WD LHL
PL810000A	i-Limb Quantum Ti Friction RHS
PL811000A	i-Limb Quantum Ti Friction LHS
PL812000A	i-Limb Quantum Ti Friction RHM
PL813000A	i-Limb Quantum Ti Friction LHM
PL814000A	i-Limb Quantum Ti Friction RHL
PL815000A	i-Limb Quantum Ti Friction LHL
i-Limb Ultra	
Codice	Descrizione
PL732000A	i-Limb UltraPT QWD RHXS
PL733000A	i-Limb UltraPT QWD LHXS
PL734000A	i-Limb UltraPT QWD RHS
PL735000A	i-Limb UltraPT QWD LHS
PL736000A	i-Limb UltraPT QWD RHM
PL737000A	i-Limb UltraPT QWD LHM
PL738000A	i-Limb UltraPT QWD RHL
PL739000A	i-Limb UltraPT QWD LHL
PL740000A	i-Limb UltraPT Flexion RHXS
PL741000A	i-Limb UltraPT Flexion LHXS
PL742000A	i-Limb UltraPT Flexion RHS
PL743000A	i-Limb UltraPT Flexion LHS
PL744000A	i-Limb UltraPT Flexion RHM
PL745000A	i-Limb UltraPT Flexion LHM
PL746000A	i-Limb UltraPT Flexion RHL
PL747000A	i-Limb UltraPT Flexion LHL
PL748000A	i-Limb UltraPT WD RHXS
PL749000A	i-Limb UltraPT WD LHXS
PL750000A	i-Limb UltraPT WD RHS
PL751000A	i-Limb UltraPT WD LHS
PL752000A	i-Limb UltraPT WD RHM
PL753000A	i-Limb UltraPT WD LHM
PL754000A	i-Limb UltraPT WD RHL
PL755000A	i-Limb UltraPT WD LHL
PL756000A	i-Limb UltraPT Friction RHXS
PL757000A	i-Limb UltraPT Friction LHXS
PL758000A	i-Limb UltraPT Friction RHS
PL759000A	i-Limb UltraPT Friction LHS
PL760000A	i-Limb UltraPT Friction RHM
PL761000A	i-Limb UltraPT Friction LHM

PL762000A	i-Limb UltraPT Friction RHL
PL763000A	i-Limb UltraPT Friction LHL
PL816000A	i-Limb UltraPT Ti QWD RHS
PL817000A	i-Limb UltraPT Ti QWD LHS
PL818000A	i-Limb UltraPT Ti QWD RHM
PL819000A	i-Limb UltraPT Ti QWD LHM
PL820000A	i-Limb UltraPT Ti QWD RHL
PL821000A	i-Limb UltraPT Ti QWD LHL
PL822000A	i-Limb UltraPT Ti Flexion RHS
PL823000A	i-Limb UltraPT Ti Flexion LHS
PL824000A	i-Limb UltraPT Ti Flexion RHM
PL825000A	i-Limb UltraPT Ti Flexion LHM
PL826000A	i-Limb UltraPT Ti Flexion RHL
PL827000A	i-Limb UltraPT Ti Flexion LHL
PL828000A	i-Limb UltraPT Ti WD RHS
PL829000A	i-Limb UltraPT Ti WD LHS
PL830000A	i-Limb UltraPT Ti WD RHM
PL831000A	i-Limb UltraPT Ti WD LHM
PL832000A	i-Limb UltraPT Ti WD RHL
PL833000A	i-Limb UltraPT Ti WD LHL
PL834000A	i-Limb UltraPT Ti Friction RHS
PL835000A	i-Limb UltraPT Ti Friction LHS
PL836000A	i-Limb UltraPT Ti Friction RHM
PL837000A	i-Limb UltraPT Ti Friction LHM
PL838000A	i-Limb UltraPT Ti Friction RHL
PL839000A	i-Limb UltraPT Ti Friction LHL

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Colocar as coberturas i-Limb Skin Active, i-Limb Skin Contour e i-Limb Skin Natural:

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i-Limb Access	
Código	Descrição
PL764000A	i-Limb Access QWD RHXS
PL765000A	i-Limb Access QWD LHXS
PL766000A	i-Limb Access QWD RHS
PL767000A	i-Limb Access QWD LHS
PL768000A	i-Limb Access QWD RHM
PL769000A	i-Limb Access QWD LHM
PL770000A	i-Limb Access QWD RHL
PL771000A	i-Limb Access QWD LHL
PL772000A	i-Limb Access Flexion RHXS
PL773000A	i-Limb Access Flexion LHXS

PL774000A	i-Limb Access Flexion RHM
PL775000A	i-Limb Access Flexion LHM
PL776000A	i-Limb Access Flexion RHL
PL777000A	i-Limb Access Flexion LHL
PL778000A	i-Limb Access WD RHXS
PL779000A	i-Limb Access WD LHXS
PL780000A	i-Limb Access WD RHS
PL781000A	i-Limb Access WD LHS
PL782000A	i-Limb Access WD RHM
PL783000A	i-Limb Access WD LHM
PL784000A	i-Limb Access WD RHL
PL785000A	i-Limb Access WD LHL
PL786000A	i-Limb Access Friction RHS
PL787000A	i-Limb Access Friction LHS
PL788000A	i-Limb Access Friction RHM
PL789000A	i-Limb Access Friction LHM
PL790000A	i-Limb Access Friction RHL
PL791000A	i-Limb Access Friction LHL
PL840000A	i-Limb Access Ti QWD RHS
PL841000A	i-Limb Access Ti QWD LHS
PL842000A	i-Limb Access Ti QWD RHM
PL843000A	i-Limb Access Ti QWD LHM
PL844000A	i-Limb Access Ti QWD RHL
PL845000A	i-Limb Access Ti QWD LHL
PL846000A	i-Limb Access Ti Flexion RHS
PL847000A	i-Limb Access Ti Flexion LHS
PL848000A	i-Limb Access Ti Flexion RHM
PL849000A	i-Limb Access Ti Flexion LHM
PL850000A	i-Limb Access Ti Flexion RHL
PL851000A	i-Limb Access Ti Flexion LHL
PL852000A	i-Limb Access Ti WD RHS
PL853000A	i-Limb Access Ti WD LHS
PL854000A	i-Limb Access Ti WD RHM
PL855000A	i-Limb Access Ti WD LHM
PL856000A	i-Limb Access Ti WD RHL
PL857000A	i-Limb Access Ti WD LHL
PL858000A	i-Limb Access Flexion RHS
PL859000A	i-Limb Access Flexion LHS
PL860000A	i-Limb Access Friction RHXS

PL861000A	i-Limb Access Friction LHXS
-----------	-----------------------------

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PL809000A	i-Limb Quantum Ti WD LHL
PL810000A	i-Limb Quantum Ti Friction RHS
PL811000A	i-Limb Quantum Ti Friction LHS
PL812000A	i-Limb Quantum Ti Friction RHM
PL813000A	i-Limb Quantum Ti Friction LHM
PL814000A	i-Limb Quantum Ti Friction RHL
PL815000A	i-Limb Quantum Ti Friction LHL
i-Limb Ultra	
Kód	Popis
PL732000A	i-Limb UltraPT QWD RHXS
PL733000A	i-Limb UltraPT QWD LHXS
PL734000A	i-Limb UltraPT QWD RHS
PL735000A	i-Limb UltraPT QWD LHS
PL736000A	i-Limb UltraPT QWD RHM
PL737000A	i-Limb UltraPT QWD LHM
PL738000A	i-Limb UltraPT QWD RHL
PL739000A	i-Limb UltraPT QWD LHL
PL740000A	i-Limb UltraPT Flexion RHXS
PL741000A	i-Limb UltraPT Flexion LHXS
PL742000A	i-Limb UltraPT Flexion RHS
PL743000A	i-Limb UltraPT Flexion LHS
PL744000A	i-Limb UltraPT Flexion RHM
PL745000A	i-Limb UltraPT Flexion LHM
PL746000A	i-Limb UltraPT Flexion RHL
PL747000A	i-Limb UltraPT Flexion LHL
PL748000A	i-Limb UltraPT WD RHXS
PL749000A	i-Limb UltraPT WD LHXS
PL750000A	i-Limb UltraPT WD RHS
PL751000A	i-Limb UltraPT WD LHS
PL752000A	i-Limb UltraPT WD RHM
PL753000A	i-Limb UltraPT WD LHM
PL754000A	i-Limb UltraPT WD RHL
PL755000A	i-Limb UltraPT WD LHL
PL756000A	i-Limb UltraPT Friction RHXS
PL757000A	i-Limb UltraPT Friction LHXS
PL758000A	i-Limb UltraPT Friction RHS
PL759000A	i-Limb UltraPT Friction LHS
PL760000A	i-Limb UltraPT Friction RHM
PL761000A	i-Limb UltraPT Friction LHM

PL762000A	i-Limb UltraPT Friction RHL
PL763000A	i-Limb UltraPT Friction LHL
PL816000A	i-Limb UltraPT Ti QWD RHS
PL817000A	i-Limb UltraPT Ti QWD LHS
PL818000A	i-Limb UltraPT Ti QWD RHM
PL819000A	i-Limb UltraPT Ti QWD LHM
PL820000A	i-Limb UltraPT Ti QWD RHL
PL821000A	i-Limb UltraPT Ti QWD LHL
PL822000A	i-Limb UltraPT Ti Flexion RHS
PL823000A	i-Limb UltraPT Ti Flexion LHS
PL824000A	i-Limb UltraPT Ti Flexion RHM
PL825000A	i-Limb UltraPT Ti Flexion LHM
PL826000A	i-Limb UltraPT Ti Flexion RHL
PL827000A	i-Limb UltraPT Ti Flexion LHL
PL828000A	i-Limb UltraPT Ti WD RHS
PL829000A	i-Limb UltraPT Ti WD LHS
PL830000A	i-Limb UltraPT Ti WD RHM
PL831000A	i-Limb UltraPT Ti WD LHM
PL832000A	i-Limb UltraPT Ti WD RHL
PL833000A	i-Limb UltraPT Ti WD LHL
PL834000A	i-Limb UltraPT Ti Friction RHS
PL835000A	i-Limb UltraPT Ti Friction LHS
PL836000A	i-Limb UltraPT Ti Friction RHM
PL837000A	i-Limb UltraPT Ti Friction LHM
PL838000A	i-Limb UltraPT Ti Friction RHL
PL839000A	i-Limb UltraPT Ti Friction LHL
i-Limb Access	
Kód	Popis
PL764000A	i-Limb Access QWD RHXS
PL765000A	i-Limb Access QWD LHXS
PL766000A	i-Limb Access QWD RHS
PL767000A	i-Limb Access QWD LHS
PL768000A	i-Limb Access QWD RHM
PL769000A	i-Limb Access QWD LHM
PL770000A	i-Limb Access QWD RHL
PL771000A	i-Limb Access QWD LHL
PL772000A	i-Limb Access Flexion RHXS
PL773000A	i-Limb Access Flexion LHXS
PL774000A	i-Limb Access Flexion RHM

PL775000A	i-Limb Access Flexion LHM
PL776000A	i-Limb Access Flexion RHL
PL777000A	i-Limb Access Flexion LHL
PL778000A	i-Limb Access WD RHXS
PL779000A	i-Limb Access WD LHXS
PL780000A	i-Limb Access WD RHS
PL781000A	i-Limb Access WD LHS
PL782000A	i-Limb Access WD RHM
PL783000A	i-Limb Access WD LHM
PL784000A	i-Limb Access WD RHL
PL785000A	i-Limb Access WD LHL
PL786000A	i-Limb Access Friction RHS
PL787000A	i-Limb Access Friction LHS
PL788000A	i-Limb Access Friction RHM
PL789000A	i-Limb Access Friction LHM
PL790000A	i-Limb Access Friction RHL
PL791000A	i-Limb Access Friction LHL
PL840000A	i-Limb Access Ti QWD RHS
PL841000A	i-Limb Access Ti QWD LHS
PL842000A	i-Limb Access Ti QWD RHM
PL843000A	i-Limb Access Ti QWD LHM
PL844000A	i-Limb Access Ti QWD RHL
PL845000A	i-Limb Access Ti QWD LHL
PL846000A	i-Limb Access Ti Flexion RHS
PL847000A	i-Limb Access Ti Flexion LHS
PL848000A	i-Limb Access Ti Flexion RHM
PL849000A	i-Limb Access Ti Flexion LHM
PL850000A	i-Limb Access Ti Flexion RHL
PL851000A	i-Limb Access Ti Flexion LHL
PL852000A	i-Limb Access Ti WD RHS
PL853000A	i-Limb Access Ti WD LHS
PL854000A	i-Limb Access Ti WD RHM
PL855000A	i-Limb Access Ti WD LHM
PL856000A	i-Limb Access Ti WD RHL
PL857000A	i-Limb Access Ti WD LHL
PL858000A	i-Limb Access Flexion RHS
PL859000A	i-Limb Access Flexion LHS
PL860000A	i-Limb Access Friction RHXS
PL861000A	i-Limb Access Friction LHXS

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L - i-Limb UltraPT Friction	PL762000A
L - i-Limb UltraPT Friction	PL763000A
S - i-Limb UltraPT Ti QWD	PL816000A
S - i-Limb UltraPT Ti QWD	PL817000A
M - i-Limb UltraPT Ti QWD	PL818000A
M - i-Limb UltraPT Ti QWD	PL819000A
L - i-Limb UltraPT Ti QWD	PL820000A
L - i-Limb UltraPT Ti QWD	PL821000A
S - i-Limb UltraPT Ti Flexion	PL822000A
S - i-Limb UltraPT Ti Flexion	PL823000A
M - i-Limb UltraPT Ti Flexion	PL824000A
M - i-Limb UltraPT Ti Flexion	PL825000A
L - i-Limb UltraPT Ti Flexion	PL826000A
L - i-Limb UltraPT Ti Flexion	PL827000A
S - i-Limb UltraPT Ti WD	PL828000A
S - i-Limb UltraPT Ti WD	PL829000A
M - i-Limb UltraPT Ti WD	PL830000A
M - i-Limb UltraPT Ti WD	PL831000A
L - i-Limb UltraPT Ti WD	PL832000A
L - i-Limb UltraPT Ti WD	PL833000A
S - i-Limb UltraPT Ti Friction	PL834000A
S - i-Limb UltraPT Ti Friction	PL835000A
M - i-Limb UltraPT Ti Friction	PL836000A
M - i-Limb UltraPT Ti Friction	PL837000A
L - i-Limb UltraPT Ti Friction	PL838000A
L - i-Limb UltraPT Ti Friction	PL839000A
i-Limb Access	
الوصف	الرمز
XS - i-Limb Access QWD	PL764000A
XS - i-Limb Access QWD	PL765000A
S - i-Limb Access QWD	PL766000A
S - i-Limb Access QWD	PL767000A
M - i-Limb Access QWD	PL768000A
M - i-Limb Access QWD	PL769000A
L - i-Limb Access QWD	PL770000A
L - i-Limb Access QWD	PL771000A
XS - i-Limb Access Flexion	PL772000A
XS - i-Limb Access Flexion	PL773000A
M - i-Limb Access Flexion	PL774000A
M - i-Limb Access Flexion	PL775000A

i-Limb Access Flexion - يد يمني مقاس L	PL776000A
i-Limb Access Flexion - يد يسرى مقاس L	PL777000A
i-Limb Access WD - يد يمني مقاس XS	PL778000A
i-Limb Access WD - يد يسرى مقاس XS	PL779000A
i-Limb Access WD - يد يمني مقاس S	PL780000A
i-Limb Access WD - يد يسرى مقاس S	PL781000A
i-Limb Access WD - يد يمني مقاس M	PL782000A
i-Limb Access WD - يد يسرى مقاس M	PL783000A
i-Limb Access WD - يد يمني مقاس L	PL784000A
i-Limb Access WD - يد يسرى مقاس L	PL785000A
i-Limb Access Friction - يد يمني مقاس S	PL786000A
i-Limb Access Friction - يد يسرى مقاس S	PL787000A
i-Limb Access Friction - يد يمني مقاس M	PL788000A
i-Limb Access Friction - يد يسرى مقاس M	PL789000A
i-Limb Access Friction - يد يمني مقاس L	PL790000A
i-Limb Access Friction - يد يسرى مقاس L	PL791000A
i-Limb Access Ti QWD - يد يمني مقاس S	PL840000A
i-Limb Access Ti QWD - يد يسرى مقاس S	PL841000A
i-Limb Access Ti QWD - يد يمني مقاس M	PL842000A
i-Limb Access Ti QWD - يد يسرى مقاس M	PL843000A
i-Limb Access Ti QWD - يد يمني مقاس L	PL844000A
i-Limb Access Ti QWD - يد يسرى مقاس L	PL845000A
i-Limb Access Ti Flexion - يد يمني مقاس S	PL846000A
i-Limb Access Ti Flexion - يد يسرى مقاس S	PL847000A
i-Limb Access Ti Flexion - يد يمني مقاس M	PL848000A
i-Limb Access Ti Flexion - يد يسرى مقاس M	PL849000A
i-Limb Access Ti Flexion - يد يمني مقاس L	PL850000A
i-Limb Access Ti Flexion - يد يسرى مقاس L	PL851000A
i-Limb Access Ti WD - يد يمني مقاس S	PL852000A
i-Limb Access Ti WD - يد يسرى مقاس S	PL853000A
i-Limb Access Ti WD - يد يمني مقاس M	PL854000A
i-Limb Access Ti WD - يد يسرى مقاس M	PL855000A
i-Limb Access Ti WD - يد يمني مقاس L	PL856000A
i-Limb Access Ti WD - يد يسرى مقاس L	PL857000A
i-Limb Access Flexion - يد يمني مقاس S	PL858000A
i-Limb Access Flexion - يد يسرى مقاس S	PL859000A
i-Limb Access Friction - يد يمني مقاس XS	PL860000A
i-Limb Access Friction - يد يسرى مقاس XS	PL861000A

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EN – Caution: Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Ossur products are used with other recommended Ossur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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