

ProCarve


DE Gebrauchsanweisung	7
EN Instructions for use	17
FR Instructions d'utilisation	27
IT Istruzioni per l'uso	37
ES Instrucciones de uso	47
NL Gebruiksaanwijzing	57
SV Bruksanvisning	67
RU Руководство по применению	77

1

7



8

9

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DE

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Table of contents

EN

1	Product description	18
1.1	Construction and Function	18
1.2	Combination possibilities	18
2	Intended use	18
2.1	Indications for use.....	18
2.2	Area of application	18
2.3	Environmental conditions	18
2.4	Service Life	19
3	Safety.....	19
3.1	Explanation of warning symbols	19
3.2	General safety instructions	19
4	Scope of delivery.....	20
5	Preparation for use	20
5.1	Optional: Replacing the socket adapter	21
5.2	Optional: conversion of the prosthetic foot	21
5.3	Bench Alignment	22
5.4	Dynamic Trial Fitting.....	22
6	Use	23
6.1	Unlocking and locking the prosthetic knee joint.....	23
6.2	Locking and unlocking the damping mechanism	24
7	Cleaning	24
8	Maintenance	24
8.1	Replacing the ski soles	25
9	Disposal	25
10	Legal information	25
10.1	Liability	25
10.2	CE Conformity.....	25
11	Technical data	25

1 Product description

INFORMATION

Last update: 2017-05-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The ProCarve is a sport prosthesis solution for skiing, snowboarding, waterskiing, wakeboarding or other sports with similar motions. The sport prosthesis solution comprises the 3R2 prosthetic knee joint and the 1E2 or 1E2=1 prosthetic foot. The prosthetic foot can also be used by itself for a TT sport prosthesis.

The prosthetic knee joint and the prosthetic foot feature adjustable shock absorbers. The function of the shock absorbers is to absorb impacts that occur when performing a sport. During flexion a pneumatic spring provides the necessary resistance, whilst hydraulics damp movements in extension. The hardness of the pneumatic spring is adjusted using a high-pressure air pump; the extension damping strength is controlled by a hydraulic valve. This enables the best settings for the chosen sport.

The 1E2=1 prosthetic foot is equipped with a polymer ring, which increases the flexion resistance further. This is practical mainly for TT fittings.

The prosthetic knee joint can be unlocked in order to make sitting down more comfortable.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

INFORMATION

The prosthetic socket must be able to withstand the higher loading of the sport.

Restrictions of the combination possibilities

Permitted combinations	
	Reference number
Lamination anchor	4R116=T, 4R119=T
Socket adapter	4R54, 4R54=10, 4R77
Double adapter	4R72=*
Tube clamp adapter	4WR95=3
Tube adapter	2WR95
Prosthetic foot	1E2, 1E2=1
Prosthetic knee joint	3R2

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

The product is not suitable for use in everyday prostheses.

2.2 Area of application

Recommended for use in sports prostheses for skiing and snowboarding sports as well as waterskiing and wakeboarding.

Approved for a body weight of **up to 100 kg**.

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -20 °C to +60°C
Allowable relative humidity: 0% to 90%, condensing
Moisture: fresh water, salt water

Unallowable environmental conditions

Mechanical vibrations or impacts (except during sports)

Perspiration, urine, acids

Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service Life

Based on ISO 10328, the prosthetic components were tested by the manufacturer for 3 million load cycles. Depending on the type of sport, this corresponds to a period of use of 3 to 5 years under high loading.

3 Safety**3.1 Explanation of warning symbols****WARNING**

Warning regarding possible serious risks of accident or injury.

**CAUTION**

Warning regarding possible risks of accident or injury.

**NOTICE**

Warning regarding possible technical damage.

3.2 General safety instructions**WARNING****Using the product**

Severe falls due to improper use of the product

- ▶ Only use the product after receiving full instruction by qualified personnel.
- ▶ Do not use the product as an everyday prosthesis.

**WARNING****Water sports and swimming with the prosthesis**

Risk of drowning due to restricted buoyancy.

- ▶ Always wear a flotation device suitable for the sport when playing water sports and swimming.
- ▶ The flotation device must be able to support your body weight together with your prosthesis.

**CAUTION****Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 18).

**CAUTION****Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

**CAUTION****Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION**Reuse on other patients and improper maintenance**

Fall due to loss of functionality or damage to the product

- ▶ Only use the product for a single patient.
- ▶ Maintain the product regularly to extend the service life of the product (see the section "Maintenance").

⚠ CAUTION**Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

3R2 ProCarve prosthetic knee joint

see fig. 2, item:	Quantity	Designation	Reference number
–	1	Instructions for use	–
–	1	Quick reference guide	647G953
1	1	Prosthetic knee joint	3R2
2	1	Knee blocking clip (black strap)	4G115

1E2, 1E2=1 ProCarve prosthetic foot

see fig. 2, item:	Quantity	Designation	Reference number
–	1	Instructions for use	–
–	1	Quick reference guide	647G953
4	1	Prosthetic foot with foot pads and ski soles	1E2*
3	1	In addition for 1E2: foot blocking clip (red strap)	4G115=1
5	1	Foot shell	4G901
6	1	High-pressure air pump	755Y68

Spare parts/accessories (not included in the scope of delivery)

Designation	Reference number
Ski sole set	2Z503=1

5 Preparation for use

⚠ CAUTION**Incorrect alignment, assembly or adjustment**

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION**Use of talcum**

Risk of injury, damage to the product due to lack of lubrication

- Do not use talcum on the product or other prosthesis components.

INFORMATION

For a waterproof prosthesis, use titanium screws and set screws.

INFORMATION

The prosthesis alignment described and the air pressure values given represent basic values. They are to be adapted to the individual requirements of the patient.

The three-dimensional alignment of the prosthesis socket and the modular components affects the function of the prosthesis. The patient can use the properties of the sports prosthesis to their fullest potential only when the product is correctly aligned.

The residual limb position, e. g. abduction, adduction and flexion, must be considered when positioning the socket connector. It is recommended to align according to the patient's everyday prosthesis.

Alignment is performed in 2 steps:

- Bench alignment
- Dynamic trial fitting

The sports prosthesis is individually adjusted during dynamic trial fitting at a sports facility. All the parameters are adjusted according to the requirements and ability of the patient during this fitting.

The 1E2* Prosthetic Foot is equipped for use in a ski binding as delivered. If the prosthetic foot is to be used in a snowboarding soft boot or another boot, the foot shell included in the scope of delivery must be attached before the bench alignment.

5.1 Optional: Replacing the socket adapter

The socket adapters on the prosthetic knee joint can be replaced. This can be done with the screws provided with the prosthetic knee joint and with a different tightening torque. Cap screws are also provided with the 4R54=10 Socket Adapter.

> Required tools and materials:

710D4 Torque Wrench, 636K13 Loctite®, 501S84=M6x18 Countersunk Head Screws

- 1) Remove the countersunk head screws from the old adapter.
- 2) Align the new adapter and position it on the prosthetic knee joint.
- 3) Apply Loctite® to the four screws and screw them in.
- 4) Tighten the posterior screws (tightening torque: **10 Nm**).
- 5) Tighten the anterior screws (tightening torque: **10 Nm**).

5.2 Optional: conversion of the prosthetic foot

Ski binding to foot shell

> Required tools and materials:

4 mm and 5 mm Allen keys, Foot Shell, 636K13 Loctite®, 501S71=M6x10 Countersunk Head Screws

- 1) Remove the 4 cap screws on the foot pad attached to the anterior side (see fig. 3, see fig. 4).
- 2) Remove the foot pad from the prosthetic foot (see fig. 5).
- 3) Remove the 4 cap screws from the foot pad attached on the posterior side and remove the foot pad.
- 4) Place the foot shell on the prosthetic foot.
- 5) Secure the 4 countersunk head screws with Loctite® and screw in the foot shell (mounting torque: **10 Nm**) (see fig. 6).

Foot shell to ski binding

> Required tools and materials:

4 mm and 5 mm Allen keys, foot pads for the ski binding, cap screws for the foot pads, 636K13 Loctite®

- 1) Remove the 4 countersunk head screws and the foot shell from the prosthetic foot (see fig. 6).
- 2) Place the anterior foot pad on the prosthetic foot (see fig. 5).
- 3) Secure the 2 longer cap screws with Loctite® and screw into the foot pad on the front side (mounting torque: **10 Nm**) (see fig. 4).

- 4) Secure the 2 shorter cap screws with Loctite® and screw into the foot pad on the side of the foot sole (mounting torque: **10 Nm**) (see fig. 3).
- 5) Attach the posterior foot pad in the same manner.

5.3 Bench Alignment

TF bench alignment

Bench alignment process (see fig. 7)			
Required materials: 710D4 torque wrench, 636K13 Loctite®, 4 mm Allen key			
①	Place the prosthetic foot on the work surface.		
②	Connect the prosthetic foot with a tube adapter and tube clamp adapter or a double adapter. Set the angle α for the chosen sport:		
	Skiing	min. 75°	max. 80°
	Snowboarding	min. 75°	max. 85°
③	INFORMATION: When positioning the alignment reference point, take into consideration the heel height of the ski boot or snowboard boot. Mount the prosthetic knee joint according to the instructions for use of the modular adapter. Position the alignment reference point (rotation axis of the prosthetic knee joint) at the height of the medial tibial plateau to floor measurement + 20 mm .		
④	Position the prosthetic socket on the prosthetic knee joint: Tuberosity reference point at the height of the tuberosity-floor measurement Set the angle β for the chosen sport and the movement style, taking the individual socket flexion into account.		
	Skiing and snowboarding	min. 120°	max. 160°
	Compensate for any height difference due to high knee flexion by using double adapters.		
⑤	Check the prosthesis alignment during the dynamic trial fitting. Make corrections if necessary.		

TT bench alignment

INFORMATION
Ottobock recommends a fitting with a thigh sleeve.

Bench alignment process (see fig. 7)			
Required materials: 710D4 torque wrench, 636K13 Loctite®, 4 mm Allen key			
①	Place the prosthetic foot on the work surface.		
②	Connect the prosthetic foot with a tube adapter and tube clamp adapter or a double adapter. Set the angle α for the chosen sport.		
	Skiing	min. 80°	max. 85°
	Snowboarding	min. 80°	max. 90°
③	INFORMATION: When positioning the alignment reference point, take into consideration the heel height of the ski boot or snowboard boot. Position the prosthetic socket on the modular adapters: Observe the MPT-floor measurement and length measurement of the everyday prosthesis Align the prosthetic socket according to the knee angle, determined by the type of sport and the movement style. Take into account the individual socket flexion.		
④	Check the prosthesis alignment during the dynamic trial fitting. Make corrections if necessary.		

5.4 Dynamic Trial Fitting

During dynamic trial fitting, the alignment and settings of the prosthesis are checked. The flexion and extension movement damping are adjusted to the patient's requirements.

Adjusting the flexion and extension damping				
Required materials: High-pressure air pump, Allen wrench (3 mm)				
①	The pneumatic springs in the shock absorber control the flexion of the prosthetic knee joint and the prosthetic foot. Their hardness is adjusted using the high-pressure air pump (see fig. 8). Unscrew the valve cap on the shock absorber. Screw the air pump on to the valve. Set the desired pressure. Note that 1 bar to 2 bar of pressure escapes when removing the air pump. Unscrew the air pump from the valve and screw the valve cap back on.			
	Prosthetic foot			
①	Skiing	TF	min. 5 bar	max. 15 bar

Adjusting the flexion and extension damping				
①	Skiing	TT	min. 20 bar	max. 25 bar
	Snowboarding	TF	min. 8 bar	max. 15 bar
		TT	min. 15 bar	max. 25 bar
	Prosthetic knee joint			
	Skiing		min. 8 bar	max. 15 bar
②	The hydraulic valves control the extension damping (see fig. 9). The adjustment influences the speed of the prosthetic knee joint extension and the plantar flexion of the prosthetic foot.			
	Set the valve to the desired value using an Allen key (3 mm).			
	Left stop (–)		Minimum damping (factory setting)	
	Right stop (+)		Maximum damping	

6 Use

⚠ CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.
- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

Leak in the product's hydraulic system

Environmental damage or skin irritation due to leaking hydraulic oil

- ▶ Before each use, check the product for leaking hydraulic oil.
- ▶ In case of leaks, do not continue using the product and have it repaired promptly.

NOTICE

Incorrect ski binding adjustment for the ProCarve

Damage to the product due to improperly adjusted trigger value (Z value)

- ▶ Have your ski binding adjusted by a ski service shop.
- ▶ Adjust the weight to trigger the ski binding on the prosthesis side to the **body weight minus 10 kg**.
- ▶ The maximum trigger torque (M_z) is **100 Nm**.

6.1 Unlocking and locking the prosthetic knee joint

The prosthetic knee joint can be unlocked in order to sit down. This enables a chairlift to be used, for example. When standing up, the lock must be engaged again manually.

Flexion (unlocking)

- > **Prerequisite:** The prosthetic knee joint is extended and not under load.
- ▶ Pull the posterior strap on the prosthetic knee joint upwards in order to unlock the prosthetic knee joint (see fig. 12).
- The prosthetic knee joint can now be flexed (see fig. 13, see fig. 14).

Extension (locking)

⚠ CAUTION

The prosthetic knee joint is not completely extended

Risk of falling due to the lock not being engaged

- ▶ Before using the prosthesis, always check whether the lock on the prosthetic knee joint has fully engaged.

> **Prerequisite:** The prosthetic knee joint is flexed.

- 1) Extend the prosthetic knee joint.
- 2) Engage the lock by pressing with the hand.
- 3) Check whether the prosthetic knee joint is engaged before using it.

6.2 Locking and unlocking the damping mechanism

Before walking, a blocking clip must be installed respectively on the damping mechanism of the prosthetic knee joint and the prosthetic foot. The blocking clips have to be removed again before skiing.

Locking the damping mechanism

> **Prerequisite:** The prosthesis is not under load.

- 1) Install the knee blocking clip (black strap) on the damping mechanism of the prosthetic knee joint (see fig. 10, see fig. 11).
- 2) Install the foot blocking clip (red strap) on the damping mechanism of the prosthetic foot.

Unlocking the damping mechanism

> **Prerequisite:** The prosthesis is not under load.

- 1) Remove the knee blocking clip (black strap) from the damping mechanism of the prosthetic knee joint (see fig. 11, see fig. 10).
- 2) Remove the foot blocking clip (red strap) from the damping mechanism of the prosthetic foot.

7 Cleaning

⚠ CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

> **Recommended cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

- 1) Clean the product with clear fresh water and a pH neutral soap.
- 2) Rinse the soap away with clear fresh water.
- 3) Dry the product with a soft cloth.
- 4) Allow to air dry in order to remove residual moisture.

8 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ Following an individual period for the patient to get accustomed to the product, check the settings of the prosthesis and, if necessary, adapt them.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the joints of the prosthesis for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

8.1 Replacing the ski soles

> Required tools and materials:

Phillips screwdriver, 501B4=4.2x16 Oval-Head Screws

- 1) Unscrew the 4 screws in the ski sole that is to be replaced (see fig. 15).
- 2) Replace the ski sole and screw it in (see fig. 16).

9 Disposal

Dispose of the product in accordance with national regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11 Technical data

Prosthetic knee joint		
Reference number	3R2	
Weight [g]	1990	
System height [mm]	241	
System height, distal [mm]	207	
System height, proximal [mm]	34	
Proximal connection	Pyramid adapter	
Distal connection	Pyramid adapter	
Max. flexion angle (unlocked) [°]	80	
Max. body weight [kg]	100	

Prosthetic foot		
Reference number	1E2	1E2=1
Weight [g]	1550	
System height		
- with attachments for ski binding [mm]	120	
- with foot shell [mm]	103	
Proximal connection	Pyramid adapter	
Max. range of motion [°]	25	
Max. body weight [kg]	100	

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Otto Bock HealthCare GmbH