
A photograph showing the lower legs and feet of two people. On the left, a woman's legs are wearing black sheer tights and black patent leather high-heeled pumps with red soles. On the right, a man's legs are wearing black knee-high socks and black lace-up oxford shoes. The background is plain white.

SIGVARIS

Wear & Care

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Uso y Cuidado

FOR EACH STEP
OF YOUR LIFE

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SIGVARIS

THANK YOU FOR CHOOSING SIGVARIS,

a global leader in graduated compression socks and hosiery. Our objective is to produce the best compression garments, which will enhance our customers' health and quality of life. SIGVARIS offers many styles of graduated compression hosiery that are both fashionable and offer all day comfort. Please take a few moments to read this information to fully enjoy the benefits of your SIGVARIS garment.

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WEAR & CARE

STOCKING APPLICATION

For best results wear your stockings daily starting as early in the day as possible. Be sure legs are dry. We recommend the use of donning gloves. Never gather the stockings or pull from the top band. A proper fit results in proper therapy.

CLOSED TOE STOCKINGS

- 1) Fold the garment down over the heel dividing the length of the stocking in half, to get your foot and heel easily in place. Remember, do not bunch or gather the stocking. *(Img. A)*
- 2) Step into the stocking and gently pull until your foot is in place. *(Img. B)*
- 3) If you are unable to get your foot in place, use the donning gloves to gently move the fabric in an upward direction. Once your heel is in place, unfold the stocking upward and gently smooth out the fabric by moving it in an upward direction. *(Img. C)*
- 4) The top band should be approximately two fingers width from the right-angle bend on calf length stockings. *(Img. D)* For thigh-highs, position the stocking comfortably on the upper thigh.

OPEN TOE STOCKING

- 5) For the open toe styles only, use the SIGVARIS Foot Slip. Never gather the stocking prior to donning. *(Img. E)*
- 6) Leave the stocking fully-extended. Pull the stocking up over the foot past the ankle. *(Img. F)*
- 7) Using the donning gloves to ease the stocking over the foot, then over the heel and then up the leg. Make sure the heel is in place first. Once your heel is in place, gently smooth out the fabric by moving it in an upward direction. Remove Foot Slip by gently pulling it out of the open toe area of the stocking. *(Img. G)*
- 8) For open toe styles, the foot band of the stocking should be placed at the base of the toes. The toes should remain free. For calf length style, the top band should be approximately one inch (two finger lengths) below bend of knee. *(Img. H)*

TO REMOVE

Remove stockings before going to bed, unless otherwise advised by a physician.

- 9) Use SIGVARIS donning gloves to make stocking removal easy. Take hold of the stockings or socks at the top band and peel downwards past the ankle. Bunching makes removal difficult. *(Img. I)*
- 10) With the flat of the hand, work the stocking over the heel and remove. *(Img. J)*

CARE RECOMMENDATIONS

Machine wash on gentle cycle using lukewarm water and a mild detergent (use SIGVARIS Washing Solution for best results). For sheer stockings, hand washing is recommended. To protect the garment during washing, we recommend using a lingerie bag.

Rinse well and air-dry. No bleach. No fabric softeners. Do not tumble dry.

INDICATIONS

8-15mmHg*: Diabetics who experience edema (see Cautions below) and others needing a low compression sock.

15-20mmHg: Tired, achy legs; Occupational or evening edema; Leg discomfort from long hours of standing or sitting; Prophylaxis of varicosities & swelling during pregnancy; Prophylaxis for legs predisposed to risk; Reduce swelling during travel; Prophylaxis of DVT during travel

18-25mmHg*: Diabetics who experience edema (see Cautions below)

20-30mmHg: Heavy, fatigued, aching legs; Mild edema in lower extremities; Mild varicosities with minimal edema; Mild varicosities during pregnancy; Minimal edema upper extremities; Post-procedure of small veins; Prevention or management of DVT

Note: Prescription Use: A prescription is recommended for stockings with the indications of 20mmHg or higher.

30-40mmHg: Moderate varicosities with mild to moderate edema (also during pregnancy); Post-procedure of larger veins to maintain treatment success; CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations); Management and treatment of DVT or Post-Thrombotic Syndrome; Superficial phlebitis; Lymphedema after decongestant therapy to maintain reduction; Orthostatic / postural hypotension; Post-fracture, post-traumatic edema

Note: Prescription Use: A prescription is recommended for stockings with the indications of 20mmHg or higher.

40mmHg or higher: Irreversible lymphedema (life-long compression); Severe Post-Thrombotic Syndrome; Severe varicosities and/or edema

CONTRAINDICATIONS

You should not wear SIGVARIS if you have any of the following medical conditions: Arterial insufficiency, intermittent claudication, ischemia, uncontrolled congestive heart failure, untreated septic phlebitis or phlegmasia coerulea dolens.

Img. A

Img. B

CAUTIONS

Skin sensitivities or allergies; Advanced neuropathy indicating sensory loss in leg or foot; Diabetes with advanced arterial compromise; Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician; Eczema, hyperdermatitis, dermatosis, dermatitis (following treatment); Incompatibility to fabric.

Please see a doctor in the event your condition worsens or any other abnormal symptoms are observed. If you experience any discomfort, remove the stocking and consult your physician.

NOTE

The manufacturer is not liable for damage/injuries which have been caused by improper handling or misuse.

* Consult your country group's website for product availability.

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Img. J

SIGVARIS CARES

Tell Us About Your Product Experience
& Be Entered to Win a Free Pair

Your feedback is greatly appreciated. Please visit <http://bit.ly/SIGcare> and take a five minute online survey about this product. Your responses will only be shared with our research and development team.

As a thank you for your time, you will be entered into our monthly drawing for a FREE pair of SIGVARIS WELL BEING OR SIGVARIS SPORTS valued up to \$80 MSRP.

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<http://bit.ly/SIGcarees>

SIGVARIS, INC.

1119 Highway 74 South
Peachtree City, GA 30269 USA
sigvarisusa.com

SIGVARIS, CORP.

8423 Chemin Dalton
Montréal, QC
H4T 1V5
sigvaris.ca

SIGVARIS SA DE CV

Periférico Sur 4600, Piso 2
Insurgentes Cuicuilco
Coyoacán, 04530
Cd de México, D.F., MX
sigvaris.mx

R670205 REV: 1015