AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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Material	Material
Aluminium (AlCuMg1 F40)	Aluminium (AlCuMg1 F40)
Polyamid (PA)	Polyamide (PA)
Polyurethan (PUR)	Polyurethane (PUR)
Aluminium (AlMg4,5Mn0,7)	Aluminium (AlMg4,5Mn0,7)
Polyoxymethylen (POM)	Polyoxymethylene (POM)
Polyethylenterephthalat (PET)	Polyethylene terephthalate (PET)
Acryl-Isocyanat (Pulverlack)	Acrylic isocyanate (Powder paint)
Edelstahl (A4)	stainless steel (A4)
Polyamid, glasfaserverstärkt (PA-GF)	Polyamide glass fiber reinforced (PA-GF)
Ethylvinylacetat (EVA)	Ethylvinylacetate (EVA)
Polyester (PES)	Polyester (PES)
Acrylat, modifiziert (Kleber)	Acrylate modified (adhesive)
Messing (CuZnx)	Brass (CuZnx)
Eisen (Fe)	Iron (Fe)
Kunstharz	Synthetic resin
Kunststofffolie	Plastic film

 **BAUERFEIND AG**
Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany
P +49 (0) 36628 66-4000
F +49 (0) 36628 66-4499
E info@bauerfeind.com
BAUERFEIND.COM


SecuTec® Genu

Knieorthese
Knee orthosis
Orthèse du genou
Ortesis de rodilla

All instructions for use:
www.bauerfeind.com/downloads

de	deutsch	2	pl	polski	17	et	eesti	32
en	english	3	cs	česky	18	lv	latviešu	33
fr	français	5	sk	slovensky	19	lt	lietuvių	34
nl	nederlands	6	hu	magyar	21	uk	українська	36
it	italiano	7	hr	hrvatski	22	he	עברית	37
es	español	9	sr	srpski	23	ja	日本語	39
pt	português	10	sl	slovenski	25	ko	한국어	40
sv	svenska	12	ro	romanian	26	ar	عربي	41
no	norsk	13	el	ελληνικά	27	zh	中文	43
fi	suomi	14	tr	türkçe	29			
da	dansk	16	ru	русский	30			

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 english

Dear Customer,

thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read through these **instructions for use** carefully. If you have any questions, please contact your doctor or medical retailer.

Intended purpose

SecuTec Genu is a medical product. It is a stabilizing orthosis that limits the range of motion in the case of complex knee injuries.

Indications

- Rupture of the anterior and / or posterior cruciate ligament (ACL / PCL)
- After ligament surgery / ligament reconstruction
- Severe and / or complex instabilities (traumatic, degenerative, e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion (e.g. after meniscus refixation or meniscal implantation)
- Patellar fracture (non-surgical and post-operative)

Risks of using this product

Important information

- The SecuTec Genu functional orthosis is a prescription product that should be worn as directed by a physician. The SecuTec Genu must be worn according to the specifications in these instructions and for the listed indications.
- Product liability cannot be assumed in the event of inappropriate use.
- No improper modifications may be made to the product. Failure to comply with this requirement may adversely affect product performance, thereby excluding any product liability.
- Should symptoms increase or should you notice any unusual changes, please contact your physician without delay.

- Combining the support with other products (e.g. compression stockings) must be discussed beforehand with your physician.
- Avoid contact with ointments, lotions, or any substances containing grease and acids.
- Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please immediately contact your medical supply retailer.
- No side effects that affect the entire body have been reported to date. This assumes the orthosis is used / fitted correctly. Any medical aids – supports and orthoses – applied externally to the body can, if tightened excessively, lead to local pressure symptoms or, rarely, constrict blood vessels or nerves.

Contraindications

No clinically significant hypersensitive reactions have been reported to date. If you have any of the following conditions, such aids should only be fitted and worn after consultation with your physician:

- 1.Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
- 2.Varicose veins (varicosis)
- 3.Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus).
- 4.Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid.

Application instructions

Size determination

Please determine the correct size in accordance with the specifications.

Procedure for the first fitting of the SecuTec Genu orthosis

Preparation (only by a trained specialist?)

- Mark the compromise axis of rotation of the knee joint (as defined by Nietert) at a 30° bend.
- The pivot point of the orthosis hinges is approximately 4 mm above an imaginary line in the center of the hinge cover.

Adjustment options for the hinge (to be performed by a trained specialist only)

- The hinge cover can be opened via a screw fastening. This also serves to fasten the movement-limiting wedges (extension – front / flexion – rear) and protects the hinge mechanism against dirt.
- Opening the screw fastening (quarter turn until the arrow points at the open lock symbol) allows the hinge protection to be removed and the desired stops to be changed.
- The desired extension and flexion wedges are selected from the enclosed frame and are used **❶** **❷**. The following limitations are possible: Extension: 0°, 10°, 20°, 30°, 45° Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Then, click the hinge cover onto the hinge, check whether it is correctly positioned and firmly attached to the hinge, and close the screw fastening **❸** **❹**.

Static fitting (only by a trained specialist)

- Perform the initial adjustment with the straps open.
- Bend the knee to be treated to a 30° angle and place the SecuTec Genu on it so that the hinge pivot point of the orthosis hinge is aligned with the compromise axis of rotation as defined by Nietert.
- Check the fit of the frame parts when they are laid flat to make sure their circumference and longitudinal direction are correct and adjust them if necessary.
- A bending iron with protective caps should be used for any extensive bending work.
- Select the required thickness of knee joint cushions and insert them (there are two thicknesses to choose from).
- If necessary, the shin cushion provided can be attached in order to redistribute pressure in the event of pressure pain in the shin region or swelling around the tibial tubercle (Tuberositas tibiae).

Putting on the SecuTec Genu

- The leg is bent at a gentle angle (approximately 30°). The position of the kneecap is then determined through palpation.
- The center of the hinge (pivot point of the orthosis hinges) is positioned at approximately the same height as the center of the kneecap.
- To treat rear instability (PCL), the (proximal) strap near the hinge can be shifted forwards (ventrally). To treat a complex instability, the enclosed sixth strap (ventral, proximal) can be used.
- Mount all the quick-fit fasteners in the frame.
- Decide on the calf strap height that provides the optimum fit for the

patient. If you mount the quick-fit fasteners for this strap in the upper part of the orthosis, thread the opposite end of the strap out of the frame and push the stopper downwards before re-fastening the Velcro closure on the strap.

- Start by tightening the calf strap, followed by the strap beneath it.
- Then tighten the bottom thigh strap at the back, followed by the strap above it. (The quick-fit fasteners are numbered so that you can easily see the order in which to close the straps **❺**.) – When closing the straps, check whether the orthosis has rotated laterally around the leg and correct this if necessary.
- Position the front strap with the protective cushion for the tibial tubercle in such a way that the strap can be closed without exerting any pressure. After all the straps have been closed, the orthosis must exert an even pressure when the muscles are tensed.
- Once the orthosis has been put on, check once more that it is in the correct position.

Taking off the SecuTec Genu

- Grip the handle of the quick-fit fastener and unhook it from the thigh frame.
- Then lift the upper orthosis frame off the thigh. During this step, the knee joint cushions stay in contact with the knee joint.
- Now remove the straps from the lower leg frame of the orthosis.
- The top and, if necessary, the bottom strap at the front can remain closed.
- After all the straps have been removed, the orthosis can be taken off the knee joint from the front.

Cleaning instructions

Never expose the orthosis to direct heat (e.g. central heating, sunlight etc.)!

- The aluminum frames of the orthosis are coated and can be cleaned using pH-neutral soap.

We recommend washing the cushions and straps by hand at 30 °C. This product has been tested by our in-house quality management system. However, should you have any complaints, please contact your local medical supply retailer.

Where to apply the orthosis

In accordance with the indications (knee). Please see intended use.

Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

Assembly and fitting instructions

The SecuTec Genu is supplied as a pre-assembled orthosis in standard sizes. The first fitting, which includes the hinge adjustments and the static fitting, must be carried out by a trained specialist.

Technical specifications / parameters

The SecuTec Genu is a functional orthosis for the knee joint. It consists of a frame, hinges, straps, quick-fit fasteners and cushions.

Notes on reuse

This product is intended for your treatment only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Genu have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Version: 2019-09

- ¹Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso
- ²A specialist is any person who is authorized, according to their state's regulations, to fit and provide instructions on the use of supports / orthoses.

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