



SensorHand Speed 8E38=8*, 8E39=8*, 8E41=8*

DE Gebrauchsanweisung	3
EN Instructions for use	20
FR Instructions d'utilisation	34
IT Istruzioni per l'uso	51
ES Instrucciones de uso	67
PT Manual de utilização	83
NL Gebruiksaanwijzing	99
SV Bruksanvisning	115
DA Brugsanvisning	130
NO Bruksanvisning	145
HU Használati utasítás	159
TR Kullanma talimatı	174
EL Οδηγίες χρήσης	189
RU Руководство по применению	203
JA 取扱説明書	220

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Date of the last update: 2021-04-15

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

Only start up the product in accordance with the information contained in the accompanying documents supplied.

1 Intended use

1.1 Medical Purpose

The 8E38=8*, 8E39=8, and 8E41=8* Ottobock SensorHand Speed are to be used **solely** for the exoprosthetic fitting of the upper extremities.

1.2 Application

The 8E38=8*, 8E39=8, and 8E41=8* Ottobock SensorHand Speed can be used for unilateral or bilateral amputees.

1.3 Application conditions

The 8E38=8*, 8E39=8, and 8E41=8* SensorHand Speed have been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

The 8E38=8*, 8E39=8, and 8E41=8* Ottobock SensorHand Speed is designed exclusively for use on one patient. Use of the product by other persons is not approved by the manufacturer.

1.4 Prosthetist qualifications

The fitting of a patient with the 8E38=8*, 8E39=8, and 8E41=8* Ottobock SensorHand Speed may only be carried out by prosthetists who have been trained and certified by Ottobock for the use of MyoBock components.

2 Safety information

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Do not use silicone spray when putting on the prosthetic glove; silicone spray may compromise secure fit of the glove. As a mounting aid the 633S2 Ottobock Procomfort Gel is recommended.
- In order to reduce the malfunctioning risk of the SensorHand Speed and of the connected components due to strong electromagnetic radiation (e.g. when near high-tension power lines, transmitters, transformers or security systems for goods in department stores), the electrode gains should be set as low as possible.
- Before interrupting or establishing electrical connections (such as when removing the hand from the prosthesis), the system must be disconnected from the power supply. For this pur-

pose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.

- If you do not achieve the expected results by adjustment and by selection of the appropriate program, please contact an Ottobock Myo-Service facility.
- Use the "Patient Information" section of this pamphlet (Section 3) to inform your patient about the function and use of the SensorHand Speed.
- The SensorHand Speed has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- Use the product only in combination with components approved by Ottobock (combination options). Ottobock shall not accept liability in the event that the product is used with components other than those Ottobock specifies.
- Any other manipulations to the prosthesis not described in these Instructions for Use are not permitted.
- The battery may only be handled by Ottobock Service Centres (batteries are not user-serviceable).
- The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7m
 - DECT cordless phones incl. base station: 0.35m
 - WiFi (routers, access points,...): 0.22m

3 Patient information

- Do not allow foreign particles or liquids to get into the SensorHand Speed mechanism.
- The SensorHand Speed should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the SensorHand Speed.
- The Auto-Grasp function which is initiated by touching the thumb sensor, requires careful handling in certain situations, i.e. when shaking hands, body hygiene, etc.
- The SensorHand Speed and SensorHand Speed components may be opened or repaired only by certified Ottobock Myo-Service technicians.
- If the SensorHand Speed will not be used for an extended period of time, it should be left in the open position to protect the sensor and control mechanism.
- Please provide your patient with the 646D165 Patient Information on the SensorHand Speed.

3.1 Driving vehicles

An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork). It is imperative that the driver is able to operate the vehicle without any risk with the SensorHand Speed turned off. Driving with the SensorHand Speed turned ON may present a risk if the SensorHand Speed inadvertently moves due to unintentional muscle contraction or other causes.

Caution! Before a SensorHand Speed with a quick-disconnect wrist can be used, the terminal device must be positioned in such a way that a slight turn of it cannot disconnect the hand from the prosthesis.

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

4 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product. The grace period is no more than one month before or three months after maintenance is due. In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover. Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty. The following components must always be sent in for maintenance and repairs: The product, battery charger and power supply. The shipping container for the loaner unit you receive must be used for sending back the components requiring inspection.

5 Scope of delivery

- 1 pc. SensorHand Speed
- 1 pc. Instructions for Use

6 Description and function

The Ottobock SensorHand Speed is a myo-electronically controlled prosthetic hand that is characterized by a high grasping speed in connection with an innovative and sensitive control concept. The Ottobock SensorHand Speed is equipped with the SUVA-Sensor² grip stabilization system "Auto-Grasp", the FlexiGrip function and a programming function.

The proportional DMC-Control (Dynamic Mode Control) allows the amputee to control grip speed and grip force proportionally to his muscle signal. If the muscle signal changes then grip speed and grip force immediately adapt to the new muscle signal.

2 Developed by Otto Bock Healthcare Products GmbH in cooperation with the Schweizerische Unfall Versicherungs Anstalt, SUVA (Swiss Insurance Agency)

The Coaxial Bushing of the 8E38=8* SensorHand Speed is marked with an orange ring, the Connecting Cable of the 8E39=8* or 8E41=8* SensorHand Speed is marked by an orange bushing.

6.1 The SUVA-sensor technology

The SUVA-Sensor Technology is integrated into the thumb (Fig. 1) and instantly senses when a gripped object changes position and begins to slip. The system then automatically and continuously, starting from the original level, increases grip force until the object is again held securely.

6.2 FlexiGrip function

The FlexiGrip function allows the user to change the position of an object within the hand without the need to open and close the hand using myo-control. The prosthetic hand follows the position shifts of an object such as a sound hand would. This gives the grip flexibility.

6.3 Programming function

In order to find the optimal adjustment you can choose between six programs. The programs are identified by coding plugs of different color (Fig. 2). Upon delivery, the SensorHand Speed is equipped with a black coding plug and preconfigured with program 1. For program selection, plug in the 757T13 MyoSelect. If you prefer coding plugs for program selection, this is also possible.

6.4 Switching the SensorHand Speed ON and OFF

By pressing the cosmetic glove you can actuate the ON/OFF switch integrated into the control electronics.

By pressing the back of the hand: Function ON (Fig. 3)

By pressing the thumb: Function OFF (Fig. 4)

7 Adjusting the electrode(s)

The best possible function of the SensorHand Speed with SUVA-Sensor Technology is achieved through optimal positioning and adjustment of the electrodes. This is easy with the 757M10/757M11=X-CHANGE Ottobock MyoBoy.

The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interferences from electric devices, the position of the electrodes should be verified and, if need be, repositioned. If the interference cannot be eliminated, please contact the Ottobock Myo-Service.

Attention!

Allow the patient to rest during the adjustment of the electrodes to avoid fatigue, and thus incorrect electrode settings (i.e. excess sensitivity).

Programs 1, 5, and 6:

Select the desired program version in the MyoSoft or PAULA program. Adjust each electrode so that the amputee can maintain the respective muscle signal for approximately two seconds at a point above the value HIGH (Fig. 7, item 1).

Program 2:

Select the desired program version in the MyoSoft or PAULA program. Adjust each electrode so that the amputee can maintain the respective muscle signal for approximately two seconds at a point above the value LOW (Fig. 7, item 2).

Program 3:

Select the AutoControl mode for the SensorHand Speed in the MyoSoft or PAULA program. Adjust each electrode so that the amputee can maintain the respective muscle signal for approximately two seconds at a point above the value ON (Fig. 8).

Program 4:

Select the VarioControl mode for the SensorHand Speed in the MyoSoft or PAULA program. Adjust each electrode so that the amputee can maintain the respective muscle signal for approximately two seconds at a point above the value HIGH (Fig. 7, item 1).

Note to Fig. 7 and 8: With programs 3 and 4 you can only see one curve!

8 Control modes

To define the desired program insert the appropriate coding plug (Fig. 2).

13E184=1 (white) 13E184=3 (green) 13E184=5 (yellow) 13E184=8 (black)
13E184=2 (red) 13E184=4 (blue) 13E184=6 (purple)

If the black coding plug (13E184=8) has been inserted into the electronic control unit of the SensorHand Speed, the 757T13 MyoSelect can be used for setting all function types. For setting the speed of the SensorHand Speed using 757T13 MyoSelect, any of the coding plugs can be inserted. To do this read the 647G131 Instructions for Use manual.

Before changing the coding plug remove inner-hand safety ring. Partially remove the inner hand. Change the coding plug (Fig. 2).

Return safety ring and inner hand.

To reset the system so that it will recognize the new coding plug, remove the battery and reinsert it. Ten program variants are available.

9 Program survey

Program 1	OPEN	CLOSE	Indication
DMC plus Sensor White coding plug	Myo-signal to the electrode	Myo-signal to the electrode	Preferred program for amputees with two strong muscle signals
Two electrodes	Proportional speed	Proportional speed	

Program 2	OPEN	CLOSE	Indication
AutoControl – LowInput Red coding plug	Myo-signal to the electrode	Myo-signal to the electrode digitally (short signal at any level)	For amputees with two weak muscle signals
Two electrodes	Proportional speed	Constant speed	
AutoControl – LowInput Red coding plug	Myo-signal to the electrode	Signal from the switch	For amputees with only 1 muscle and weak muscle signal
One electrode & one optional Ottobock Switch	Proportional speed	Constant speed	

AutoControl – LowInput Red coding plug	Hand opens as long as opening side of the switch is activated	Signal from closing side of switch: Hand closes	For patients with too weak or no muscle signal
Any Ottobock switch	Constant speed	Constant speed	

Program 3	OPEN	CLOSE	Indication
AutoControl Green coding plug	Quick, sustained myo-signal to the electrode	Very slow muscle relaxation: Hand stays open: Quick muscle relaxation: Hand closes	For amputees with only 1 muscle and very weak muscle signal
One electrode	Constant speed	Constant speed	
AutoControl Green coding plug	Hand opens as long as switch is activated	Hand closes automatically as soon as switch is released	For amputees with weak or no muscle signals
Any Ottobock switch	Constant speed	Constant speed	

Program 4	OPEN	CLOSE	Indication
VarioControl Blue coding plug	Speed and strength of muscle tension to the electrode	Speed and strength of muscle relaxation to electrode	For amputees with 1 muscle and strong muscle signal or tendency to co-contraction
One electrode	Proportional speed	Proportional speed	
VarioControl Blue coding plug	Speed and strength of pull to linear transducer	Speed of release of the linear transducer	For amputees with weak or no muscle signal
One linear transducer	Proportional speed	Proportional speed	

Program 5	OPEN	CLOSE	Indication
VarioDual Yellow coding plug	Speed and strength of muscle tension to first electrode	Speed and strength of muscle relaxation to first electrode Proportional speed.	Control for amputees with 2 strong muscle signals
Two electrodes	Proportional speed	Grip force proportional to strength of muscle signal to second electrode.	

Program 6	OPEN	CLOSE	Indication
DMC plus Sensor Purple coding plug	Myo-signal to electrode	Myo-signal to electrode	Program for amputees with 2 strong muscle signals
Two electrodes	Proportional speed	Proportional speed	SUVA Sensor and FlexiGrip can be switched off.

All power forces (e.g. 100 N ~ 10 Kp) are values without obligation which are given to show the functional principle of the various modes.

9.1 Program 1: **DMC plus Sensor**
 White coding plug
 Control with 2 electrodes

This control corresponds to the DMC-Control but is additionally provided with the “SUVA-Sensor” grip stabilization system. Level of grip speed and grip force is determined by signal strength (resulting from muscle tension).

Open: Proportionally controlled by the opening electrode.

Close: Proportionally controlled by the closing electrode.

Example 1: A weak muscle signal generates the smallest grip force (10N) for grasping an object. If the Sensor recognizes that the object is slipping it automatically increases grip force up to 50 % higher (15N) than the initial grip strength. FlexiGrip becomes active at 20N. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

Example 2: A stronger muscle signal generates a higher grip force and the automatic grip force increases to a maximum of 100N to prevent the object from slipping. FlexiGrip becomes active at 130N. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

Program 1: DMC plus Sensor (white coding plug)				
OPEN	CLOSE	Initial Grip Strength	Auto-Grasp Grip Force Increase	Flexi-Grip Function
Myo-signal to the electrode	Myo-signal to the electrode	proport.: 0 N – 100 N	proportional: Up to 50% higher than initial grip strength	dependent on initial grip force, each slightly above the max. grip force
Speed: proportional 15 mm/s to 300 mm/s	Speed: proportional 15 mm/s to 300 mm/s		e.g. Initial grip strength 10 N = increase of grip force up to max. 15 N	FlexiGrip effective: min. at 20 N max. at 130 N

The user can always independently overcome the AutoGrasp feature with a stronger muscle signal and increase the grip force, if necessary, up to a maximum of 100N.

Attention: A short “OPENING” impuls will stop the AutoGrasp as well as the FlexiGrip function.

9.2 Program 2: **AutoControl – LowInput**
 Red coding plug
 Control with 2 electrodes,
 1 electrode and 1 switch, or
 1 switch

The hand closes at maximum speed and grips the object with the minimal grip force (10 N).

If the Sensor senses the object slipping, it automatically increases grip force up to the maximum (100 N).

FlexiGrip becomes active at 130 N. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

Open: proportionally controlled by the opening electrode.

Close: at highest speed by short muscle signal at any level over the ON-Threshold to closing electrode.

Program 2: 2 Electrodes AutoControl – LowInput (red Coding Plug)				
OPEN	CLOSE	Initial Grip Strength	AutoGrasp Force Increase	FlexiGrip Function
Myo-signal to the electrode Speed: proportional 15 mm/s to 300 mm/s	Myo-signal to the electrode digitally (short signal at any level) Speed: constant 300 mm/s	10 N	up to max. 100 N	effective at maximal grip force

9.2.1 Control with 1 electrode and 1 switch

Open: proportionally controlled by the opening electrode.

Close: at highest speed by short activation of switch

Attention: Plug the red/white cable of the 13E99 Connecting Cable to closing side of the 9E169 Coaxial Plug (Fig. 5: (1) Battery Connecting Cable, (2) Electrode, (3) red/white wire of 13E99 Connecting Cable)

Program 2: Electrode and Switch AutoControl – LowInput (red coding plug)				
OPEN	CLOSE	Initial Grip Strength	Autograsp Force Increase	FlexiGrip Function
Myo-signal to the electrode: Hand opens Speed: proportional 15 mm/s to 300 mm/s	Signal to switch: Hand closes Speed: constant 300 mm/s	10 N	up to max. 100 N	effective at maximal grip force

9.2.2 Control with 1 Switch:

This program may only be used with any Ottobock MyoBock Switch.

Open: at maximal speed; as long as the open-side of the switch is activated. The Hand stays open.

Close: at maximal speed by activating the close-side of the Switch.

Program 2: Switch AutoControl – LowInput (red Coding Plug)				
OPEN	CLOSE	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip Function
Hand opens as long as open-side of switch is activated. Speed: constant 300 mm/s	Signal to close-side of switch: Hand closes Speed: constant 300 mm/s	10 N	up to max. 100 N	effective at maximal grip force

9.3 Program 3:

AutoControl

Green coding plug

Control with 1 electrode, 1 linear transducer,
or 1 switch

The Hand closes at maximal speed and grips an object with lowest grip force (10 N).
If the sensor senses a change in position it automatically causes the grip to continuously increase up to the grip force necessary (max. 10N). FlexiGrip becomes effective at 130 N. Once the grasp is secure the SensorHand Speed reverts to the initial grip force.

9.3.1 Control with 1 electrode:

- Open:

at maximum speed with a quick, constant muscle signal to the electrode.
- Close:

at maximum speed by quick **relaxation** of muscle
- Stop:

by very slow **muscle relaxation** to the electrode; Hand stops in an opened position.
- Example 1:

After opening the Hand the amputee relaxes the muscle very slowly. The opening position remains unchanged.
- Example 2:

After opening the Hand the amputee relaxes the muscle at the highest speed. The Hand closes automatically at maximum speed and begins to grip the object with 10N grip force.

If the sensor senses a change in the position of the object, the Hand – if necessary – will automatically readjust up to the maximum grip force (100 N).

Program 3: 1 Electrode AutoControl (green coding plug)				
OPEN	CLOSE	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip FunCtion
Quick constant Myo-Signal to the electrode	Very slow muscle relaxation to the electrode: Hand stops in opened position Quick muscle relaxation to electrode: Hand closes	10 N	up to max. 100 N	effective at maximal grip force
Speed: constant 300 mm/s	Speed: constant 300 mm/s			

9.3.2 Control with 1 linear transducer:

- Open:

at maximum speed by fast pull on linear transducer
- Close:

at maximum speed by fast release on linear transducer
- Stop:

by very slow release on linear transducer. Hand stops in an open position.
- Example 1:

After opening the amputee very slowly releases pull on linear transducer. The opening position remains unchanged.
- Example 2:

After opening the amputee release pull on linear transducer with highest possible speed.

The hand automatically closes with highest speed and starts to grasp object with 10 N grip force.

If the sensor senses a change in the position of the object, the Hand – if necessary – will automatically readjust up to the maximum grip force (100 N).

Program 3: 1 linear transducer AutoControl (green coding plug)				
Open	Close	Initial Grip Strength	AutoGrasp Grip Force	Flexi-Grip Function
High speed pull on linear transducer	Very slow pull release on linear transducer Hand stops in open position	10 N	up to max. 100 N	effective at maximal grip force
	Quick pull release on linear transducer Hand closes			
Speed: constant 300 mm/s	Speed: constant 300 mm/s			

9.3.3 Control with 1 switch:

Open: at maximum speed, as long as switch is activated

Close: after releasing the switch the Hand automatically closes at maximum speed and begins to grasp the object with 10 N grip force.

Attention: Plug the red/white cable of the 13E99 Connecting Cable to closing side of the 9E169 Coaxial Plug (Fig. 6: (1) Battery Connecting Cable, (2) Electrode, (3) red/white wire of 13E99 Connecting Cable)

Program 3: 1 Switch AutoControl (green coding plug)				
OPEN	CLOSE	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip Function
Opens as long as switch is activated.	Closes automatically as soon as switch is released	10 N	up to max. 100 N	effective at maximal grip force
Speed: constant 300 mm/s	Speed: constant 300 mm/s			

9.4 Program 4:

VarioControl

Blue coding plug

Control with 1 electrode or 1 linear transducer

With this program, the opening speed is determined by the level and speed of the muscle signal. The closing speed is determined by the reduction of the muscle tension.

FlexiGrip is effective at maximal grip force. Once the grasp is secure the SensorHand Speed reverts to the initial grip force.

Open: Proportional. Opening speed is determined by level and speed of **muscle tension**.

Close: Proportional. Closing speed is determined by speed and strength of **muscle relaxation**. This also affects the level of the maximal AutoGrasp grip force increase.

Stop: Very slow **muscle relaxation** to the electrode, Hand stops in an opened position.

Example 1: After opening the amputee very slowly relaxes the muscle. Closing is simultaneous to the length of muscle relaxation time. The object is grasped with low force (10 N). There is no automatic grip force increase.

Example 2: After opening the amputee relaxes the muscle at maximum speed. The Hand closes at maximum speed und begins to grasp the object with 10 N. If the Sensor senses a change in position, it automatically readjusts – if necessary – up to the maximum grip force (100 N).

Program 4: 1-electrode control VarioControl (blue coding plug)				
OPEN	CLOSE	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip Function
Speed and level of muscle tension to electrode	Speed and level of muscle relaxation to electrode	10 N	None at low to medium Closing speed	Effective at 15 N with low to medium Closing speed
Speed: proportional 15 mm/s to 300 mm/s	Speed: proportional 15 mm/s to 300 mm/s	10 N	At medium to high Closing speed up to max. 100 N	Effective at 130 N with medium to high Closing speed

9.4.1 Control with 1 linear transducer

Open: proportional. Opening speed is determined by speed and force of pull on linear transducer.

Close: proportional. Closing speed is determined by release speed of pull on linear transducer. This also determines the level of the maximal grip force increase.

Stop: very slow pull release on linear transducer. Hand remains in opened position.

Example 1: After opening the amputee slowly releases pull on linear transducer.

Closing takes place analog to time of muscle relaxation with slow speed. The object is grasped with low force (10 N). There is no automatic grip force increase.

Example 2: After opening the amputee releases the pull on the linear transducer at highest possible speed. The hand closes with highest possible speed and starts to grasp an object with 10 N grip force. If the Sensor senses a change in position, it automatically readjusts – if necessary – up to maximum grip force (100 N).

Program 4: 1 linear transducer VarioControl (blue coding plug)				
Open	Close	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip Function
Speed and force of pull on linear transducer	Release speed on linear transducer	10 N	At low to medium closing speed: none	At low to medium closing speed: Effective from 15 N
Speed: proportional 15 mm/s to 300 mm/s	Speed: proportional 15 mm/s to 300 mm/s	10 N	At medium to high closing speeds: up to max. 100 N	At medium to high closing speeds effective at maximal grip force

9.5 Program 5:

VarioDual
Yellow coding plug
Control with 2 electrodes

With this program, the opening speed is determined by speed and level of muscle tension. Until reaching a minimal grip force of approx. 10 N closing speed is dependent on muscle relaxation speed. Grip force is determined by the following or the simultaneous muscle signal to the second electrode.

FlexiGrip is effective, depending on the initial grip strength, at slightly over the maximal AutoGrasp Grip Force. Once the grasp is secure, the SensorHand Speed reverts to the initial grip force.

Electrode 1:

- Open: Proportional. The opening speed is determined by speed and level of **muscle tension**.
- Close: Proportional. The closing speed is determined by speed and level of **muscle relaxation**. Grip force is approx. 10 N.
- Stop: with very slow **muscle relaxation** to the electrode; Hand stops in an opened position

Electrode 2:

- Gripping: **Grip force build-up is determined by the level of the muscle signal to the second electrode.** Maximal grip force is approx. 100 N.
- Example 1: After opening, the amputee relaxes the muscle. Closing takes place proportionally to the muscle relaxation. The object is gripped with minimal grip force (10 N). FlexiGrip is effective at 20 N. Once the object is securely grasped, the SensorHand Speed reverts to the minimal grip force.
- Example 2: After gripping, as described in Example 1, the object must be grasped with a higher grip force. The amputee must create a muscle signal to the second electrode. Proportionally a grip force between 10 N and 100 N can be built-up. If the object starts to slip the grip force is increased up to 50%. FlexiGrip is effective at approx. the double value of the initial grip strength up to 130 N. Once the object is securely grasped, the SensorHand Speed reverts to the initial grip force.

Program 5: 2 Electrodes VarioDual (yellow coding plug)				
Open	Close	Initial Grip Strength	AutoGrasp Grip Force Increase	FlexiGrip Function
Speed and level of muscle tension to the 1st electrode	Close: Speed and level of muscle relaxation to the 1st electrode	10 N	No AutoGrasp	effective at 20 N
Speed: proportional 15 mm/s to 300 mm/s	Speed: proportional 15 mm/s to 300 mm/s Grip force build-up: Grip force dependent on level of muscle signal to 2nd electrode. Grip force: proportional 10 N to 130 N	Proportional: 10 N to 100 N	proportional: up to max. 1,5-fold of initial grip strength	dependent on initial grip strength, effective at slightly over the maximal grip force increase min. from 20 N max. from 130 N

9.6 Programm 6:

DMC plus Sensor can be switched off **Purple coding plug** **Control with 2 electrodes**

This control mode corresponds to program 1; the SUVA-Sensor as well as the FlexiGrip can be temporarily deactivated.

9.6.1 Switching the SUVA-Sensor and the FlexiGrip Function On and Off

In order to grasp soft and yielding objects such as soft foamed plastics you can manually deactivate the SUVA-Sensor. Open the SensorHand Speed as far as it will go and keep it open by muscle signal. Simultaneously apply some pressure to the SUVA-Sensor (press to a hard object). A short vibrating signal will confirm the switch-off. In order to switch the SUVA-Sensor on, follow the same procedure. Two short vibrating signals confirm the activation of the SUVA-Sensor*.

Note:

Please note that with a switched-off SUVA-Sensorik the grip force will not be automatically adjusted and that grasped objects may slip from the hand.

After inserting a battery, vibrating signals inform the amputee about the current mode:

single vibration	Sensor is switched off
two vibrations	Sensor is switched on

10 Handling the batteries

Only a fully charged battery should be used to operate the SensorHand Speed. Always keep a charged, back-up battery available as a replacement. An intelligent "battery manager" alerts the patient to low voltage conditions: The hand will become slower and the gripping force will decrease. At the same time the battery is protected against harmful total discharge.

We recommend operating the SensorHand Speed with the EnergyPack (757B20) to fully profit for a long time from the full capability.

With restricted capacity, the SensorHand Speed can optionally be operated with the EnergyPack (757B21), the X-ChangePack (757B15), or the MyoEnergy Integral (757B25=*, 757B35=*).

For detailed instructions on how to handle the batteries, please refer to the instructions for use enclosed with the batteries.

Attention:

As soon as the SensorHand Speed recognises a fully charged EnergyPack (757B20, 757B21) or MyoEnergy Integral (757B25=*, 757B35=*), it automatically switches to Li-Ion battery technology. Once this has happened, only this technology should be used to operate the system. If an X-ChangePack (757B15) is employed after the control has switched to Li-Ion battery technology, the battery can not be used to its full capacity. The MyoSelect (757T13) can be used, however, to reset the SensorHand Speed control for use with X-ChangePack.

For operational security and reliability, only 757B20, 757B21, 757B25=*, 757B35=* or 757B15 Batteries are to be used.

Upon delivery the SensorHand Speed is (factory default) set for NiMH batteries.

11 Technical data

Static current	2 mA
Operating temperature range	0–70 °C (32–158 °F)
Opening width	100 mm
Proportional speed	15–300 mm/s
Proportional grip force	0–approx. 100 N
Power supply	EnergyPack 757B20 (7.2 V)
Service life of the terminal device	5 years
Battery service life	2 years

With limitations:

EnergyPack	757B21 (7.2 V)
MyoEnergy Integral	757B25=*, 757B35=* (7.4 V)
X-ChangePack, Interchangeable Battery NiMH	757B15 (6 V)

Ambient conditions

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-con- densing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-con- densing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-con- densing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-con- densing

12 Symbols used

MD Medical device

13 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer assumes no liability for damage caused by component combinations (use of third-party products) and applications which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

14 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Ensure that the product is not damaged and that no liquid penetrates into the system components.

Then, dry the prosthesis component with a soft cloth.

15 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

16 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Vienna · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com