










Instructions for Use

MIAMI TLSO 456

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	18
NO	Bruksanvisning	21
DA	Brugsanvisning	24
SV	Bruksanvisning	27
EL	Οδηγίες Χρήσης	30
FI	Käyttöohjeet	34
NL	Gebruiksaanwijzing	37
PT	Instruções de Utilização	40
PL	Instrukcja użytkowania	43
CS	Návod k použití	46
TR	Kullanım Talimatları	49
RU	Инструкция по использованию	52
JA	取扱説明書	55
ZH	中文说明书	58
KO	사용 설명서	61

ENGLISH

PRODUCT DESCRIPTION

Miami TLSO 456 is a reusable rigid thoraco-lumbo-sacral spine orthosis implemented as a modular compression belt. The device is designed to support and provide motion restriction to a patient's spine during routine daily use within the thoracolumbar and lumbar spine.

The device is a non-invasive, adjustable brace for single patient use. Device length can be adjusted to patient morphology requirements. A pulley system is used to tighten the belt and belt panels to the patient.

List of Parts (Fig. 1)

- a. Chest Buckle
- b. D-Ring
- c. Overlap closure/belt ends
- d. Shoulder straps
- e. PTE
- f. Belt
- g. Strap anchor

INDICATIONS FOR USE

Gross thoracic and lumbar spine immobilization in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative spinal pathologies
- Post-surgical stabilization
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

CONTRAINDICATIONS

Do not use with unstable, displaced fractures.

INTENDED USE

Patients requiring a gross thoracic and lumbar spine immobilization in all three planes. Sizing is appropriate for most adults; the product is non-gender specific.

GENERAL SAFETY INSTRUCTIONS

Warning:

The Miami TLSO with posterior thoracic extensions is a supportive brace only and is not intended or guaranteed to prevent spinal injury. It should be applied by a licensed healthcare practitioner and as directed by a physician or other qualified medical authority. Össur accepts no responsibility for any injury which might occur while using the orthosis.

Precaution:

Do not overtighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system.

If brace causes discomfort, irritation, or is damaged, contact your healthcare provider for further instructions and discontinue use.

This brace is single patient use only.

The braces are not designed to be worn while bathing.

INSTRUCTIONS FOR FITTING

MIAMI TLSO BRACE (PTE) APPLICATION

1. Log roll the patient onto their side. Omit this step if fitting in standing position. If standing, patient must stand with back straight and chest out.
2. Position the posterior panel centered over the spine, with the bottom of the posterior panel at approximately the sacrococcygeal joint (**Fig. 2**).
3. Log roll the patient back into the supine position. Omit this step if fitting in standing position.
4. Wrap the lumbar belt around the waist and through the Össur Fit Tool (**Fig. 3**). Adjust the length of the belt arms to position tool centrally. Secure the hook and loop tabs to desired sizing (**Fig. 4**).
5. Remove Össur Fit Tool from the belt.
Note: If an Össur Fit Tool is not available, use sizing guide (**Fig. 5**) to estimate best patient fit based on waist circumference. Simply fold or trim the belt arms at the suggested location. If the patient's waist exceeds 50" (127 cm), an extension belt can be added to each belt arm to reach a maximum belt circumference of 70" (178 cm).
6. Place belt ends into the Overlap Closures (**Fig. 6**). To obtain the lowest profile fit, the belt arm ends can be trimmed if later adjustments are not required.
Note: Once belt ends are secured in Overlap Closures, the Overlap Closures will serve as belt ends.
Note: Panels can be customized by heat-molding, grinding, or trimming, to achieve correct patient fit.
7. Close the belt arms around the waist.
8. Remove strap anchor from belt, place the shoulder straps over the shoulder, and re-fasten the strap anchor onto the belt (**Fig. 7**).
Note: Place the strap anchor on a smooth area of the belt as posterior to the patient as possible without introducing axilla discomfort.
9. Slide fingers through the holes in the two compression system handles and pull until the belt is at the appropriate tightness (**Fig. 8**). Adjust the compression system cord length if necessary. (**Fig. 9**)
10. If shoulder pad position needs to be improved, log roll the patient, remove alligator tabs, trim straps, replace alligator tabs, and log roll the patient back into the supine position.
11. Pull webbing to tighten shoulder straps.
12. Adjust the chest strap position by sliding the d-ring (**Fig. 10**).
13. Fasten chest buckle and tighten chest strap.
14. Upon sitting, standing or ambulation, the brace may require fine tuning for comfort, fit and stabilization. This may include adjusting or tightening straps, compression system, and/or rigid panels.

BRACE REMOVAL

1. Detach chest strap buckle
2. Detach compression system handles and loosen. Affix them onto the belt after loosening.
3. Detach the overlap closure and remove brace.
4. Before reapplying the brace and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart.

PANEL ADJUSTMENTS

The posterior panels can be modified as necessary to optimize patient fit and comfort. Simply remove the panels from their fabric sleeve and adjust the rigid panel with a heat gun and/or trimming device.

MAINTENANCE & CLEANING

The Miami TLSO system is designed for minimum maintenance and care. The padded covers may be removed from the rigid panels for washing.

- Hand-wash using mild soap and rinse thoroughly
- Air dry flat
- Do not machine wash or place in tumble dryer

PRODUCT MATERIAL LIST

Nylon, glass-filled nylon, polyamide, spandex, thermoplastic polyurethane, polyurethane, ester polyurethane, polyester, polyurethane foam, bamboo, acrylonitrile butadiene styrene, steel, stainless steel, ethylene vinyl acetate, glass, polyethylene, low density polyethylene, aluminum.

ENVIRONMENTAL CONDITIONS

The device can be used in normal conditions (home or health care), both indoors and outdoors.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

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EN – Caution: This product has been designed and tested as a single use product, and is not intended for multiple use. In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure. If any problems occur with the use of this product, immediately contact your medical professional.

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