







Instructions for Use



TOTAL KNEE® 1900


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	12
FR	Notice d'utilisation	19
ES	Instrucciones para el uso	26
IT	Istruzioni per l'uso	33
DA	Brugsanvisning	46
SV	Bruksanvisning	53
EL	Οδηγίες χρήσης	59
NL	Gebruiksaanwijzing	73
PT	Instruções de Utilização	80
PL	Instrukcja użytkowania	87
TR	Kullanım Talimatları	94
RU	Инструкция по использованию	100
JA	取扱説明書	107
ZH	中文说明书	113
KO	사용 설명서	119


8


9

10


11


12

A

B

DESCRIPTION

The device is a polycentric knee joint with a Geometric Locking System, and a unique swing phase control through the use of friction within an elastic polymer medium.

Proximal adapter options are available, and the device has a Distal Tube Connector.

Interchangeable Bumpers (**Fig. 7A**) allow for knee flexion to be adjusted to suit the patient's preference and activity levels. The adjustable Extension Promoter (**7D**) assists in limiting heel rise and moving the device into a fully extended position.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for low to moderate impact use, e.g., walking.

The weight limit for the device is 100 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Proximal Adapter Options

The following proximal adapter options are available for the device (**Fig. 1**):

- 4-Prong Socket Adapter (1A),
- Female Pyramid Adapter (1B),
- Euro 4-Hole Adapter (1C), or
- Male Pyramid Adapter (1D).

Caution: Only use the specified Total Knee proximal adapters, please refer to Össur catalogue for more information.

1. Fully engage the threads of the selected proximal adapter.
Caution: Adapters from other manufacturers may not fully engage all the threads and compromise the strength of the prosthesis.
2. Align the Pinch Bolt (1E) of the Male Pyramid Adapter (1D), 4-Prong Socket Adapter (1A), and the Euro 4-hole Adapter (1C) perpendicular to the anterior/posterior axis of the device, except for the Female Pyramid Adapter (1B), which should be angled 45°.
Caution: Maximum connection strength will not be established unless the adapter is correctly aligned.
3. For secure fastening to the device, tighten the adapter's Pinch Bolt according to its instructions for use.
Caution: Do not tighten the proximal adapter against the housing of the device.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

Caution: Contact between socket/adapter and proximal 2/3 portion of the device's back link must be avoided when fully flexing the knee (**Fig. 2**). Contact in this area can damage the device. Acceptable contact is between the socket and the distal 1/3 portion of the device's back link.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through pivot axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through pivot axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 4**).

- Proximal Adapter Pinch Bolt (**Fig. 1E**): Please refer to the respective proximal adapter's instructions for use.
- Tube Clamp Screw (**Fig. 5A**): 15 Nm

Bumper Selection (Fig. 6,7)

The selection of a bumper determines the amount of stance flexion.

1. Carefully insert a flat head screwdriver or a similar flat tool under the bottom of the Bumper Cover (7C) and remove it (**Fig. 6**).
2. Remove the Bumper (7A) by inserting a small and flat hand tool underneath it and lift upwards. Leave the underlying Shim (7B) in place.
3. Select a Bumper to install based on which Total Knee device is being used, and on the patient's weight:

Weight kg	Weight lbs	Color	Type
35–50	80–110	Green	Soft
50–70	110–155	Black	Regular
70–90	155–200	Red	Firm
90–100	200–220	Yellow	Extra Firm

4. Re-install the Bumper Cover (7C) by hooking the groove under the chassis front and pushing the two pins into the chassis holes.

Influencing Factors on Bumper Selection:

- Active patients may require a stiffer bumper.
- Patients who are accustomed to wearing a prosthesis without a stance-flexion feature may initially prefer a stiffer Bumper. After a trial period of two weeks, selecting a softer Bumper is recommended. Increasing the softness of the Bumper allows the patient to get familiarized with the device and facilitates the transition to using its stance-flexion feature.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Ensure the patient can activate and deactivate the geometric lock.

Swing Control

The device is designed to provide smooth and fluid friction swing control through the mechanical manipulation of a unique polymer medium. Additionally, the swing control can be enhanced by using the Extension Promoter and optional Extension Assist.

Caution: For safety, please make initial assessments and adjustments with the patient standing between parallel bars.

Polymer Friction Adjustment (Fig. 8)

1. Loosen the Set Screw (8A) by a ½ turn. Do not completely remove the Set Screw from the right side of the device's housing unit.

2. On the anterior side of the device, turn the Friction Adjustment Screw (8B) clockwise to increase the resistance to flexion and extension. Do not exceed $\frac{1}{4}$ of a clockwise turn, which is the maximum setting for the polymer friction (8C).
3. After obtaining the desired friction setting, tighten the Set Screw to secure the polymer friction setting.

Caution: Do not overtighten the Set Screw. Apply 2 Nm (17 in-lbs) of torque.

Bumper adjustment

As needed, try different Bumpers, and select the one that best suits the need of the patient.

The final bumper selection might differ from the initial selection based on the patient's weight (please refer to the previous section: Bumper Selection).

The torque required to release the Geometric Lock and the stance flexion motion:

- increases when using a soft Bumper, and
- decreases when using a firm Bumper.

Depending on which foot is selected, combine:

- softer bumper with soft heels, and
- firmer bumper with firm heels.

Bumper height adjustment

A thin plastic plate, or a Shim (**Fig. 7B**), is placed underneath the Bumper (7A) to adjust its height. There are two types of Shims:

- Grey: 0.38 mm thick, and
- Clear: 0.76 mm thick (the device is shipped with the Clear Shim installed).

Adding Shims changes the release point of the Geometric Lock that allows for knee flexion during the swing phase. Adjusting the height of the Bumper affects knee stability, along with ease of walking and sitting:

- Higher Bumper: facilitates the release into swing flexion and decreases stance flexion.
- Lower Bumper: increases the device's stability and stance flexion, while delaying the release into swing flexion.

Gradually increase the height of the bumper until the patient feels premature release of the knee. Reduce the number of Shims until the patient finds a suitable balance between device stability and at which point the knee releases into swing flexion.

Caution: Do not install more Shims than the two different types combined. Adding more Shims can deactivate the Geometric Lock, and negatively affect the stability of the knee during stance.

Note: If no Shims are used, the device will exhibit some flexion/extension play during stance.

Extension Promoter

The Extension Promoter reduces excessive heel rise and promotes faster knee extension. It is delivered in the minimum setting.

- Flex the knee to access the Extension Promoter Adjustment Screw (7D).
- To reduce excessive heel rise, especially when walking fast, turn the Extension Promoter Adjustment Screw clockwise.

Extension Assist

The Extension Assist is an optional feature that facilitates full knee extension. Use of the Extension Assist is only recommended when the patient is unable to extend the knee to engage the Geometric Lock. Three springs of different stiffness are included with the Extension assist:

- Yellow: regular (comes pre-installed with the unit),
- Blue: firm, and
- Red: extra firm.

Installation of the Extension Assist (Fig. 9)

1. Remove the Bumper Cover (A), Bumper (B), and Shims (C).
2. Install the Spring Cup (D) into the Hole (E) inside the Distal Tube Connector.
3. Firmly seat the Spring Cup by gently tapping it with a rubber mallet (or a similar hand tool).
4. Place the Ferrule (F) over the Cable (G).
5. Place the Spring (H) over the Cable and Ferrule (try the yellow spring first).
6. Insert the Cable and Spring into the Spring Cup, with the ball end of the Cable first.
7. Slide ball end of the Cable through the Hole (E) in the Tube Connector and insert the ball into the Slot on the device's Bottom Link (I) by using needle-nose pliers or a small screwdriver.
8. Re-install the selected Shims, Bumper, and Bumper Cover.

Follow-up check

- Reassess adjustments after a couple of weeks.
- It may be necessary to change the Shim and Bumper adjustments after the patient has gained some experience with the device.
- Over time, the patient's increased device experience, confidence during use, and proficiency may warrant a change to a firmer Bumper and/or an increase in Shim thickness.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

The stance flexion Bumper must be replaced if cracked or worn, depending on the patient's activity level.

Check for damage, excessive wear and dirt.

If the device is subjected to excessive moisture or a corrosive environment, it is recommended to clean and lubricate the device frequently.

In a slightly flexed position, three sets of roller bearings will be visible (**Fig. 10A**). Place few drops of general purpose oil or sewing machine oil on the roller bearings. Move the knee several times and wipe off excess oil with a soft cloth. No other parts need external lubrication. Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

SUPPLEMENTARY DEVICES

Cosmetic Options (Fig. 11,12)

The device is intended to protect the prosthesis and to restore the appearance of a missing limb.

- A discontinuous cosmetic solution that does not interfere with knee function is recommended.
- Continuous cosmetic solutions are also a possibility.

Following cosmetic solutions are available, please refer to Össur catalogue for further information:

- Shin Fairing, Black (11A),
- Knee Cap, Large, Beige (11B),
- Knee Cap, Small, Black (11C),
- Cosmetic Foam Cover, Regular (12A), and
- Cosmetic Foam Cover, Large (12B).

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3–5 years of use.

ISO 10328 - P5 - 100 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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