X-ACT FIT SHOULDER

- ESTABILIZADOR DE HOMBRO X-ACT FIT
- X-ACT FIT-SCHULTER
- X-ACT FIT SHOULDER
- ORTHÈSE D'ÉPAULE X-ACT FIT
- X-ACT FIT-SCHOUDERBRACE
- SUORTE DE OMBRO X-ACT FIT
- RAMENNÍ ORTÉZA X-ACT FIT
- RAMENNÁ ORTÉZA X-ACT FIT
- X-ACT FIT-SKULDER
- X-ACT FIT AXEL
- X-ACT OLKAPÄÄ KUNTOON
- X-ACT SUNN SKULDER
- ORTEZA X-ACT FIT SHOULDER
- ORTEZĂ PENTRU UMĂR X-ACT FIT

ENGLISH..... 2

ESPAÑOL 4

DEUTSCH..... 6

ITALIANO..... 8

FRANÇAIS..... 10

NEDERLANDS..... 12

PORTUGUÊS..... 14

ČEŠTINA..... 16

SLOVENČINA..... 18

DANSK.....20

SVENSKA.....22

SUOMI..... 24

NORSK.....26

POLISH28

ROMANIAN30

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











DJO, LLC
5919 Sea Otter Place,
Suite 200
Carlsbad, CA 92010
U.S.A.

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS: The X-Act Fit shoulder brace provides support and immobilization of the shoulder, elbow, wrist and hand following surgery and other shoulder conditions. Rotator cuff repairs; Bankart and capsular shift procedures; Soft-tissue strains; Shoulder dislocation/subluxation; Non or minimally displaced proximal humerus fractures; Pectoral tendon rupture and distal biceps rupture; AC-Joint instabilities; Proximal biceps tenodesis, SLAP-Repairs; ORIF clavicle fractures; Humeral shaft fractures; Scapular fractures; SC joint separation, UCL recon of the elbow and other elbow fracture indications; Non- and Post-surgical rehabilitation.

CONTRAINDICATIONS: None

WARNING AND CAUTIONS:

- Device is supportive only and does not guarantee or intend to prevent injury.
- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not use over open wounds.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- DJO, LLC accepts no responsibility for any injury which might occur while using DJO, LLC products.
- Do not use this device if it was damaged and/or packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

PRE-APPLICATION INFORMATION:

SIZING:

1. Measure patient's waist circumference 1" (2-5cm) above belly button.
2. The WAIST BELT is packaged for the largest configuration and will fit up to 58" waist.
3. To adjust the WAIST BELT for smaller sizes, remove the Y-Tab closure on end of belt and trim WAIST BELT to fit patient. Reattach the Y-Tab closure. Each trim line is 2"(5cm) apart.

PREPARING THE BRACE:

1. Prepare the brace by opening the WAIST BELT CLOSURE MITT.

APPLICATION INFORMATION:

1. FORM THE WAIST BAND :

1. Shape the malleable waist band to fit patient's torso.

NOTE: Do not remove plastic sleeve until step 3.

2. POSITION WAIST BAND & SIZE WAIST BELT:

1. (Right Shoulder shown) Place the malleable WAIST BAND on the affected side of body approximately 1" (2.5 cm) above the belly button. (Fig C)
2. Wrap the WAIST BELT around the torso and attach the BELT CLOSURE STRAP to FRONT PANEL of WAIST BAND. (Fig D)
3. Tighten WAIST BELT using the BELT ADJUSTMENT STRAP as needed for secure fit. (Fig E) The ARM REST BAR should be pointing forward and in a neutral position.

NOTE: Excess WAIST BELT may be trimmed as needed. Remove the Y-tab on the WAIST BELT, trim belt and re-apply Y-tab. Secure WAIST BELT to WAIST BAND.

4. Movable HIP PAD may be positioned on inside of WAIST BELT on unaffected side to provide additional comfort & support.

3. ADJUSTING FOR ABDUCTION / ER / IR POSITIONS:

The DonJoy X-Act Fit Shoulder is packaged and configured for a neutral abduction position.

UP TO 45 DEGREES OF ABDUCTION OR ER/IR

1. Remove Plastic Sleeve on ARM REST BAR.
2. Shape the malleable ARM REST BAR into the desired arm position for up to 45 Abduction or for External/Internal Rotation. Secure the AB STRUT under the ARM REST BAR. (Fig F)

45-TO-90-DEGREE ABDUCTION

1. Remove Plastic Sleeve on ARM REST BAR.
2. Locate the button on the EZ AXIS HINGE plate on the WAIST BAND. Push in the BUTTON and pull away from WAIST BAR to rotate the ARM REST BAR so it's at a 90-degree angle to WAIST BAND. Push down on EZ AXIS HINGE plate and make sure there is an AUDIBLE CLICK locking the ARM REST BAR in place.

NOTE: ARM REST BAR should be pointing down toward patient's hip. (Fig G)

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3. Form the ARM REST BAR to desired level of Abduction.
4. Remove the CINCH STRAP around the AB BAR. Detach AB EXTENSION and position along AB STRUT for optimal support of ARM REST BAR. Tightly secure the CINCH STRAP around the AB STRUT & EXTENSION. The hook on the top of the AB EXTENSION should be securely attached to the ARM REST BAR. (Fig H). The AB STRUT 'lip' may be formed if desired, to be flush with AB Extension.

NOTE: For some patients, the WAIST BAND may need to be placed slightly higher on patient's torso to obtain 90-degree abduction.

4. SIZING AND APPLICATION OF ARM SLING:

APPLY ARM SLING

1. Remove forearm, thumb, and bicep straps from the ARM ENVELOPE. Open QUICK RELEASE BUCKLES at front of ARM ENVELOPE. Place forearm into ARM ENVELOPE as far back as possible.

ADJUST ARM SLING LENGTH

1. To adjust ARM ENVELOPE for a smaller size, remove the hook closure tabs with QUICK RELEASE BUCKLES from front of ARM ENVELOPE, fold the end of the ARM ENVELOPE inward and adhere the internal hook and loop tabs at the desired length. (Fig I).
2. Re-attach the hook closure tabs with QUICK RELEASE BUCKLES to either side of the ARM ENVELOPE (as far forward as possible), making sure the QUICK RELEASE BUCKLE is positioned above top of ARM ENVELOPE. (Fig J)
3. Place arm in ARM ENVELOPE and secure THUMB STRAP between the thumb and forefinger at front of sling. Secure FOREARM and BICEP STRAPS at desired locations.
4. Place shoulder strap over the opposite shoulder. Connect shoulder strap QUICK RELEASE BUCKLES to its mate on front of ARM ENVELOPE. (Fig K). Trim and tighten straps as needed.

NOTE: Final trimming of straps should be completed AFTER abduction bar is positioned.

5. ATTACH SLING TO ARM REST BAR:

1. Position ARM SLING ENVELOPE along the hook and loop strip of the ARM REST BAR with elbow a few inches behind the bend on the ARM REST BAR. Make sure the ARM ENVELOPE is firmly attached to the ARM REST BAR(Fig L).
2. Adjust and trim shoulder strap, front and back, as needed

OPTIONAL: SHOULDER RING APPLICATION

1. Remove the standard shoulder strap from the ARM ENVELOPE.
2. Lengthen FRONT and BACK SHOULDER STRAPS to their longest length for first time fitting. Attach the posterior strap to the O-ring on back of the ARM SLING ENVELOPE.
3. Open the SHOULDER RING and wrap under unaffected arm, then over top of shoulder. Secure SHOULDER RING with hook and loop closure (Fig. M). Make sure the hook and loop closure of SHOULDER RING is away from patient with WIDEST ENDS over top of shoulder. Trim and/or adjust the SHOULDER RING circumference as needed.

POSITIONING OF SWIVEL D-RING STRAP TABS

1. Adjust the front and back SWIVEL D-RING STRAP TABS on SHOULDER RING so they are approximately 1"- 2" from top of shoulder (Fig. N). Placement of the SWIVEL D-RING STRAP TABS on the shoulder ring helps to maintain straps and arm in proper position.
2. Secure SHOULDER STRAP with Quick Release Buckles to front of sling.
3. Adjust and Trim excess SHOULDER STRAP as needed.

CLEANING INSTRUCTIONS: Separate the rigid components from the soft goods prior to cleaning.

Hand wash in water (86°F/30°C) using soap. Air dry.

MATERIAL CONTENTS: Nylon, Aluminum, Stainless Steel, Polyurethane Foam, Polyurethane.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

RX ONLY.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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