

# ANKLIZER® II PUMP WALKER • ANKLIZER® II PUMP WALKER - LOW TOP

## Application Instructions

0814 1151 - 0814 1155 • 0814 1161 - 0814 1165



**Step 1.** Open hook and loop closure straps and remove from D-rings. Open liner. Place the foot and leg inside the liner, making sure the heel is firmly back against the liner. If you find there is extra space, pads may be inserted at either side of the ankle or at the heel.



**Step 2.** Close the liner over toes and top of foot. Close the liner around the shin. Tightly secure the hook and loop closures. Fold down the circular tab over the front of the ankle.

An additional padded strap is attached for additional padding at the front of the ankle.

Position the heel of the foot against the rocker boot. Brace uprights should fall midline over the ankle.



Feed closure straps through D-rings and fasten closure straps around uprights and liner starting at the toes and working upward, alternating the direction until all straps are secured. Strap adjustment may be necessary as swelling patterns change throughout the rehabilitation process.



Upright

**Step 3.** Determine the position of the uprights with respect to the liner. If the upright falls on top of an overlap, the liner should be opened to allow the upright to lie underneath permitting removal of the brace without disturbing the position of the upright.



**Step 4.** Slide the sleeves from the uprights and press the uprights firmly against the liner.



**Step 5.** Fasten closure straps around uprights and liner by alternating the direction of the top three straps until they are secured.

**Note: Hi-Top has 3 straps for the upper part of the walker and the Low-Top has 1 strap for the upper part of the walker.**



**Step 6.** Inflate air chamber by first turning knob clockwise until knob offers gentle resistance; then pressing and releasing the air bulb repeatedly until the liner feels stable and supportive on the leg. There should be equal pressure on both sides of the leg, from calf to ankle. **Do not over-inflate chamber.**

Air can be released from the air chamber by slowly turning the knob counter-clockwise, and halting air release by turning the knob clockwise until the knob offers gentle resistance. To re-inflate, repeat beginning of **Step 6**. When prescribed by physician, remove leg from product by releasing closure straps and opening liner.



HI-TOP



LOW-TOP

Reorder Number		Size	Female Shoe Size	Male Shoe Size	Unit
HI-TOP	LOW-TOP				
0814 1151	0814 1161	X-Small	4 1/2	3 1/2	Each
0814 1152	0814 1162	Small	5 - 8 1/2	4 - 7	Each
0814 1153	0814 1163	Medium	9 - 11	7 1/2 - 10	Each
0814 1154	0814 1164	Large	11 1/2 - 13	10 1/2 - 12	Each
0814 1155	0814 1165	X-Large	13 1/2 Plus	12 1/2	Each

### BIRD & CRONIN INC.

Quality Medical Products To Help Improve Your Life  
1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail sales@birdcronin.com • www.birdcronin.com

©2013 Bird & Cronin, Inc. All Rights Reserved 10001151