

SPORTSMAN'S CHOICE™ TENNIS ELBOW SUPPORT

Application Instructions



1. With your free hand, position the tennis elbow support below the elbow as shown in the picture.
2. Slip hook and loop closure strap through D-ring and fasten securely.
3. Re-adjust as needed for comfort.
4. If pain increases or persists, discontinue use and contact your healthcare provider.

HAND WASH COLD • AIR DRY

INDICATIONS: Soft tissue injuries where support and heat are prescribed.

CONTRAINDICATIONS: If rash or blisters occur, apply talcum powder or stockinette sleeve under the support. If this does not clear up the problem, discontinue use. Supports SHOULD NOT be worn while sleeping or by persons with circulatory problems.

WASHING INSTRUCTIONS: Rinse the support in mild soapy or clear, cool water after every use. AIR DRY. DO NOT use solvents or chlorine bleach and DO NOT put in dryer.

RUBBER MANUFACTURERS WARNING: Contact of this material with the skin may cause irritation. Persons with known susceptibility to dermatitis or with broken skin should avoid contact with this material. Anyone who develops skin irritation from contact with this material should discontinue its use immediately and consult a physician.

WARRANTY: Covers defects in workmanship or materials for a period of three (3) weeks from purchase date. DOES NOT accept any liability for any damage or injury rising from defective materials or workmanship or problems covered from CONTRAINDICATIONS.



Single Patient Use

BIRD & CRONIN INC.

Quality Medical Products To Help Improve Your Life
1200 Trapp Road • Eagan, Minnesota USA 55121
651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934
e-mail sales@birdcronin.com • www.birdcronin.com

©2012 Bird & Cronin, Inc. All Rights Reserved - 10007426A