

Included in package:



1 Fabrication Dummy

Manufactured by:



419 N. Curtis Rd., ID 83706 USA
(208) 429-0026
www.coyotedesign.com



CD103FD

Fabrication Dummy

For CD103 Air-Lock or CD117 Easy Off Lock
CD103D Deep Air-Lock or CD117D Easy Off Deep lock



2-Year warranty against manufacturer defects, excessive wear or breakage

Made in U.S.A.

External Prosthetic Components



CD103FD.revA.02082016



EC REP

Advena Ltd
Pure Offices, Plato Close
Tachbrook Park
Warwick, CV34 6WE, UK

Included in package:



1 Fabrication Dummy

Manufactured by:



419 N. Curtis Rd., Boise, ID 83706 USA
(208) 429-0026
www.coyotedesign.com



CD103FD

Fabrication Dummy

For CD103 Air-Lock or CD117 Easy Off Lock
CD103D Deep Air-Lock or CD117D Easy Off Deep lock



2-Year warranty against manufacturer defects, excessive wear or breakage

Made in U.S.A.

External Prosthetic Components



CD103FD.revA.02082016



EC REP

Advena Ltd
Pure Offices, Plato Close
Tachbrook Park
Warwick, CV34 6WE, UK

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103FD



Put CD103FD fabrication dummy over anchor.



Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



Pressing in around fabrication dummy will help remove any air pockets.



If needed use parachute cord or rope to pull plastic tight to dummies edge.



Let plastic cool before cutting fabrication dummy free.



Cut tooling free with cast saw, knife or utility knife.



Remove tooling from anchor.



Sand flexible inner down so anchor is almost 100% exposed



Pull a PVA bag over the anchor and cast. Glue your lock to anchor and start laminating.

Flexible Inner Socket with Coyote Designs Fabrication Dummy CD103FD



Put CD103FD fabrication dummy over anchor.



Pull preferred flexible plastic over socket and installed fabrication dummy using a vacuum system.



Pressing in around fabrication dummy will help remove any air pockets.



If needed use parachute cord or rope to pull plastic tight to dummies edge.



Let plastic cool before cutting fabrication dummy free.



Cut tooling free with cast saw, knife or utility knife.



Remove tooling from anchor.



Sand flexible inner down so anchor is almost 100% exposed



Pull a PVA bag over the anchor and cast. Glue your lock to anchor and start laminating.