smooth curves ahead

*curvas suaves à frente • In Zukunft läuft alles glatt • curvas suaves a la vista
الانحناء بسهولة • courbes douces devant • curvas suaves à vista
впереди плавные кривые • 前部平滑曲线*

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
- 3mm Hex Key
- 4mm Hex Key
- FootHorn

MOUNTING

Use only high quality endoskeletal components.

ASSEMBLY AND DISASSEMBLY (for sock replacement)


Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

HYDRAULIC RANGE

The Odyssey_{K2} foot has 12° of hydraulic motion. The foot is designed to provide 1° hydraulic dorsiflexion from the neutral standing position. This is intended to prevent midstance drift; the sensation of user instability (lacking toe support) during standing.


The foot was developed utilizing a flexible composite base. As a result, a typical user will experience an additional range of dynamic motion during ambulation.

Note: Excessive angular adjustment will affect the hydraulic range of the foot. After making an alignment change, ensure that the user retains 1° of hydraulic dorsiflexion.

BENCH ALIGNMENT

The weight line should divide the foot at 1/3 heel lever to 2/3 toe lever. It will fall just anterior of the pyramid.


STATIC ALIGNMENT

Using a 3mm hex key for adjustment, the hydraulic valves should be set at minimum resistance. Have the user stand comfortably and evaluate the heel-toe balance of the foot. They will sense the hydraulic movement of the ankle but should not feel like they are falling forward or backward. Use alignment to position the foot at the point where they feel the most balanced.

SYMPTOM	ALIGNMENT CHANGE
Falling Backward	Shift foot posterior relative to the socket
Falling Forward	Shift foot anterior relative to the socket

DYNAMIC ALIGNMENT

Have the user begin by walking on level ground, to evaluate the heel-toe resistance and gait timing. Using a 3mm hex key, adjust for planterflexion resistance first, then dorsiflexion. Finalize dynamic alignment by observing the user walking on an inclining-declining surface (ramp). Make further adjustments to the resistance valves as necessary.

Dorsiflexion resistance affects the user's gait through Midstance, as the body travels over the foot.

Planterflexion resistance affects the user's gait from Heel Strike to Foot Flat.

DESIRED RESULT	VALVE ADJUSTMENT
Firmer Planterflexion	Turn P-valve clockwise (more resistance)
Softer Planterflexion	Turn P-valve counterclockwise (less resistance)
Firmer Dorsiflexion	Turn D-valve clockwise (more resistance)
Softer Dorsiflexion	Turn D-valve counterclockwise (less resistance)

⚠️ PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

ADDITIONAL CONSIDERATIONS

Have the user practice standing up from a seated position in order to acclimate to the motion of the ankle. Use caution when driving. Make sure the user is comfortable with the motion of the ankle if using the Odyssey_{K2} as their driving foot.

GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.

WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder may effect the function of the CPI Sock and lead to noise.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

WARRANTY INSPECTION / MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Hydraulic Assembly
- CPI Sock
- Composites and Adapters
- Foot Shell

WARRANTY INSPECTION SCHEDULE FOR ODYSSEY_{K2} : 6 MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365


College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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The Odyssey_{K2} is covered by the following patents:

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MADE IN THE USA

©2018 College Park Industries, Inc. All rights reserved. College Park, OdysseyK2, Intelliweave, and Gait Matching are registered trademarks of College Park Industries, Inc.

college park

800.728.7950 | 586.294.7950 | college-park.com

COLLEGE PARK INDUSTRIES, INC
27955 College Park Dr. Warren, MI 48088 USA

EMERGO EUROPE
Prinsessegracht 20, 2514 AP The Hague, Netherlands

**Australian Sponsor
EMERGO AUSTRALIA**

Level 20, Tower II, Darling Park, 201 Sussex Street,
Sydney, NSW 2000 Australia

430 INS OK2 TIS 180316