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## USE

This knee is intended exclusively for lower limb prosthesis.

It is specifically recommended for any patient (from low activity level to active) up to 125 kg (carried load included).

## FEATURES

In extension this joint system with linkages puts the instantaneous center of rotation backward the load line. This ensures a good stability in stance phase. Every flexion and extension speed is adjustable independently (See page 7: Fig. 3).

The extension assist strength is adjustable through an extension spring cylinder (See page 7 : Fig.3) *Note: a « short extension spring cylinder » kit is available (ref. 1P11067) and allows to shorten the total height of the knee by 22 mm. This kit includes a short cylinder, an adapted spring and assembling instructions.*

This knee weighs about 1050gr for the 1P120 and 1P121 models; and about 1080gr for the 1P120-KD and 1P121-KD models.


The knee is designed to allow flexion up to 145°. However flexion can be limited by the socket dimensions or the cosmetic cover.

## ASSEMBLING (Fig. 1a et 1b)


Connections with the other prosthetic components

This knee is available in four different models :

- 1P120 (grey) and 1P121 (orange): upper connection by a male pyramid, convenient for TF amputation or hip disarticulation. Lower connection for Ø34 tubular components 1G01-P6 and 1D52-P6 (See Fig.1b for tightening)
- 1P120-KD (grey) and 1P121-KD (orange): upper connection by a M36x1,5 screw, convenient for long stump or knee disarticulation. Lower connection for Ø34 tubular components 1G01-P6 and 1D52-P6 (See Fig.1b for tightening)

**FINISHING:** One-part cosmetic cover 1G16 (complete or partial)

**Inform your patient:** The warranty is not applicable to damages due to a bad assembly, an inadequate alignment, a very dusty environment without adapted protection, or any inappropriate use.

Using a cosmetic cover is highly recommended to prevent your fingers or clothing from being jammed in the joint. Avoid exposing the knee to situations that may cause the corrosion of the metallic parts (fresh water, sea water, acid, ...). In case of fall, shock or abnormal behaviour of the joint, contact your CPO for knee control.


Fig. 2

**STABILITY ADJUSTMENT** (Fig 2)**To increase dynamism**

- Use the different adjustment screws of the connectors to bring the load line closer to the knee joint axis (near 0), but not over it, and/or bring the knee joint axis forwards.
- Check the alignment with an appropriate tool (plumb line or laser).

Note : This type of alignment decreases stability, which eases knee flexion during walking.

**To increase stability**

- Use the different adjustment screws of the connectors to bring the load line away from the knee joint axis (up to 10 mm) and/or bring the knee joint axis backwards.
- Check the alignment with an appropriate tool (plumb line or laser).

Note : This type of alignment increases stability, which makes the patient feeling safer during walking.

Note: The KD model of this knee (connection with M36x1.5 thread) does not allow the stability to be modified. The socket should be made accordingly.

The HYDEAL II knee joint comes adjusted for medium walking pace (knobs **E** and **F** in medium position and assist cylinder **R** screwed home unless three turns).

To adapt the initial adjustments to the patient's walking pace, adjust knobs **E** and **F** using the hex wrench delivered with the knee. We recommend you to make the adjustments in the following order :

**1: FLEXION** >> screw on knob **F** to decrease flexion angle, or screw out knob **F** to increase flexion angle (*adjustment range : from one turn on to one turn out*)

**2: EXTENSION** >> Screw on knob **E** to decrease the impact at the end of extension (*Adjustment range : from one turn on to one turn out*)

**3: EXTENSION ASSIST:** If the optimal adjustment of the patient's walking pace cannot be reached >> use assist cylinder **R** to modify the strength of the extension assist : screw on to increase the strength of the assist (quicker return) and screw out to decrease the strength of the assist (slower return).


Fig. 3

When all the adjustments are finalized, check that the knee can reach full extension with low walking speeds (risk of falling).

This hydraulic knee is designed mainly to control swing phase. For example the knee may be damaged if the hydraulic system bears the body weight when resistance is adjusted at a high level during flexion.

## MAXIMAL FLEXION LIMIT (Fig.4)

In full flexion, socket **E** should never touch linkages **A** and **B**. The socket should bear on extension limit **J** in order to protect the knee.


Fig. 4

## MAINTENANCE

Check at least once a year the knee functions and make the necessary maintenance operations if needed.

Spare parts :

- extension limit : **EE015** (x1) >> to be sticked with XC051 Proteor adhesive
- maximal flexion limit : **EY001** (x1) >> self-adhesive

In order not to damage the knee, do not use talcum powder to eliminate possible noises due to cosmetic cover friction, but a silicon spray . Talcum powder damages the mechanical parts, which may lead to function defects and risk of falling for the patient.  
**The application of talcum powder will render all claims against PROTEOR null and void..**

Except from knobs **E** and **F** (Fig. 3) and tightening screws **V1** and **V2** (Fig. 1b) of the Ø34 tube, screwing on or loosening any other screws of the knee is strictly forbidden.

Never grease the knee roll bearings because this could damage them very soon.

**RECYCLING** This product includes items made of various materials : elastomers, plastics, aluminum, steel and brass. It also contains oil. They must be recycled in compliance with the laws in force.

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