

Wound Care & Compression

Recommendations for venous leg ulcers: Clinical practice guidelines of the Society for Vascular Surgery and the American Venous Forum.

Patients with venous leg ulcer healing. Guideline 5.1

- We recommend compression therapy over no compression therapy to increase venous leg ulcer healing rate (GRADE -1: Level of evidence - a)

Patients with venous leg ulcer healed and recurrence. Guideline 5.2

- We suggest compression therapy to decrease the risk of ulcers recurrence. (GRADE 2: Level of evidence - b)

Patients with treatment of leg ulcers. Guideline 5.3

- We suggest the use of multi component compression bandage over single – component bandages for the treatment of venous leg ulcers. (Grade 2: Level of evidence - b)

Patients with arterial insufficiencies. Guideline 5.4

- We suggest VLU underlying arterial disease; do not use compression bandages of stockings if ABI is .5 or less.

Primary Venous Disease (CEAP C3-C4) due to reflux. Guideline 8.1

- Recommend compression 20-30 mmHg knee high or thigh high. (Grade 2: Level evidence - c)

HCPCS Codes

- Inelastic Wraps (A6545)
- Compression Stockings (A5631- 30-40 mmHg), (A6532- 40-50 mmHg)
- * Juzo Compression Wrap is PDAC approved



Visit www.juzousa.com or www.juzo.ca for more information.

Juzo USA
P.O. Box 1088
Cuyahoga Falls, OH 44223
1 888-255-1300
www.juzousa.com
support@juzousa.com

Juzo Canada Ltd.
1100 Burloak Dr. Ste. 300
Burlington, ON L7L 6B2
1 888-255-1300
www.juzo.ca
support@juzo.ca

® Juzo, Fibersoft, Keep Moving & Freedom in Motion are registered trademarks of Julius Zorn, Inc.
© 2016 Julius Zorn, Inc.

Juzo Compression Solutions

For Vascular & Wound Care



Look for Juzo. Find Comfort.




Improving the Patient’s Quality of Life

Compression therapy is the standard of care for the management of wound care. Compression therapy products help manage leg ulcer healing rates and decrease the risk of ulcer recurrence. They are used for management of venous hypertension, lymphedema and venous ulcers to reduce edema and aid the return of venous blood to the heart.

Compression therapy has shown to improve the quality of life for patients when used correctly. Symptoms of pain, swelling, skin discoloration, cosmetic appearance of the leg, increased activity levels, improved sleep, and decrease in depression symptoms have been documented improvements with the use of compression therapy¹.

The use of compression may also decrease the cost associated with venous disease. A patient with a healed ulcer using compression may save an estimated \$17,000².

1. American Venous Forum Handbook.
2. Korn P, Patel ST, Heller JA, Deitch JS, Krishnasastry KV, Bush HL, et al. Why insurers should reimburse for compression stockings in patients with chronic venous stasis. J Vasc Surg 2002;35:950-7.



MANAGEMENT OF VENOUS STASIS ULCERS

After wound cleaning, a short-stretch compression bandage can be applied. This provides an optimal method of managing the circumference fluctuations that arise early in the treatment course.



Juzo® Compression Wrap

The Juzo Compression Wrap is a short-stretch garment that can be employed as an alternative to conventional multi-layer bandages. Thanks to its short stretch properties, it delivers a very high working pressure and a low resting pressure. The fastener makes it easy to apply, and easy to adjust for fluctuations in body circumference. The special notches and minimal overlying straps resist gapping and bulk. The calf wrap can be worn in combination with a foot wrap.



Advantages of Juzo® Compression Wrap

- High working pressure thanks to the short-stretch material
- Adjustable to the individual wearer when managing edema
- Faster application compared to multi-layer bandages
- Package includes a stockinette and compressive anklets
- The item is reversible (black or beige) to suit the individual patient’s preference
- Made in the U.S.A.

PLEASE NOTE:
If there are any open wounds in the region being treated, these wound areas must first be covered with a suitable wound dressing or bandage before applying the compression device.

ONGOING MANAGEMENT

As soon as edema has been brought under control, the patient can commence wearing a Juzo Ulcer Pro compression stocking system. Targeted compression therapy, combined with more exercise, can lead to good treatment outcomes for affected patients.



Juzo Ulcer Pro

The Juzo UlcerPro™ is a two layered system that can be applied over top of the wound dressings, while helping to hold them in place and provide therapeutic graduated compression.

The Ulcer Pro dual layer system includes a pair of open toe, dual stretch, knee high compression stockings and two mild compression liners for greater ease in donning over the wound dressing.

Kit available in 30-40 or 40-50 mmHg compression.



Juzo Ulcer Liner

The Juzo Ulcer Liner offers mild compression. It can be worn by itself or used as an under layer for a Juzo compression stocking.

Juzo Silver Ulcer Liner

Juzo also offers a Juzo Silver Ulcer Liner, knitted with X-Static™ fiber and offers mild compression. It can be worn by itself or used as an under layer for a Juzo compression stocking.



Juzo Soft Silver Stockings

Made with X-Static™ fiber which inhibits the growth of bacteria on the surface of the fabric providing anti-odor and antimicrobial performance for the life of the product.

- ANTIMICROBIAL- Protects fabric from the growth of bacteria and fungi.
- ANTI-ODOR- Reduces odor and keeps garment fresher.
- REGULATES TEMPERATURE- Keeps you cooler in the summer and warmer in the winter.
- PERMANENT- Does not wash out. Lasts for the life of the product.

MANAGEMENT OF RECURRENCES

Once the ulcer has healed, compression stockings can be used to maintain the therapy results over the long term. These are selected according to the degree of edema formation – if the leg has minor swelling, we recommend our Juzo Soft or Basic line. For more extreme swelling, use the Juzo Dynamic Line.



Juzo Dual Stretch

Juzo DualStretch™ compression stocking is designed to be therapeutic and provide greater ease of donning.

The foot and ankle area offers two-way stretch for the perfect fit and comfort while the one-way stretch zone helps with putting on and taking off the garment.

The Juzo Dual Stretch can be used as a base layer under the Juzo Compression Wrap System, as an alternative foot compression solution, worn as the top layer for the Juzo Ulcer Pro or worn alone as a compression garment.



Juzo Soft

Juzo’s Soft Stockings deliver the softest relief for legs and are easy to put on and off. The opaque knit leaves the garment virtually inconspicuous, without sacrificing comfort or looks. Available in open or closed toe and 15-20, 20-30 & 30-40 mmHg compression. Juzo Soft is available in standard and seasonal colors.



Juzo Dynamic

The Juzo Dynamic is an all purpose, durable seamless compression stocking, providing the best containment for swelling. All of the attributes, styles and features provide many compression options and solutions. The Juzo Dynamic is available in 20-30, 30-40 and 40-50 mmHg and open or closed toe. Also available in MAX sizing.



For value line compression stockings, use Juzo Basic.