Instructions for use

Cheetah
(Flex-Sprint III)

Life Without Limitations®

.....	3
ⒺⒺ Instructions for use.....	5
ⒺⒺ Gebrauchsanweisung.....	7
ⒺⒺ Guide de fabrication.....	9
ⒺⒺ Instrucciones para el uso.....	11
ⒺⒺ Istruzioni per l'uso.....	13
ⒺⒺ Brugsanvisning.....	15
ⒺⒺ Bruksanvisning.....	17
ⒺⒺ Gebruiksaanwijzing.....	19
ⒺⒺ Instruções de Utilização.....	21
ⒺⒺ 取扱説明書.....	23

EN - Caution: Össur products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when Össur products and components are used with other recommended Össur or authorized components. If un-usual movement or product wear is detected in a structural part of a prosthesis at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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




Cheetah™ (Flex- Sprint III)

Cheetah System

- Custom foot for track & field sports
- Preferred transtibial sprinting foot
- Clearance for hardware
Minimum: 250mm (10") for Pylon Adapter (D)
Maximum: 32mm (1 1/4") for lamination adapter (C)
- Maximum weight 147 kg (325lbs) for lamination adapter
- One adult foot size

Connector Options (A)

Pylon Connector FSX50003 (1-6), FSX50004 (7-9)

- Used for clearance greater than 250mm (10")
- Adapter is laminated direct to the back of the socket protruding from the apex of the distal socket to allow attachment of the foot.

Bolts go through adapter first then through foot module. Apply Loctite 410 to all threads and torque to following specifications:

Category	Bolt Size	Torque
M 1-6	5/16"	33 Nm (24 ft-lbs)
L 7-9	3/8"	45 Nm (33 ft-lbs)

Lamination Connector (B)

FSX50001 (1-6), FSX50002 (7-9)

- Used for clearance less than 250mm (10")
- Adapter is laminated direct to the back of the socket

Bolts go through adapter first then through foot module. Apply Loctite 410 to all threads and torque to following specifications:

Category	Bolt Size	Torque
M 1-6	5/16"	33 Nm (24 ft-lbs)
L 7-9	3/8"	45 Nm (33 ft-lbs)

The Cheetah™ may also be connected by a direct lamination.

Friction Pads

Always use between T and foot module. If alignment wedges are used, place friction pads between each layer.

Bench Alignment

Sagittal plane (C+D)

- Load line bisects socket and falls through the ground contact point or approximately 5cm (2") posterior.
- Determine the appropriate socket flexion according to amputee requirements.

Coronal plane

- Introduce appropriate socket abduction/adduction (using a jig)
- 5°-7° additional external rotation
- Achieve proper height approximately 25mm (1") longer than the sound side with shoe to compensate for vertical compression and plantar flexion

Note: The recommendations for bench alignment are used as a starting point. Definitive alignment will be achieved during dynamic assessment.

Cutting Foot Module

- Before cutting check for enough clearance.

- Measure for length and add 50mm (2") to compensate for measurement error. The final height is approximately 25mm (1") taller than the sound side.
- Clamp foot module in vice. Use a sharp fine tooth hacksaw blade.
- Sand away any sharp edges.
- Secure foot to using fitting clamps or synthetic plaster bandages.
- Cut again if necessary.
- Lowest hole position is in 30 cm (12") when pylon is vertical.

Drilling Foot Module Attachment Holes

- Back-up foot module when drilling using wood block or discarded carbon fiber section and clamp securely.
- Use drill guide template to locate holes.
- Use sharp drill bits. Suggestion: Bullet nose or solid carbide tipped drill bits.
- Drill size:
Cat 1-6 8,4mm (21/64")
Cat 7-9 11mm (7/16")
- Drill speed 1200-1800 RPM.
- Feed drill slowly with light intermittent pressure - do not force.
- Foot is now ready for assembly.

Dynamic Alignment

Physical ability and length of residual limb will result in different alignment. Best performance will be achieved with best achievable symmetry of gait. Each athlete will have her/ his own optimal alignment, which relates to the individual's running style. Optimal alignment is critical in order to reduce energy expenditure and eliminate excessive muscle and ligament strain while running.

- Check for best symmetry of gait
- Stride length can be affected by A-P position socket to foot
- Vertical deflection is in the average between 50mm (2") to 75mm (3")

Trouble shooting

Exclude reasons related to socket fit. Try to achieve symmetry in running as close as possible. The symmetry in running should be assessed at maximum speed.

Further recommendations:

For track & field events shoes usually feature spikes:

The sole of the spare shoe is then attached to the prosthetic foot in the following way:

- Cut off the sole area where the spikes are attached (1).
- Remove the insole
- Bond leather to the back side of the prepared part.
- Sand the surfaces (2).
- Bond assembly to the ground contact area (3).

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Össur Iberia S.A.