

DE	Gebrauchsanweisung	3
EN	Instructions for Use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização	28
NL	Gebruiksaanwijzing	33
SV	Bruksanvisning	38
DA	Bruksanvisning	42
NO	Bruksanvisning	47
FI	Käyttöohje	52
PL	Instrukcja użytkowania	57
HU	Használati utasítás	62
CS	Návod k použití	67
RO	Instrucțiuni de utilizare	71
HR	Upute za uporabu	76
TR	Kullanma talimatı	81
EL	Οδηγίες χρήσης	86
RU	Руководство по применению	91
JA	取扱説明書	96
ZH	使用说明书	100

4R84=D
4R84=D-62

4R72=D
4R72=D-62
4R75=D-70

①

■	●	▲
DE Einzelteile	Einzelteile-Pack	Mindestmenge
EN Single components	Single-Component Pack	Minimum order quantity
FR Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
IT Singoli componenti	Imballo dei singoli componenti	Quantità minima
ES Componentes	Kit componentes	Cantidad mínima
PT Componentes	Conjunto de peças de reposo	Quantidade mínima
NL Onderdelen	Onderdelenpakket	Minimikvantitet
SV Separata delar	Service-Set	Minimum aantal
DA Komponenter	Komponentpakke	Mindste mængde
NO Enkle komponenter	Enkeltkomponentpakke	Minste ordremenge
FI Yksittäisosat	Yksittäisosapakkaus	Vähimmäismäärä
PL Podzespoły	Zestaw naprawczy	Ilość minimalna
HU Alkatrészek	Alkatrész-csomag	Minimum mennyiség
CS Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
RO Componente individuale	Set componente individuale	Cantitate minimă
HR Pojedinačni dijelovi	Pakiranje pojedinačnih dijelova	Minimalna količina
TR Müferit parçalar	Müferit parça ambalajı	Minimum miktar
EL Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
RU Отдельные детали	Комплект деталей	Минимальное количество
JA 1 個から発注いただける部品	セットで発注いただく部品	複数入パックで発注いただく部品
ZH 单个部件	维修组件	最低起订量

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CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Date of the last update: 2013-08-14

- Please read this document carefully.
- Follow the safety instructions.

The scope of delivery is shown on the front page.

1 Single Components (fig. 1)

1.1 Single Component ■

—

1.2 Single-Component Pack ●

—

1.3 Minimum Order Quantity ▲

(1) 506G3=M8x12-V Set Screw

2 Description

2.1 Intended Use

The structural parts shown on the front page are to be **solely** used for the prosthetic fitting of the lower limb.

2.2 Field of Application

Field of application according to the **Ottobock MOBIS Mobility System**:

Approved for a **patient weight of up to 150 kg (330 lbs)**.

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water, and acids.

Using a medical device in such environmental conditions will render all claims against Otto Bock HealthCare GmbH null and void.

Please also inform your patients!

2.3 Function

The 4R84=D and 4R84=D-62 Connection Adapters provide connection between pyramid receiver and Ø 30 mm tube clamp.

The 4R72=D and 4R72=D-62 Connection Adapters provide connection between pyramid adapter and Ø 30 mm tube clamp.

The 4R75=D-70 Connection Adapter provides connection between pyramid adapter and Ø 34 mm tube clamp.

3 Technical Data

	4R84*		4R72*		4R75
	*=D	*=D-62	*=D	*=D-62	=D-70
Weight	65 g	145 g	70 g	150 g	170 g
System height min.	–	20 mm	–	67 mm	76 mm
System height max.	–	48 mm	–	96 mm	106 mm
System height	19 mm	–	66 mm	–	–
Diameter	30 mm	30 mm	30 mm	30 mm	34 mm
Material	Titanium	Stainless steel	Titanium	Stainless steel	Stainless steel
Patient weight max.	150 kg	150 kg	150 kg	150 kg	150 kg
Mobility grade	–	–	–	–	–

*Stainless steel

4 Handling

4.1 Assembly

The Connection Adapters 4R84=D-62, 4R72=D-62 and 4R75=D-70 must be shortened according to the measurements of the amputee.

INFORMATION

The section surface must be perpendicular to the adapter longitudinal axis. Carefully debur the tube cutting edge inside and outside.

Loosen the cap screw of the tube clamp adapter by two turns and insert the tube 5 mm deep into the tube clamp adapter while slightly rotating the tube. Next, place the tube clamp adapter with the tube inserted vertically on the worktop and insert the tube up to the bottom of the clamp.

4.2 Adjustments and Final Assembly

CAUTION

Risk of falling due to changed prosthesis geometry. Dirt on the contact surfaces can cause the adapters to turn while using the prosthesis.

Prior to final assembly, the contact surfaces of the adapters must be cleaned with 634A3 Acetone!

Static corrections are possible at any time during alignment and trial walking as well as after completion of the prosthesis using the four 506G3 Adjustment Screws. To exchange a module or during disassembly the previous position may be retained by loosening only the two most deeply screwed in, adjacent set screws.

Replace set screws that appear too long (to prevent foam cover damage) or too short (stability) with the appropriate screw length.

Adjustment screws are available in the following sizes.

Article No.	Length
506G3=M8x12-V	12 mm
506G3=M8x14	14 mm
506G3=M8x16	16 mm

Use 710D4 Torque Wrench to tighten the adjustment screws.

CAUTION

Risk of falling as a result of loosened set screws. To prevent this risk, proceed as follows:

To tighten the set screws: Use 710D4 Torque Wrench.

Trial fitting: Tighten the 506G3 Set Screws clockwise to **10 Nm**.

Completion: Secure the 506G3 Set Screws with 636K13 Loctite® and pre-tighten clockwise to **10 Nm**. Then complete assembly by tightening to **15 Nm**.

4.3 Maintenance Instructions

INFORMATION

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Disposal

5.1 Disposal Information

INFORMATION

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

6 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

7 Warranty

Ottobock warrants this product/these products from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Ottobock distribution company for more information.

8 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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