






28R140N Dorso Arexa

DE	Gebrauchsanweisung	5
EN	Instructions for use	10
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização	34
NL	Gebruiksaanwijzing	40
SV	Bruksanvisning	46
DA	Brugsanvisning	51
NO	Bruksanvisning	57
FI	Käyttöohje	62
PL	Instrukcja użytkowania	68
CS	Návod k použití	74
BG	Инструкция за употреба	79
TR	Kullanma talimatı	85
RU	Руководство по применению	91
JA	取扱説明書	97
ZH	使用说明书	103

Material	Rostfreier Stahl / stainless steel, PE Schaumstoff / PE foam, Frotteebezüge / Terry Cloth Covers Kunststoff / plastic: PA, TPU, PP
-----------------	--

INFORMATION

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INFORMATION

Last update: 2015-09-23

► Please read this document carefully.

► Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 28R140N Dorso Arexa hyperextension brace.

2 Product Description

2.1 Components (see Fig. 1)

Item	Component	Item	Component
1	Sternal support, padded	8	Textile cover
2	Base plate, padded	9	Waistband, two-piece
3	Ratchet closure	10	Screws for waistband
4	Belt guiding element	11	Torx wrench
5	Extension element	12	Clip rivet for pads
6	Symphysis support, padded	13	Waistband pad
7	Back support, padded	–	Sternal support, short (to order if needed, article no.: 29R331=S)

2.2 Design

The ventral section consists of a base plate (2) with adaptable extension elements (5), two ratchet closures (3), a sternal support (1) and a symphysis support (6). The supports (1, 6) are continuously height-adjustable. The supports with the tubes are attached to the base plate by special clamping elements. The sternal support is also available in a short version (option).

The flexible, two-piece waistband (9) is attached to the back support (7). The waistband slides into the ratchet closures (3) and engages audibly. The base plate and waistband as well as the sternal, symphysis and back supports are padded.

3 Application

3.1 Indications for use

The modular brace is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

The brace must be used according to the indications and following the instructions of a physician.

3.2 Indications

Stable traumatic vertebral body fractures in the area of the lower thoracic spine and upper lumbar spine (also osteoporotic) without neurological damage (Th10 to L2).

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support/brace will be applied; sensory and circulatory disorders in the body region being treated.

3.4 Effects

The brace straightens the spinal column in the sagittal plane and stabilizes it. Based on the three-point principle with the back support (7), sternal, and symphysis supports (1, 6), the brace causes a lordosis, thereby relieving the anterior vertebral bodies.

4 Safety

4.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury

NOTICE Warning regarding possible technical damage

4.2 General safety instructions

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety instructions

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema, or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.

► **Inform the patient.**

⚠ CAUTION

Contact with heat, embers, or fire

Injuries (e.g., burns) caused by melting of the material

- Keep the product away from open flame, embers, and other sources of heat.

► **Inform the patient.**

NOTICE

Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e. g. sand)

Damage and premature product wear

- Following any contact with the substances identified above, promptly clean the product in accordance with the chapter "Cleaning".

► **Inform the patient.**

NOTICE

Improper molding or application

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- The product may only be molded by qualified personnel.
- Do not make any improper changes to the product.
- Always apply the product according to the information in the instructions.

► **Inform the patient.**

5 Handling

INFORMATION

- The daily duration of use and period of application are generally determined by the physician.
- The initial fitting and application of the product must be carried out by qualified personnel.
- Instruct the patient in the handling and care of the product.
- Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size selection

The brace is available in a universal size. It can be used for fitting of waist-band circumference up to 150 cm/59 inches.

5.2 Fitting and application

- > The patient is lying or standing.
- 1) Loosen the screws in the base plate (see fig. 2).
- 2) Taking into account the maximum extension length marked on the tubes, adjust the height of the sternal and symphysis supports to the patient (see fig. 3).
- 3) Tighten the screws in the base plate.
- 4) If adaptation to the anatomical body contours of the patient is required, shape the tubes between the base plate and the supports manually over a suitable edge (see fig. 4). **Functional risks are associated with, for example, repeated cold forming and the resulting reduction in breaking elongation (strain hardening).**
- 5) If required, remove the same number of extension elements not needed to adjust the circumference from both sides and then reattach the belt guiding elements (see fig. 5). Shorten the base plate pad accordingly (see fig. 6).
- 6) Secure the two-piece waistband to the back support in the optimum position for the patient using countersunk head screws (see fig. 7).
- 7) Position the back support dorsally and the base plate ventrally on the patient.
- 8) Insert the waistband symmetrically into the ratchet closures (see fig. 8). Fine-tuning is done by using the ratchet closures (see fig. 9).
- 9) Check the proper fit of the brace and correct if required.
- 10) Adjust the length of the waistband pad (see fig. 10) and attach it with clip rivets (see fig. 11).
- 11) Shorten the ends of the waistband protruding from the ratchet closures, leaving some in reserve, and round them so they are free of burrs (see fig. 12).
- 12) To open the brace, push one of the buttons on the ratchet closure (see fig. 13).

Verification

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the brace is applied properly and fits correctly.
- ▶ **Inform the patient.**

Before handing over the product to the patient, the following points must be verified after applying the brace to the patient:

- Pressure, placement, and fit of the three supports (1, 6, 7).
- Tension and proper positioning of the waistband (9) and base plate (2).
- Fit and hyperextension effect in the sitting position.

Instructing the patient

⚠ CAUTION

Loosening the ratchet closures

Loss of the support function, injuries due to wear or unintentional release of the ratchet closures

- ▶ At each application, ensure that the ratchet closures audibly engage in the waistband. The waistband must not slip out under load.
- ▶ Practice the application and removal of the product with the patient.
- ▶ In case of damage or noticeable changes in the product, the patient should consult a specialist workshop.

INFORMATION

Continuous, direct skin contact with the product can lead to skin irritation. Only wear the product over clothing (e.g. a T-shirt or undershirt) in daily use. For limited times (e.g. water therapy), the product can also be worn directly on intact skin.

5.3 Cleaning

The brace is not damaged by water, which means that personal hygiene (e.g. showering) and therapy in water that contains chlorine or salt are, in principle, possible when wearing the brace. Subsequently, the brace must be cleaned according to the following instructions.

- 1) Remove the textile covers from the supports and wash at 30 °C in the washing machine.
- 2) Clean the brace with a damp cloth and neutral soap.
- 3) Rinse the brace with water and dry it with a cloth.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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Otto Bock HealthCare GmbH