


4R220 Dynamic Vacuum System


DE	Gebrauchsanweisung	8
EN	Instructions for use	20
FR	Instructions d'utilisation	31
IT	Istruzioni per l'uso	44
ES	Instrucciones de uso	56
PT	Manual de utilização	68
NL	Gebruiksaanwijzing	80
SV	Bruksanvisning	92
DA	Brugsanvisning	103
NO	Bruksanvisning	114
FI	Käyttöohje	126
PL	Instrukcja użytkowania	137
HU	Használati utasítás	149
CS	Návod k použití	161
RO	Instrucțiuni de utilizare	172
HR	Upute za uporabu	184
SL	Navodila za uporabo	195
SK	Návod na používanie	206
BG	Инструкция за употреба	218
TR	Kullanma talimatı	230
EL	Οδηγίες χρήσης	241
RU	Руководство по применению	255
ZH	使用说明书	268
KO	사용 설명서	277



INFORMATION

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INFORMATION

Last update: 2016-07-01

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 4R220* Dynamic Vacuum System is a system for providing an active vacuum. It consists of a component with an integrated piston pump, a liner and a sealing sleeve.

The Dynamic Vacuum System uses the pistoning movements between the residual limb and prosthetic socket to generate the vacuum. The piston of the pump is equipped with magnets. A metallic counter-piece is installed on the liner. Air is suctioned out of the space between the liner and prosthetic socket into the cylinder during the swing phase, and this air is pressed out during the stance phase. The result is a permanent vacuum, with a level that adapts to the user's activity level.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

Combination possibilities	
Designation	Reference number
Liner	6Y94
Socket adapter	4R22, 4R23, 4R37, 4R51, 4R54, 4R55, 4R73=A, 4R73=D, 4R74, 4R77, 4R95

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Product approved for max. body weight

- The maximum approved body weight is specified in the technical data (see Page 31).

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Excessive strain on the product
Risk of injury due to breakage of load-bearing components
► Use the product according to the specified field of application (see Page 21).

⚠ CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product for a single patient.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION**Product contains a magnet**

Influence on or damage to devices or items due to a strong magnetic field

- Keep the product away from devices and items that are sensitive to magnetic fields.

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

4R220 Dynamic Vacuum System		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Cylinder body	–
1	Piston	–
1	Stop ring	4X339
4	Snap bushing	5X163
4	Countersunk head screw	501S128=M6x22
4		501S128=M6x25
1	Plaster screw	503S3
1	Valve bushing	–
1	Duckbill valve	–
1	Piston dummy	–
1	Silicone dummy	–
1	Valve dummy	–
1	Spacer plate	4X314
1	Mounting wrench	4X338
1	Polylub GLY 801 special grease	633F30=2

Spare parts (not included in scope of delivery)		
Quantity	Designation	Reference number
Single-component pack, consisting of:		4X320
1	Piston	–
1	Stop ring	–
Single-component pack, consisting of:		4X322
1	Valve bushing	–
2	Duckbill valve	–
Dummy set, consisting of:		4X326

Spare parts (not included in scope of delivery)		
Quantity	Designation	Reference number
1	Piston dummy	–
1	Silicone dummy	–
1	Valve dummy	–
1	Plaster screw	–

5 Preparation for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

INFORMATION

Not all of the materials listed may be available in your country.

- If the materials are not available, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Preparing for Socket Fabrication

INFORMATION

The Dynamic Vacuum System can be used with a total surface weight-bearing socket or a specific weight-bearing socket. Observe the following instructions:

- Maintain the residual limb length with liner unchanged while taking the plaster cast.
- Do not shorten the plaster model during modelling.

- > **Required materials:** Piston dummy, silicone dummy, short or long valve dummy, 99B25 nylon stockinette, 99B81=70X19X5 PVA bag, string, 633W8 wax, 636K8* Plastaband
- 1) Distally flatten the plaster model so that the piston dummy can be positioned in the longitudinal residual limb axis (see fig. 2). Position the piston dummy on the plaster model and loosely secure it with a plaster screw.
 - 2) Pull a nylon stockinette over the plaster model.
 - 3) Tie off the nylon stockinette between the piston dummy and plaster model and trim off the excess nylon stockinette (see fig. 3).
 - 4) Turn the plaster screw all the way in (see fig. 4).

5) **INFORMATION: No PVA bag is required for fabricating a check socket.**

Soak the shorter PVA bag and pull it over the plaster model.

6) Use string to tie off the PVA bag in the undercut of the piston dummy (see fig. 5).

7) Cut off the excess PVA bag with a scalpel.

8) Position the silicone dummy on the piston dummy. Check that there is no excess PVA bag protruding (see fig. 6).

9) **INFORMATION: Do not close the opening in the cylinder body from the inside. Screw in the valve dummy after positioning the cylinder body.**

Position the cylinder body on the piston dummy. Perform the mediolateral alignment of the exhaust opening depending on the prosthetic alignment.

10) Insulate the heads of the cap screws with wax.

11) Fill the heads of the cap screws with Plastaband (see fig. 7).

12) Screw the valve dummy into the exhaust opening until the O-ring closes the opening (see fig. 8).

13) Close the slot in the valve dummy with Plastaband.

5.2 Fabricating the Prosthetic Socket

5.2.1 Optional: fabricating a check socket

> **Recommended materials:** ThermoLyn (e.g. 616T52* ThermoLyn rigid, 616T83* ThermoLyn clear), 616S5* BetaSil

1) Prepare for fabrication of the prosthetic socket (see Page 24).

2) Attach 2 strips of Plastaband lengthwise over the cylinder on both the anterior and the posterior side. This makes it easier to remove the cylinder after the trial fitting (see fig. 9).

3) Attach a strip of Plastaband circularly over the proximal edge of the cylinder to seal the prosthetic socket (see fig. 9).

4) Complete the vacuum forming process.

5) Finish the prosthetic socket (see Page 28).

6) **If the system leaks:** seal the transition from the cylinder body to the inside of the prosthetic socket with BetaSil.

5.2.2 Laminating a Transtibial Socket

INFORMATION

The reinforcement described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the reinforcement.

- > **Required materials:** 623T3=8 or 623T3=10 perlon stockinette, 81A1=8 or 81A1=10 nylon stockinette, 616B1=25x* carbon fibre cloth strap, 616G2 carbon UD stockinette, 616G15 woven carbon fibre stockinette, 99B81=100X19X5 PVA bag, 617H119 Orthocryl lamination resin 80:20 PRO, string, 627B40 polyethylene adhesive tape
- 1) Prepare for fabrication of the prosthetic socket (see Page 24).
- 2) Cut 3 pieces of perlon stockinette (each twice the length of the plaster model).
- 3) Pull a perlon stockinette over the plaster model to the edge. Tie off the second half of the perlon stockinette and fold it over the plaster model (see fig. 10).
- 4) Pull the 2 other pieces of perlon stockinette over the plaster model the same way.
- 5) Expose the valve dummy using a scalpel (see fig. 11).
- 6) Wrap string around the valve dummy and tie the perlon stockinette off tightly in the distal groove of the cylinder body (see fig. 12).
- 7) **INFORMATION: Tie this off tightly so the system is airtight after the lamination process.**
Tie off the 6 layers of erlon stockinette in the undercut of the cylinder body so they are up against the cylinder body (see fig. 13).
- 8) Position a layer of carbon fibre tape from the medial condyle over the cylinder body to the lateral condyle (see fig. 14).
- 9) Position a layer of carbon fibre tape from the MPT (mid-patella-tendon) point over the cylinder body to the opening for the back of the knee (see fig. 15).
- 10) Position a layer of carbon fibre tape circularly at the height of the MPT (mid-patella-tendon) point (see fig. 16).
- 11) Pull a nylon stockinette over the plaster model.
- 12) Expose the valve dummy.
- 13) Wrap string around the valve dummy and tie the nylon stockinette off tightly in the distal groove of the cylinder body (see fig. 17).
- 14) Tie the nylon stockinette off tightly in the undercut of the cylinder body with perlon string (see fig. 18).
- 15) Cut off a piece of woven carbon fibre stockinette (1.3 times the length of the plaster model).
- 16) Pull the woven carbon fibre stockinette over the plaster model to the edge (see fig. 19).
- 17) Tie off the excess woven carbon fibre stockinette at the distal end and pull it over the plaster model (see fig. 20).

- 18) Tie off the woven carbon fibre stockinette tightly in the distal groove of the cylinder body.
- 19) Tie off the woven carbon fibre stockinette tightly in the undercut of the cylinder body using perlon string.
- 20) Expose the valve dummy.
- 21) Pull a nylon stockinette over the plaster model.

Prosthetic socket with carbon finish

- 1) Cut a piece of carbon UD stockinette (twice the length of the plaster model).
- 2) Pull the carbon UD stockinette over the plaster model up to the edge.
- 3) Tie off the excess carbon UD stockinette at the distal end and fold it over the plaster model.
- 4) Pull a nylon stockinette over the plaster model.
- 5) Soak the longer PVA bag and pull it over the plaster model.
- 6) Perform the casting process with Orthocryl. Approx. **30 %** more Orthocryl is required due to the 6 layers of perlon stockinette.
- 7) **Once the lamination resin is sufficiently distributed:** displace excess lamination resin in the area of the cylinder body to the distal end by wrapping with polyethylene adhesive tape. Wrap tightly so the system becomes airtight.
- 8) Allow lamination resin to cure.
- 9) Finish the prosthetic socket (see Page 28).

Prosthetic socket without carbon finish

- 1) Cut a piece of carbon UD stockinette (1.5 times the length of the plaster model).
- 2) Pull the carbon UD stockinette over the distal half of the plaster model (see fig. 21).
- 3) Tie off the excess carbon UD stockinette at the distal end and fold it over the plaster model (see fig. 22).
- 4) Tie off the carbon UD stockinette in the undercut of the cylinder body (see fig. 23).
- 5) Cut a piece of perlon stockinette (twice the length of the plaster model).
- 6) Pull the perlon stockinette over the plaster model to the edge. Tie off the second half of the perlon stockinette and fold it over the plaster model (see fig. 24).
- 7) Soak the longer PVA bag and pull it over the plaster model (see fig. 25).
- 8) Perform the casting process with Orthocryl. Approx. **30 %** more Orthocryl is required due to the 6 layers of perlon stockinette.

- 9) **Once the lamination resin is sufficiently distributed:** displace excess lamination resin in the area of the cylinder body to the distal end by wrapping with polyethylene adhesive tape. Wrap tightly so the system becomes airtight.
- 10) Allow lamination resin to cure.
- 11) Finish the prosthetic socket (see Page 28).

5.3 Finishing the Prosthetic Socket

- > **Required materials:** mounting wrench, snap bushings, 617H46 silicone bonding agent, 633F30=2 Polyub GLY 801
- 1) Mark and trim the contour of the prosthetic socket.
 - 2) Sand free the valve dummy and remove the Plastaband.
 - 3) Unscrew the valve dummy and retain it. After fabricating the check socket, the valve dummy is still required for the fabrication of the final prosthetic socket.
 - 4) Remove the prosthetic socket from the plaster model.
 - 5) Remove the piston dummy and silicone dummy from the plaster model and insert into the cylinder so it is closed (see fig. 26).
 - 6) Close the opening in the piston dummy with Plastaband.
 - 7) Screw in the valve dummy (see fig. 27).
 - 8) Sand the contour of the prosthetic socket.
 - 9) Sand the distal end of the socket down to the heads of the cap screws (see fig. 28).
 - 10) Remove the cap screws.
 - 11) Sand the distal end of the socket flat and check that its surface is even.
 - 12) Remove any dirt from the 4 openings in the distal end of the socket. This ensures that the snap bushings will lie flat against the cylinder body.
 - 13) **INFORMATION: The snap bushings must protrude approx. 1 mm from the laminate. This ensures forces are transferred to the Dynamic Vacuum System.**
Apply silicone bonding agent to the snap bushings and insert them into the openings with the groove end first (see fig. 29).
 - 14) Unscrew the valve dummy.
 - 15) Remove the piston dummy and silicone dummy from the prosthetic socket.
 - 16) Wipe out the cylinder space with a lint-free cloth.
 - 17) Grease the cylinder space and O-rings of the valve bushing with the supplied Polyub GLY 801.

18) **NOTICE! Do not allow the duckbill valve to come into contact with grease.**

Insert the duckbill valve into the valve bushing with the pointed end first (see fig. 30).

19) Screw the valve bushing into the prosthetic socket (torque value: **3 Nm**).

20) Insert the piston into the cylinder to the stop (see fig. 31).

21) Insert the stop ring into the thread of the cylinder and tighten it with the mounting wrench (see fig. 32).

5.4 Mounting the socket adapter

The connection to the distal prosthetic components is established by a socket adapter. A spacer plate is positioned on the snap bushings in the distal end of the socket when mounting the socket adapter. To ensure forces are transferred to the Dynamic Vacuum System, the spacer plate must not lie against the laminate.

> **Required tools and materials:**

710D4 torque wrench, 636K13 Loctite®, spacer plate, 501S128=M6x22 or 501S128=M6x25 countersunk head screws

1) Set the spacer plate onto the prosthetic socket (see fig. 33).

2) **If the pyramid receiver or the pyramid adapter is rotatable:** Place the pressure plate on the adapter.

3) Apply the adapter to the prosthetic socket (see fig. 34).

4) **Optional:** Align the pyramid receiver or the pyramid adapter.

5) **NOTICE! Only use the screws specified here.**

Select suitable countersunk head screws.

Fixed pyramid adapter/pyramid receiver: 501S128=M6x22

Swivelling pyramid adapter/pyramid receiver: 501S128=M6x25

6) Secure the screws with Loctite®.

7) Screw in the 2 posterior countersunk head screws and tighten them (torque value: **12 Nm**).

8) Screw in the 2 anterior countersunk head screws and tighten them (torque value: **12 Nm**).

6 Use

6.1 Donning the prosthesis

1) Apply the liner to the residual limb so that the distal end is aligned with the longitudinal axis of the residual limb. Unroll the liner onto the residual limb so that it has no wrinkles, without shifting the soft tissue.

- 2) Slide into the prosthetic socket with the liner. The magnets in the piston of the Dynamic Vacuum System connect to the metallic counter-piece in the liner.
- 3) Unroll the sealing sleeve over the prosthetic socket onto the thigh.

6.2 Taking off the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Pull the residual limb and liner out of the prosthetic socket.

7 Cleaning

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Dirt in the prosthetic socket or on the liner can impair the function of the Dynamic Vacuum System.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ To lubricate the cylinder or piston seal, use **Polylub GLY 801**.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R220
Weight [g]	210
System height [mm]	37
Material	Aluminium
Max. body weight [kg]	150

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Otto Bock HealthCare GmbH