


50S10 Malleo Immobil Walker high	
50S11 Malleo Immobil Walker low	
50S12 Malleo Immobil Air Walker high	
50S14 Malleo Immobil Air Walker low	
50S15 Malleo Immobil ROM Walker	
50S17 Malleo Immobil Achilles Walker	

DE	Gebrauchsanweisung	6
EN	Instructions for Use	12
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	33
PT	Manual de utilização	40
NL	Gebruiksaanwijzing	46
SV	Bruksanvisning	53
DA	Brugsanvisning	59
NO	Bruksanvisning	66
FI	Käyttöohje	72
PL	Instrukcja użytkowania	79
HU	Használati utasítás	85
CS	Návod k použití	92
RO	Instrucțiuni de utilizare	98
HR	Upute za uporabu	106
SL	Navodila za uporabo	112
BG	Инструкция за употреба	118
TR	Kullanma talimatı	125
EL	Οδηγίες χρήσης	132
RU	Руководство по применению	139
ZH	使用说明书	147

Malleo Immobil	Art-Nr.	Max. Körpergewicht / User weight limit	Schuhgröße / Shoe Size		
			EU	US	
				Männer / Men's	Frauen / Women's
Walker high	50S10-1=S	100 kg 220 lbs	31-35	1-3½	2-4½
	50S10-1=M		36-42	4-9	5-10
	50S10-1=L		43-45	9½-11	10½-12
	50S10-1=XL		46-50	11½-15	12½-13
Walker low	50S11-1=S	125 kg 275 lbs	31-35	1-3½	2-4½
	50S11-1=M		36-42	4-9	5-10
	50S11-1=L		43-45	9½-11	10½-12
	50S11-1=XL		46-50	11½-15	12½-13
Air Walker high	50S12-1=S	100 kg 220 lbs	31-35	1-3½	2-4½
	50S12-1=M		36-42	4-9	5-10
	50S12-1=L		43-45	9½-11	10½-12
	50S12-1=XL		46-50	11½-15	12½-13
Air Walker low	50S14-1=S	125 kg 275 lbs	31-35	1-3½	2-4½
	50S14-1=M		36-42	4-9	5-10
	50S14-1=L		43-45	9½-11	10½-12
	50S14-1=XL		46-50	11½-15	12½-13
ROM Walker	50S15-1=S	100 kg 220 lbs	31-35	1-3½	2-4½
	50S15-1=M		36-42	4-9	5-10
	50S15-1=L		43-45	9½-11	10½-12
	50S15-1=XL		46-50	11½-15	12½-13
Achilles Walker	50S17-1=S	100 kg 220 lbs	31-35	1-3½	2-4½
	50S17-1=M		36-42	4-9	5-10
	50S17-1=L		43-45	9½-11	10½-12
	50S17-1=XL		46-50	11½-15	12½-13

Material	Liner: PA, PU Schaumstoff/PU Foam Walker: Glasfaserverstärkter Kunststoff/glass-fibre reinforced plastic
-----------------	---

INFORMATION

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INFORMATION

Date of the last update: 2014-04-23

- Please read this document carefully.
- Follow the safety instructions.

The instructions for use provide you with important information for adjusting and applying the 50S10 Malleo Immobil Walker high, 50S11 Malleo Immobil Walker low, 50S12 Malleo Immobil Air Walker high, 50S14 Malleo Immobil Air Walker low, 50S15 Malleo Immobil ROM Walker, and 50S17 Malleo Immobil Achilles Walker immobilisation orthoses.

2 Application

2.1 Application

The orthosis is intended **exclusively** for orthotic fittings of the foot/ankle joint area and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

2.2 Indications

50S10 Malleo Immobil Walker high

- Postoperative, conservative treatment of stable forefoot, metatarsus and/or ankle fractures
- Postoperative, conservative treatment of distal fibula fractures
- Conservative treatment of forefoot, metatarsal and tarsal injuries
- Ligament, soft tissue and tendon injuries

50S11 Malleo Immobil Walker low

- Postoperative, conservative treatment of stable forefoot, metatarsus and/or ankle fractures
- Soft tissue injuries of the foot
- Exostosis treatment for hallux valgus

50S12 Malleo Immobil Air Walker high

- Postoperative, conservative treatment of stable forefoot, metatarsus and/or ankle fractures
- Postoperative, conservative treatment of distal fibula fractures
- Soft tissue injuries of the foot
- Severe ankle sprains

50S14 Malleo Immobil Air Walker low

- Postoperative, conservative treatment of stable forefoot, metatarsus and/or ankle fractures
- Soft tissue injuries of the foot
- Exostosis treatment for hallux valgus

50S15 Malleo Immobil ROM Walker

- Postoperative, conservative treatment of stable forefoot, metatarsus and/or ankle fractures
- Postoperative, conservative treatment of distal fibula fractures
- Postoperative, conservative treatment of Achilles tendon ruptures
- Conservative treatment of forefoot, metatarsal and tarsal injuries
- Ligament, soft tissue and tendon injuries

50S17 Malleo Immobil Achilles Walker

- Postoperative, conservative treatment of Achilles tendon ruptures

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases / injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the support / orthosis will be applied; sensory and circulatory disorders in the body region being treated and in the area of the legs.

2.4 Effects

The orthosis serves to immobilise or limit the movement (only in the case of the ROM Walker model) of the foot/ankle joint area.

3 Safety

3.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury

NOTICE Warning regarding possible technical damage

3.2 General Safety Instructions

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- Observe the safety information in these instructions for use.

- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product is for single patient use only.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C or above +60 °C (e.g. sauna, excessive sunlight, drying on a radiator).
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the physician treating you.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table)

4.2 Adaptation

- 1) Open all the hook-and-loop closures on the orthosis and remove the liner (see Fig. 1 and see Fig. 2).

INFORMATION: For the Air Walker high, Air Walker low and Achilles Walker models, the top hook-and-loop closure is removed from the Walker together with the liner.

- 2) **Only Achilles Walker and, optionally, Walker high and Air Walker high models:** Position the heel wedge (29S35) in the liner.
- 3) Position the foot in the liner (see Fig. 3).
- 4) Close all the hook-and-loop closures on the liner (see Fig. 4 and see Fig. 5).
- 5) Position the foot with the liner in the Walker (see Fig. 6). Check the position of the heel whilst doing this.

INFORMATION: The foot angle should be 90° (except models with the 29S35 heel wedge).

- 6) Press on the lateral guides medially and laterally (see Fig. 7).
→ The liner is secured in the Walker.
- 7) Close all hook-and-loop closures on the Walker one after the other, starting at the bottom (see Fig. 8 and see Fig. 9).

Setting the Joint Angles (Only ROM Walker Model)

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.
- ▶ **Inform the patient.**

The joints on the ROM Walker can be adjusted in 7.5° increments from 45° dorsal flexion to 45° plantar flexion. Take the following steps to enable plantar or dorsal flexion movement:

- 1) Move the Plexiglas disc so that the window is over the pin (see Fig. 10).
- 2) Remove the pins from the joint (see Fig. 11).
- 3) Set the same angle medially and laterally.
- 4) Reinsert the pins (see Fig. 12).

Alternatively, the joints can be completely fixated:

- ▶ Tighten the Allen head screw (see Fig. 13, item 1) with the enclosed Allen key.

4.3 Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- 1) Open all hook-and-loop closures on the orthosis.
- 2) Position the foot in the orthosis. Check the position of the heel whilst doing this (see Fig. 14).
- 3) Close all hook-and-loop closures on the liner (see Fig. 15).
- 4) Close all hook-and-loop closures on the Walker one after the other, starting at the bottom (see Fig. 16 and see Fig. 17).

Increasing compression (only Air Walker high, Air Walker low and Achilles Walker models)

⚠ CAUTION

Pumping up the air cushion too much

Risk of local pressure and constriction of blood vessels and nerves due to excessively tight application

- ▶ Only fill the air cushion with enough air to impede relative movements between the body part and the orthosis.
- ▶ **Inform the patient.**

- 1) To inflate, turn the rotary valve to the right [Inflate] (see Fig. 18, item 1).
- 2) Press the pump as often as required to achieve the desired compression (see Fig. 18, item 2).
- 3) **Optional:** to deflate, turn the rotary valve to the left [Deflate] (see Fig. 18, item 1).

4.4 Cleaning

Clean the orthosis regularly:

Textile component:

- 1) Remove the textile component from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile component in warm water at 30 °C with a standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. sunshine, stove or radiator).

Composite component/heel wedge:

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat (e.g. sunshine, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The

manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH