ENGLISH

MD Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional. The device must be fitted and adjusted by a healthcare professional. For more detailed device information, see **Table 1**.

INDICATIONS FOR USE AND TARGET PATIENT POPULATION

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The weight limit for the device is 136kg. The device is for low to moderate impact use, e.g., walking.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional. **Warning:** Only use Össur replacement parts for the device for ensured durability The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Rotation Adjustments:

- Loosen two neighbouring set screws on female adapter.
- Set adapter to desired rotation.
- Secure set screws with medium strength threadlocker and torque.

Warning: For secure fastening use medium strength threadlocker and apply rated torque twice, see **Table 1**.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

ENVIRONMENTAL CONDITIONS

The device is Weatherproof. A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed. No contact with salt water or chlorinated water is allowed. Dry with a cloth after contact with fresh water or humidity. Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

SPECIFICATIONS

Table 1 - Device Specifications

Part #	Description	Material	Tightening Torque	
			Mount Screws	Set Screws
M625400	4-hole adapter with pyramid rotatable	Stainless Steel	10 Nm	None
M626400	4-hole adapter with pyramid receiver rotatable			15 Nm

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.


LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

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