















Instructions for Use

ICELOCK[®]
544 SERIES

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	19
NO	Bruksanvisning	22
DA	Brugsanvisning	25
SV	Bruksanvisning	28
EL	Οδηγίες χρήσης	31
FI	Käyttöohjeet	35
NL	Gebruiksaanwijzing	38
PT	Instruções de Utilização	41
JA	取扱説明書	45
ZH	中文说明书	48

INTRODUCTION

Thank you for choosing a high quality component from Össur. We can assure you that you have made the right choice. Please carefully review and follow these instructions before fabricating the prosthesis or accessory products. Remember to inform the patient concerning proper handling and care of the prosthesis and its subcomponents, as set forth in these instructions. This is important to achieve patient comfort and maximum product life.

The users of Icelock 544 are clinical specialists and end user. This device should only be used with a prescription from a physician that includes the necessary treatment information.

INDICATION FOR USE

The Icelock 544 Series components are specifically designed to support negative pressure suspension of prostheses to liners for Lower extremity amputees. Weight Limit: 136 kg.

ISO 10328 - P6 - 136 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

CONTRAINDICATIONS FOR USE

Össur does not recommend prescribing Icelock 544 for amputees that are not capable to use the device correctly according to instructions, for example due to lack of hand dexterity or cognitive disability.

PRODUCT SELECTION AND FITTING

TYPE OF VALVE PLATE

Icelock 544 offers a number of different options regarding connectors, check valve configurations and release buttons.

The following chart shows the different options available:

Item	Check valve to socket (Sleeve only)	Extended Air Chamber	Release button	Expulsion port to air (bypass)	Integrated filter for pump	Barb connector to external Pump
L-544511			Yes		Yes	Yes
L-544513		Yes	Yes		Yes	Yes
L-544521			Yes	Yes		
L-544512	Yes	Yes			Yes	Yes
L-544541			Yes		Yes	
L-544520	Yes			Yes	Yes	
L-544540					Yes	

L-544511 is specifically designed for Unity® vacuum pumps, while other configurations are more suitable in other scenarios, L-544521 suits well with Iceross Seal-In liners without elevated vacuum since it has a bypass port and a Release button to let air in. Consult with Össur Customer Service if in any doubt with product selection.

CAUTION !

Never use an autoexpulsion(check) valve without a functional release valve when combined with Iceross Seal-In® or other solutions where liner is sealed against the inside wall of a socket. Some means of allowing air back in the socket must be in place. Never operate this product if a change in its operation is observed. In such cases the user should consult the clinical specialist.

Icelock 544 is a single patient device.

DAILY CARE

Proper product care is important for maintaining product lifetime and ensure correct function. Make sure not to let aggressive cleaners, solvents, acids, salt water, perspiration or any other harmful substance come in contact with Icelock 544 Series components.

PATIENT CONTROLS

Icelock 544 series components equipped with a release button require the patient to press the release button when doffing the prosthesis in case of an internal seal, like Iceross Seal-In products. Compressing the release button upon donning will also facilitate faster and more complete engagement and is strongly recommended when using Iceross Seal-In liners. Capability to follow this process is essential.

STORAGE

Icelock L-544 Series can be stored as long as it is not corroded and seals are soft and pliable. Cool, dark place is preferred for long term storage.

FABRICATION

Icelock L-544 Series can either be set up on a test socket for temporary test fit purposes, or laminated into a definitive socket.

Test Socket Procedure

1. Using a designated 2-component PU adhesive putty, bond the L-544400 (**Figure 1**) assembly to desired location on test socket (**Figure 2**).
2. Wrap test socket and adapter firmly with glass reinforced casting tape, enabling test fitting of the prosthesis (**Figure 3**).
3. Drill through housing dummy into cast using a sharp 3mm drill bit. Take care to align drill with centerline (**Figure 4**).
4. Remove housing dummy and replace with Housing, there is one seal on each end of Housing (**Figure 5**).
5. Complete assembly applying Loctite 243 to 4-hole screws and torque to 10Nm (**Figure 6**).

Definitive Socket Procedure

1. Drape a layer of 3-4mm PETG or equivalent rigid thermoplastic sheet over cast positive (**Figure 7**).
2. Buff distal area. Using a designated 2-component PU adhesive putty,

bond the L-544400 assembly to desired location directly on thermoplastic layer. Take care to fill the center area for a bubble free bond between plastic dummy and thermoplastic layer **(Figure 9 & 10)**.

3. Tie reinforcement material in groove on socket adapter and do a second tie between socket and adapter. Decide amount of reinforcement according to conventional practices taking into account user weight and activity level **(Figure 11)**.
4. Apply PVA bag and impregnate laminate with conventional methods **(Figure 12)**.
5. Tie around lamination dummy with electrical tape to remove resin distal to socket adapter. This will produce clean edge and minimize finishing work **(Figure 12)**.
6. Once resin has set, drill through housing dummy into cast using a sharp 3mm drill bit. Take care to align drill with centerline **(Figure 13)**. Remove excess resin, PVA bag and plastic dummy.
7. Remove/Break plaster from inside of socket.
8. Insert and hand tighten Housing, there is one seal on each end of Housing. Make sure smaller seal sits flat against bottom of hole when inserting Housing **(Figure 14)**, larger seal sits between Housing and Valve Plate **(Figure 15)**.
9. Complete assembly applying Loctite 243 to 4-hole screws and torque to 10Nm, make sure threads are clean and provide full engagement of screws without resistance. Use screws supplied with kit only.

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

For more information please visit www.ossur.com.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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