










Instructions for Use

ICEROSS SEAL-IN[®] X LINERS

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	21
IT	Istruzioni per l'uso	26
NO	Bruksanvisning	31
DA	Brugsanvisning	36
SV	Bruksanvisning	41
EL	Οδηγίες Χρήσης	46
FI	Käyttöohjeet	52
NL	Gebruiksaanwijzing	57
PT	Instruções de Utilização	62
PL	Instrukcja użytkowania	67
CS	Návod k použití	72
TR	Kullanım Talimatları	77
RU	Инструкция по использованию	82
JA	取扱説明書	87
ZH	中文说明书	92
KO	사용 설명서	96

ENGLISH

Iceross Seal-In X Liners are referred to as the **device** in the following document. This document provides important information on the indications for use, fit, and handling of the device. This document is intended for a certified prosthetist and the user of the device.

DEVICE DESCRIPTION

The device is a silicone liner available in two variants:

- Transtibial, 3 mm / 6 mm profiles;
- Transfemoral, standard / conical.

The device should be paired with a suitable **Iceross Seal-In X Seal**. Together they form a suspension system where a distal vacuum chamber is created below the level of the seal that reliably suspends the residual limb in the prosthesis.

INDICATIONS FOR USE

The device is a non-invasive, single patient, reusable prosthetic device designed for daily use.

- Transtibial: Low to high impact levels.
- Transfemoral: Low to extreme impact levels.

The device is to be used exclusively for the exoprosthetic fitting of amputations of the lower limb, for example due to:

- Trauma
- Vascular Disease
- Cancer
- Congenital Defects

The device can be used by persons from the aforementioned groups regardless of age and level of hand dexterity.

NB: Users should have an appropriate level of physical fitness to use an external prosthesis for load transfer while ambulating.

CONTRAINDICATIONS

Not known.

SAFETY

Warnings

- Misuse of the device may result in potential loss of suspension.
- The device should not be used by individuals with the following conditions:
 - Extremely short residual limb (**refer to tables 1 and 2**).
 - Extreme residual limb volume fluctuations.
 - Cognitive limitations hindering safe device usage.

Residual limb length is measured as follows;

- Transtibial: From mid-patella tendon to distal end.
- Transfemoral: From perineum to distal end, with the tissue hanging down.

Minimum residual limb length is defined as follows;

TABLE 1

TRANSTIBAL	
Liner Size	Min. Length (cm)
16-23,5	11
25-30	12
32-36	13
40	14

TABLE 2

TRANSFEMORAL	
Liner Size	Min. Length (cm)
23,5-28	16
30-34	17
36-40	18
45	19
50	20
55	21

Precautions

- In cases of functional change or loss, the user should immediately cease using the device, and contact their practitioner. This may be the result of device damage or changes in the residual limb.
- The device should be used in combination with a suitable **Iceross Seal-In X Seal**. If the device is used as a cushion liner without a seal, a sleeve must be used to maintain suspension.

SIZING

Precaution: Selecting the correct device size is essential. Selecting a size which is too small/tight may result in pistoning, numbness, or blistering. Selecting a size which is too large/loose, may result in increased perspiration and movement of the limb inside the device, which may lead to blisters and rashes. If any of the above symptoms, or any other indication of improper fit are experienced, the user should contact their practitioner immediately.

To determine the correct device size:

- Measure the circumference of the residual limb 4cm from the distal end.
- Choose the corresponding device size (**Figure 1**). If the measurement taken does not have a corresponding device size, round the measurement down.

NB: Never choose a device size which exceeds the measurement, i.e. do **NOT** round up.

E.g.: If the residual limb measures 24,5cm at the indicated location, choose device size 23,5 **NOT** size 25.

Trimming

The device may be trimmed to increase range of motion, or for user preference (**Figure 2**). The **Iceross Curve Master** is specifically designed to produce a smooth, rounded edge when trimming, minimizing the risk of skin irritation and device tears.

Precaution: Do not cut the device below the proximal trim-lines of the socket. Excessive trimming may reduce the suction between the device and limb, and compromise suspension.

DONNING

Precautions

- Failure to adhere to these instructions may compromise comfort and/or suspension.
- Do NOT apply lotion to the residual limb immediately prior to donning the device. To prevent device damage, the residual limb should be clean and dry.
- Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.
- Ensure the inside of the device is clean, dry, and free from any foreign objects that may cause skin irritation.

1. Orientation:

- Transtibial device should be oriented with the logo on the front of the residual limb (**Figure 3**).
- Transfemoral device requires no specific orientation.

2. Grip the top of the device from the inside, and slide it over the hand until the inner surface is fully exposed, taking care not to damage it with fingernails (**Figure 4**).

3. After exposing as much of the distal end of the device as possible, position it against the distal end of the residual limb (**Figure 5**).

4. With light compression, roll the device all the way upwards onto the limb (**Figures 6 and 7**), taking care not to damage it with fingernails. Do NOT tug or pull.

5. Check that no air pockets are present.

SOCKET CONSTRUCTION GUIDELINES

Recommended socket materials to ensure long-term device function:

- All PETG materials
- Laminated sockets
- Thermoplastic inner socket materials with smooth surface
- EVA inner socket materials

The inner surface of the socket must be airtight and as smooth as possible, irrespective of socket material.

To achieve the smoothest surface possible on the plaster positive model, it is recommended to use fine abrasives for final smoothing, followed by silicone spray.

Use of soft/smooth socket lining material will prolong the device sealing function.

An expulsion valve with release function, such as the Icelock 500 series, is essential for proper device suspension functioning. The socket valve must be positioned as distally as possible, below the **Iceross Seal-In X Seal**.

Socket Donning

Precautions

- Do NOT use alcohol spray, as this will cause the low-friction coating used on **Iceross Seal-In X Seals** to break down, making correct socket donning difficult, and reducing device functionality.
- Do NOT use creams or lotions with the device, as this will cause the **Iceross Seal-In X Seal** to soften and expand over time, affecting its ability to be held in tension against the device.

NB: An *Iceross Seal-In X Seal* ring must be donned before donning the socket.

With the device and seal ring correctly configured, simply step into the socket, expelling all air until fully engaged into the distal end (**see Figures 8 & 9**).

Socket Fit

Achieving the intended function of the *Iceross Seal-In System* is dependent on the fit of the hard socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- User can don the socket easily without the use of lubricant aids.
- Full contact between the device lips and the socket inner wall.
- No air pockets develop when the user loads the prosthesis.
- Good proximal fit.
- No pistoning in the socket.
- No visible movement when pulling on the prosthesis.
- Good rotational stability.
- Full movement possible without compromising distal vacuum.

Precaution: If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

DAILY CARE

NB: Daily cleaning of the residual limb is essential. Use of a mild, pH balanced, 100% fragrance- and dye-free liquid soap is recommended. If skin is dry, apply a pH balanced, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Perspiration

Precaution: If device use results in excessive perspiration such that stability or suspension are compromised, it is recommended to remove the device and dry off both the device and the residual limb.

Cleaning

It is essential to take proper care of the device, and wash it daily following use. The device can be used immediately after cleaning and drying.

- Remove the device, turn inside out (**Figure 4**).
- Wash with pH balanced, 100% fragrance- and dye-free liquid soap.
- The device can also be machine washed (40 °C / warm) with a mild detergent. Fabric softeners, bleaches, and other products/cleaning solutions may result in device damage and should NOT be used.
- After washing, rinse the device thoroughly with hot water (**Figure 10**), and pat dry on both sides with a lint-free cloth.

NB: Always return the device to its neutral state (with the textile cover facing out) immediately after cleaning. Do NOT store inside-out as device damage may result.

Inspection

Precautions

- The device should always be checked for wear and damage. Users are advised to contact their practitioner immediately if any damage is evident, as it may compromise device effectiveness.
- Both the practitioner and user must ensure the device is NOT exposed to glass, carbon fibers, or other foreign particles. Such substances can become embedded in the device material causing skin irritation. Washing alone may not be sufficient to eliminate the problem. The device should be returned to the practitioner if it is inadvertently exposed to foreign substances or chemicals.
- Users should be aware that other common household or bath products may cause or contribute to skin irritation, e.g. soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol-sprays.

STORAGE

Store the device with the textile cover facing out in a cool and dry place. For extended periods of storage, place the device in a clean plastic bag.

WARRANTY

Össur offers a 6-month (from invoice date) limited warranty on the device.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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