

Instructions for Use

SEAL-IN X SEALS


		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	24
NO	Bruksanvisning	29
DA	Brugsanvisning	33
SV	Bruksanvisning	38
FI	Käyttöohjeet	43
EL	Οδηγίες Χρήσης	47
NL	Gebruiksaanwijzing	52
PT	Instruções de Utilização	57
PL	Instrukcja użytkowania	62
CS	Návod k použití	67
TR	Kullanım Talimatları	72
RU	Инструкция по использованию	77
JA	取扱説明書	82
ZH	中文说明书	86
KO	사용 설명서	90


1






2

3

Iceross Seal-In X Seals are referred to as **the device** in the following document. This document provides important information on the indications for use, fit, and handling of the device. This document is intended for a certified prosthetist and the user of the device.

DEVICE DESCRIPTION

The device consists of a hypobaric membrane and textile donning aid. The unique anti-friction coating makes donning and doffing a prosthetic socket easier, without use of lubricant aids. When the device is correctly positioned onto a donned **Iceross Seal-In X Liner**, it will be held in tension against the silicone bands on the exterior surface of the liner, then compressed against the interior surface of the prosthetic socket. This creates a distal vacuum chamber that reliably suspends the residual limb in the prosthesis.

The device is offered in 3 variants:

- **Iceross Seal-In X Grip**
- **Iceross Seal-In X Volume**
- **Iceross Seal-In X Classic**

INDICATIONS FOR USE

The device is a non-invasive, single patient, reusable prosthetic device designed for daily use.

In combination with an **Iceross Seal-In X Liner**:

- Transtibial liner variant supports low to high impact levels.
- Transfemoral liner variant supports low to extreme impact levels.

The device is to be used exclusively for the exoprosthetic fitting of amputations of the lower limb, for example due to:

- Trauma
- Vascular Disease
- Cancer
- Congenital Defects

The device can be used by individuals from the aforementioned groups regardless of age and level of hand dexterity.

NB: Users should have an appropriate level of physical fitness to use an external prosthesis for load transfer while ambulating.

CONTRAINDICATIONS

Not known.

SAFETY

Warning: A strong cautionary statement regarding possible dangers related to the intended use or foreseeable misuse of a device.

Precaution: A cautionary statement requiring action to be taken in advance of device use, or in the event of malfunction or change in device performance. Intended to ensure safe and effective use of device, and to protect the user against greater risk of harm and/or possible failure of the device, when used as intended.

Warning

- The device should not be used by individuals with the following conditions:
 - Extremely short residual limb (**refer to tables 1 and 2**).
 - Extreme residual limb volume fluctuations.
 - Cognitive limitations hindering safe device usage.

Residual limb length is measured as follows;

- Transtibial: From mid-patella tendon to distal end.
- Transfemoral: From perineum to distal end, with the tissue hanging down.

Minimum residual limb length is defined as follows;

TABLE 1

TRANSTIBIAL	
Liner Size	Min. Length (cm)
16-23,5	11
25-30	12
32-36	13
40	14

TABLE 2

TRANSFEMORAL	
Liner Size	Min. Length (cm)
23,5-28	16
30-34	17
36-40	18
45	19
50	20
55	21

Precautions

- In cases of functional change or loss, the user should immediately cease using the device, and contact their practitioner. This may be the result of device damage or changes in the residual limb.
- The device should be used in combination with a suitable **Iceross Seal-In X Liner**.

SIZING

Precaution: Selecting the correct device size is essential. Selecting a size which is too small/tight may result in excessive pressure on the residual limb, which can affect residual limb health. It may also result in over-stretching and permanent deformation of the device, affecting device performance.

One of the following methods may be used to determine correct device size:

Method 1 – Seal Ring Kits

The preferred method of determining correct device size is by using an **Iceross Seal-In X Kit** (please refer to catalog for further details).

1. Place a device from the **Iceross Seal-In X Kit** at appropriate location on top of liner, considering the recommended distance between lowest socket edge and uppermost seal lip (**Figure 1**).
2. Assess for any deformation of the device and residual limb. The device should sit securely in place without excessive compression (**Figure 2**).
3. If the device does not fit at the appropriate location, select another size or relocate if appropriate, to obtain proper tension.

Method 2 - Measurement

Legend for figure 3:

X: Seal size

Y: Circumference (cm)

A: Tight fit

B: Loose fit

Measure the circumference (in centimeters) at the appropriate location, directly on the skin.

Refer to the graph in **Figure 3** to estimate correct device size.

Examples:

- Circumference of 40 cm corresponds to device sizes 44 or 47.
- Circumference of 26 cm corresponds to device sizes 28 or 30.

DONNING

1. The device must be clean, and free from dirt prior to use. If necessary, rinse with water, and dry with a lint-free cloth.
2. Use fabric end of device to pull on, and adjust to the desired location (**Figures 4 & 5**).

Precautions:

- Silicone portion of device must not extend proximal to the socket brim (**Figure 6**).
- Device must be aligned with the silicone bands on the **Iceross Seal-In X Liner** to achieve a proper seal (**Figure 6**).
- Device must not be worn directly on an exposed residual limb.
- Do NOT use creams, lotions, or other lubricants with the device, as this will cause it to soften and expand over time, affecting its ability to be held in tension on the **Iceross Seal-In X Liner**.

NB: The uppermost seal lip of the device should be no closer to the socket brim than indicated in **Figure 1**.

SOCKET CONSTRUCTION GUIDELINES

Recommended socket materials to ensure long-term device function:

- All PETG materials
- Laminated sockets
- Thermoplastic inner socket materials with smooth surface
- EVA inner socket materials

The inner surface of the socket must be airtight and as smooth as possible, irrespective of socket material.

To achieve as the smoothest surface possible on the plaster positive model, it is recommended to use fine abrasives for final smoothing, followed by silicone spray.

Use of soft/smooth socket lining material will prolong device sealing function.

An expulsion valve with release function, such as in the **Icelock 500** series, is essential for proper device suspension functioning. The socket valve must be positioned as distally as possible, below the device.

Socket Donning

Precautions

- Do NOT use alcohol spray, as this will cause the low-friction coating of the device to break down, making correct socket donning difficult, and reducing device functionality.

- Do NOT use creams or lotions with the device, as this will cause it to soften and expand over time, affecting its ability to be held in tension against the **Iceross Seal-In X Liner**.

With the **Iceross Seal-In X Liner** and device correctly configured, simply step into the socket, expelling all air until fully engaged into the distal end (**Figure 7**).

Socket Fit

Achieving the intended function of the **Iceross Seal-In System** is dependent on the fit of the hard socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- User can don the socket easily without the use of lubricant aids.
- Full contact between the device lips and the socket inner wall.
- No air pockets develop when the user loads the prosthesis.
- Good proximal fit.
- No pistoning in the socket.
- No visible movement when pulling on the prosthesis.
- Good rotational stability.
- Full movement possible without compromising distal vacuum.

Precaution: If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

Seal-In Socks

Users who experience residual limb volume fluctuations and resulting compromised suspension, are advised to use **Iceross Seal-In Socks**. These should be placed below the seal lip (**Figures 1A & 8**).

Precautions

- If volume fluctuation cannot be adequately compensated for by this means, the socket will need to be remade to suitably reduced dimensions.
- The **Iceross Seal-In Socks** should never be placed between the **Iceross Seal-In X Liner** and the device (**Figure 8**).

DAILY CARE

NB: The residual limb should be inspected daily by the user, to verify that skin/limb health is satisfactory, and continued use of the device is safe.

Cleaning

Rinse the silicone portion of the device with water daily, and dry with a lint-free cloth. This will prevent the build-up of dirt and other residue, which may affect the seal-in properties and hygiene of the device over time.

NB: Do NOT use soap or any other cleaning aids with the device.

Inspection

Precaution: The device should ALWAYS be checked for holes, other damage, or wear. Users are advised to contact their practitioner immediately if damage is evident, as it may compromise device effectiveness.

WARRANTY

Össur offers a 6-month (from invoice date) limited warranty on the device.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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