

Never apply a B-Cool® 2.0 gel pack(s) directly onto the skin.

An additional protective barrier such as a B-Cool® Wrap or cloth towel is recommended to protect against frostbite to the affected area.

In the unlikely event there is exposure to the eyes, immediately flush with water and consult with a healthcare professional as a precaution.

If the B-COOL® 2.0 gel pack is punctured, dispose of it and re-order a replacement B-COOL® 2.0 gel pack to continue treatment.

COLD THERAPY APPLICATION INSTRUCTIONS

- 1. Place the B-COOL® 2.0 gel pack(s) into the freezer for at least 4 hours before using.**
- 2. Place in the B-COOL® 2.0 Wraps or cloth towel, then apply to affected area.**
- 3. Remove pack if it becomes uncomfortable. Always follow your physician's or healthcare professional's specific treatment instructions.**

