

B-COOL™ Hip Wrap



1. Remove wrap and insert gel packs into pockets.

Note: Place gel packs into freezer until frozen

2. Place wrap over affected hip.
3. Attach with tab and secure shorter strap around the thigh to hold in place.
4. Attach with tab and secure longer strap around the waist.
5. Adjust straps accordingly for comfort.
6. Never apply a frozen or heated pack directly onto the skin.
7. Always follow your physician's or healthcare provider's specific cold therapy treatment.

| Reorder No. | Dimensions | Description |
|-------------|--------------------------------------|------------------------------------|
| 0814 3333 | 11 3/4" x 12 3/4" (29.85 x 32.39 cm) | Wrap Only |
| 0814 3331 | 11 3/4" x 12 3/4" (29.85 x 32.39 cm) | Wrap with B-Cool 6" x 10" Gel 2/PK |
| 0814 3332 | 11 3/4" x 12 3/4" (29.85 x 32.39 cm) | Wrap with B-Cool 2.0 2/PK |

BIRD & CRONIN INC.

Designers and Manufacturers of Quality Medical Products

1200 Trapp Road • Eagan, Minnesota USA 55121

800/328-1095 • Fax 800/279-7934 • 651/683-1111 • Fax 651/688-9855

email sales@birdcronin.com • www.birdcronin.com

© 2011 Bird & Cronin, Inc. All Rights Reserved 10003333cm