

BAKER NIGHT SPLINT

Application Instructions

Reorder No.	Size
0814 4802	Small
0814 4803	Medium
0814 4804	Large
0814 4805	X Large
0814 4806	XX-Large

1. Flex knee when applying.
2. Attach the instep strap first to hold heel into splint.
3. Extend knee to check tightness of instep strap.
4. Next adjust remaining straps and close buckles.
5. For additional stretching place toe wedge under liner.
6. Prone position with pillows under the shins is acceptable.

CAUTIONS

Over stretching, especially during extended wear (several hours) may produce temporary numbness and tingling. This is usually a rare situation and can be avoided by maintaining the stretch arc. If symptoms such as above occur, loosen the stretch straps or remove the splint. Symptoms should resolve within several minutes. If symptoms persist or condition worsens, discontinue use of the splint and contact your physician.

WARNING

***DO NOT WALK ON PRODUCT.
REMOVE BEFORE WALKING.***

CARE

Liner, Pads, Straps: Hand wash with mild detergent. Air Dry.

Shell: Clean with damp cloth and mild detergent.
Wipe Dry.

BIRD & CRONIN INC.

Designers and Manufacturers of Quality Medical Products

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail sales@birdcronin.com • www.birdcronin.com